---

**AB-TPOK-3**  
**(P002046)**

---

DiaSorin

---

---

English.....	p. 1
Italiano.....	p. 9
Français.....	p. 17
Deutsch.....	p. 26
Español.....	p. 35
Português.....	p. 44
Magyar.....	p. 52
Czesky.....	p. 60
Ελληνικά.....	p. 68

**TPO ANTIBODIES IMMUNORADIOMETRIC ASSAY KIT**  
**Procedure for quantitative determination of**  
**autoantibodies to human thyroid peroxidase (IgG to TPO)**  
**in human serum or plasma samples**

For in vitro use only

**1. INTRODUCTION**

Patients affected by autoimmune thyroid diseases (Hashimoto's disease, primary hypothyroidism, idiopathic myxoedema and Graves-Basedow disease) circulate autoantibodies that react with a variety of human thyroid antigens, such as thyroglobulin (hTg), thyroid microsomal antigen (M), second antigen of the colloid and TSH receptor.

Assay of thyroglobulin (anti-hTg) and microsomal (anti-M) autoantibodies is presently the method most commonly used in clinical practice for aetiopathogenic diagnosis of autoimmune thyroid diseases. Such antibodies, especially microsomal antibodies, are present at elevated levels in the majority of patients affected by Graves-Basedow disease and in almost all patients affected by Hashimoto's disease or idiopathic myxoedema. Conversely, thyroglobulin and microsomal antibodies are never found elevated in other thyroid diseases, such as sporadic or endemic non-toxic goitre and malignant or benign neoplasms. Thyroglobulin and microsomal antibodies are generally absent in normal subjects; however, microsomal antibodies may be observed in females over 50 with no clinical signs of thyroid disease.

Microsomal antibodies were originally detected by immunofluorescence and complement fixation, but are currently titrated by passive haemagglutination or radio- and enzyme immunoassays. The reliability of these assays is proportional to the antigen purification. Whilst hTg is well identified and purified since long, microsomal antigen structure is unknown and it has not yet been purified. Therefore, microsomal antigen preparations currently used for autoantibodies detection are forcefully contaminated by other antigens, such as hTg. This fact has made it necessary to add excess hTg to the reaction mixture when assaying microsomal antibodies in order to avoid false positive results caused by the interference of thyroglobulin antibodies.

Recently, thyroid peroxidase (TPO) has been identified as the main and possibly single autoantigenic component of microsomal antigen at biochemical, immunologic and molecular level. TPO plays a key role in thyroid hormone synthesis, by catalyzing iodide oxidation and its incorporation into hTg thyrosyl residues.

The availability of anti-TPO monoclonal antibodies has made it possible to purify TPO by affinity chromatography and therefore to develop immunoassays specific for TPO autoantibodies. These methods have the obvious advantages of using a purified antigen and of requiring no pre-adsorption of the specimen to be assayed with hTg.

Recent studies have demonstrated that immunoassays of TPO are not only clinically reliable, but also more specific than microsomal antibody assays employing a non-purified antigen.

**2. PRINCIPLE OF THE ASSAY**

The assay is an immunoradiometric (IRMA) method. In the first assay stage, IgG to TPO contained in calibrators or samples binds to solid-phase TPO. In the second assay stage, the <sup>125</sup>I-tracer (<sup>125</sup>I-labelled protein A) is added, which binds to the solid-phase TPO-IgG complex. After a second incubation, the amount of <sup>125</sup>I-tracer bound to the solid phase is proportional to the concentration of IgG to TPO present in calibrators or samples to be assayed. At the end of each incubation, the unbound material is removed by aspiration and washing. The method adopted for B/F separation is based on the use of coated tubes, where recombinant TPO is fixed on the tube walls.

The assay of IgG to TPO employs protein A as tracer, which is a bacterial cell wall protein derived from *Staphylococcus aureus*, capable of binding the Fc fragment of IgG molecules.

### 3. REAGENTS PROVIDED IN THE KIT

Coated tubes	100
<sup>125</sup> I-tracer	1 vial
Anti-TPO calibrators	5 vials
Control serum	1 vial
Sample diluent / Zero calibrator	1 bottle
Incubation buffer	1 bottle
Wash buffer	2 bottles
Number of tubes	100

**STORAGE:** Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming.

Reagents from different batches must not be mixed.

#### 3.1. Coated tubes

The inner surface of each tube is coated with recombinant thyroid peroxidase (baculovirus).

Before use, bring the coated tubes to room temperature prior to opening the box, to avoid condensation of humidity.

Securely reseal the box containing unused tubes. Do not mix different batches of coated tubes.

#### 3.2. <sup>125</sup>I-tracer (red): ready-to-use reagent

The vial contains 21 mL protein A labelled with <sup>125</sup>I, BSA, phosphate-citrate buffer, preservatives and an inert red dye. Radioactivity is 444 kBq (12 µCi) or less on the calibration date.

#### 3.3. Anti-TPO calibrators: ready-to-use reagent

Each vial contains 0.5 mL of prediluted human serum containing anti-TPO antibodies, calf serum, BSA, PBS buffer and preservatives. The calibrator concentrations are the following: 10 - 30 - 100 - 300 - 1000 U/mL. The kit calibrators are referenced to the NIBSC 66/387 anti-TMS Reference Preparation. *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

*Because the calibrators have been prediluted 1:51 and values assigned accordingly, the U/mL concentrations of the unknown specimens can be determined directly by interpolation from the calibration curve. A multiplying factor is necessary only when the specimen dilution is greater than 1:51.*

#### 3.4. Control serum: lyophilized reagent

The vial contains prediluted human serum containing anti-TPO antibodies, calf serum, BSA, PBS buffer and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Reconstitute the vial contents with 1 mL distilled water. The resulting solution is stable at 2-8°C until the kit expiry date. The control serum *should not be further diluted*. Read the value directly from the calibration curve.

#### 3.5. Sample diluent / Zero calibrator: ready-to-use reagent

The bottle contains 50 mL PBS buffer, calf serum, BSA and preservatives. The reagent is used as zero calibrator and to dilute samples.

*The reagent is common to AB-TPOK-3 and AB-HTGK-3 kits.*

### 3.6. Incubation buffer (blue): ready-to-use reagent

The bottle contains 52 mL TRIS buffer solution, BSA, detergents, preservatives and an inert blue dye.

### 3.7. Wash buffer: reagent in solution (10x)

Each bottle contains 50 mL 0.5% Triton X-100 and saline solution.

Dilute the contents of each bottle to 500 mL with deionized water. The resulting solution is stable at 2-8°C until the kit expiry date. The reagent is used to rinse coated tubes.

## 4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Glassware.
- Disposable polystyrene tubes.
- Micropipettes with disposable tips (10, 50, 200, 500, 1000 µL) (10, 50 µL: trueness ± 3%, precision 2%; 200, 500, 1000 µL: trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Horizontal shaker capable of achieving a shaking speed of 300-350 rpm.
- Device for dispensing and aspiration of wash buffer capable of delivering 2-3 mL per wash cycle for two wash cycles.
- Gamma counter suitable for counting <sup>125</sup>I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

## 5. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulants citrate, EDTA and heparin have been tested and may be used with this assay. Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

### Sample dilution

All samples should be diluted 1:51 with sample diluent (3.5) before assaying. Dispense 10 µL sample and 500 µL sample diluent into polystyrene tubes and mix with a Vortex.

If high levels of IgG to TPO are expected, an additional dilution with the sample diluent supplied in the kit should be performed. *Calibrators and control serum are ready to use and must not be diluted.*

## 6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each calibrator and sample.

Dispense reagents *in the bottom of coated tubes*. Operate according to the following scheme:

reagents	tubes	Calibrators 0-5	Diluted samples
Calibrators		50 µL	–
Diluted samples		–	50 µL
Incubation buffer		200 µL	200 µL

- **Mix** the contents of tubes and **incubate for 2 hours at room temperature** while continuously shaking (300-350 rpm).
- Carefully **aspirate** the incubation mixture and **wash** twice with 2 mL wash buffer.

Be sure that the aspirator tip touches the bottom of the coated tube so that all the liquid is removed. Failure to remove adhering solution adequately may result in poor reproducibility and spurious results. No trace of dye should still be visible.

- **Dispense 200 µL tracer** into all tubes. Prepare two non-coated tubes for total activity computation containing only 200 µL tracer and set them aside until counting.
- **Incubate for 2 hours at room temperature** while continuously shaking (300-350 rpm).
- Carefully **aspirate** tracer and **wash** twice with 2 mL wash buffer. Operate as above.
- **Measure the radioactivity** of tubes.

## 7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Compute the B/T ratio for each calibrator and unknown sample as follows:

$$B/T\% = \frac{\text{calibrator or sample mean counts}}{\text{total activity mean counts}} \times 100$$


Plot in log-log coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of IgG to TPO concentration expressed as U/mL on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1). Directly from the calibration curve, read the IgG to TPO concentration of each sample expressed as U/mL. If the sample was diluted beyond the initial 1:51 factor, the antibody concentration value derived from the diluted specimen must be multiplied by the additional dilution factor. Keep in mind that the calibrators are already diluted 1:51.

### Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/T x 100
Total activity	113,962	–
Zero calibrator/sample diluent	43	0.04
10 U/mL	954	0.84
30 U/mL	2,495	2.19
100 U/mL	8,308	7.29
300 U/mL	21,395	18.77
1000 U/mL	65,080	57.11
Sample	6,914	6.07

By interpolation from the calibration curve, the sample is found to contain 90.3 U/mL anti- TPO IgG.

**Fig. 1**

### Interpretation of results

The first calibrator (10 U/mL) represents the system's cut-off.

Samples with IgG to TPO values greater than or equal to the cut-off value should be graded positive. Samples with IgG to TPO values less than the cut-off value should be graded negative.

*However, samples with IgG to TPO ranging within  $\pm 20\%$  of the cut-off value should be graded disputable. Each laboratory should establish its own reference ranges.*

### 8. EXPECTED VALUES

During the clinical investigation 359 subjects were taken into consideration: 200 euthyroid subjects, 98 hypothyroid patients affected by Hashimoto's autoimmune thyroiditis, and 61 hyperthyroid patients affected by Graves' disease. In this population 200 subjects were graded negative and 159 subjects were graded positive for anti-TPO IgG with a reference test. The negative population encompassed subjects from all the three groups.

With DiaSorin kit 195 subjects were graded negative (anti-TPO IgG levels below 10 U/mL) and 157 subjects were graded positive (anti-TPO IgG levels above 10 U/mL).

Concordance of results was 98% (352/359).

However, the values reported above are merely indicative: each laboratory should establish its own reference ranges.

### 9. SPECIFIC PERFORMANCE CHARACTERISTICS

#### 9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

**Interference.** Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by anticoagulants (citrate, EDTA, heparin), slight haemolysis (up to 20 mg/dL haemoglobin), lipaemia (up to 10000 mg/dL triglycerides), or bilirubinaemia (up to 50 mg/dL bilirubin).

**Cross-reactions.** No cross-reactivity is observed due to the presence of thyroglobulin autoantibodies up to a concentration of one million IU/mL.

## 9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 0.8 U/mL at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator (sample diluent), that is, two standard deviations above zero.

## 9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	25	25	25
Mean (IU/mL)	41.04	121.13	244.00
Standard deviation	1.70	7.40	6.00
Coefficient of variation (%)	4.10	6.10	2.40

Reproducibility	D	E	F	G
Number of determinations	25	25	25	25
Mean (U/mL)	21.4	71.5	146.5	680.5
Standard deviation	1.9	5.6	7.0	43.9
Coefficient of variation (%)	9.1	7.8	4.8	6.5

## 9.4. Trueness

The assay trueness has been checked by the dilution test.

**Dilution test.** Six sera with high anti-TPO antibodies concentration were tested after serially diluting with sample diluent. Examples are given in the following table.

Dilution	Expected concentration, U/mL	Measured concentration, U/mL	% Recovery
neat	–	437.6	–
1:2	218.8	196.8	89.9
1:4	109.4	109.4	100.0
1:8	54.7	54.9	100.4
1:16	27.3	25.1	91.9
1:32	13.7	13.6	99.3
neat	–	207.9	–
1:2	104.0	108.4	104.3
1:4	52.0	55.5	106.8
1:8	26.0	22.9	88.1
1:16	13.0	12.2	93.9
1:32	6.5	5.7	87.7

## 9.5. High-dose saturation effect

Whenever samples containing extremely high antibody concentrations are tested in a two-step sandwich method, the saturation effect can mimic concentrations lower than real. However, a well-optimized method excludes grossly underestimated results (as it may occur in one-step sandwich methods).

Concentrations of IgG to TPO up to 150,000 U/mL were tested in this assay and it was observed that such levels produce analytical signals that are consistently above the top calibrator concentration (saturation curve). For correct quantification, samples containing IgG to TPO levels greater than that of the top calibrator should be diluted with sample diluent and retested. The results must then be multiplied by the additional dilution factor beyond the initial 1:51 to obtain the IgG to TPO levels of the neat specimens.

## 10. LIMITATIONS OF THE PROCEDURE

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration and washing are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration of incubation mixture and rinsing of tubes
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

## 11. WARNINGS AND PRECAUTIONS

Test components contain sodium azide as a preservative. Because sodium azide may form explosive lead or copper azide in plumbing, it is recommended that drains be thoroughly flushed with water after disposal of solutions containing sodium azide (Council Directive 99/45/EC).

R 22 – Harmful if swallowed.

R 31 – Contact with acids liberates toxic gas.

S 28 – After contact with skin, wash immediately with plenty of water.

S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

## 12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

## 13. BASIC RULES OF RADIATION SAFETY

### Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 12.5 µCi (462 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the

U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

### SCHEME OF THE ASSAY

- 1 - PREPARE REAGENTS. DILUTE SAMPLES 1:51.
- 2 - IDENTIFY COATED TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS ACCORDING TO THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

REAGENTS \ TUBES	T	CAL 0-5	DILUTED SAMPLES
CALIBRATORS	—	50 $\mu$ L	—
DILUTED SAMPLES	—	—	50 $\mu$ L
INCUBATION BUFFER	—	200 $\mu$ L	200 $\mu$ L

- 4 - INCUBATE FOR TWO HOURS AT ROOM TEMPERATURE WHILE SHAKING.
- 5 - ASPIRATE THE INCUBATION MIXTURE AND RINSE TWICE WITH 2 mL WASH BUFFER.
- 6 - DISPENSE 200  $\mu$ L TRACER INTO ALL TUBES.
- 7 - INCUBATE FOR TWO HOURS AT ROOM TEMPERATURE WHILE SHAKING.
- 8 - ASPIRATE TRACER AND RINSE TWICE WITH 2 mL WASH BUFFER.
- 9 - MEASURE THE RADIOACTIVITY OF TUBES.

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B. CZARNOCKA et al.

Purification of the human thyroid peroxidase and its identification as the microsomal antigen involved in autoimmune thyroid disease.  
*FEBS*, **190** (1) : 147 (1985).

H. ENGLER et al.

Assessment of anti-thyroglobulin and anti-microsomal autoantibodies in patients with autoimmune thyroid disease: comparison of hemagglutination assay, enzyme-linked immunoassay and radioligand assay.  
*Chim. Clin. Acta*, **179** : 251 (1989).

M. GREINER et al.

A modified ROC analysis for the selection of cut-off values and the definition of intermediate results of serodiagnostic tests.  
*J. Immunol. Meth.*, **185** : 123 (1995).

J. GUSTAFSSON et al.

Purification of the human thyroid peroxidase using ion-exchange liquid chromatography.  
*J. Clin. Lab. Immunol.*, **23** : 57 (1987).

Y. KOHNO et al.

Autoantibodies to thyroid peroxidase in patients with chronic thyroiditis: effect of antibody binding on enzyme activities.  
*J. Exp. Immunol.*, **65** : 534 (1986).

T. KOTAMI

Detection of autoantibodies to thyroid peroxidase in autoimmune thyroid disease by micro-ELISA and immunoblotting.  
*J. Clin. Endocrinol. Metab.*, **62** (5) : 928 (1986).

F. LIBERT et al.

Thyropoxidase, an auto-antigen with a mosaic structure made of nuclear and mitochondrial gene modules.  
*EMBO J.*, **6** (13) : 4193 (1987).

M. LUDGATE, G. VASSART

Review – The molecular genetics of three thyroid autoantigens: thyroglobulin, thyroid peroxidase and the thyrotropin receptor.  
*Autoimmunity*, **7** : 201-217 (1990).

S. MARIOTTI et al.

Comparison of serum thyroid microsomal and thyroid peroxidase autoantibodies in thyroid diseases.  
*J. Clin. Endocrinol. Metab.*, **65** (5) : 987 (1987).

S. MARIOTTI et al.

The thyroid microsomal-microvillar autoantigen.  
*Folia Allergol. Immunol. Clin.*, **34** (1) : 1 (1987).

A. PINCHERA et al.

Significance of thyroid autoantibodies in autoimmune thyroid diseases.  
In: *Autoimmunity and the Thyroid*, P.G. Wolfish, R.J. Wall, R. Volpi eds., Academic Press, Orlando, p. 139 (1985).

L. PORTMANN et al.

Anti-thyroid peroxidase antibody in patients with autoimmune thyroid disease: possible identity with anti-microsomal antibody.  
*J. Clin. Endocrinol. Metab.*, **61** (5) : 1001 (1985).

S.H. ROMAN et al.

Enzyme-linked immunosorbent microassay and hemagglutination compared for detection of thyroglobulin and thyroid microsomal autoantibodies.  
*Clin. Chem.*, **30** (2) : 246 (1984).

S. SHUMAN

Thyroid antigens and autoimmunity.  
*Adv. Immunol.*, **14** : 85 (1971).

R. VOLPI


Autoimmune thyroid disease. A perspective.  
*Mol. Biol. Med.*, **3** : 25 (1986).

A.P. WEETMAN, A.M. MCGREGOR

Autoimmune thyroid disease: development in our understanding.  
*Endocrinol. Rev.*, **5** : 309 (1984)

DiaSorin Inc.  
1951 Northwestern Avenue  
P.O. Box 285  
Stillwater, MN 55082-0285 U.S.A.  
Phone: 651-439-9710  
FAX: 651-351-5669

DiaSorin S.p.A.,  
13040 Saluggia (VC) Italy

APM0130

34695 7/10