INSIK-5
(P2796)

DiaSorin

English	p. 1
Italiano.....	p. 9
Français	p. 17
Español	p. 26

INSULIN RADIOIMMUNOASSAY KIT

Procedure for quantitative determination of insulin in human serum or plasma samples

For in vitro diagnostic use only

1. INTRODUCTION

Insulin is a polypeptide hormone of molecular weight 6000 daltons, composed of two peptide chains, A and B, joined by two cross-linked disulphide bonds and synthesized by the beta cells of the islets of Langerhans of the pancreas.

Insulin influences most of the metabolic functions of the body. Its best known action is to lower the blood glucose concentration by increasing the rate at which glucose is converted to glycogen in the liver and muscles and to fat in adipose tissue, by stimulating the rate of glucose metabolism and by depressing gluconeogenesis.

Insulin stimulates the synthesis of proteins, DNA and RNA in cells generally, and promotes the uptake of amino acids and their incorporation into muscle protein. It increases the uptake of glucose in adipose tissue and its conversion into fat and inhibits lipolysis.

Insulin primary action is on the cell membrane, where it probably facilitates the transport of glucose and amino acids into the cells. At the same time it may activate intracellular enzymes such as glycogen synthetase, concerned with glycogen synthesis.

2. PRINCIPLE OF THE ASSAY

The principle of the assay is based on the competition between labelled insulin and insulin contained in calibrators or samples to be assayed for a fixed and limited number of antibody binding sites. After the incubation, the amount of labelled insulin bound to the antibody is inversely related to the concentration of unlabelled insulin present in calibrators or samples. The method adopted for B/F separation is based on the use of a precipitating reagent, in which a second antibody is pre-precipitated and in excess.

3. REAGENTS PROVIDED IN THE KIT

¹²⁵ I-labelled insulin	1 vial
Insulin antiserum	1 vial
Insulin calibrators	6 vials
Control serum	1 vial
Precipitating reagent	2 bottles
Number of tubes	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming.

Reagents from different batches must not be mixed.

3.1. ¹²⁵I-labelled insulin (red): lyophilized reagent

The vial contains porcine insulin labelled with ¹²⁵I, BSA, phosphate buffer, preservatives and an inert red dye. Radioactivity is 51 kBq (1.37 µCi) or less on the calibration date.

Reconstitute the vial contents with 10 mL distilled water. The resulting solution is stable for four weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.2. Insulin antiserum (blue): lyophilized reagent

The vial contains antiserum raised in guinea-pigs, BSA, phosphate buffer, preservatives and an inert blue dye. Reconstitute the vial contents with 10 mL distilled water. The resulting solution is stable for four weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.3. Insulin calibrators: lyophilized reagent

The vials contain increasing amounts of human insulin, human serum and preservatives. The calibrators are referenced to WHO 66/304 international standard. *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

Reconstitute the contents of the zero calibrator vial with 2 mL distilled water and the contents of the 1-5 calibrator vials with 1 mL distilled water. The resulting solutions contain 0 - 10 - 25 - 50 - 100 - 200 µIU/mL insulin respectively and are stable for four weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.4. Control serum: lyophilized reagent

The vial contains human insulin, human serum and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Reconstitute the vial contents with 1 mL distilled water. The resulting solution is stable for four weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.5. Precipitating reagent: ready-to-use reagent

Each bottle contains 53 mL polyethylene glycol, TRIS buffer, antibody to guinea-pig IgG raised in goats, non-specific guinea-pig IgG and preservatives.

4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Disposable polystyrene tubes.
- Micropipettes with disposable tips (100, 1000 µL) (trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Multisample centrifuge capable of achieving 1500-2000 x g*.
- Water aspirator (optional).
- Gamma counter suitable for counting ¹²⁵I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

5. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulant EDTA has been tested and may be used with this assay. Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

If insulin levels greater than 200 µIU/mL are expected, the samples should be diluted with the zero calibrator.

6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time.

Perform all assay steps in the order given and without any appreciable delays between the steps.

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

A clean, disposable tip should be used for dispensing each calibrator and sample.

- Operate according to the following scheme:

reagents \ tubes	Total activity	Calibrators 0-5	Samples
Calibrators	-	100 μ L	-
Samples	-	-	100 μ L
Tracer	100 μ L	100 μ L	100 μ L
Antiserum	-	100 μ L	100 μ L

- **Mix** the contents of tubes with a Vortex and **incubate for 1.5 hours at room temperature** (procedure A) or **overnight** (18-24 hours) **at 2-8°C** (procedure B).
- Let the bottle of **precipitating reagent reach room temperature** and **mix** well by repeatedly tilting end over end.
- **Dispense 1 mL precipitating reagent** into all tubes (except for total activity).
- **Mix** the contents of tubes with a Vortex and **let stand for 15 min at room temperature**.
- **Centrifuge** all tubes for 15 min at 1500-2000 x g*.
- **Discard the supernatant** by aspiration or decantation. When aspirating, avoid disturbing the precipitate. When decanting, let the tubes drain mouth down on blotting paper.
- **Measure the radioactivity** of precipitate.

7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Evaluate the binding ability as follows:

$$(B/T)_0\% = \frac{\text{zero calibrator mean counts}}{\text{total activity mean counts}} \times 100$$

Express the mean counts for each calibrator and unknown sample as a percentage of zero calibrator mean counts:

$$B/B_0\% = \frac{\text{calibrator or sample mean counts}}{\text{zero calibrator mean counts}} \times 100$$

Plot in linear-linear or semilog coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of insulin concentration expressed as μ U/mL on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1). Directly from the calibration curve, read the insulin concentration of each sample expressed as μ U/mL. If the sample was diluted, the insulin concentration value derived from the diluted specimen must be multiplied by the dilution factor.

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/Bo x 100
Total activity	18,622	–
Zero calibrator	9,738	100
10 $\mu\text{IU/mL}$	8,287	85.1
25 $\mu\text{IU/mL}$	6,213	63.8
50 $\mu\text{IU/mL}$	4,460	45.8
100 $\mu\text{IU/mL}$	3,009	30.9
200 $\mu\text{IU/mL}$	2,220	22.8
Sample	4,966	51.0

By interpolation from the calibration curve, the sample is found to contain 35 $\mu\text{IU/mL}$ insulin.


Fig. 1

8. EXPECTED VALUES

Insulin levels normally range from 5 to 25 $\mu\text{IU/mL}$. However, the range given is merely indicative: each laboratory should establish its own reference ranges.

9. SPECIFIC PERFORMANCE CHARACTERISTICS

9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by EDTA, haemolysis (up to 1000 mg/dL haemoglobin), lipaemia (up to 500 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin) or one freeze-thaw cycle of samples. Insulin values determined in plasma samples obtained with heparin or citrate are significantly lower than those determined in serum samples.

Cross-reactions. The percentage of cross-reactions, calculated according to Abraham, shows the specificity of the antiserum used.

- Human insulin	100%
- Bovine insulin	100%
- Porcine insulin	100%
- Rat insulin	100%
- Porcine proinsulin	28.0%
- Bovine glucagone	0.46%
- Porcine glucagone	0.03%
- Porcine C-peptide	< 0.01%
- Human C-peptide	< 0.01%

9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 3.0 $\mu\text{IU}/\text{mL}$ at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator, that is, two standard deviations below zero.

9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	10	10	10
Mean ($\mu\text{IU}/\text{mL}$)	24.1	73.6	130.8
Standard deviation	1.6	7.8	7.2
Coefficient of variation (%)	6.6	10.6	5.5

Reproducibility	A	B	C
Number of determinations	20	20	20
Mean ($\mu\text{IU}/\text{mL}$)	22.6	79.6	135.5
Standard deviation	1.4	8.6	13.1
Coefficient of variation (%)	6.2	10.8	9.7

9.4. Trueness

The assay trueness has been checked by the dilution and recovery tests.

Dilution test. One serum with high insulin concentration was tested after serially diluting with the zero calibrator.

Dilution	Expected concentration, $\mu\text{U/mL}$	Measured concentration, $\mu\text{U/mL}$	% Recovery
neat	–	93.0	–
1:2	46.5	45.1	97.0
1:4	23.3	21.8	93.9
1:8	11.6	11.2	96.0
1:16	5.8	5.4	92.2

Recovery test. Two sera containing insulin were tested as such and after mixing with increasing amounts of insulin.

Added concentration, $\mu\text{U/mL}$	Expected concentration, $\mu\text{U/mL}$	Measured concentration, $\mu\text{U/mL}$	% Recovery
–	–	32.2	–
114	145.0	160.0	110.5
57	88.0	87.0	98.9
27	58.0	61.0	105.2
15	46.0	48.0	105.7
–	–	57.8	–
114	169.0	168.0	99.5
57	112.0	105.0	94.0
27	82.0	80.0	97.7
15	70.0	65.0	93.4

10. LIMITATIONS OF THE PROCEDURE

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration or decanting are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration or decanting of incubation mixture
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

11. WARNINGS AND PRECAUTIONS

Test components contain sodium azide as a preservative. Because sodium azide may form explosive lead or copper azide in plumbing, it is recommended that drains be thoroughly flushed with water after disposal of solutions containing sodium azide (Council Directive 99/45/EC).

- R 22 – Harmful if swallowed.
- R 31 – Contact with acids liberates toxic gas.
- S 28 – After contact with skin, wash immediately with plenty of water.
- S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Sodium merthiolate (Council Directive 99/45/EC):

- R 20/21/22 – Harmful by inhalation, in contact with skin and if swallowed.
- R 33 – Danger of cumulative effects.
- S 28 – After contact with skin, wash immediately with plenty of water.
- S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

13. BASIC RULES OF RADIATION SAFETY

Reagents Containing Iodine-125

This kit contains radioactive material that does not exceed 1.5 µCi (54 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

SCHEME OF THE ASSAY

- 1 - RECONSTITUTE REAGENTS.
- 2 - IDENTIFY ASSAY TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS BY THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

REAGENTS \ TUBES	T	CAL 0-5	SAMPLES
CALIBRATORS	-	100 μ L	-
SAMPLES	-	-	100 μ L
TRACER	100 μ L	100 μ L	100 μ L
ANTISERUM	-	100 μ L	100 μ L

- 4 - INCUBATE FOR 1.5 HOURS AT ROOM TEMPERATURE (PROCEDURE A) OR OVERNIGHT AT 2-8°C (PROCEDURE B).
- 5 - ADD 1 mL PRECIPITATING REAGENT TO ALL TUBES (except for total activity) AND THOROUGHLY MIX.
- 6 - LET STAND FOR 15 MIN AT ROOM TEMPERATURE.
- 7 - CENTRIFUGE FOR 15 MIN AT 1500-2000 x g*.
- 8 - DISCARD SUPERNATANT.
- 9 - MEASURE THE RADIOACTIVITY OF PRECIPITATE

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- S. BABA et al.
Metabolism of insulin and proinsulin.
Excerpta Med. I.C.S., **468** : 270 (1979).
- R.A. DE FRONZO, J.D. TOBIN, R. ANDRES
Glucose tramp technique: a method for quantifying insulin secretion and resistance.
Am. J. Physiol., **237** : E214-E223 (1979).
- E. FERRANNINI, A. PILO
Pattern of insulin delivery after intravenous glucose injection in man and its relation to plasma glucose disappearance.
J. Clin. Invest., **64** : 243 (1979).
- E. FERRANNINI et al.
Independent stimulation of glucose metabolism and Na⁺-K⁺ exchange by insulin in the human forearm.
Am. J. Physiol., **255** : E953-E958 (1988).
- A.M. FOURNIER et al.
Blood pressure, insulin, and glycemia in nondiabetic subjects.
Am. J. Med., **80** : 861-864 (1986).
- R.O.B. GANS et al.
Renal and cardiovascular effects of exogenous insulin in healthy volunteers.
Clin. Sci., **80** : 219-225 (1991).
- R.O.B. GANS et al.
The effect of angiotensin I-converting enzyme inhibition on insulin action in healthy volunteers.
Eur. J. Clin. Invest., **21** : 527-533 (1991).
- L.C. GIUDICE et al.
Insulin-like growth factors and their binding proteins in the term and preterm human fetus and neonate with normal and extremes of intrauterine growth.
J. Clin. Endocrinol. Metab., **80** : 1548-1555 (1995).
- G. GRAGNOLI, A.M. SIGNORINI, S. TANGANELLI
Plasma levels of insulin, C-peptide and glucagon in liver cirrhosis.
J. Endocrinol. Invest., **4** : 1 (1981).
- K. KOSAKA, Y. AKAUMA, R. HOGURA, T. KUZUYA
Plasma insulin response during the 100 g glucose tolerance test in diabetic patients.
Excerpta Med. I.C.S., **468** : 283 (1979).
- H. KUZUKA et al.
Dissociation of peripheral insulin and C-peptide responses in diabetes mellitus.
Excerpta Med. I.C.S., **468** : 217 (1979).
- M. LAAKSO, S.V. EDELMAN, G. BRECHTEL, A.D. BARON
Decreased effect of insulin to stimulate skeletal muscle blood flow in obese man. A novel method for insulin resistance.
J. Clin. Invest., **85** : 1844-1852 (1990).
- G. MARCHESINI et al.
Insulin and glucagon levels in liver cirrhosis. Relationship with plasma amino-acid imbalance of chronic hepatic encephalopathy.
Dig. Dis. Sci., **24** : 594 (1979).
- M. MODAN et al.
Hyperinsulinemia. A link between hypertension, obesity and glucose intolerance.
J. Clin. Invest., **75** : 809-817 (1985).
- S. MÜLLNER, H. NEUBAUER, W. KÖNIG
A radioimmunoassay for the determination of insulins from several animal species, insulin derivatives and insulin precursors in both their native and denatured state.
J. Immunol. Meth., **140** : 211-218 (1991).
- T. POLLARE, H. LITHELL, C. BERNE
A comparison of the effects of hydrochlorothiazide and captopril on glucose and lipid metabolism in patients with hypertension.
N. Engl. J. Med., **321** : 868-873 (1989).
- D. PORTE, J.B. HALTER
The endocrine pancreas and diabetes mellitus.
In: Textbook of Endocrinology, R.H. Williams ed., W.B. Saunders, Philadelphia, p. 715 (1971).
- A. P. ROCCHINI et al.
Insulin and renal sodium retention in obese adolescents.
Hypertension, **14** : 367-374 (1989).

J.I. STARR et al.

Insulin, proinsulin and C-peptide.

In: Methods of Hormone Radioimmunoassay, B.M. Jaffe, H.R. Behrman eds., Academic Press, New York, p. 613 (1979).

K.W. TAYLOR

Insulin.

In: Hormone Assays and Their Clinical Application, J.A. Loraine, E.T. Bell eds., Churchill Livingstone, Edinburgh, p. 333 (1976).

R. TREVISAN et al.

Role of insulin and atrial natriuretic peptide in sodium retention in insulin-treated IDDM patients during isotonic volume expansion.

Diabetes, **39** : 289-298 (1990).

M. TROVATI et al.

Insulin influences the renin-angiotensin-aldosterone system in humans.

Metabolism, **38** : 501-503 (1989).

M.A. WILSON, L.E.M. MILES

Radioimmunoassay of insulin.

In: Handbook of Radioimmunoassay, G.E. Abraham ed., M. Dekker Inc., New York, p. 275 (1977).

R.L. YOUNG, R.J. FUCHS, M.J. WALTJEN, M.T. COOPERMAN


Glucose insulin response to oral glucose in a healthy obese population.

Diabetes, **28** : 208 (1979).

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710
FAX: 651-351-5669

APM0134

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

34672 7/10