VIP

¹²⁵I RIA Kit

For the quantitative determination of vasoactive
intestinal peptide (VIP) in plasma

Instruction Manual

Manuel d'Instructions

Testanleitung

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REF: 39125

DiaSorin

Stillwater, Minnesota 55082-0285, U.S.A.

TABLE OF CONTENTS

English	Page 1
Français	13
Deutsch	26
Italiano.....	40
Português.....	53

VASOACTIVE INTESTINAL PEPTIDE (VIP) RADIOIMMUNOASSAY

1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE

This kit is intended for the quantitative determination of vasoactive intestinal peptide (VIP) in plasma by radioimmunoassay (RIA).

2. SUMMARY AND EXPLANATION

VIP is a 28 amino acid compound originally extracted from the small intestine of the pig.^{1,2} Its amino acid sequence, compared to several related peptides, is as follows:

	1	2	3	4	5	6	7	8	9	10	11	12	13	14
VIP	His-	Ser-	Asp-	Ala-	Val-	Phe-	Thr-	Asp-	Asn-	Tyr-	Tyr-	Arg-	Leu-	Arg
Secretin	His-	Ser-	Asp-	Gly-	Thr-	Phe-	Thr-	Ser-	Glu-	Leu-	Ser-	Arg-	Leu-	Arg
Glucagon	His-	Ser-	Gln-	Gly-	Thr-	Phe-	Thr-	Ser-	Asp-	Tyr-	Ser-	Lys-	Tyr-	Leu
GIP	Tyr-	Ala-	Gln-	Gly-	Thr-	Phe-	Ile-	Ser-	Asp-	Tyr-	Ser-	Lys-	Ala-	Met
	15	16	17	18	19	20	21	22	23	24	25	26	27	28
VIP	Lys-	Gln-	Met-	Ala-	Val-	Lys-	Lys-	Tyr-	Leu-	Asn-	Ser-	Ile-	Leu-	Asn (NH ₂)
Secretin	Asp-	Ser-	Ala-	Arg-	Leu-	Gln-	Asp-	Leu-	Leu-	Gln-	Gly-	Leu-	Val	(NH ₂)
Glucagon	Asp-	Ser-	Arg-	Arg-	Ala-	Gln-	Asp-	Phe-	Val-	Gln-	Trp-	Leu-	Met-	Asn-Thr
GIP	Asp-	Lys-	Ile-	Arg-	Gln-	Gln-	Asp-	Phe-	Val-	Asn-	Trp-	Leu-	Leu-	Ala-etc.

There are similarities between the amino-terminal region of VIP and the amino-terminal regions of secretin, glucagon, and gastric inhibitory polypeptide (GIP). The carboxy-terminal region (i.e., residues 15-28), however, is quite different.

VIP has widespread occurrence throughout the gastrointestinal tract, where Bloom and Polak find tissue concentrations of VIP to be higher than any other gut peptide.³ Peripheral effects of VIP in the gut include: increased alkaline juice flow from the pancreas¹ inhibition of gastric acid production^{4,5} stimulation of insulin release⁶ hyperglycemia as a result of glycogenolysis⁷ stimulation of small intestinal juice flow and increased cyclic adenosine monophosphate (cAMP) content of the intestine.^{8,9}

VIP is also found in high concentration throughout the central nervous system, where it has strong vasodilator and hypotensive effects on the cardiovascular system, and in the respiratory system where it influences bronchodilation and augments ventilation.¹⁰ VIP is located in both neurons and endocrine cells, although more predominantly in the former, and appears to be released through neuronal stimulation.¹⁰⁻¹⁴

VIP measurement is used to research certain types of pancreatic islet tumors which secrete excessive amounts of the peptide. This condition has been designated as the Verner-Morrison syndrome,¹⁵ and is characterized by refractory watery diarrhea, hypokalemia, and achlorhydria (WDHA syndrome). Individuals with this syndrome have greatly elevated plasma VIP levels, often thousands of times the normal concentration.¹⁶⁻¹⁸ Another form of diarrhea, which manifests symptoms very similar to those of the Verner-Morrison syndrome, is the result of islet hyperplasia. Unlike the Verner-Morrison syndrome, however, this condition is characterized by normal VIP levels; these cases have been designated by Bloom and Polak as the pseudo-Verner-Morrison syndrome.^{3,19} Tumors other than those of the islet which produce elevated VIP levels include: ganglioneuroblastoma, bronchiogenic carcinoma, pheochromocytoma, medullary thyroid carcinoma, and retroperitoneal histiocytoma.^{15,21-23}

3. PRINCIPLE OF THE ASSAY

VIP is measured by a radioimmunoassay using delayed tracer addition to increase sensitivity. Sample and rabbit anti-VIP are added, followed by a 24 hour incubation at 2-8°C. ¹²⁵I VIP is then added, followed by a second incubation for 24 hours at 2-8°C. Pre-precipitated carrier, second antibody, and polyethylene glycol are added in a single step. The assay can be centrifuged and decanted after a minimum of 2 hours incubation at 2-8°C.

4. REAGENTS PROVIDED IN THE KIT

VIP Citrate Buffer	1 vial/2 mL
VIP Calibrator 0	1 vial/5 mL
VIP Calibrators (1-5)	5 vials/2 mL
VIP Antiserum	1 vial/12.5 mL
¹²⁵ I VIP	1 vial/12.5 mL
Precipitating Complex (GAR-PPT)	1 vial/35 mL
VIP Controls (Level 1 and 2)	2 vials/2 mL
Number of tests	125

STORAGE: Upon receipt, and prior to reconstitution, store all reagents at 2-8°C. After reconstitution store all reagents at -15° or lower until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 VIP Citrate Buffer: lyophilized reagent

BSA-citrate buffer contains thimerosal as a preservative. Reconstitute the vial with 2 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.2 VIP Calibrator 0: lyophilized reagent

Citrated plasma is diluted in BSA-citrate buffer containing thimerosal and other stabilizers. Reconstitute the vial with 5 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.3 VIP Calibrators: lyophilized reagent

VIP calibrators, at nominal concentrations ranging from 25-400 pg/mL, are prediluted in VIP-free plasma and BSA-citrate buffer. Exact concentration values are assigned according to each lot. Reconstitute each vial with 2.0 mL of purified water, mix and allow them to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using. This kit is standardized with Synthetic VIP. The calibrator demonstrates commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

4.4 VIP Antiserum: lyophilized reagent

Rabbit anti-VIP serum is diluted in BSA-borate buffer containing thimerosal. Reconstitute the vial with 12.5 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.5 ¹²⁵I VIP: lyophilized reagent

VIP is labeled with iodine-125 and diluted in BSA-citrate-EDTA buffer containing thimerosal and other stabilizers. Reconstitute the vial with 12.5 mL of purified water, mix and allow it to stand for 15-20 minutes until the contents are completely dissolved; mix thoroughly before using.

4.6 Precipitating Complex (GAR-PPT): lyophilized reagent

Normal rabbit serum, pre-precipitated with goat anti-rabbit serum and polyethylene glycol (PEG), is diluted in BSA-borate buffer containing thimerosal and other preservatives. Reconstitute each vial with 35 mL of purified water; mix thoroughly until the suspension appears homogeneous and then allow it to stand for a minimum of 30 minutes at room temperature with occasional mixing.

4.7 VIP Control (level 1 and level 2): lyophilized reagent

Defibrinated human plasma is spiked, if necessary, with the appropriate amount of VIP to obtain a concentration within a specified range. Sodium azide is added as a preservative. Reconstitute each vial with 2.0 mL of purified water, mix and allow them to stand for 15-20 minutes until the contents are completely dissolved.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV), or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING THIMEROSAL

Some reagents in this kit contain thimerosal which contains a mercury compound. Disposal of elemental mercury, inorganic mercury, mercury oxides, and mercury compounds should be done in strict compliance with all local, state, and federal regulations.

WARNING: This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 2 μ Ci (74 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.

5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

7. COLLECTION AND PREPARATION OF SPECIMEN

Two hundred microliters in duplicate, of EDTA plasma containing Aprotinin, are required for the VIP assay. Plasma samples may be assayed directly, with no extraction required.

The use of plasma containing EDTA and Aprotinin is necessary in the initial sample collection for VIP.

Collect blood in a 5 mL or 10 mL evacuated glass tube using EDTA as an anticoagulant. Aprotinin should be added immediately to the collection tube to prevent degradation of the VIP in the specimen at a concentration of 2,500 KIU/5 mL of whole blood. Immediately centrifuge the samples for 15 minutes at $760 \times g^{**}$ to obtain hemolysis-free plasma. Separate the serum from the cells and place in storage tubes. If the specimen is not used immediately, store it at -15°C or lower. Specimens should not be repeatedly frozen and thawed.

All plastics, glassware or other materials coming into contact with the specimen should be entirely free of any contamination.

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Temperature controlled centrifuge to accommodate 12 x 75 mm tubes.
- 8.3 Gamma scintillation counter capable of counting iodine-125.
- 8.4 Vortex.
- 8.5 Suggested pipetting devices:
 - a. Micropipettors calibrated to deliver 100 μL and 200 μL .
 - b. Repeating dispensers calibrated to deliver 100 μL , 500 μL , and 1 mL.
- 8.6 0.85% saline.
- 8.7 Purified water.

9. ASSAY PROCEDURE

CAUTION: SET UP ASSAY ON **CRUSHED ICE**. VIP IS LABILE IN PLASMA STORED AT ROOM TEMPERATURE; THEREFORE, THIS PRECAUTION IS NECESSARY TO OBTAIN VALID RESULTS.

- 9.1 Reconstitute lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 25°C . Gently mix all reagents before using.

$$** g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.2 Set up labeled 12 x 75 mm disposable glass tubes in duplicate according to the Scheme of the Assay, on the back page.
- 9.3 Add reagents as follows:
- a. **Total count tubes**
Set aside until step 5
 - b. **Nonspecific binding tubes (NSB)**
200 μ L of Calibrator 0
100 μ L of VIP citrate buffer
 - c. **Calibrator 0**
200 μ L of Calibrator 0
100 μ L of VIP antiserum
 - d. **VIP calibrators (A-E)**
200 μ L of VIP calibrator
100 μ L of VIP antiserum
 - e. **Controls and unknown samples**
200 μ L of plasma
100 μ L of VIP antiserum
- 9.4 Vortex the tubes gently without foaming and incubate them for 24 hours (\pm 2 hours) at 2-8°C.
- 9.5 Add 100 μ L of 125I VIP to all tubes.
- 9.6 Vortex the tubes gently without foaming and incubate them for 24 hours (\pm 2 hours) at 2-8°C.
- 9.7 Add 1 mL of 0.85% saline to all tubes except the total count tubes.
- 9.8 Vigorously mix the GAR-PPT; add 500 μ L to all tubes except the total count tubes.
- 9.9 Vortex the tubes gently without foaming and incubate them for 2 hours (\pm 15 minutes) at 2-8°C.
- 9.10 Centrifuge the tubes using 760 x g* for 20 minutes at 20-25°C.
- 9.11 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a minimum time of 2 minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 9.12 Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for 60 seconds or longer. (see Limitations of the Procedure section).

10. PROCEDURAL COMMENTS

- 10.1 Assay all samples in duplicate to ensure confidence in values obtained.
- 10.2 Add each aliquot of reagent to the lower third of assay tube to ensure complete mixture of reagents.
- 10.3 Some manufacturers' disposable tubes yield elevated nonspecific bindings.
- 10.4 If you choose to aspirate the supernatant from the precipitate, be careful not to disturb the precipitate.
- 10.5 Any value greater than the highest calibrator must be diluted with calibrator 0 and assayed again. Correct result using the appropriate dilution factor.
- 10.6 When performing the centrifugation step, the temperature must be controlled so that it does not exceed 25°C.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

10.7 To completely monitor the consistent performance of an RIA there are additional factors which may be checked. DiaSorin suggests a check of the following parameters to assure consistent kit performance.

a. Total Counts

b. Maximum Binding

Average counts per minute (CPM) of calibrator 0 Tubes/Average CPM of Total Count Tubes.

c. Nonspecific Binding

Average CPM of NSB Tubes/Average CPM of Total Count Tubes.

d. Slope of Calibrator Curve

For example, monitor the 80, 50 and 20% suppression points of the calibrator line.

11. QUALITY CONTROL

Each laboratory should include at least two kit controls in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included with this kit. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

12. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either a linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is %B/B₀ versus log concentration.

12.1 Calculate the average CPM for each calibrator, control and unknown sample.

12.2 Subtract the average CPM of the NSB tubes from all counts.

12.3 Divide the corrected CPM of each calibrator, control, or sample by the corrected CPM of the calibrator 0.

$$\% B/B_0 = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of Calibrator 0} - \text{CPM of NSB}} \times 100$$

12.4 Using two cycle semi-log or log-logit graph paper, plot percent B/B₀ for the VIP calibrators (vertical axis) versus the concentration (horizontal axis).

12.5 Draw a best fit line through the points.

12.6 Interpolate the levels of VIP in the unknown samples from the plot.

12.7 If an unknown sample has been diluted, correct for the appropriate dilution factor.

12.8 Calculate maximum binding by dividing CPM of calibrator 0 by the average total counts obtained in the total count tubes.

TABLE III
DiaSorin VIP RIA Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent (B/T)	Percent Bound (B/B ₀)	Conc. (pg/mL)
Total Count	10,928	10,930				
	10,931					
NSB	808	932		8.5		
	1,055					
Calibrator 0	4,603	4,709	3,777	43.1	100.0	
	4,834					
Calibrators (pg/mL)						
A (19)	4,034	4,046	3,114		82.5	
	4,058					
B (30)	3,542	3,558	2,626		69.5	
	3,573					
C (54)	2,768	2,767	1,835		48.8	
	2,766					
D (140)	1,948	1,988	1,056		28.0	
	2,028					
E (290)	1,297	1,352	420		11.1	
	1,408					
Unknown Samples						
1	2,051	2,060	1,128		29.9	126
	2,070					
2	3,210	3,209	2,277		60.3	40
	3,207					

Typical sample data and a calibrator curve are shown in TABLE III and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

VIP SAMPLE CALIBRATOR CURVE


FIGURE 1

The DiaSorin QC lab uses a smoothed spline curve fit.

13. LIMITATIONS OF THE PROCEDURE

- 13.1 If the initial concentration of the unknown sample is greater than the highest calibrator, dilute with calibrator 0 only.
- 13.2 Specimens should not be repeatedly frozen and thawed.
- 13.3 Counting times should be sufficient to prevent statistical error (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).

14. EXPECTED VALUES

CLINICAL INTERPRETATION OF VALUES

Forty-five normal fasting persons (mean age 28) were sampled for plasma VIP determination. In 21 males, the value was 42.8 ± 11.9 pg/mL VIP (mean \pm standard deviation), and in 24 females, 42.6 ± 7.2 pg/mL. Since there is no significant difference between the values for males and females, the overall average is 42.7 ± 9.5 pg/mL in fasting adults. A tentative two standard deviation range would be 23-63 pg/mL. This normal range may be broader if persons from a hospital population without gastrointestinal disease are used as a control group. Others report the normal range to go up to 200 pg/mL. The normal range in the assay at Ohio State University is undetectable to 170 pg/mL, with an average of 62 ± 44 pg/mL. Since the Ohio State assay report values below 50 as undetectable, their assay is not quite as sensitive as the DiaSorin procedure.

A series of doubled-blinded clinical specimens,* including several VIPomas and ganglioneuroblastomas, were obtained from Dr. Thomas O'Dorisio at Ohio State University. Later, he revealed his results on these same samples; these are shown in TABLE IV. The data demonstrate many extremely elevated VIP values, which are characteristic of VIP-producing tumors, and show that the DiaSorin VIP kit does indeed measure elevated values in humans and compares well with an established research laboratory. The regression equation for the first 22 specimens (excluding the 4 normals with undetectable levels reported by Ohio State) is $y = 0.78x + 32$, where y is the DiaSorin value and x is the Ohio State value. The correlation coefficient is 0.94.

The control samples used in the Ohio State laboratory were also measured with the DiaSorin kit. Their low control (Ohio State range <50-100 pg/mL) had a mean value of 53 pg/mL and the high control (Ohio State range 450-1014 pg/mL) was found to be 1077 pg/mL in 3 separate analyses with the DiaSorin kit. Several of these Ohio State samples from normal persons had values ranging from 65-176 pg/mL on the DiaSorin assay. The broad range of pathological samples in WDHA syndrome is typical of the literature and quite similar to the situation of gastrin levels in gastrinoma.

* These specimens had been stored for an unknown amount of time in Dr. O'Doriso's laboratory. The number of freeze and thaws between their analysis and our analysis is unknown.

TABLE IV
Results of Double-Blind Split Samples With
Dr. Thomas O'Doriso at Ohio State University
values: picograms/ML

Case	VIP DiaSorin	VIP Ohio State	Comments
I. Ganglioneuroblastoma			
1	1888	1700	symptomatic WDHA
1A	248	135	asymptomatic
1B	464	430	symptomatic WDHA
2	3776	>2000	
2A	4672	1700	
2B	352	225	post resection with diarrhea
II. VIPoma			
1	368*	960	
1A	184	290	post resection
2	1152	1200	metastatic
3	304	290	metastatic after resection
4	208	230	post resection
III. Watery diarrhea with no documented diagnosis			
1	11250	>14500	
2	368	185	
3	928	880	neg. exploration
IV. Watery diarrhea with other tumors			
1	192	390	hypernephroma
2	2360	1500	smallcell carcinoma
3	208	330	hypernephroma with vasodilation
4	272	360	pheochromocytoma
V. Normal patients			
1	64	<50	
2	96	<50	
3	104	<50	
4	96	<50	
5	176	110	
VI. Unknown			
1	720	1150	
2	512	1200	
3	320	170	

* DiaSorin assay was done after several freeze and thaws.

** Regression analysis excluding <50 pairs; $y = 0.79x + 335$; $r = 0.94$;
 $y = \text{DiaSorin VIP}$; $x = \text{Ohio State VIP}$

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Precision

Intra Assay Variation (values = 5)

Sample Number	Mean Value	S.D	%C.V.
Low	62.6	1.2	1.9
High	87.8	2.6	3.0

Inter Assay Variation (values = 73)

Sample Number	Mean Value	S.D	%C.V.
Low	56.3	5.4	9.58
High	80.4	5.7	7.07

15.2 TRUENESS: THE ASSAY TRUENESS HAS BEEN VERIFIED BY THE DILUTION TEST AND THE RECOVERY TEST.

Linearity (Parallelism)

Serial Dilution Study of five Patient Samples (Values = pg/mL)

Sample Number	1/8	1/16	1/32
1	1584	1696	1888
2	3040	4192	4672
3	>3200	3440	3776
	1/10	1/15	1/20
4	2630	2130	2360
	1/100	1/150	1/200
5	7200	11250	10400

Recovery (Accuracy) (Values = pg/mL)

TABLE V

Recoveries of VIP at Several Level Approximating the Normal Range from Fresh Plasma Containing EDTA and Aprotinin (Each point is the mean of duplicate analyses)

Plasma Sample	VIP Added (pg/mL)	VIP Found (pg/mL)	Recovery %	VIP Plasma Sample	VIP Added (pg/mL)	VIP Found (pg/mL)	Recovery %
1	0	14	--	3	0	14	--
	25	38	97		25	33	85
	50	50	78		50	48	75
	100	85	75		100	86	75
2	0	2	--	4	0	14	--
	25	23	85		25	38	90
	50	48	92		50	50	75
	100	80	78		100	80	68
			Mean Recovery	25	89%		
				50	80%		
				100	74%		

15.3 Analytical Sensitivity

When defined as the apparent concentration at 3 standard deviations from the counts at maximum binding, the minimum detectable amount is 1.0 pg/tube (5 pg/mL).

15.4 Analytical Specificity

Comparison of the cross-reactivity of the VIP antibody was made with the following peptides:

Peptides	% Cross-reactivity
Somatostatin	<0.1
pMotilin	<0.1
pGastrin Inhibitory Polypeptide (GIP)	<0.1
Glucagon	<0.1
pGastrin Releasing Peptide (GRP)	<0.1
Gastrin ¹⁷ I	<0.1
Gastrin ³⁴ I	<0.1
pSecretin	<0.1
hPancreatic Polypeptide (PPP)	<0.1
hProinsulin	<0.1
pInsulin	<0.1
Cholecystokinin (CCK8)	<0.1
Substance P	<0.1
Neurotensin	<0.1
Methionine Enkephalin	<0.1
Leucine Enkephalin	<0.1
Bombesin	<0.1

REFER TO LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

1. CAUTION: Set up assay on crushed ice. VIP is labile in plasma stored at room temperature; therefore, this precaution is necessary to obtain valid results
2. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow reagents to reach above 25°C. Identify tubes in duplicate.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0	Cal 1-5	Controls and unknown samples
VIP Citrate Buffer	-	100 µL	-	-	-
Calibrator (0)	-	200 µL	200 µL	-	-
Calibrators (1-5)	-	-	-	200 µL	-
Controls and Unknown Samples	-	-	-	-	200 µL
Antiserum	-	-	100 µL	100 µL	100 µL

4. Cover the tubes with parafilm and vortex gently.
5. Incubate for 24 hours +/- 2 hours at 2-8°C.
6. Dispense 100 µL of the tracer into all wells.
7. Cover and vortex gently.
8. Incubate for 24 hours +/- 2 hours at 2-8°C.
9. Dispense 1 mL of 0.85% saline into all wells except the total count tubes.
10. Dispense 500 µL of the GAR-PPT into all wells except the total count tubes.
11. Cover the tubes and vortex gently.
12. Incubate for 2 hours +/- 15 minutes at 2-8°C.
13. Centrifuge using 760 x g* for 20 minutes.
14. Decant the supernatants.
15. Count each tube in a gamma counter for 60 seconds or longer.

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$$*g = (1118 \times 10^8) (\text{raio em cm}) (\text{rpm})^2$$

REFERENCES/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFIA/REFERÊNCIAS

1. Said, S.I. and V. Mutt, "Isolation from Porcine-Intestinal Wall of a Vasoactive Octa-cosapeptide Related to Secretin and to Glucagon," **European Journal of Biochemistry**, 28:199, (1972).
2. Said, S.I. and V. Mutt, "Polypeptide with Broad Biological Activity: Isolation from Small Intestine," **Science**, 169:1217, (1970).
3. Bloom, S.R. and J.M. Polak, "VIP Measurement in Distinguishing Verner-Morrison Syndrome and Pseudo Verner-Morrison Syndrome," **Clinical Endocrinology**, 5:223, (1976).
4. Barbezat, G.O. and M.I. Grossman, "Intestinal Secretion: Stimulation by Peptides," **Science**, 174:422, (1971).
5. Makhlof, G.M. and S.I. Said, "The Effect of Vasoactive Intestinal Peptide (VIP) on Digestive and Hormonal Function," In **Gastrointestinal Hormones**, Ed. J.C. Thompson, Texas University Press, pp. 599-610, (1975).
6. Kerins, C. and S.I. Said, "Hyperglycemic and Glycogenolytic Effects of Vasoactive Intestinal Polypeptide," **Proceedings of the Society for Experimental Biology and Medicine**, 142:1014, (1972).
7. Schebalin, M., S.I. Said and G.M. Makhlof, "Stimulation of Insulin and Glucagon Secretion by Vasoactive Intestinal Peptide," **American Journal of Physiology**, 232: 197, (1977).
8. Schwartz, C.J., D.V. Kimberg, H.E. Sheerin, M. Field and S.I. Said, "Vasoactive Intestinal Peptide Stimulation of Adenylate Cyclase and Active Electrolyte Secretion in Intestinal Mucosa," **Journal of Clinical Investigation**, 54:536, (1974).
9. Krejs, G.J., J.S. Fordtran, S.R. Bloom, J. Fahrenkrug and O.B. Schaffalitzky De Muckadell, "Effect of VIP Infusion on Water and Ion Transport in the Human Jejunum," **Gastroenterology**, 78:722, (1980).
10. Said, S.I., "Vasoactive Intestinal Polypeptide (VIP): Current Status," In **Gastrointestinal Hormones**, Ed. J.C. Thompson, Texas University Press, pp. 591-597, (1975).
11. Bryant, M.G., S.R. Bloom, J.M. Polak, R.H. Albuquerque, I. Modlin and A.G.E. Pearse, "Possible Dual Role for Vasoactive Intestinal Peptide as Gastrointestinal Hormone and Neurotransmitter Substance," **Lancet**, 1:991, (1976).
12. Domschke, S., W. Domschke, S.R. Bloom, P. Mitznegg, S.J. Mitchell, G. Lux and U. Strunz, "Vasoactive Intestinal Peptide in Man: Pharmacokinetics, Metabolic and Circulatory Effects," **Gut**, 19:1049, (1978).
13. Fahrenkrug, J., "Vasoactive Intestinal Polypeptide: Measurement, Distribution and Putative Neurotransmitter Function," **Digestion**, 19:149, (1979).
14. Uddman, R., J. Fahrenkrug, L. Malm, J. Alumets, R. Hakanson and F. Sundler, "Neuronal VIP in Salivary Glands: Distribution and Release," **Acta. Physiol. Scand.**, 110:31, (1980).
15. Verner, J.V. and A.B. Morrison, "Islet Cell Tumor and a Syndrome of Refractory Watery Diarrhea and Hypokalemia," **American Journal of Medicine**, 25:374, (1958).
16. Bloom, S.R. and J.M. Polak, "Glucagonomas, VIPomas and Somatostatinomas," **Clinics in Endocrinology and Metabolism**, 9(2):285, (1980).
17. Modlin, I.M. and S.R. Bloom, "VIPomas and the Watery Diarrhea Syndrome," **South African Medical Journal**, 54:53, (1978).
18. Yamaguchi, K., K. Abe, S. Miyakawa, S. Ohnami, M. Sakagami and N. Yanaihara, "The Presence of Macromolecular Vasoactive Intestinal Polypeptide (VIP) in VIP-Producing Tumors," **Gastroenterology**, 79(4), (1980).

19. Kidd, G.S., M. Donowitz, T. O'Dorisio, S. Cataland and F. Newman, "Mild Chronic Watery Diarrhea-Hypokalemia Syndrome Associated with Pancreatic Islet Cell Hyperplasia," **The American Journal of Medicine**, 66:883, (1979).
20. O'Dorisio, M.S., T.M. O'Dorisio, S. Cataland and S.P. Balcerzak, "Vasoactive Intestinal Polypeptide as a Biochemical Marker of Polymorphonuclear Leukocytes," **Journal of Laboratory Clinical Medicine**, 96(4):666, (1980).
21. Kaplan, S.J., C.T. Holbrook, H.G. McDaniel, W.L. Buntain and W.M. Crist, "Vasoactive Intestinal Peptide Secreting Tumors of Childhood," **American Journal Diseases of Child.**, 134:21, (1980).
22. Iida, Y., O. Nose, H. Kai, A. Okada, T. Mori, P. Lee, K. Kakudo and N. Yanaihara, "Watery Diarrhea with a Vasoactive Intestinal Peptide-producing Ganglioneuroblastoma," **Archives of Disease in Childhood**, 55:929, (1980).
23. Sagor, G.R., M.G. Bryant, M.A. Ghatei, R.M. Kirk and S.R. Bloom, "Release of Vasoactive Intestinal Peptide in the Dumping Syndrome," **British Medical Journal**, 282:507, (1981).

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso
	Maximum temperature	Température maximale	Höchsttemperatur	Temperatura massima
IVD	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnostica in vitro.
LOT	Lot No.	No. de lot	Chargen-Nr.	Lotto n°
Ab	Antiserum	Antisérum	Antiserum	Antisiero
Ab PEG	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reagente precipitante
Ag ¹²⁵ I	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	¹²⁵ I-markiertes Antigen	Tracciatore: antigene etichettato con ¹²⁵ I
BUF	Buffer	Tampon	Puffer	Tampone
CAL	Calibrator	Étalon	Kalibrator	Calibratore
CONTROL	Control serum	Sérum de contrôle	Kontrollserum	Siero di controllo
	Radioactive	Radioactif	Radioaktiv	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo


SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482
In the United Kingdom Call: +44(0) 1344 401 430 FAX: +44(0) 7884 050812
10403
34675 7/10

PRINTED IN U.S.A.

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DiaSorin Inc.