**CYCLO-Trac® SP-Whole Blood
Radioimmunoassay for Cyclosporine**

For the quantitative measurement of cyclosporine
(Cyclosporin A, Cyclosporin, Sandimmune®) in whole blood

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1. INTENDED USE

FOR IN VITRO DIAGNOSTIC USE.

The CYCLO-Trac[®] SP Whole Blood radioimmunoassay (RIA) kit is intended for the quantitative measurement of cyclosporine (Cyclosporin A, Cyclosporin, Sandimmune[®]) in whole blood as an aid in the management of heart, kidney, and liver transplant patients receiving cyclosporine.

CAUTION: U.S. FEDERAL LAW RESTRICTS THIS DEVICE TO SALE AND DISTRIBUTION BY, OR ON THE ORDER OF, A PHYSICIAN, OR TO A CLINICAL LABORATORY; AND USE IS RESTRICTED TO, BY, OR ON THE ORDER OF, A PHYSICIAN. THE CYCLO-TRAC SP RIA KIT IS FOR USE ONLY WITH CYCLO-SPORINE.

2. SUMMARY AND EXPLANATION

Cyclosporine (Cyclosporin A, Cyclosporin, Ciclosporin, Sandimmune) is a cyclic undecapeptide of fungal origin with powerful immunosuppressive properties.^{1, 2} While its mechanism of action is still under investigation, cyclosporine appears to affect the metabolism of both T-helper lymphocytes and T-suppressor lymphocytes, which impairs the immune system.³⁻⁵

Though this drug has demonstrated utility in preventing the rejection of new organ transplants,⁶⁻⁹ a number of adverse effects,^{10,11} the most serious of which are nephro- and hepa-toxicity associated with high blood levels of this drug,¹²⁻¹⁶ have been reported. Moreover, absorption and metabolism of cyclosporine is highly variable from patient to patient.¹⁷⁻¹⁹ Thus, dosages must be individualized and monitored to achieve optimal immunosuppression and minimize adverse side effects.^{18,19} Nearly all of the cyclosporine present in whole blood is bound to red blood cells and plasma proteins.²⁰⁻²² This distribution is affected by changes in temperature and hematocrit.^{21,22} To ensure consistent results, a standardized protocol for blood collection is necessary.¹⁹

Two methods are currently available to monitor cyclosporine levels in patients: high performance liquid chromatography (HPLC) and radioimmunoassay (RIA). First generation RIA's employed polyclonal antibodies which exhibit high cross-reactivity with circulating metabolites.²³ Cyclosporine values generated by these RIA's invariably exceeded those generated by HPLC.²⁴ Second generation RIA's are based on monoclonal antibody technology.²⁵ Both specific and nonspecific monoclonal antibodies have been developed. The specific antibody measures only the parent compound of cyclosporine and generates values similar to those obtained by HPLC. The nonspecific antibody measures both parent compound and metabolites and produces results similar to those of first generation RIA's.²⁶

3. PRINCIPLE OF THE ASSAY

The DiaSorin CYCLO-Trac SP Whole Blood RIA employs a specific monoclonal antibody that measures only cyclosporine and demonstrates no significant cross-reactivity with its metabolites. A methanol extraction step is required for the calibrators, controls, and samples prior to assay. The assay itself consists of a single incubation step and uses a double antibody phase separation technique. To simplify the procedure, the first antibody (mouse monoclonal specific to Cyclosporin A), and the second antibody (donkey anti-mouse) are premixed with a carrier (normal mouse serum) to produce a single reagent called Anti-CYCLO-Trac SP ImmunoSep.

To determine cyclosporine levels with this kit, whole blood samples, calibrators, and controls are first extracted with methanol. The methanol extracts are then combined with Iodine-125 labeled cyclosporine tracer and Anti-CYCLO-Trac SP ImmunoSep. Following a 1 hour incubation, the tubes are centrifuged, decanted, and then counted. The amount of radio-activity remaining in the pellet is inversely proportional to the concentration of cyclosporine found in the sample.

4. REAGENTS PROVIDED IN THE KIT

CYCLO-Trac SP NSB ImmunoSep (BLUE)	1 vial/9 mL
CYCLO-Trac SP Calibrators (0-6)	1 vial/2 mL, 5 vials/1 mL
CYCLO-Trac SP ImmunoSep	2 vials/100 mL
CYCLO-Trac SP Tracer (RED)	1 vial/20 mL
CYCLO-Trac SP Controls	2 vials/1 mL
Number of tests	200

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent at 2-8°C until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

Reagents from different batches must not be mixed.

4.1 CYCLO-Trac SP NSB ImmunoSep (blue): ready to use reagent

Donkey anti-mouse serum/normal mouse serum premix is diluted in phosphate-buffered saline containing 0.08% sodium azide and a surfactant. Blue dye is added.

Shake vigorously immediately before use.

4.2 CYCLO-Trac SP Calibrators (0-6): ready to use reagent

Cyclosporine-free serum contains 0.1% sodium azide and a surfactant. Six cyclosporine calibrators, at nominal concentrations ranging from 0-1200 ng/mL, are prediluted in human serum containing 0.1% sodium azide and a surfactant. Exact concentration values are assigned according to each lot and appear on vial labels. The DiaSorin CYCLO-Trac SP calibrator has been calibrated against the USP Cyclosporine Reference Standard. Any comparison with other products or procedures should be done with this reference standard. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

4.3 Anti-CYCLO-Trac SP ImmunoSep: ready to use reagent

Monoclonal mouse anti-cyclosporine antibody and donkey anti-mouse/normal mouse serum premix is diluted in phosphate-buffered saline containing 0.08% sodium azide and a surfactant. **Shake vigorously immediately before use.**

4.4 ¹²⁵I CYCLO-Trac SP (red): ready to use reagent

A derivative of cyclosporine is labeled with Iodine-125 and diluted in citrate-EDTA buffer which contains 0.05% sodium azide, human serum and a surfactant. Red dye is added.

4.5 CYCLO-Trac SP Controls (level 1 and level 2): ready to use reagent

Human serum is spiked with the appropriate amount of cyclosporine A to obtain a concentration within a specified range. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. 0.1% sodium azide and other stabilizers are added.

5. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 9 μCi (333 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of abnormal particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.
- 6.5 Poor reproducibility between duplicates.

7. SPECIMEN COLLECTION AND PREPARATION

One hundred microliters of EDTA whole blood are required for the CYCLO-Trac SP whole blood extraction procedure.

NOTE: Heparinized samples should not be used.^{19, 27}

COLLECTION AND STORAGE OF WHOLE BLOOD

Collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. EDTA (7.2 mg/ 5 mL blood) should be used as an anticoagulant. No further additives or preservatives are required to maintain the integrity of the sample. If specimens are not used immediately, store at 2-8°C or lower for up to 7 days. For prolonged storage (up to 3 months) freeze samples at -20°C or lower.^{28, 29}

Whole blood samples that have been stored at 4°C, or have not been immediately processed, should be thoroughly mixed prior to taking samples for testing in the CYCLO-Trac SP whole blood kit.

8. EQUIPMENT AND MATERIAL REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Temperature controlled centrifuge to accommodate 12 x 75 mm tubes.
- 8.3 Gamma scintillation counter capable of counting Iodine-125.
- 8.4 Vortex.
- 8.5 Pipetting devices:
 - a. 400 µL, 100 µL and 50 µL positive displacement pipettes.
 - b. Repeating dispensers calibrated to deliver 1.0 mL and 100 µL.
- 8.6 Methanol (HPLC grade).
- 8.7 Caps for extraction tubes (Scientific Products No. T1226-32).

9. ASSAY PROCEDURE

- 9.1 Label one 12 x 75 glass test tube for each calibrator, control and unknown sample.
- 9.2 Dispense 400 µL of methanol (HPLC grade) into each tube using a positive displacement pipette.* Cover tubes with parafilm if samples are not yet ready to be added.
- 9.3 Add the calibrators, the controls and the whole blood samples into the methanol using a 100 µL positive displacement pipette.*

NOTE: Slowly dispense each sample with the pipette tip placed slightly under the surface of the methanol.
- 9.4 Immediately after the sample is added, cap each tube and vortex it for 15 seconds.
- 9.5 Centrifuge the capped tubes using 1600 x g** at 20-25°C for 5 minutes.
- 9.6 Take care not to allow the methanol extracts to evaporate before pipetting into the assay tubes.
- 9.7 Set-up labeled 12 x 75 mm test tubes according to the protocol Scheme of the Assay (back cover).
- 9.8 Allow reagents to equilibrate to room temperature (20-25°C) before using in the assay.

* Alternately, 200 µL of sample may be extracted with 800 µL of methanol.

** $g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$

- 9.9** Add reagents as follows:
- a. Total count tubes**
100 µL of ¹²⁵I CYCLO-Trac SP (red)
 - b. Nonspecific binding tubes (NSB)**
50 µL of extracted 0 calibrator
100 µL of ¹²⁵I CYCLO-Trac SP (red)
1.0 mL of thoroughly mixed CYCLO-Trac SP NSB ImmunoSep (blue)
 - c. CYCLO-Trac SP 0 calibrator**
50 µL of extracted CYCLO-Trac SP 0 Calibrator
100 µL of ¹²⁵I CYCLO-Trac SP (red)
 - d. CYCLO-Trac SP calibrators (A-E)**
50 µL of extracted CYCLO-Trac calibrator
100 µL of ¹²⁵I CYCLO-Trac SP (red)
 - e. Control and unknown samples**
50 µL of extracted sample
100 µL of ¹²⁵I CYCLO-Trac SP (red)
- 9.10** Shake the bottle of Anti-CYCLO-Trac SP ImmunoSep vigorously; immediately add 1.0 mL of well-mixed reagent to all tubes except the total count and NSB tubes.
- 9.11** Vortex the tubes vigorously and incubate them for 1 hour at 20-25°C.
- 9.12** Centrifuge the tubes using 1600 x g* for 20 minutes at 20-25°C.
- 9.13** Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a maximum time of 2 minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning tubes upright.
- 9.14** Using a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy (see Limitations of Procedure section).

10. PROCEDURAL COMMENTS

- 10.1** In order to obtain maximum accuracy and precision, the methanol extraction steps must be carefully followed; positive displacement pipettes (SMI) must be used to dispense both the methanol and the samples. Positive displacement pipettes are also required for transferring 50 µL of methanol extract to each assay tube.
- 10.2** Add each aliquot of reagent to the lower third of assay tube to ensure complete mixture of reagents.
- 10.3** Some manufacturers' disposable tubes may yield elevated nonspecific bindings or cause increased pellet slippage.
- 10.4** ImmunoSep reagents must be vigorously mixed immediately before use.
- 10.5** If a sample reads greater than the highest calibrator, re-extract the sample using a blood: methanol ratio greater than 1:4 (1:9 or 1:19 are recommended). Multiply the assay results by the appropriate correction factor.
- 10.6** If the tubes cannot be decanted within 5 minutes after the centrifuge has stopped, it is recommended to re-centrifuge the tubes.
- 10.7** Run the assay immediately after the extraction step.

* $g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$

11. QUALITY CONTROL

Each laboratory should include at least two control samples in every assay to ensure the validity of each assay's results. If you are using unassayed controls and have not established a mean and standard deviation, a mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included with this kit. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. It is not recommended that kit reagents be interchanged with the exception of ¹²⁵I CYCLO-Trac SP tracer.

In order for a laboratory to completely monitor the consistent performance of an RIA assay, there are additional factors which must be checked. DiaSorin suggests a regular check of the following parameters to assure consistent kit performance.

a. Total Counts

b. Maximum Binding

Counts per minute (CPM) of 0 Calibrator Tube/Average CPM of Total Count Tubes.

c. Nonspecific Binding

CPM of NSB Tube/Average CPM of Total Count Tubes.

d. Slope of Calibrator Curve

For example, monitor the 80, 50 and 20% points of the calibrator line for run-to-run reproducibility.

12. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIA's. Each is based on obtaining a calibration curve by plotting the extent of binding against the concentrations of the calibration calibrators. The curve for the CYCLO-Trac SP assay is obtained by plotting % B/B₀ vs. log concentration.

12.1 Calculate the average CPM for each calibrator, control and patient sample.

12.2 Subtract the average CPM of the NSB tubes from all counts to obtain the corrected CPM.

12.3 Divide the corrected CPM of each calibrator, control or patient sample by the corrected CPM of the 0 calibrator.

$$\% B/B_0 = \frac{\text{CPM of calibrator or patient sample} - \text{CPM of NSB}}{\text{CPM of 0 calibrator} - \text{CPM of NSB}} \times 100$$

12.4 Using three cycle semi-log graph paper, plot % B/B₀ of the calibrators (Y-axis) against their corresponding concentration (X-axis).

12.5 Draw a best-fit curve through the points.

12.6 Interpolate cyclosporine concentrations of the patient samples from the calibrator curve.

In practice, most laboratories will employ a computer based data reduction system. DiaSorin recommends using either a spline or 4-parameter logistic (4PL) curve fitting program.

NOTE: DO NOT use a linear concentration scale (X-axis) and or a logit B/B₀ scale (Y-axis).

REDUCTION DATA

The DiaSorin QC lab uses a smoothed spline curve fit.

TABLE I
DiaSorin CYCLO-Trac SP Whole Blood Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent Bound (B/T)	Percent (B/B ₀)	Conc. (ng/mL)
Total Count	51045	50892				
	50738					
NSB	1221	1328		2.6		
	1434					
0 Calibrator	25323	25432	24104	50.0		
	25541					
Calibrators (ng/mL)						
A (20)	22781	22806	21478		89.1	
	22830					
B (61)	19850	19624	18296		75.9	
	19397					
C (154)	14843	14911	13583		56.4	
	14978					
D (379)	9506	9570	8242		34.2	
	9633					
E (1099)	4893	4858	3530		14.6	
	4823					
Unknown Samples						
1	16225	16122	14794		61.4	126
	16019					
2	9764	9570	8242		34.2	379
	9375					

Typical sample data and a calibrator curve are shown in TABLE I and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

CYCLO-Trac[®] SP WHOLE BLOOD SAMPLE CALIBRATOR CURVE


FIGURE 1

13. LIMITATIONS OF THE PROCEDURE

- 13.1 Specimens should not be repeatedly frozen and thawed.
- 13.2 Counting times should be sufficient to prevent statistical error.
- 13.3 Monoclonal antibodies such as OKT3 have been evaluated and shown not to interfere in the assay.

14. EXPECTED VALUES

Therapeutic blood levels for cyclosporine in kidney, heart, and liver transplant patients have not been well established. This has in part been due to the varied treatment protocols for transplantation practiced, the methods available for measuring cyclosporine and the concomitant administration of other immunosuppressive medications with cyclosporine.^{19, 30, 31}

A clinical evaluation using the CYCLO-Trac SP Whole Blood kit was conducted at the University Hospitals of London, Ontario and Albany, New York, in kidney, heart, and liver transplant patients. These evaluations measured temporal trough cyclosporine blood levels in transplant patients. Trough blood values were observed post transplant from 4 days to 90 days in kidney (n = 10), heart (n = 10), and liver (n = 10) transplant patients at the University Hospital, London Ontario. In addition, a similar evaluation was conducted with kidney transplant patients (n = 18) at University Hospital, Albany, New York. At Albany, trough blood levels were observed from transplant to 199 days post transplant.

TABLE II shows the mean cyclosporine blood values and significant events (rejection or nephrotoxicity) observed in transplant patients from the University Hospital, London, Ontario. At this center, patients' blood levels were maintained at 400 to 650 ng/mL cyclosporine depending upon transplant type and clinical condition of the patients.

TABLE II

Temporal cyclosporine whole blood values in kidney, heart and liver transplant patients as measured by the CYCLO-Trac SP Whole Blood Kit.

Post Operative	Kidney ^a				Heart ^b				Liver ^c					
	N	Event ^d	Mean	S.D.	N	Event ^d	Mean	S.D.	N	Event ^d	Mean	S.D.		
0	-	-	-	-	-	-	-	-	-	-	-	-		
4	8	2	438	269	9	-	161	115	-	-	-	-		
7	10	3	392	138	10	-	326	176	9	-	372	111		
11	9	6	428	207	10	-	563	119	10	-	496	136		
14	10	5	387	219	10	-	598	162	10	1	620	231		
21	10	3	487	283	10	5	505	185	10	-	478	109		
30	9	3	416	163	9	4	450	102	10	1	608	217		
60	9	1	579	474	9	3	482	171	8	1	437	66		
90	9	1	252	105	8	1	433	173	8	-	347	102		
Total events		24					13					3		

^a Cyclosporine: 3.75 mg/kg to 5 mg/kg I.V. changing to 12 mg/kg in two divided doses as soon as possible, prednisone 2-5 mg/kg with taper to 20 mg/day, methylprednisolone 1-2 mg/kg. If no kidney function, cyclosporine was discontinued and OKT3 given 5 mg/kg for 10 days.

^b Cyclosporine begun 3rd post operative day 5 mg/kg in two divided doses, prednisone 15 mg daily taper gradually, azathioprine in addition as indicated, ALG 50 mg/kg given once and continued only as indicated, OKT3 5 mg/kg I.V. for 7 days as indicated.

^c Cyclosporine given I.V. initially with oral dose started at day 5 at 4 mg/kg and maintained at desired blood level, prednisone 15 mg daily decreasing gradually, azathioprine 2-2.5 mg/kg and discontinued when platelets begin to fall, methylprednisolone 500 mg intraoperatively and 15 mg I.V. for 10 days, OKT3 5 mg/kg for 10 days.

^d Nephrotoxicity or rejection.

TABLE III shows the mean cyclosporine blood values and corresponding creatinine levels in 18 kidney transplant patients from the University Hospital, Albany, New York.³² At Albany, patients were maintained at 100 to 300 ng/mL trough cyclosporine level.

TABLE III

Temporal cyclosporine whole blood values in 18 kidney transplant patients as measured by the CYCLO-Trac SP Whole Blood Kit.

Week	Post Operative ^a Day	N ^b	Cyclosporine (ng/mL)		N	Creatinine Levels (mg/dL)			
			Mean	S.D.		Mean	S.D.	MIN	MAX
1	0-7	12	151	128	143	4.9	4.0	0.7	18.0
2	8-14	16	185	96	125	2.7	2.2	1.1	8.2
3	15-21	19	173	94	94	2.7	1.9	1.1	9.4
4	22-28	17	198	146	76	3.0	1.9	1.1	8.5
5	29-35	22	210	168	70	2.9	1.7	1.2	8.2
6	36-42	17	164	63	68	2.5	1.1	1.1	5.6
7	43-49	14	175	53	55	2.4	1.1	1.1	5.4
8	50-56	17	198	99	42	2.1	1.0	1.2	5.4
9	57-63	16	176	54	40	3.1	2.7	1.0	10.9
10	64-70	4	199	59	34	3.5	3.3	0.9	11.7
11	71-77	18	146	73	27	3.0	3.1	0.9	11.1
12	78-84	5	142	76	18	3.5	3.3	1.1	11.2
13	85-91	1	304	-	7	4.2	0.3	3.6	4.6
14	92-98	1	317	-	10	3.4	1.1	1.7	4.8
15	99-105	2	121	10	9	4.1	1.6	2.0	5.5
16	106-112	1	49	-	5	5.4	2.5	1.7	8.7
17	113-119	1	79	-	3	1.6	0.1	1.6	1.7

^a Patients treated with cyclosporine 2-5 mg/kg (at twelve hour intervals), azathioprine 4 mg/kg preoperatively and tapered to 0.75 to 1.0 mg/kg, and methylprednisolone given as a 300 mg bolus and then adjusted to 20 mg/day by 30 days post transplant and 10 mg/day by 60 days post transplant.

^b Individual cyclosporine determinations.

The results of these studies demonstrate that cyclosporine blood levels were maintained at desired blood levels for each transplant type as indicated for the transplant protocols utilized by each center.

Because cyclosporine can produce different blood levels in different patients, interpretation of results and adjustment of cyclosporine dosage must take into consideration post transplant period, other medications administered and relevant clinical information concerning the patient.

DRUG INTERACTIONS

The following is reported by Sandoz LTD., Basle Switzerland for the drug Sandimmune® (cyclosporine):

Cyclosporine is extensively metabolized by the liver. Therefore, its whole blood level may be influenced by drugs that affect the hepatic microsomal enzymes, particularly the cytochrome P-450 system. Substances known to inhibit these enzymes will decrease hepatic metabolism and increase cyclosporine whole blood levels. Drugs that have been reported to increase the whole blood concentration of cyclosporine include danazol, ketoconazole, erythromycin, nicardipine, diltiazem, methyl-testosterone, and oral contraceptives. Substances that are inducers of cytochrome P-450 activity will increase hepatic metabolism and decrease cyclosporine whole blood levels.

Drugs that have been reported to decrease the whole blood concentration of cyclosporine include rifampin, phenytoin, phenobarbital, I.V. trimethoprim with sulfamethoxazole, and carbamazepine. Monitoring of blood levels of cyclosporine and appropriate dosage adjustment are essential when these drugs are used concomitantly.

It has been found that prednisolone clearance is reduced in patients treated with cyclosporine and that whole blood levels of cyclosporine, as measured by RIA, increase following the administration of high dose methylprednisolone.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Precision

Within-assay variation (values = ng/mL)

	Mean Value	S.D.	%C.V.	N
LOW	46	4.9	10.7	25
MEDIUM	186	6.0	3.2	25
HIGH	625	19.2	3.1	25

Between-assay variation results from multi-center proficiency study using whole blood samples.

	Mean Cyclosporin (ng/mL)	S.D.	%C.V.	N
Center 1*	91	15.0	16.5	5
	203	27.9	13.7	5
	543	53.7	9.9	5
Center 2	71	6.7	9.5	5
	173	12.9	7.58	5
	466	41.5	8.9	5
Center 3	78	5.4	6.8	5
	195	18.4	9.4	5
	500	50.2	10.0	5
Center 4	87	6.9	7.9	5
	200	13.5	6.8	5
	505	42.5	8.4	5
Center 5	78	3.6	4.6	5
	192	10.2	5.3	5
	497	23.9	4.8	5

* Indicates between assay variation seen in a laboratory just beginning to run radio-immunoassays.

15.2 TRUENESS: THE ASSAY TRUENESS HAS BEEN CHECKED BY THE LINEARITY TEST AND THE RECOVERY TEST.

Linearity (Parallelism)

Average recovery for 30 different serially diluted patient samples (concentration range 132-1716 ng/mL).

$$\% \text{ Recovery} = 100 \times \frac{(\text{dilution factor} \times \text{value of diluted sample})}{(\text{value of neat sample})}$$

Dilution	1/2	1/4	1/8
% Recovery	96.1	93.3	95.9
S.D.	13.7	19.2	22.6
% C.V.	14.3	20.6	23.6

Accuracy

Recovery Study: Average recoveries at 3 levels of added cyclosporine for 5 samples, 6 replicates each (n=30).

	Added Cyclosporine (ng/mL)	% Recovery	S.D.	% C.V
Level 1	95	119	17	14
Level 2	212	109	11	10
Level 3	467	97	7	7

15.3 Analytical Sensitivity (Limits of Detection)

When defined as the apparent concentration at three standard deviations from the counts at maximum binding, the mean sensitivity of the CYCLO-Trac SP Whole Blood Assay is 8.7 ng/mL (n=20).

15.4 Analytical Specificity

Average cross-reactivity of cyclosporine metabolites in the CYCLO-Trac SP Whole Blood assay. Cross-reactivity is defined as:

$$100 \times \frac{(\text{Cyclosporine concentration at 50\% B/B}_0)}{(\text{Metabolite concentration at 50\% B/B}_0)}$$

Cyclosporine Metabolite	% Cross-reactivity
1	1.7
8	<0.1
13	0.1
17	0.7
18	0.2
21	0.8

15.5 Interference

The following substances were spiked into aliquots of a whole blood pool containing cyclosporine (160-250 ng/mL). None of these substances interfered with cyclosporine measurement when present in concentrations listed below.

Substance	Concentration
Cholesterol	100-1000 mg/dL
Triglycerides	300-1000 mg/dL
Bilirubin	5-30 mg/dL
Azathioprine	50-200 ng/mL
Cimetidine	500-1000 ng/mL
Nifedipine	100-200 ng/mL
6-methyl-prednisolone	20-60 ng/mL
OKT3	1500 ng/mL

The following list of coadministered drugs have been tested at levels up to 100 µg/mL and found not to interfere with the assay:

Acetaminophen	Lidocaine
N-Acetylprocainamide	Metoclopramide
Cephalosporine	Nafcillin
Chloramphenicol	Penicillin
Digoxin	Propranolol
Disopyramide	Rifampin
Ethosuximide	Salicylic Acid
Gentamycin Sulfate	Tobramycin
Kanamycin	Verapamil

Comparison of DiaSorin CYCLO-Trac SP Whole Blood Assay to HPLC

Study	Transplant Type	Line Equation	Regression Coefficient	n
Pharmacokinetic	N/A	$Y = 0.9121X + 47$	$r = 0.9719$	373
University of Pennsylvania	Kidney	$Y = 0.8022X + 34$	$r = 0.9585$	82
	Heart	$Y = 0.9196X + 50$	$r = 0.9223$	42
	Liver	$Y = 0.9572X + 48$	$r = 0.9067$	40
Albany Medical Center	Kidney	$Y = 1.038X + 27$	$r = 0.9124$	187
University of Western Ontario	Kidney	$Y = 1.3074X + 29$	$r = 0.9703$	73
	Heart	$Y = 1.4672X + 16$	$r = 0.9696$	75
	Liver	$Y = 1.3575X - 12$	$r = 0.9067$	63

Correlation data presented represents the results of studies done in a number of laboratories using independently validated HPLC procedures. No comparison should be made to this data, and any comparison of methods should be done only after careful standardization of the methods has been achieved.

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1. Tutschka, P.J., "Cyclosporin A as a Tool for Experimental and Clinical Transplantation," **Biomedicine**, 36:341, (1982).
2. Tutschka, P.J., "Cyclosporin A-Immunobiological Aspects," **Hepatology**, 2 (5):645,(1982).
3. Britton, S. and R. Palacios, "Cyclosporin A-Usefulness, Risks and Mechanism of Action," **Immunological Review**, 65:5, (1982).
4. Keown, P.A., G.L. Essery, C.R. Stiller, N.R. Sinclair, R. Mullen and R.A. Ulan, "Mechanisms of Immunosuppression by Cyclosporin," **Transplantation Proceedings**, 13 (1):386, (1981).
5. White, D.J.G., "Cyclosporin A Clinical Pharmacology and Therapeutic Potential," **Drugs**, 24:332, (1982).
6. Merion, R.M., D.J. White, S. Thiru, D.B. Evans, and R.Y. Calne, "Cyclosporine: Five Year's Experience in Cadaveric Renal Transplantation," **New England Journal of Medicine**, 310 (3), 148-54, (1984).
7. Bieber, C., et. al., **Heart Transplantation**, (3), 81, (1983).
8. Busuttill, R.W., L.I. Goldstein, G.M. Danovich, M.E. Ament and L.D. Memsic, "Liver Transplantation Today," **Annals of Internal Medicine**, 104 (3), 337-389, (1986).
9. Cohen, D.J., "Cyclosporine: Clinical Applications and Monitoring Methods," **Laboratory Management**, pp. 29-34, (1985).
10. Sells, R.A., **Transplantation Proc.**, 15, 2495, (1983).
11. Albrechtsen, D., G. Soedal, K.J. Berg, H. Bondevik, I.B. Brekke, P. Fauchald, A. Jakobsen, B.R. Opedal, H.E. Rugstad, E. Thorsby, et. al., "Cyclosporine in Living Related and Cadaveric Renal Transplantation," **Transplantation Proc.**, 16 (5), 1911-4, (1984).
12. Atkinson, K., et. al., **Transplantation Proc.**, 15, 2761, (1983).
13. Hows, J.M., P.M. Chipping, S. Fairhead, J. Smith, A. Baughan and E.C. Gordon-Smith, "Nephrotoxicity in Bone Marrow Transplant Recipients Treated with Cyclosporin A," **Br. J. Haematol**, 54 (1), 69-78, (1983).
14. Kahan, B.D., C.A. Wideman, M. Reid, S. Gibbons, M. Jarowenko, S. Flechner and C.T. Van Buren, "The Value of Serial Serum Trough Cyclosporine Levels in Human Renal Transplantation," **Transplantation Proc.**, 16 (5), 1195-9, (1984).
15. Keown, P.A. and et. al., "In: Transplantation and Clinical Immunology," **Excerpta Medica**, Vol. 13, J.L. Toursine (Ed.), pp. 34, (1981).
16. Bowers, L.D. and D.M. Canafax, "Cyclosporine: Experience with Therapeutic Monitoring," **Therapeutic Drug Monitoring**, 6 (2), 142-7, (1984). See also D.J. Freeman, *Am. R. Cell*, **Phys. Surg. Can.**, 17, 301, (1984).
17. Giacherio, D. and T. Annesley, "Cyclosporin: Monitoring the Dosage of a New Immuno-suppressant," **MediLab**, April/May, pp. 20-25, (1985).
18. Kahan, B.D., et. al., **Transplantation Proc.**, 15, 446, (1983).
19. "National Academy of Clinical Biochemistry/American Association for Clinical Chemistry Task Force on Cyclosporine Monitoring," **Clinical Chemistry**, 33 (7), 1269-88, (1987).
20. Yatscoff, R.W., D.N. Rush and J.R. Jeffery, "Effects of Sample Preparation on Con-centrations of Cyclosporin A Measured in Plasma," **Clinical Chemistry**, 30 (11): 1812, (1984).
21. Annesly, T., K. Matz, D. Giacherio and C. Feldkamp, "Cyclosporine Distribution in Blood: Concomitant Effect of Hematocrit, Concentration, Metabolites and Temperature," **Journal of Clinical Immunoassay**, 10 (4), 226-231, (1987).

22. Annesly, T.M., D.A. Giacherio and C.S. Feldkamp, "The Concentration Related Distribution of Cyclosporine in Blood," **The Journal of Clinical Immunoassay**, 9 (1), Spring (1986).
23. Donatsch, P., E. Abisch, M. Homberger, R. Traber, M. Trapp and R. Voges, "A Radioimmunoassay to Measure Cyclosporin A in Plasma and Serum Samples," **Journal of Immunoassay**, 2 (1), 19-32, (1981).
24. Lensmeyer, G.L. and B.L. Fields, **Clinical Chemistry**, 31, 196, (1985).
25. Ball, P.E., H. Munzer, H.P. Keller, E. Abisch and J. Rosenthaler, "Specific ³H Radioimmunoassay with a Monoclonal Antibody for Monitoring Cyclosporine in Blood," **Clinical Chemistry**, 34:257-261, (1988).
26. Holt, D.W., A. Johnston, J.T. Marsden, L. Vernillet, P.A. Keown, T.G. Rosano, L.M. Shaw and J. Rosenthaler, "Monoclonal Antibodies for Radioimmunoassay of Cyclosporine: a Multicenter Comparison of Their Performance with the Sandoz Polyclonal Radioimmunoassay Kit," **Clinical Chemistry**, 34 (6):1091-6, (1988).
27. Potter, J.M. and H. Self, "Cyclosporin A: Variation in Whole Blood Levels Related to In Vitro Anticoagulant Usage," **Therapeutic Drug Monitoring**, 8 (1): 122-123, (1986).
28. Hows, J.M., and J.M. Smith, "In Vitro Stability of Cyclosporin," **Journal of Clinical Pathology**, 36:720-721, (1983).
29. Rudy, J.L. and J. Ruthledge, "Preparation and Optimal Storage of a Whole Blood Control Pool for Cyclosporin A," **Clinical Chemistry**, 32:2204-5, (1986).
30. Shaw L, "Advances in Cyclosporine Pharmacology, Measurement and Therapeutic Monitoring," **Clinical Chemistry**, 35 (7):1299-1308, (1989).
31. Moger, T., G. Post, S. Steviuff and C. Anderson, "Cyclosporine Nephrotoxicity is Minimized by Adjusting Dosage on the Basis of Drug Concentration in Blood," **Mayo Clinic Proc.**, 63:241-274, (1988).
32. Rosano T.G., B.M. Freed, M.A. Pell and N. Lempert, "Cyclosporine Metabolites in Human Blood and Renal Tissue," **Transplant Proceedings**, 18 p(6) Suppl. 5:35-40, (1986).

Sandimmunne® is a registered trademark of Sandoz LTD., Basle, Switzerland.
 CYCLO-Trac® is a registered trademark of DiaSorin Inc.

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
IVD	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
LOT	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
Ab	Antiserum	Antisérum	Antiserum	Antisero	Antisiero
Ag ¹²⁵ I	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Trazador: antígeno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
BUF	Buffer	Tampon	Puffer	Tampón	Tampone
CAL	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
CONTROL	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Radioactive	Radioactif	Radioaktiv	Radiactivo	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482	
In the United Kingdom Call: +44(0) 1344 401 430 FAX: +44(0) 7884 050812	
11598	34671 7/10

PRINTED IN U.S.A.