**GammaDab<sup>®</sup> Gastrin**  
**<sup>125</sup>I RIA Kit**

For the quantitative determination of  
gastrin levels in serum

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### 1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

The GammaDab<sup>®</sup> [<sup>125</sup>I] Gastrin Radioimmunoassay Kit is to be used for the quantitative determination of gastrin levels in serum.

### 2. SUMMARY AND EXPLANATION

Gastrin stimulates the stomach to secrete gastric acid more effectively than other gastric hormones.<sup>1</sup> Endocrine G cells produce, store and release gastrin within the pyloric and upper duodenal mucosa. Big gastrin (G-34), a linear peptide with thirty-four amino acids, and little gastrin (G-17), with seventeen amino acids, are the most prevalent forms of circulating gastrin.<sup>2</sup> Both polypeptides include one tyrosine residue, which exists either in the nonsulfated (G-17-I, G-34-I) or sulfated state (G-17-II, G-34-II). Other gastrin molecules including gastrin-13 (I and II) also circulate in the human bloodstream, and differ in biological activity, molecular size and N-terminal peptide chain length.<sup>3-6</sup> Determination of serum gastrin levels can aid in the diagnosis of various diseases, including peptic ulcer, pernicious anemia and Zollinger-Ellison (ZE) syndrome.<sup>1</sup> Gastrinomas, or tumors secreting excessive gastrin, can occur along the duodenum or among pancreatic non-beta islet cells, causing ZE syndrome.<sup>7,8</sup> Treating ZE patients with H<sub>2</sub> receptor antagonists, notably cimetidine, has achieved considerable success.<sup>9,10</sup>

ZE syndrome patients generally have serum gastrin values from 300-250,000 pg/mL.<sup>1</sup> Peptic ulcer and pernicious anemia patients may also present elevated serum gastrin values. One report has suggested that many ZE syndrome patients remain undetected and are treated simply as peptic ulcer patients.<sup>8</sup>

Three stimulation tests may aid in distinguishing ZE syndrome from other diseases. Serum gastrin levels are measured before and after:

- 2.1 a protein meal
- 2.2 calcium infusion of 12 to 20 mg/kg for 2-3 hours
- 2.3 secretin injection of 1-3 units

TABLE IV in Expected Values section summarizes typical responses.

### 3. PRINCIPLES OF THE ASSAY

The GammaDab [<sup>125</sup>I] Gastrin Radioimmunoassay Kit, based on competitive binding principles of radioimmunoassay, utilizes a precipitating antiserum reagent to separate antibody-bound tracer from unbound tracer. Nonradioactive gastrin from serum samples, gastrin calibrators or controls competes with a constant amount of tracer for binding sites on gastrin antibody, which is held at a limiting concentration. The amount of tracer which will bind to the antibody is inversely proportional to the amount of nonradioactive gastrin present in the assay tube.

The precipitating antiserum reagent containing second antibody in polymer solution, is used to separate the antibody-bound tracer from unbound tracer by immunoprecipitation. The assay tubes are centrifuged and supernatants decanted. The antibody-bound tracer, which is in the precipitate, is counted in a gamma counter. A calibrator curve is constructed and the gastrin concentration in the samples is interpolated from the calibrator curve.

### 4. REAGENTS PROVIDED IN THE KIT

Gastrin Antiserum	1 vial / 13 mL
Gastrin Calibrators	6 vials / 1 mL
Gastrin Control	1 vial / 1 mL
Gastrin Precipitating Reagent	1 vial / 130 mL
Gastrin Tracer	1 vial / 13 mL
Number of tests	125

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent at 2-8, or as mentioned below, until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

Reagents from different batches must not be mixed

**4.1 [125I] Gastrin Tracer: ready to use reagent**

Each vial contains approximately 2 µCi tracer (<1 µg/mL Gastrin) in 13 mL of tris buffer and human serum albumin with 0.1% sodium azide as a preservative.

**4.2 Rabbit Anti-Gastrin Serum: ready to use reagent**

Each vial contains 13 mL of rabbit anti-gastrin serum (<1 µg/mL) in normal rabbit serum and tris buffer with 0.1% sodium azide as a preservative.

**4.3 Gastrin Precipitating Antiserum Reagent: ready to use reagent**

Each bottle contains 130 mL of sheep anti-rabbit serum (<20%) in a water soluble polymer and tris buffer with 0.1% sodium azide as a preservative.

**4.4 Gastrin Control: ready to use reagent**

Each vial contains 1 mL of gastrin in tris buffer and human serum albumin with 0.1% sodium azide as a preservative. The nominal value is 300 ± 60 pg/mL.

**4.5 Gastrin Calibrators: ready to use reagent**

Each vial contains 1 mL of gastrin in tris buffer and human serum albumin with 0.1% sodium azide as a preservative. Calibrators are provided at the following concentrations: 0, 40, 100, 200, 500 & 1000 pg/mL, respectively. Store at 2-8°C for up to one month. For longer storage, the calibrators should be kept frozen. The DiaSorin Gastrin Calibrators are calibrated using the current production master lot. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

**5. WARNINGS AND PRECAUTIONS**

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

**REAGENTS CONTAINING HUMAN SOURCE MATERIAL**

**Treat as potentially infectious.**

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus, hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4<sup>th</sup> ed., May 1999 or current edition.

**REAGENTS CONTAINING SODIUM AZIDE**

**CAUTION:** Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, U.S. 1976.

**European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

## REAGENTS CONTAINING IODINE-125

Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation there from, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

## 6. COLLECTION AND PREPARATION OF SERUM

Serum samples should be drawn from patients after overnight fasting, preferably twelve hours or more, since circulating gastrin levels rise after eating.

Collect a venous blood sample aseptically using an evacuated glass tube without additives. Within one hour of collection, separate the serum fraction after centrifugation for use in the assay. Lipemic samples may interfere with assay results.

Serum samples may be stored up to eight hours at 2-8°C, or below -20°C up to two weeks.<sup>12</sup> Avoid freeze-thaw cycles and self-defrosting freezers. Immediately prior to assay, thaw samples at room temperature and mix well before pipetting. Do not thaw samples in a 37°C water bath.

## 7. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 7.1 Test tube rack.
- 7.2 Polypropylene test tubes (12 x 75 mm).
- 7.3 Precision pipet (100 µL).
- 7.4 Semi-automatic pipet (1.0 mL).
- 7.5 Vortex mixer.
- 7.6 Magnetic mixer and spin bar (optional).
- 7.7 Centrifuge (2-8°C).
- 7.8 Waste beaker.
- 7.9 Plastic-backed absorbent paper.
- 7.10 Semi-logarithmic graph paper (available upon request).

## 8. ASSAY PROCEDURE

- 8.1 Allow all reagents to reach ambient temperature and mix well without foaming prior to use.
- 8.2 Label polypropylene test tubes in duplicate according to the following scheme. Total Counts [T<sub>1</sub>,T<sub>2</sub>] and B<sub>0</sub> tubes [1,2] may be required for certain data reduction and quality control programs but may be omitted if the calibrator curve is plotted on semi-logarithmic graph paper.

Tube No.	Contents of Tubes		Code Letter
T <sub>1</sub> ,T <sub>2</sub>	Total Counts (Tracer)		
1,2	Gastrin Blank, B <sub>0</sub>	0 pg/mL	A
3,4	Gastrin Calibrator	40 pg/mL	B
5,6	Gastrin Calibrator	100 pg/mL	C
7,8	Gastrin Calibrator	200 pg/mL	D
9,10	Gastrin Calibrator	500 pg/mL	E
11,12	Gastrin Calibrator	1000 pg/mL	F
13,14	Gastrin Control		
15,16	Patient "X" Serum Sample		
17,18	Patient "Y" Serum Sample		

- 8.3 Add to the appropriate duplicate tubes:
- 100 µL of each gastrin blank and calibrator.
  - 100 µL of gastrin control.
  - 100 µL of each patient sample.
- 8.4 Add 100 µL of gastrin tracer to all tubes including Total Counts. Cap the Total Count tubes and set aside.
- 8.5 Add 100 µL of rabbit anti-gastrin serum to all tubes except Total Counts. Mix by vortexing at low speed.
- 8.6 Incubate tubes sixty minutes at room temperature.
- 8.7 Add 1.0 mL of precipitating antiserum reagent to all tubes except Total Counts. Mix by vortexing at low speed.  
TO INSURE A HOMOGENEOUS MIXTURE, STIR THE PRECIPITATING ANTI-SERUM REAGENT CONTINUOUSLY WITH A MAGNETIC MIXER AND SPIN BAR OR HAND SWIRL EVERY 10-15 PIPETTINGS.
- 8.8 Incubate tubes ten minutes at room temperature.
- 8.9 Centrifuge all tubes except Total Counts at 2-8°C for ten minutes at a minimum relative centrifugal force (RCF, TABLE I) of 1500 x g.
- 8.10 Decant each tube except Total Counts into a waste beaker. Tap the rim of each tube vigorously onto absorbent paper to remove any residual supernatant.  
FAILURE TO REMOVE ADHERING SOLUTION ADEQUATELY FOLLOWING DECONTATION MAY RESULT IN POOR REPLICATION AND SPURIOUS VALUES.
- 8.11 Count all tubes, including Total Counts, for one minute in a gamma counter with the window suitably adjusted for iodine-125.
- 8.12 Calculate results. Refer to Results section.

**TABLE I**  
**Computing Relative Centrifugal Force (RCF)\***

Speed of Centrifuge (rpm)	RADIUS (inches)								
	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0
3000	767	895	1022	1150	1278	1406	1534	1661	1789
3500	1044	1218	1392	1566	1740	1914	2087	2261	2435

\*RCF VALUE (x g) =  $R \times N^2 \times 284 \times 10^{-7}$

R = radius in inches

N = speed of centrifuge measured in revolutions per minute (rpm)

## 9. QUALITY CONTROL

Each laboratory should utilize controls at several levels to monitor assay performance. The controls should be treated as unknowns. Quality control charts should be maintained to follow the performance of the controls. Appropriate statistical methods should be used to evaluate trends. Acceptable performance limits should be ascertained for the individual laboratory.

## 10. CALCULATION OF RESULTS

- 10.1 Record the counts per minute (CPM) bound for each tube.
- 10.2 Plot the CPM bound for the gastrin calibrators (vertical axis) versus gastrin concentration (horizontal axis) on semi-logarithmic graph paper.
- 10.3 Draw the best fitting curve. Refer to FIGURE 1 and TABLE III for a typical calibrator curve and data.
- 10.4 Locate the CPM bound for each tube which corresponds to each control and patient sample on the vertical axis and follow a horizontal line intersecting the calibrator curve. At point of intersection, read gastrin concentration (pg/mL) from the horizontal axis.

## 11. S.I. UNITS

Gastrin serum calibrators, control and patient values can be converted directly from picograms per milliliter (pg/mL) to nanograms per liter (ng/L), of S.I. units. Refer to TABLE II for values of gastrin calibrators listed in both units.

**TABLE II**  
**Concentration of Calibrators**


Cat. No.	Calibrator Levels	S.I. Units
CA-2946	0 pg/mL	0 ng/L
CA-2947	40 pg/mL	40 ng/L
CA-2948	100 pg/mL	100 ng/L
CA-2949	200 pg/mL	200 ng/L
CA-2950	500 pg/mL	500 ng/L
CA-2951	1000 pg/mL	1000 ng/L

**TABLE III**

**Recording the Data  
Do not use to calculate unknowns**

Tube No.	Contents of Tubes	CPM Bound	Gastrin Concentration (pg/mL)
T1	Total Counts (Tracer)	25320	—
T2	Total Counts (Tracer)	25005	—
1	Gastrin Blank, B <sub>0</sub> 0 pg/mL	17121	—
2	Gastrin Blank, B <sub>0</sub> 0 pg/mL	17129	—
3	Gastrin Calibrator      40 pg/mL	15125	—
4	Gastrin Calibrator      40 pg/mL	15287	—
5	Gastrin Calibrator      100 pg/mL	12933	—
6	Gastrin Calibrator      100 pg/mL	12692	—
7	Gastrin Calibrator      200 pg/mL	10056	—
8	Gastrin Calibrator      200 pg/mL	10282	—
9	Gastrin Calibrator      500 pg/mL	6554	—
10	Gastrin Calibrator      500 pg/mL	6649	—
11	Gastrin Calibrator      1000 pg/mL	4387	—
12	Gastrin Calibrator      1000 pg/mL	4353	—
13	Gastrin Control	8510	304
14	Gastrin Control	8492	305
			Av. 305
15	Patient "X" Serum Sample	14142	64
16	Patient "X" Serum Sample	13989	67
			Av. 65
17	Patient "Y" Serum Sample	9879	217
18	Patient "Y" Serum Sample	9864	218
			Av. 217

**FIGURE 1**

**12. LIMITATIONS OF THE PROCEDURE**

- 12.1 Lipemic samples may interfere with assay results.
- 12.2 Remove any visible particles by centrifugation before assaying.
- 12.3 The user should note that erroneous results may be caused by improperly storing or handling patient samples.
- 12.4 Failure to obtain appropriate gastrin values for controls may indicate imprecise manipulations, improper handling or deterioration of reagents.

### 13. EXPECTED VALUES

When fasting serum gastrin levels are elevated, the stimulation tests outlined in TABLE IV can aid in differentiating between patients with Zollinger-Ellison syndrome and those with duodenal ulcer.<sup>1,11</sup>

**TABLE IV**  
**Stimulation Tests for Zollinger-Ellison Syndrome**  
**Effect on Serum Gastrin Level**

	High Protein Meal	Intravenous Infusion of Calcium	Secretin Injection
<b>Normal</b>	large increase 2-3 times	small increase	small decrease
<b>Zollinger-Ellison Syndrome</b>	no change	large increase 2-3 times	large increase 3 times
<b>Duodenal Ulcer</b>	large increase 5-10 times	small increase	small decrease

**Reference range:**

Less than 108 pg/mL

[Less than 108 ng/L]

This range may be used as a guideline until each laboratory has established its own reference range. This reference range has been established based upon a study of 12-hour fasting specimens in a population of 99 normal patients.

### 14. SPECIFIC PERFORMANCE CHARACTERISTICS

#### 14.1 Analytical Sensitivity

The sensitivity of the calibrator curve is defined as the smallest single value which can be distinguished from zero. This value was calculated from the 95% confidence limits (2 standard deviations) for thirty replicates at the zero point of the calibrator curve. The calculated sensitivity is 6.0 pg/mL.

#### 14.2 Precision

The intra-run precision was determined from the mean of twenty simultaneous assays per sample. The inter-run precision was determined from the mean of the average of duplicates for twenty separate runs.

Precision Sample	Number of Assays	Standard Mean (pg/mL)	Intra-Run Deviation (pg/mL)	Coefficient of Variation (%)
Serum Pool A	20	64	2.6	4.1
Serum Pool B	20	239	7.9	3.3
Serum Pool C	20	526	15.8	3.0

Inter-Run Precision Sample	Number of Separate Runs	Mean (pg/mL)	Standard Deviation (pg/mL)	Coefficient of Variation (%)
Serum Pool A	20	64	3.1	4.8
Serum Pool B	20	236	8.3	3.5
Serum Pool C	20	484	30.6	6.3

**14.3 Trueness: The assay trueness has been checked by the recovery assay.**

**Recovery from Pooled Sera**

Spiked samples were prepared by adding aliquots of a known high gastrin stock solution to aliquots of a previously analyzed serum pool. Each spiked sample was assayed in quadruplicate. Percent recovery is calculated as (recovered/expected) x 100.

Contribution from Spike (pg/mL)	Contribution from Pool (pg/mL)	Expected Value (pg/mL)	Recovered Value (pg/mL)	Percent (%) Recovery
837	36	873	934	107
605	42	647	688	106
324	50	374	383	102
216	53	269	267	99
53	57	110	102	93

**14.4 Avidity**

The calculated affinity constant of the gastrin antiserum is approximately  $9 \times 10^{10}$  liters/mole.

**14.5 Analytical Specificity**

Data on the cross-reactivity of the antiserum used in this kit are shown expressed as the ratio of gastrin concentration to the cross-reacting compound concentration at 50% inhibition of maximum binding.

Compound	% Cross-reactivity
Gastrin G-17-I	100
Gastrin G-17-II	64
Gastrin G-34-I	22
Gastrin G-13-I	29
CCK Peptide	0.3
Secretin	0.3

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






1. RIA for Physicians: An Aid in Differential Diagnosis, Vol. 1, No. 4, ed. Travis, **Scientific Newsletters**, Inc., Anaheim, CA, (1976).
2. Bentley, P.H., G.W. Kenner and R.C. Sheppard, "Human Gastrin: Isolation, Structure and Synthesis," **Nature**, 209:583, (1966).
3. Rehfeld, J.F., et. al., "Gastrin Heterogeneity in Serum and Tissue: A Progress Report," **Gastrointestinal Hormones**, Thompson (ed.), U. Texas Press, Austin, TX, 43, (1975).
4. Walsh, J.H. and M.I. Grossman, "Medical Progress: Gastrin (First of Two Parts)," **N. Eng. J. Med.**, 292:1324 & 1377, (1975).
5. Yalow, R.S., "Heterogeneity of Peptide Hormones with Relation to Gastrin," **Gastrointestinal Hormones**, Thompson (ed.), U. Texas Press, Austin, TX, 25, (1975).
6. Rehfeld, J.F., "REVIEW: Gastrins in Serum, A review of Gastrin Radioimmunoanalysis and the Discovery of Gastrin Heterogeneity in Serum," **Scand. J. Gastroent.**, 8:577, (1973).
7. McGuigan, J.E. and W.L. Trudeau, "Immunochemical Measurement of Elevated Levels of Gastrin in the Serum of Patients with Pancreatic Tumors of the Zollinger-Ellison Variety," **N. Eng. J. Med.**, 278:1308, (1968).
8. Isenberg, J.I., et. al., "Progress in Gastroenterology: Zollinger-Ellison Syndrome," **Gastroenterology**, 65:140, (1973).
9. Siepler, J.K., et. al., "H<sub>2</sub> Receptor Antagonists," **J. Hosp. Pharm.**, 35:141, (1978).
10. Stage, J.G., F. Stadil and K. Fischerman, "New Aspects in the Treatment of the Zollinger-Ellison Syndrome," **Cimetidine**, Crentzfeldt (ed.), Excerpta Medica, Amsterdam & Oxford, 137, (1978).
11. Stadil, F., "Chapter 30: Gastrinomas," **Gastrointestinal Hormones**, Jerzy Glass (ed.), Raven Press, N.Y., 729, (1980).
12. Rayford, P.L., K. Hejtmancik and J.C. Thompson, "Radioimmunoassay: Gastrointestinal Hormones," **World J. Surg.**, 3, 423, (1979).

### SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Antiserum	Antisérum	Antiserum	Antisuero	Antisiero
	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reactivo precipitante	Reagente precipitante
	Tracer: antigen labelled with <sup>125</sup> I	Traceur : antigène marqué à l'iode <sup>125</sup>	Tracer: <sup>125</sup> I-markiertes Antigen	Trazador: antigéno etiquetado con <sup>125</sup> I	Tracciatore: antigene etichettato con <sup>125</sup> I
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Buffer	Tampon	Puffer	Tampón	Tampone
	Radioactive	Radioactif	Radioaktiv	Radiactivo	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,  
13040 Saluggia (VC) Italy

Phone: 651-439-9710  
FAX: 651-351-5669

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482	
In the United Kingdom Call: 44 118 9364200 FAX: 44 118 9792061	
12752	27678 12/03

PRINTED IN U.S.A.

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DiaSorin Inc.