---

**HGH-CTK/ irma**  
(P2244)

---

DiaSorin

---

---

English.....	p. 1
Italiano.....	p. 8
Français.....	p. 16
Deutsch.....	p. 25
Español.....	p. 33
Português.....	p. 41
Magyar.....	p. 49
Ελληνικά.....	p. 57

## **hGH IMMUNORADIOMETRIC ASSAY KIT**

### **Procedure for quantitative determination of human growth hormone (hGH) in human serum or plasma samples**

For in vitro use only

#### **1. INTRODUCTION**

Human growth hormone (hGH) is the hormone most abundantly secreted by the pituitary. During childhood and adolescence it causes body statural growth, and throughout life it affects major metabolic processes, in such a way that growth is encouraged. It is not, however, essential to health in adult life.

hGH is synthesized by the alpha cells of the anterior pituitary. It is composed of a single chain of 191 amino acid residues with a molecular weight of 21,500 daltons and two disulphide bonds.

hGH circulates in different molecular forms, called Big-Big GH, Big GH and Little GH. Little GH is the most abundant form and probably the original monomeric form, while the other two seem to be polymeric forms.

hGH does not have a specific target organ, but receptors have been demonstrated on human liver cells and circulating lymphocytes. The metabolic effects of hGH can be generally divided into two main categories: *a)* anti-insulin direct action (lipolytic and causing hyperglycaemia); *b)* insulin-like indirect action (anti-lipolytic and mitogenic). The latter action causes body statural growth which is the classic action of hGH. This is exerted through somatomedins, relatively low-molecular weight hormones (M.W. 7,000 daltons) of chemical structure very similar to that of proinsulin.

At least two somatomedins have been identified, one with insulin-like properties, the other one more directly involved in body statural growth. Unlike all other polypeptide hormones, somatomedins mostly circulate bound to carrier proteins which prolong their half-lives, therefore making their concentrations relatively stable.

#### **2. PRINCIPLE OF THE ASSAY**

The assay is a non-competitive immunoradiometric (IRMA) method, based on the use of antibody-coated tubes, employing two mouse monoclonal antibodies directed against two different epitopes on the hGH molecule. Monoclonal IgG to hGH is coated on the tube walls and labelled monoclonal IgG to hGH is used as tracer.

The assay features a single incubation during which hGH contained in calibrators or samples allows the tracer to bind to the solid phase. After the incubation, the amount of labelled antibody bound to the solid phase is proportional to the concentration of hGH present in calibrators or samples. At the end of incubation, the unbound material is removed by aspiration and washing.

#### **3. REAGENTS PROVIDED IN THE KIT**

Coated tubes	100
<sup>125</sup> I-tracer	1 vial
hGH calibrators	7 vials
Control serum	1 vial
Wash buffer	2 bottles
Number of tubes	100

*STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.*

*Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.*

*When reconstituting the contents of the vials, mix gently to avoid foaming.*

*Reagents from different batches must not be mixed.*

##### **3.1. Coated tubes**

The inner surface of each tube is coated with biotinylated IgG to hGH (mouse monoclonal) directed against an hGH epitope that is different from the one against which the IgG used for the tracer is raised.

*Before use, bring the coated tubes to room temperature prior to opening the box, to avoid condensation of humidity. Securely reseal the box containing unused tubes. Do not mix different batches of coated tubes.*

**3.2. <sup>125</sup>I-tracer (red):** ready-to-use reagent

The vial contains 11 mL IgG to hGH (mouse monoclonal) labelled with <sup>125</sup>I, BSA, horse serum, mouse IgG, phosphate buffer, preservatives and an inert red dye. Radioactivity is 333 kBq (9 µCi) or less on the calibration date.

**3.3. hGH calibrators:** ready-to-use reagent

The vials contain increasing amounts of hGH, human serum, phosphate buffer and preservatives. hGH concentrations are as follows: 0 - 1.5 - 4.5 - 9 - 24 - 60 - 150 µIU/mL (0 - 0.5 - 1.5 - 3 - 8 - 20 - 50 ng/mL). The zero calibrator vial contains 3 mL; the 1-6 calibrator vials contain 0.5 mL. *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

The kit calibrators are referenced to the 2nd International Standard NIBSC 98/574, Recombinant DNA-derived hGH (1 ng = 3 µIU).

**3.4. Control serum:** lyophilized reagent

The vial contains hGH, human serum and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. Reconstitute the vial contents with 1 mL distilled water. Store the resulting solution for two weeks at 2-8°C or in deep-frozen aliquots (-20°C or below) for extended storage.

**3.5. Wash buffer:** reagent in solution (10x)

Each bottle contains 50 mL 0.5% Triton X-100 and saline solution.

Dilute the contents of each bottle to 500 mL with deionized water. The resulting solution is stable at 2-8°C until the kit expiry date. The reagent is used to rinse coated tubes.

**4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED**

- Distilled or deionized water.
- Glassware.
- Micropipettes with disposable tips (50, 100, 1000 µL) (50 µL: trueness ± 3%, precision 2%; 100, 1000 µL: trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Horizontal shaker capable of achieving a shaking speed of 300-350 rpm (optional).
- Device for dispensing and aspiration of wash buffer capable of delivering 2-3 mL per wash cycle for two wash cycles.
- Gamma counter suitable for counting <sup>125</sup>I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

**5. SPECIMEN COLLECTION AND PREPARATION**

Either human serum or plasma may be used. The anticoagulants citrate, EDTA and heparin have been tested and may be used with this assay. Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

*If hGH levels greater than 150 µIU/mL (50 ng/mL) are expected, the samples should be diluted with the zero calibrator.*

**6. ASSAY PROCEDURE**

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time.

Perform all assay steps in the order given and without any appreciable delays between the steps. The reagents must be dispensed as quickly as possible (max 30 min).

A clean, disposable tip should be used for dispensing each calibrator and sample.

- Dispense reagents in the bottom of coated tubes. Operate according to the following scheme:

	tubes	Calibrators 0-6	Samples
reagents			
Calibrators		50 µL	–
Samples		–	50 µL
Tracer		100 µL	100 µL

- Prepare two non-coated tubes for total activity computation containing only 100 µL tracer and set them aside until counting.
- **Mix** the contents of tubes with a Vortex and **incubate** according to the procedure selected:
  - **90 min at room temperature**, while continuously shaking (300-350 rpm) (procedure A)
  - **overnight (18-22 hours) at room temperature** (procedure B).
- Carefully **aspirate** the incubation mixture and **rinse twice** with 2 mL **wash buffer**.  
*Be sure that the aspirator tip touches the bottom of the coated tube so that all the liquid is removed. Failure to remove adhering solution adequately may result in poor reproducibility and spurious results. No trace of red dye should still be visible.*
- **Measure the radioactivity** of tubes.

## 7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Compute the B/T ratio for each calibrator and unknown sample as follows:

$$B/T\% = \frac{\text{calibrator or sample mean counts}}{\text{total activity mean counts}} \times 100$$

Plot in log-log coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of hGH concentration expressed as µIU/mL or ng/mL on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1). Directly from the calibration curve, read the hGH concentration of each sample expressed as µIU/mL or ng/mL. If the sample was diluted, the hGH concentration value derived from the diluted specimen must be multiplied by the dilution factor.

### Calculation example (procedure A)

*The following data must only be considered an example and should not be employed instead of the data obtained by the user.*

Description	cpm	B/T%
Total activity	112,199	100
Zero calibrator	23	0.02
1.5 µIU/mL - 0.5 ng/mL	891	0.79
4.5 µIU/mL - 1.5 ng/mL	2,861	2.55
9 µIU/mL - 3 ng/mL	5,598	4.99
24 µIU/mL - 8 ng/mL	13,308	11.86
60 µIU/mL - 20 ng/mL	29,546	26.33
150 µIU/mL - 50 ng/mL	59,715	53.22
Sample	3,636	3.24

By interpolation from the calibration curve the sample is found to contain 5.76 µIU/mL (1.92 ng/mL) hGH.


Fig. 1

## 8. EXPECTED VALUES

Given the variability of hGH secretion within a single individual and between different individuals, as reported in the literature, the basal test is considered of limited diagnostic significance and assessments based on the response of hGH to challenge tests are to be preferred.

It should be remembered, however, that the basal reference values reported below are to be considered merely indicative. Each laboratory should establish its own assessment methods and reference ranges.

Subjects	hGH, $\mu$ IU/mL	hGH, ng/mL
Prepubertal age	0.0 - 4.0	0.0 - 1.3
Pubertal age: males	0.0 - 10.0	0.0 - 3.3
females	0.0 - 30.0	0.0 - 10.0
Adult age: males	0.0 - 14.0	0.0 - 4.7
females	0.0 - 32.0	0.0 - 10.7

## 9. SPECIFIC PERFORMANCE CHARACTERISTICS

### 9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

**Interference.** Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by anticoagulants (citrate, EDTA, heparin), haemolysis (up to 200 mg/dL haemoglobin), lipaemia (up to 500 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin) or one freeze-thaw cycle of samples.

**Cross-reactions.** The kit recognizes the two monomeric variant forms of hGH (M.W. 20,000 and 22,000 daltons). Specificity was evaluated as the apparent hGH concentration obtained when testing high concentrations of potentially interfering hormones.

Hormone	Concentration	Apparent hGH level
hPRL	100 mIU/mL 3rd IS 84/500	3 $\mu$ IU/mL (1 ng/mL)
hCG	100 IU/mL 3rd IS 75/537	undetectable
hTSH	1 IU/mL 2nd IRP 80/558	undetectable
hLH	2 IU/mL 1st IRP 68/40	undetectable
hFSH	2 IU/mL 2nd IRP 78/549	undetectable
hCS	3 $\mu$ g/mL	undetectable

### 9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 0.12  $\mu$ IU/mL (0.04 ng/mL) at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator, that is, two standard deviations above zero.

### 9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	10	10	10
Mean ( $\mu$ IU/mL)	36.80	16.40	3.42
Standard deviation	0.68	0.30	0.12
Coefficient of variation (%)	1.9	1.9	3.9

Reproducibility	A	B	C
Number of determinations	10	10	10
Mean ( $\mu$ IU/mL)	34.28	15.82	4.02
Standard deviation	0.86	0.58	0.16
Coefficient of variation (%)	2.5	3.7	4.1

### 9.4. Trueness

The assay trueness has been checked by the dilution and recovery tests.

**Dilution test.** Two sera with high hGH concentration were tested after serially diluting with the zero calibrator.

Dilution	Expected concentration, $\mu$ IU/mL	Measured concentration, $\mu$ IU/mL	% Recovery
neat	–	77.40	–
1:2	38.70	37.40	96.6
1:4	19.35	18.60	96.1
1:8	9.68	9.40	97.2
1:16	4.84	4.80	99.2
1:32	2.42	2.40	99.2
neat	–	44.40	–
1:2	22.20	23.00	103.6
1:4	11.10	10.80	97.3
1:8	5.55	5.40	97.3
1:16	2.78	2.80	100.9
1:32	1.39	1.40	100.9

**Recovery test.** Two sera containing hGH were tested as such and after mixing with increasing amounts of hGH.

Added concentration, $\mu\text{IU/mL}$	Expected concentration, $\mu\text{IU/mL}$	Measured concentration, $\mu\text{IU/mL}$	% Recovery
–	–	33.40	–
50.00	81.74	78.80	96.4
25.00	56.74	54.40	95.9
12.50	44.24	44.60	100.8
6.50	38.24	39.80	104.1
2.10	33.84	34.60	102.3
–	–	13.00	–
50.00	62.36	62.50	100.2
25.00	37.36	38.40	102.8
12.50	24.86	25.00	100.6
6.50	18.86	19.20	101.9
2.10	14.46	14.60	101.0

### 9.5. High-dose hook effect

Whenever samples containing extremely high antigen concentrations are tested in a one- step *sandwich* method, the hook effect can mimic concentrations lower than real. For correct quantification, samples containing hGH levels greater than that of the top calibrator should be diluted with the zero calibrator and retested. The results must then be multiplied by the dilution factor to obtain the levels of the neat specimens.

The kit has been designed in such a way that doses up to 7.5 mIU/mL (2.50  $\mu\text{g/mL}$ ) hGH produce an analytical signal which is above that of the top calibrator.

### 10. LIMITATIONS OF THE PROCEDURE

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results. A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration and washing are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration of incubation mixture and rinsing of tubes
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

### 11. WARNINGS AND PRECAUTIONS

Test components contain sodium azide as a preservative. Because sodium azide may form explosive lead or copper azide in plumbing, it is recommended that drains be thoroughly flushed with water after disposal of solutions containing sodium azide (Council Directive 99/45/EC).

R 22 – Harmful if swallowed.

R 31 – Contact with acids liberates toxic gas.

S 28 – After contact with skin, wash immediately with plenty of water.

S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

## 12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an overkill approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

## 13. BASIC RULES OF RADIATION SAFETY

### Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 8.85  $\mu\text{Ci}$  (328 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

## SCHEME OF THE ASSAY

- 1 - RECONSTITUTE THE CONTROL SERUM.
- 2 - IDENTIFY COATED TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS ACCORDING TO THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

	TUBES	CAL 0-6	SAMPLES
REAGENTS			
CALIBRATORS		50 $\mu\text{L}$	–
SAMPLES		–	50 $\mu\text{L}$
TRACER		100 $\mu\text{L}$	100 $\mu\text{L}$

- 4 - INCUBATE FOR:
  - . 90 MIN AT ROOM TEMPERATURE WHILE SHAKING (Procedure A)
  - . OVERNIGHT AT ROOM TEMPERATURE (Procedure B).
- 5 - ASPIRATE THE INCUBATION MIXTURE AND RINSE TWICE WITH 2 mL WASH BUFFER.
- 6 - MEASURE THE RADIOACTIVITY OF TUBES.

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- R.L. ABRAMS et al.  
Hypothalamic regulation of growth hormone secretion.  
J. Clin. Invest., **43** : 1242 (1964).
- G. BANFI, M. MARINELLI, M. PONTILLO, P. BONINI  
Standardisation with synthetic 22-kDa monomer human growth hormone reduces discrepancies between two monoclonal immunoradiometric assay kits.  
Clin. Chem., **38** (10) : 2107-2110 (1992).
- R. BOUILLON, P. DE MOOR  
Heterogeneity of human growth hormone in serum.  
Ann. Endocr., **35** : 606 (1974).
- A.C. CELNIKER, A.B. CHEN, R.M. WERT Jr., B.M. SHERMAN  
Variability in the quantitation of circulating growth hormone using commercial immunoassays.  
J. Clin. Endocrinol. Metab., **68** : 469-475 (1989).
- W.H. DAUGHADAY  
The adenohipophys. In: «Textbook of endocrinology», R.H. Williams ed., W.B. Saunders, Philadelphia, p. 73 (1981).
- E. GHIGO et al.  
Growth hormone-releasing activity of hexarelin, a new synthetic hexapeptide, after intravenous, subcutaneous, intranasal, and oral administration in man.  
J. Clin. Endocrinol. Metab., **78** : 693-698 (1994).
- W.M. HUNTER  
Growth hormone. In: «Hormone assays and their clinical application», J.A. Loraine, E.T. Bell eds., Churchill Livingstone, Edinburgh, p. 221 (1976).
- U.J. LEWIS et al.  
Human growth hormone: a complex of proteins.  
Rec. Progr. Horm. Res., **36** : 477 (1980).
- C.H. LI  
Human growth hormone: perspectives on its chemistry and physiology.  
In: «Advances in human growth hormone research», S. Raiti ed., U.S. Department of Health, Education and Welfare Publ., DHEW publ. no. (NIH) 74-612, p. 321 (1974).
- G.M. MOLINATTI et al.  
Radioimmunoassay of human growth hormone.  
J. Nucl. Biol. Med., **13** : 26 (1968).
- G.T. PEAKE, J. MORRIS, M.T. BUCKMAN  
Growth hormone. In: «Methods of hormone radioimmunoassay», B.M. Jaffe, H.R. Behrman eds., Academic Press, New York, p. 223 (1979).
- E. SCHAFF-BASS et al.  
Advances in diagnosis and treatment of short stature with special reference to the role of growth hormone.  
J. Pediatr., **104** : 801 (1984).
- B.T. TRUDD  
Growth, growth hormone and the somatomedins: a historical perspective and current concepts.  
Ann. Clin. Biochem., **28** : 542-555 (1991).

DiaSorin Inc.  
1951 Northwestern Avenue  
P.O. Box 285  
Stillwater, MN 55082-0285 U.S.A.  
Phone: 651-439-9710  
FAX: 651-351-5669

APM0133

DiaSorin S.p.A.,  
13040 Saluggia (VC) Italy

34663 7/10