TRYPsik
(P2573)

DiaSorin

English	p. 1
Italiano.....	p. 10
Français	p. 19
Deutsch	p. 28
Español	p. 37

TRYPSIN RADIOIMMUNOASSAY KIT
Procedure for quantitative determination of
trypsin-like immunoreactivity (TLI)
in human serum or plasma samples

1. INTRODUCTION

Assayable species

The molecular species actually measured in a radioimmunoassay system employing trypsin as immunogen, calibrator and material to be labelled was initially believed to be trypsin. Further experiments brought about the eventual identification of the assayable species as cathodic trypsinogen, i.e. the inactive precursor of the enzyme involved in RIA setting-up. It must be noted that circulating trypsin, if any, should escape RIA detection due to formation of non-immunoreactive complexes with serum inhibitors (i.e. α_2 -macroglobulin).

The expression of analytical data as trypsin equivalent concentration (i.e. ng/mL trypsin or more precisely cathodic trypsin) has been used to date by all the researchers concerned, including the authors who first recognized the analyte as trypsinogen. This way of data expressing is maintained as being fully acceptable in terms of analytical self-consistency. Fig. 1 shows the linear correlation of the responses obtained using cathodic trypsinogen and trypsin as calibrator, proving the feasibility of a conversion between the two unit systems. The species measured is usually referred to as trypsin-like immunoreactivity (TLI).

The use of trypsin RIA to assay trypsinogen is justified by the fact that the enzyme can only be effectively purified by Trasylol affinity-chromatography.


Fig. 1: Check for the immunological equivalence between cathodic trypsin and trypsinogen. The data refer to the mean \pm S.D. recorded in 5 independent experiments.

Biochemical features

Trypsinogen is a zymogen, which is synthesized in the exocrine acinar cells of the pancreas and accumulated in the nuclei of the pancreatic cells. Two forms of trypsinogen exist, giving origin to two distinct trypsin isoenzymes. These species are referred to as *cathodic* and *anodic* trypsinogen and trypsin, according to the respective electrophoretic behaviour.

Trypsinogen has a molecular weight of 24,000 daltons and it is formed by a single chain of 229 amino acid residues. When zymogens are secreted into the gastrointestinal tract, they are converted into their active forms by the selective hydrolytic cleavage of one or more specific peptide bonds in the zymogen molecule. Trypsinogen is converted into active trypsin in the intestinal lumen by removal of a hexapeptide from the N-terminal end through the action of the proteolytic enzyme enterokinase.

Clinical relevance

The clinical relevance of a serum TLI assay lies in the fact that, being synthesized only by the pancreatic acinar cells, it is considered a specific indicator of pancreatic damage, unlike enzymes alpha-amylase and lipase usually employed in pancreatic diagnosis. This remark is supported by the fact that totally pancreatectomized patients present zero TLI levels, but alpha-amylase values within the normal range.

As it is impossible to measure the enzymatic activity of TLI in serum, because of the presence of many inhibitors, the availability of a reliable RIA system for the determination of TLI in serum is of great interest.

In **acute pancreatitis**, TLI levels reach a maximum value of five to ten times higher than the upper normal limit and remain elevated for four to five days after the onset of abdominal pain. No overlap of these elevated levels with normal levels was ever observed, provided that patients affected by severe chronic renal failure (decreased TLI excretion) are not considered. Therefore, when interpreting high levels, the renal status of the patients should be evaluated.

Up to now no serum test has been available to diagnose **chronic pancreatitis**. As pancreatic acini atrophy occurs in this disease, about 50% of patients show a trypsin level below the lower normal limit.

The TLI assay has been proposed as a screening test for the early diagnosis of **cystic fibrosis in newborns**. TLI levels are significantly raised in early stages of cystic fibrosis and fall to subnormal levels as the disease progresses. Abnormally low TLI values have been reported in patients with **diabetes**.

Also, patients affected by **pancreatic carcinoma** seem to present abnormal TLI levels.

2. PRINCIPLE OF THE ASSAY

The principle of the assay is based on the competition between labelled TLI and TLI contained in calibrators or samples to be assayed for a fixed and limited number of antibody binding sites. After the incubation, the amount of labelled TLI bound to the antibody is inversely related to the concentration of unlabelled TLI present in calibrators or samples. The method adopted for B/F separation is based on the use of a precipitating reagent in which a second antibody is pre-precipitated and in excess.

3. REAGENTS PROVIDED IN THE KIT

¹²⁵ I-labelled cathodic trypsin	1 vial
Cathodic trypsin antiserum	1 vial
Zero calibrator	1 vial
Cathodic trypsin calibrators	6 vials
Control serum	1 vial
Precipitating reagent	2 bottles
Number of tubes	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming.

Reagents from different batches must not be mixed.

3.1. ¹²⁵I-labelled cathodic trypsin: lyophilized reagent

The vial contains cathodic trypsin labelled with ¹²⁵I, BSA, phosphate buffer and preservatives. Radioactivity is 57 kBq (1.55 µCi) or less on the calibration date.

Reconstitute the vial contents with 10 mL distilled water. The resulting solution is stable for three weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.2. Cathodic trypsin antiserum: lyophilized reagent

The vial contains antiserum raised in rabbits to human cathodic trypsin, BSA, phosphate buffer and preservatives. Reconstitute the vial contents with 10 mL distilled water. The resulting solution is stable for three weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

3.3. Zero calibrator: ready-to-use reagent

The vial contains 10 mL bovine serum and preservatives.

3.4. Cathodic trypsin calibrators: lyophilized reagent

The vials contain increasing amounts of cathodic trypsin, BSA and preservatives.

Reconstitute the contents of the 1-6 calibrator vials with 1 mL zero calibrator. The resulting solutions contain 5 - 10 - 30 - 80 - 200 - 400 ng/mL trypsin and are stable for three weeks at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage. *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

3.5. Control serum: lyophilized reagent

The vial contains cathodic trypsin, human serum and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Reconstitute the vial contents with 1 mL distilled water. Store the resulting solution in deep-frozen aliquots (-20°C or below).

3.6. Precipitating reagent: ready-to-use reagent

Each bottle contains 53 mL polyethylene glycol, TRIS buffer, antibody to rabbit IgG raised in goats, non-specific rabbit IgG and preservatives.

4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Disposable polystyrene tubes.
- Micropipettes with disposable tips (100, 1000 µL) (trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Multisample centrifuge capable of achieving 1500-2000 x g*.
- Water aspirator (optional).
- Gamma counter suitable for counting ¹²⁵I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

5. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulant EDTA has been tested and may be used with this assay (see §9.1 for heparin and citrate; do not use oxalate). Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles. *Decomplemented serum must not be used.*

6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each calibrator and sample.

- Operate according to the following scheme:

reagents \ tubes	Total activity	Calibrators 0-6	Samples
Calibrators	-	100 µL	-
Samples	-	-	100 µL
Tracer	100 µL	100 µL	100 µL
Antiserum	-	100 µL	100 µL

- **Mix** the contents of tubes with a Vortex and **incubate for one hour at room temperature**.
If low values of TLI are expected (e.g. chronic pancreatitis), sensitivity can be enhanced by adding the tracer 30 min after the other reagents kept at room temperature (30 min preincubation + one hour incubation at room temperature). Conversely, if high values of TLI are expected (e.g. acute pancreatitis), it is suggested to dilute the samples 1:2 with the zero calibrator (50 µL sample + 50 µL zero calibrator).
- Let the bottle of **precipitating reagent reach room temperature** and **mix** well by repeatedly tilting end over end.
- **Dispense 1 mL precipitating reagent** into all the tubes (except for total activity).
- **Mix** the contents of tubes with a Vortex and **let stand for 15 min at room temperature**.
- **Centrifuge** all tubes for 15 min at 1500-2000 x g*.
- **Discard the supernatant** by aspiration or decantation. When aspirating, avoid disturbing the precipitate. When decanting, let the tubes drain mouth down on blotting paper.
- **Measure the radioactivity** of precipitate.

7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Evaluate the binding ability as follows:

$$(B/T)_0\% = \frac{\text{zero calibrator mean counts}}{\text{total activity mean counts}} \times 100$$

Express the mean counts for each calibrator and unknown sample as a percentage of zero calibrator mean counts:

$$B/B_0\% = \frac{\text{calibrator or sample mean counts}}{\text{zero calibrator mean counts}} \times 100$$

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

Plot in linear-linear or semilog coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of TLI concentration expressed as ng/mL on the abscissa (x axis). A calibration curve is thus obtained (Fig. 2). Directly from the calibration curve, read the TLI concentration of each sample expressed as ng/mL. If the sample was diluted, the TLI concentration value derived from the diluted specimen must be multiplied by the dilution factor.

Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/Bo x 100
Total activity	20,540	—
Zero calibrator	11,235	100
5 ng/mL	9,830	87.5
10 ng/mL	8,898	79.2
30 ng/mL	6,348	56.5
80 ng/mL	3,876	34.5
200 ng/mL	2,326	20.7
400 ng/mL	1,626	14.5
Sample	6,680	59.5

By interpolation from the calibration curve, the sample is found to contain 26 ng/mL TLI.


Fig. 2

8. EXPECTED VALUES

The ranges given below are merely indicative: each laboratory should establish its own reference ranges.

Subjects	TLI, ng/mL	
	Mean \pm S.D.	Range
Normal controls	32.0 \pm 10.9	10 - 57
Chronic pancreatitis	16.5 \pm 12.9	< 1.5 - 47
Acute pancreatitis	350 \pm 260	92 - 850
Total pancreatectomy	< 1.5	–

9. SPECIFIC PERFORMANCE CHARACTERISTICS

9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by EDTA, haemolysis (up to 1000 mg/dL haemoglobin), lipaemia (up to 500 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin) or one freeze-thaw cycle of samples. TLI values determined in plasma samples obtained with heparin or citrate are about 10% lower than those determined in serum samples.

Cross-reactions. The RIA system was characterized for specificity to trypsin-related substances. Complete discrimination from bovine trypsin, negligible interference with anodic species and large cross-reactivity of cathodic trypsinogen are apparent.

9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 2.5 ng/mL at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator, that is, two standard deviations below zero.

9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	20	20	20
Mean (ng/mL)	10.10	45.10	97.60
Standard deviation	0.55	2.50	7.10
Coefficient of variation (%)	5.40	5.50	7.30

Reproducibility	A	B	C
Number of determinations	10	10	10
Mean (ng/mL)	10.70	46.20	101.20
Standard deviation	0.75	3.50	9.70
Coefficient of variation (%)	7.00	7.60	9.60

9.4. Trueness

The assay trueness has been checked by the dilution test. The recovery test has little meaning in heterologous systems like this (calibrator and analyte are different). In this particular case, trypsin calibrator could not be recovered as it forms non-immunoreactive complexes with the inhibitory systems present in the unknown sample.

Dilution test. One serum with high trypsin concentration was tested after serially diluting with the zero calibrator.

Dilution	Expected concentration, ng/mL	Measured concentration, ng/mL	% Recovery
neat	–	226.9	–
1:2	113.5	119.8	105.6
1:4	56.7	61.7	108.8
1:8	28.4	34.7	122.3
1:16	14.2	19.2	135.4
1:32	7.1	9.4	132.6

10. LIMITATIONS OF THE PROCEDURE

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration or decanting are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration or decanting of incubation mixture
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

11. WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

13. BASIC RULES OF RADIATION SAFETY

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 1.5 μCi (55 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

SCHEME OF THE ASSAY

- 1 - RECONSTITUTE THE REAGENTS.
- 2 - IDENTIFY ASSAY TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS BY THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

REAGENTS	TUBES	T	CAL 0-6	SAMPLES
CALIBRATORS		-	100 µL	-
SAMPLES		-	-	100 µL
TRACER		100 µL	100 µL	100 µL
ANTISERUM		-	100 µL	100 µL

- 4 - INCUBATE FOR ONE HOUR AT ROOM TEMPERATURE.
- 5 - ADD 1 mL PRECIPITATING REAGENT TO ALL THE TUBES (except for total activity) AND MIX THOROUGHLY.
- 6 - LET STAND FOR 15 MIN AT ROOM TEMPERATURE.
- 7 - CENTRIFUGE FOR 15 MIN AT 1500-2000 x g*.
- 8 - DISCARD SUPERNATANT.
- 9 - MEASURE THE RADIOACTIVITY OF PRECIPITATE.

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REFERENCES/BIBLIOGRAFIA/BIBLOGRAPHIE/LITERATUR/BIBLIOGRAFIA

T.E. ADRIAN et al.

Plasma trypsin in chronic pancreatitis and pancreatic adenocarcinoma.
Clin. Chim. Acta, 97 : 205 (1979).

A. ANDRIULLI et al.

Clinical evaluation of serum cationic trypsin-like immunoreactivity in patients with different pancreatic diseases.
Ital. J. Gastroenterol., 12 : 260 (1980).

A. ANDRIULLI et al.

Circulating trypsin-like immunoreactivity in chronic pancreatitis.
Dig. Dis. Sci., 26 : 532 (1981).

J.R. CROSSLEY, R.B. ELLIOT, P.A. SMITH

Neonatal screening for cystic fibrosis, using immunoreactive trypsin assay in dried blood spots.
Clin. Chim. Acta, 113 : 111 (1981).

P. DANDONA, E. ELIAS, A.G. BECKETT

Serum trypsin concentrations in diabetes mellitus.
Brit. Med. J., 2 : 1125 (1978).

D.R. GAMBLE, A. MOFFAT, V. MARKS

Serum immunoreactive trypsin concentrations in infectious and non-infectious illnesses and in juvenile diabetes.
J. Clin. Path., 32 : 897 (1979).

M.C. GEOKAS et al.

Determination of human pancreatic cathodic trypsinogen in serum by radioimmunoassay.
Am. J. Physiol., 236 : E77 (1979).

A.F. HEELEY

Screening for cystic fibrosis in the newborn using a dried blood disc I.R.T. assay.
The Proceedings of the First U.K. Symposium on Neonatal Screening for Cystic Fibrosis, held at the Royal Society of Medicine, London, on 5th November 1980, p. 21.

L.T. KIRBY et al.

Use of a dried blood spot in immunoreactive trypsin assay for detection of cystic fibrosis.
Clin. Chem., 27 (5) : 678 (1981).

G. LAKE-BAKAAR, S. MCKAVANAGH, J.A. SUMMERFIELD

Urinary immunoreactive trypsin excretion: a non-invasive screening test for pancreatic cancer.
Lancet, ii : 878 (1979).

G. LAKE-BAKAAR

Measurement of trypsin in duodenal juice by radioimmunoassay.
Gut, 21 : 402 (1980).

R. MALVANO

Radioimmunoassay of circulating trypsin-like immunoreactivity: methodological aspects.
The Proceedings of the First U.K. Symposium on Neonatal Screening for Cystic Fibrosis, held at the Royal Society of Medicine, London, on 5th November 1980, p. 6.

S. RECCHIA et al.


A comparative evaluation of two commercial RIA kits for trypsin-like immunoreactivity.
J. Nucl. Med. All. Sci., 24 (3-4) : 197 (1980).

Serum Immunoreactive Trypsin

Reports from the International Symposium Advances in the Diagnosis and in the Therapy of Pancreatic Diseases, G. Labò, P. Vezzadini eds.
Scand. J. Gastroenterol., 15 : Suppl. 62 (1980).

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710
FAX: 651-351-5669

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

APM0137

34664 7/10