TESTO-CTK
(P3093)

DiaSorin

English	p. 1
Italiano.....	p. 9
Français	p. 17
Deutsch.....	p. 26
Español	p. 35
Ελληνικά.....	p. 44


**TESTOSTERONE DIRECT
RADIOIMMUNOASSAY KIT**

**Procedure for quantitative determination of testosterone
in human serum or plasma samples**

For in vitro use only

1. INTRODUCTION

Testosterone, the main androgen produced by the Leydig cells of the testes, is a steroid hormone of molecular weight 288.4 daltons, having the following chemical formula:

Though mainly synthesized in the testes, testosterone is also produced – in much smaller amounts – in the adrenal cortex and in the ovaries. This is because a common meta-bolic pathway is followed when a steroid hormone is synthesized.

Testosterone is transported in the plasma by a beta-globulin, called testosterone-binding globulin. It is estimated that about 97-99% of the circulating testosterone is bound. The remainder, present as free component, is assumed to be the metabolically active portion of the total testosterone level.

To perform its biological activities, testosterone appears to require conversion into an active and more potent metabolite, dihydrotestosterone, and this can occur both in the circulation and in target tissues.

At the time of puberty, pituitary gonadotrophins increase, stimulating testicular maturation: LH acts directly on the interstitial mesenchymal elements, causing them to secrete testosterone and estrogens, as they develop into mature Leydig cells. FSH acts on the seminiferous tubules to induce and maintain normal spermatogenesis.

At the onset of puberty, testosterone stimulates the development of male sex characteristics (such as external genitalia, accessory sex organs, hair growth, linear growth, voice timbre, psyche, muscle and bone tissue mass) which it will maintain throughout life.

The age at which maximum testosterone levels are observed varies from 15-20 to 25 years according to the theories. Testosterone is secreted episodically throughout the day, so that multiple peaks can be observed.

2. PRINCIPLE OF THE ASSAY

The principle of the assay is based on the competition between labelled testosterone and testosterone contained in calibrators or samples to be assayed for a fixed and limited number of antibody binding sites. After the incubation, the amount of labelled testosterone bound to the antibody on the tube walls is inversely related to the concentration of unlabelled testosterone present in calibrators or samples. The method adopted for B/F separation is based on the use of antibody-coated tubes, where the antibody is coated on the tube walls.

3. REAGENTS PROVIDED IN THE KIT

Coated tubes	100
¹²⁵ I-labelled testosterone	1 bottle
Testosterone calibrators	7 vials
Control serum	1 vial
Number of tubes	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

3.1. Coated tubes

The inner surface of each tube is coated with testosterone antiserum raised in rabbits.

Before use, bring the coated tubes to room temperature prior to opening the box, to avoid condensation of humidity.

Securely reseal the box containing unused tubes. Do not mix different batches of coated tubes.

3.2. ¹²⁵I-labelled testosterone (red): ready-to-use reagent

The bottle contains 52 mL testosterone labelled with ¹²⁵I, BSA, phosphate-citrate buffer, preservatives and an inert red dye. Radioactivity is 59 kBq (1.6 µCi) or less on the calibration date.

3.3. Testosterone calibrators: ready-to-use reagent

Each vial contains steroid-free serum, preservatives and testosterone at the following concentrations: 0 - 0.25 - 0.5 - 1.0 - 2.5 - 5.0 - 10.0 ng/mL (0 - 0.87 - 1.73 - 3.47 - 8.7 - 17.3 - 34.7 nmol/L).

The zero calibrator vial contains 1 mL in human serum; the 1-6 calibrator vials contain 0.5 mL in newborn calf serum.

As no international standard preparation is currently available, the kit calibrators are referenced to an internal reference preparation (99% purity by HPLC). *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

3.4. Control serum: lyophilized reagent

The vial contains testosterone, human serum and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Reconstitute the vial contents with 1 mL distilled water. The resulting solution is stable for one week at 2-8°C. Store in deep-frozen aliquots (-20°C or below) for extended storage.

4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Micropipettes with disposable tips (50, 500 µL) (50 µL: trueness ± 3%, precision 2%; 500 µL: trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Thermostatically-controlled heating block or water bath capable of maintaining 37° ± 1°C.
- Device for aspiration of incubation mixture.
- Gamma counter suitable for counting ¹²⁵I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

5. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulants EDTA and heparin have been tested and may be used with this assay. Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

If testosterone levels greater than 10 ng/mL are expected, the samples should be diluted with the zero calibrator.

6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time. Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each calibrator and sample.

- Dispense reagents *in the bottom of coated tubes*. Operate according to the following scheme:

reagents	tubes	Calibrators 0-6	Samples
Calibrators		50 µL	–
Samples		–	50 µL
Tracer		500 µL	500 µL

- **Mix** the contents of tubes with a Vortex and **incubate for 3 hours at 37°C**.
- Carefully **aspirate** the incubation mixture. *Be sure that the aspirator tip touches the bottom of the coated tube so that all the liquid is removed. Failure to remove adhering solution adequately may result in poor reproducibility and spurious results. No trace of dye should still be visible.*
- **Measure the radioactivity** of tubes.

7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Compute the B/Bo ratio for each calibrator and unknown sample as follows:

$$B/Bo\% = \frac{\text{calibrator or sample mean counts}}{\text{zero calibrator mean counts}} \times 100$$

Plot in semilog or linear-linear coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of testosterone concentration expressed as ng/mL or nmol/L on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1). Directly from the calibration curve, read the testosterone concentration of each sample expressed as ng/mL or nmol/L. If the sample was diluted, the testosterone concentration value derived from the diluted specimen must be multiplied by the dilution factor.

Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/Bo x 100
Zero calibrator	10,559	100
0.25 ng/mL - 0.87 nmol/L	7,816	74.0
0.5 ng/mL - 1.73 nmol/L	6,999	66.3
1.0 ng/mL - 3.47 nmol/L	5,537	52.4
2.5 ng/mL - 8.7 nmol/L	4,184	39.6
5.0 ng/mL - 17.3 nmol/L	3,347	31.7
10.0 ng/mL - 34.7 nmol/L	2,350	22.3
Sample	5,376	50.9

By interpolation from the calibration curve, the sample is found to contain 1.2 ng/mL (4.3 nmol/L) testosterone.

B/Bo x 100


Fig. 1

8. EXPECTED VALUES

The ranges given below are merely indicative: each laboratory should establish its own reference ranges.

Normal males	2.8 - 9.0 ng/mL	10 - 30 nmol/L
Normal females	0.1 - 1.0 ng/mL	0.3 - 3.5 nmol/L

Conversion to nmol testosterone/L is possible using the following formula:
 $\text{nmol testosterone/L} = \text{ng testosterone/mL} \times 3.47$.

9. SPECIFIC PERFORMANCE CHARACTERISTICS

9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by anticoagulants (EDTA, heparin), haemolysis (up to 200 mg/dL haemoglobin), lipaemia (up to 500 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin). Testosterone values determined in plasma samples obtained with citrate are about 30% lower than those determined in serum samples.

Cross-reactions. The percentage of cross-reactions, calculated according to Abraham, shows the specificity of the antiserum used.

- Testosterone	100%
- 5-dihydrotestosterone	6.9%
- Androstenedione	1.1%
- 19-norethisterone	0.33%
- Danazol	0.15%
- Dehydroepiandrosterone	$4.5 \times 10^{-2}\%$
- 5 alpha-androstan, 3 alpha, 17 beta diol	$2.5 \times 10^{-2}\%$
- Testosterone glucuronide	$2.0 \times 10^{-2}\%$
- Cortisol	} $< 1.0 \times 10^{-2}\%$
- Corticosterone	
- Cholesterol	
- Deoxycorticosterone	
- 17-epitestosterone	
- Testosterone sulphate	
- Androsterone	
- 3 beta-oxyandrostan, 5 beta, 17 one	
- Progesterone	
- Estrone	
- Estradiol, Ethynylestradiol	
- Estriol	
- Dexamethasone	

9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 0.02 ng/mL (0.07 nmol/L) at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator, that is, two standard deviations below zero.

9.3. Precision

Different sample pools, containing different concentrations of analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	14	14	14
Mean (ng/mL)	1.13	3.61	7.27
Standard deviation	0.091	0.158	0.262
Coefficient of variation (%)	8.06	4.38	3.61

Reproducibility	D	E	F	G
Number of determinations	10	10	10	10
Mean (ng/mL)	0.419	2.164	4.650	9.342
Standard deviation	0.032	0.135	0.339	0.661
Coefficient of variation (%)	7.58	6.23	7.28	7.08

9.4. Trueness

The assay trueness has been checked by the dilution and recovery tests.

Dilution test. Two sera with high testosterone concentration were tested after serially diluting with the zero calibrator.

Dilution	Expected concentration, ng/mL	Measured concentration, ng/mL	% Recovery
neat	–	5.38	–
1:2	2.69	2.60	96.7
1:4	1.34	1.30	96.7
1:8	0.67	0.67	99.7
1:16	0.34	0.34	101.2
1:32	0.17	0.18	107.1
neat	–	5.13	–
1:2	2.57	2.79	108.8
1:4	1.28	1.32	102.9
1:8	0.64	0.68	106.0
1:16	0.32	0.31	96.7
1:32	0.16	0.15	93.6

Recovery test. Two sera with low testosterone concentration were tested as such and after mixing with increasing amounts of testosterone.

Added concentration, ng/mL	Expected concentration, ng/mL	Measured concentration, ng/mL	% Recovery
–	–	0.30	–
3.88	4.18	4.10	98.1
2.46	2.76	2.50	90.6
1.11	1.41	1.50	106.4
0.69	0.99	1.00	101.0
0.39	0.69	0.70	101.4
–	–	4.84	–
3.88	8.72	8.90	102.1
2.46	7.30	7.40	101.4
1.11	5.95	6.40	107.6
0.69	5.53	5.70	103.1
0.39	5.23	5.30	101.3

10. LIMITATIONS OF THE PROCEDURE

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration of incubation mixture
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

11. WARNINGS AND PRECAUTIONS

Test components contain sodium azide as a preservative. Because sodium azide may form explosive lead or copper azide in plumbing, it is recommended that drains be thoroughly flushed with water after disposal of solutions containing sodium azide (Council Directive 99/45/EC).

R 22 – Harmful if swallowed.

R 31 – Contact with acids liberates toxic gas.

S 28 – After contact with skin, wash immediately with plenty of water.

S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Sodium merthiolate (Council Directive 99/45/EC):

R 20/21/22 – Harmful by inhalation, in contact with skin and if swallowed.

R 33 – Danger of cumulative effects.

S 28 – After contact with skin, wash immediately with plenty of water.

S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

13. BASIC RULES OF RADIATION SAFETY

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 1.7 μCi (63 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. It's receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state of which the Commission has entered into an agreement for the

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

SCHEME OF THE ASSAY

- 1 - RECONSTITUTE THE CONTROL SERUM.
- 2 - IDENTIFY COATED TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS ACCORDING TO THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

	TUBES	CAL 0-6	SAMPLES
REAGENTS			
CALIBRATORS		50 μL	-
SAMPLES		-	50 μL
TRACER		500 μL	500 μL

- 4 - INCUBATE FOR 3 HOURS AT 37°C.
- 5 - ASPIRATE THE INCUBATION MIXTURE.
- 6 - MEASURE THE RADIOACTIVITY OF TUBES.

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F.J. AULETTA, B.V. CALDWELL, G.L. HAMILTON

Androgens: testosterone and dihydrotestosterone.

In: «Methods of hormone radioimmunoassay», B.M. Jaffe, H.R. Behrman eds., Academic Press, New York, p. 715 (1979).

C.W. BARDIN, C.A. PAULSEN

The testes

In: «Textbook of endocrinology», R.H. Williams ed., W.B. Saunders, Philadelphia, p. 293 (1981).

W.P. COLLINS et al.

Radioimmunoassay of plasma testosterone.

J. Steroid Biochem., **3** : 333 (1972).

D.C. CUMMING, S.R. WALL

Non-sex hormone-binding globulin-bound testosterone as a marker for hyperandrogenism.

J. Clin. Endocrinol. Metab., **61** (5) : 873 (1985).

P. DEWIS, D.C. ANDERSON

A practical approach to the hirsute patient.

Practitioner, **226** (2) : 223 (1982).

M.L. DUFAU et al.

Radioimmunoassay of plasma testosterone.

Clin. Chim. Acta, **37** : 109 (1972).

S. FURUYAMA, D.M. MAYES, C.A. NUGENT

A radioimmunoassay for plasma testosterone.

Steroids, **16** : 415 (1972).

S.P. GLICKMAN, R.L. ROSENFELD, R.M. BERGENSTAL, J. HELKE

Multiple androgenic abnormalities, including elevated free testosterone, in hyperprolactinemic women.

J. Clin. Endocrinol. Metab., **55** (2) : 251 (1982).

D.M. HALPIN, J.M. BURRIN, G.F. JOPLIN

Serum testosterone levels in women with Cushing's disease.

Acta Endocrinol., **122** (1) : 71 (1990).

A.A.A. ISMAIL, J.A. LORAINE

Testosterone.

In: «Hormone assays and their clinical application», J.A. Loraine, E.T. Bell eds., Churchill Livingstone, Edinburgh, p. 580 (1976).

W.W. LUTHOLD et al.

Serum testosterone fractions in women: normal and abnormal clinical states.

Metabolism, **42** (5) : 638 (1993).

E. NEACSU, L. SIMIONESCU, A. ZIMEL, A. CARAGHEORGHEOPOLO

The development of a radioimmunoassay system for testosterone (T) and dihydrotestosterone (DHT). Part 3: The preparation of radioiodinated testosterone.

Endocrinologie, **28** (3-4) : 127 (1990).

R.S. RITTMASER, D.M. ARAB, L. LEHMAN

Dose-response effect of depot leuprolide acetate on serum androgens in hirsute women.

Fertil. Steril., **65** (5) : 912 (1996).

G.L. ROYER et al.

Relationship between age and levels of total, free, and bound testosterone in healthy subjects.

Curr. Ther. Res., **35** (3) : 345 (1984).

E.A. RUCHHOFT, K.E. ELKIND-HIRSCH, R. MALINAK

Pituitary function is altered during the same cycle in women with polycystic ovary syndrome treated with continuous or cyclic oral contraceptives or a gonadotropin-releasing hormone agonist.

Fertil. Steril., **66** (1) : 54 (1996).

L. SIMIONESCU, E. NEACSU, A. ZIMEL, A. CARAGHEORGHEOPOL

The development of a radioimmunoassay system for testosterone (T) and dihydrotestosterone (DHT). Part 2: The preparation of antisera to testosterone.

Endocrinologie, **28** (3-4) : 107 (1990).

K.B. SINGH, D.K. MAHAJAN, J. WORTSMAN

Effect of obesity on the clinical and hormonal characteristics of the polycystic ovary syndrome.

J. Reprod. Med., **39** (10) : 805 (1994).

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710
FAX: 651-351-5669

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

APM0136

34669 7/10