DiaSorin Inc.
 1951 Northwestern Avenue, P.O. Box 285,
 Stillwater, MN 55082-0285 U.S.A.
 Phone: 651-439-9710 - Fax: 651-351-5669

IRMA-mat[®] Ferritin

Instructions for Use

REF 324.511

100 tests

English

IVD

For professional use only!

Intended Use

In vitro assay for the quantitative determination of ferritin in human serum.

Summary and Explanation of the Test

Ferritin is a high-molecular iron storage protein (MW approx. 450,000) which can incorporate varying amounts of iron in the form of iron-(III)-hydroxide/phosphate (max. 4,500 iron ions/molecule) (9). The iron represents approximately 25 % of the total molecular weight.

Ferritin can mainly be found in the cytoplasm of the reticuloendothelial cells, in liver cells and to a small extent in the precursors of the red cells in the bone marrow (8). It stores iron in a biologically mobilisable form.

The determination of serum ferritin has its significance in the diagnosis and monitoring of iron deficiency. Reduced ferritin serum levels indicate iron deficiency; iron overload or idiopathic haemochromatosis are associated with elevated ferritin serum concentrations (3).

Elevated ferritin levels have also been found in a variety of malignancies, especially in Hodgkin- or Non-Hodgkin lymphoma, acute or chronic leukaemia as well as in breast or bronchial carcinoma (1, 2, 4, 5, 6, 7).

Principles of the Procedure

Two-site immunoradiometric assay (sandwich principle) using two highly specific monoclonal antibodies. Antibody-coated polystyrene tubes serve as solid phase. The tracer antibody and the coated antibody react simultaneously with the ferritin present in patient samples or standards. Excess tracer is removed by a washing step and the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

CONTENTS

Determinations	100
¹²⁵ I-anti-ferritin, monoclonal (mouse), red radioactive content (kBq / μ Ci)	32 mL < 456 /12
7 standards A-G (0.6 mL) in human serum albumin (The exact concentration is indicated on each vial label)	1 set
Diluent (0 ng/mL) in human serum albumin	11 mL
Test tubes, coated with anti-ferritin, monoclonal (mouse)	2 x 50
Control serum	0.6 mL
Quality control report	1

Material Required but not Provided

- Micropipets (25 μ L, 300 μ L) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively an appropriate automated analyser system, if available
- 0.9% NaCl solution
- Uncoated polystyrene tubes for the dilution of sera and controls

Warnings and Precautions for Users

Reagents Containing Human Source Material

TREAT AS POTENTIALLY INFECTIOUS.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 12 μ Ci (456 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to “Decontamination of Laboratory Sink Drains to Remove Azide Salts,” in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly (avoid foam formation).

Storage of Reagents

Store all reagents at 2–8°C until the expiry date printed on the package label.

- **Keep upright for storage.**
- **Keep away from direct light.**

Sample Collection and Storage

- Collect samples using standard procedures.
- Sample material: serum; disturbances with plasma are not known.
- Storage at 2 - 8°C: 24 h
- For longer storage periods: freeze to below -20°C.
- Freezing and thawing samples once does not affect the test results.
- Stored samples should be thoroughly mixed prior to use (Vortex mixer).
- Do not use samples which are agglutinated, lipemic, hemolysed, icteric or contaminated.

Interfering Substances

No interference with test results is seen by concentrations of bilirubin < 0.125 mg/mL, haemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The single components of each kit are carefully matched. In case of exchange or mixture of any components from different lots the manufacturer does not guarantee reliable results. See bottom of the kit for the lot numbers of all components.
- Strictly adhere to the pipetting step sequence.
- The measuring time at the gamma counter must be adjusted for counting at least 1 minute.
- This kit must not be used after the expiry date printed on the package label.
- Observe quality control guidelines for medical laboratories.
- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that standards and samples are assayed in duplicate.

If values are expected out of the standard range (> 1000 ng/mL), samples should be further diluted with the diluent (e.g. by factors 10, 100, 1000).

Alternatively to the manual performance and evaluation an appropriate automated analyser system can be used on responsibility of the laboratory.

1. Pipet 25 μL of each standard, control or patient sample onto the bottom of the corresponding coated tube.
2. Add 300 μL ^{125}I -anti-ferritin, vortex.
3. Incubate for 2 hours (\pm 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes three times with 2 mL of 0.9% NaCl solution.
6. Measure radioactivity (cpm) in all tubes.

25 μL	Pipet standard, control or patient sample
300 μL	Add ^{125}I -anti-ferritin
	Mix
2 h (\pm 5 min)	Incubate at room temperature (18–25°C) on a shaker*
	Aspirate
3 x 2 mL	Wash with 0.9% NaCl solution
1 min	Measure (gamma scintillation counter)

* Keep shaking conditions constant!

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

Calculation of Results

The standard curve is established manually as follows:

1. Determine the mean CPM for each pair of duplicate tubes.
2. Divide the mean CPM of each standard (B) by the mean CPM of the highest standard (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding ($\%B/B_{\text{max}}$) for each standard.
3. On semi-log paper, plot the relative binding of each standard ($\%B/B_{\text{max}}$) on the Y-axis versus the corresponding concentrations on the X-axis.
4. Read sample concentrations (ng/mL) directly off the standard curve by their corresponding relative binding ($\%B/B_{\text{max}}$).

Samples with counts higher than the 1000 ng/mL standard must be diluted with the kit diluent and re-assayed. For the final concentration of such samples, the appropriate dilution factor has to be considered. An example of a standard curve is given in the QC-certificate. This curve must not be used for the calculation of unknown samples.

The instrumental calculation of radio-immunological measuring values is based on a spline approximation.

Quality Control

Observe quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory. This control should be simultaneously tested with each test run and treated like patient samples.

The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Expected Values

Normal range (evaluated in serum of healthy blood donors)
(5 - 95th percentile)

Men (n = 99)	15	–	332 ng/mL
Women (premenopausal; n = 73)	6	–	85 ng/mL

Ferritin levels are age-dependent and are, in women, severely affected by the reproductive status (iron depletion due to menstruation) so that each laboratory should establish its own reference range.

Limitations of the Procedure

Elevated ferritin serum levels have been found in cases of liver diseases, acute and chronic inflammation and different types of infection. Therefore, ferritin serum levels have to be interpreted in context with the clinical picture and other diagnostic procedures.

All tests, in which antigen is incubated together with labelled antibodies and immobilised antibodies in a liquid phase, bear the risk that samples containing extremely high concentrations of the antigen will give measuring values below those of the highest standard. In case of the IRMA-mat® Ferritin, this phenomenon is observed at concentrations exceeding 25,000 ng/mL. If such values are suspected, measurement should be repeated after further dilution (e.g. by factors 10, 100, 1000) of the specimen.

There is but limited comparability with Ferritin values obtained with kits from other manufacturers because of the tissue-specificity (liver, spleen) of the antigens used in these kits.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may give falsely elevated or decreased values. Although HAMA-neutralising agents are added, extremely high HAMA serum concentrations may occasionally influence results. These samples should not be used for the IRMA-mat® Ferritin assay.

Analytical Data

Calibration

The standard used in the IRMA-mat® Ferritin assay is orientated using the reference standard NIBSC 94/572.

Measuring range

IRMA-mat® Ferritin allows to measure concentrations between 1 and 1000 ng/mL.

High-dose hook

No high-dose hook effect was observed for ferritin concentrations up to 25,000 ng/mL.

Precision

Intra-assay Mean value (ng/mL)	CV (%)	n=	Inter-assay Mean value (ng/mL)	CV (%)	n=
25.3	4.3	9	34.5	9.0	12
162.0	3.0	9	233.8	8.0	12
263.0	3.0	9	332.2	8.0	12

Analytical sensitivity

The detection limit of the IRMA-mat® Ferritin assay is < 1.0 ng ferritin/mL. This limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest ferritin concentration that can be differentiated from zero with statistical significance.

Specificity

The following cross-reactivities were found:

human liver ferritin	100 %
human spleen ferritin	60 %
human heart ferritin	< 10 %

Linearity upon Dilution

A patient's serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Dilution	Measured value (ng/mL)	Expected value (ng/mL)	Recovery (%)
1	708.3	728.9	97.2
1 : 1.25	611.8	581.9	105.1
1 : 2.5	282.4	287.8	98.1
1 : 5	136.6	140.7	97.0

Recovery

A patient's serum with low ferritin content was spiked with different amounts of ferritin and then measured.

Original concentration: 33.5 ng/mL

Measured value (ng/mL)	Expected value (ng/mL)	Recovery (%)
232.0	232.0	100.0
215.3	192.3	111.9
111.4	112.9	98.7
72.7	73.2	99.4

DiaSorin Inc.
1951 Northwestern Avenue, P.O. Box 285,
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710 - Fax: 651-351-5669

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