ACTH IRMA

Immunoradiometric Assay

For the quantitative determination of adrenocorticotropic hormone (ACTH) in human EDTA plasma

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1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE ONLY.

This kit is intended for the quantitative determination of human adrenocorticotrophic hormone (ACTH) in EDTA plasma. The measurement of plasma ACTH is used in the diagnosis of adrenal dysfunctions.

2. SUMMARY AND EXPLANATION

The measurement of ACTH in human plasma has been of great clinical interest for many years. In 1948 Sayers et. al. described a method for determining ACTH based on its bio-activity.¹ Since that time radioimmunoassays (RIA) measuring ACTH were developed but were based on a primary antibody which did not measure biological activity. The DiaSorin ACTH IRMA is designed to detect whole molecule ACTH (1-39). ACTH IRMA is a more sensitive method for the detection of ACTH in plasma than RIA, allowing for clearer differentiation of normal and ACTH-suppressed levels.²

ACTH is a secretory product of the pituitary adrenocorticotrophic cells within the anterior pituitary gland.³ It is synthesized as part of a large 265 amino acid precursor molecule POMC (proopiomelanocortin).⁵ ACTH is a 39 amino acid hormone with the following structure:

NH₂-Ser-Tyr-Ser-Met-Glu-His-Phe-Arg-Trp-Gly-Lys-Pro-Val-Gly-Lys-Lys-Arg-Arg-Pro-Val-Lys-Val-Tyr-Pro-Asn-Gly-Ala-Glu-Asp-Glu-Ser-Ala-Glu-Ala-Phe-Pro-Leu-Glu-Phe-OH.

The sequence 1-13 is the same as alpha melanocyte stimulating hormone. Sequence 1-26 is considered the active portion of the molecule and is required for full biological activity.⁴

The ACTH molecule is a basic protein and labile to the action of proteolytic enzymes. It functions by increasing the synthesis and release of all adrenal corticosteroids from the adrenal cortex including cortisol, aldosterone, and the adrenal androgens. The major factors effecting ACTH release include CRH (corticotropin releasing hormone), free cortisol concentration in plasma, stress and the sleep-wake cycle. Control of ACTH concentration is mediated through a negative feedback mechanism and the central nervous system. ACTH rises in response to reductions in plasma cortisol levels. Conversely, as cortisol rises, ACTH release is inhibited at the level of the pituitary.³

ACTH is secreted by the pituitary gland in a pulsatile manner but follows an overriding diurnal rhythm. Levels of ACTH in the normal individual are typically highest in the early morning (6-8 A.M.) and lowest late in the day. Because of this, it is customary to perform ACTH blood draws early in the day. Conditions that have been demonstrated to upset the normal diurnal rhythm include Cushing's syndrome and ectopic ACTH syndrome as well as physiologic stress such as surgery.⁶⁻⁸

ACTH results are used in conjunction with other diagnostic assays to evaluate patients with various diseases involving the adrenal cortex. Cushing's syndrome is a metabolic disorder resulting from an increased production of cortisol by the adrenal gland. The diagnosis of Cushing's syndrome depends upon the demonstration of increased cortisol production and the lack of dexamethasone suppression. Once the diagnosis has been established, further testing is performed to determine its cause.^{3,9,10} This syndrome is commonly caused by adrenal hyperplasia, secondary to pituitary ACTH overproduction. The overproduction of cortisol is due to either hypersecretion of pituitary ACTH or ectopic production of ACTH by nonendocrine tumors.

The term "Cushing's Disease" is traditionally used to define an individual who has an ACTH-producing pituitary tumor; however, it is used by some clinicians to describe anyone who has hypersecretion of pituitary ACTH regardless of the presence of a tumor.

When Cushing's syndrome is caused by excessive secretion of pituitary ACTH, elevation of ACTH may occur. The site of ACTH production in corticotroph microadenomas has been localized using the technique of bilateral inferior petrosal sinus catheterization.^{11,12}

Ectopic ACTH syndrome is another possible cause of Cushing's syndrome. In these patients, non-endocrine malignant tumors produce high levels of ACTH and POMC resulting in very high cortisol production which is not subject to feedback inhibition at any level.^{13,14}

When Cushing's syndrome is due to a primary adrenal abnormality, such as an adrenal tumor, the adrenal gland acts independently of ACTH and pituitary secretion is suppressed. These patients have very low or undetectable ACTH concentrations and no response is seen with dexamethasone suppression or metyrapone stimulation tests.³

Adrenal insufficiency, is a condition which leads to inadequate cortisol production. It can be due to a primary inability of the adrenal to produce sufficient quantities of cortisol or secondary failure due to insufficient production of ACTH. The basal ACTH concentration of patients with primary adrenocortical insufficiency (Addison's disease) will be markedly elevated. When these patients are administered endogenous ACTH, the adrenals are unresponsive. Secondary adrenal insufficiency is caused by pituitary ACTH deficiency. Plasma ACTH levels may help to distinguish between primary and secondary adrenal insufficiency, since they are elevated in the former and reduced or absent in the latter.³

Nelson's syndrome is a condition caused by the development of large pituitary tumors in patients who have undergone bilateral adrenalectomy for Cushing's syndrome. ACTH levels are typically followed in these patients to help monitor clinical course and response to therapeutic regimens.^{15,16}

3. PRINCIPLE OF THE ASSAY

The DiaSorin ACTH IRMA is a highly specific and sensitive immunoradiometric assay for the measurement of human ACTH 1-39. This is accomplished using two ACTH antibodies, each one specific for different regions of the ACTH molecule. The iodinated tracer used in this system contains a purified polyclonal goat antibody specific for ACTH 26-39 and an iodine-125 labeled monoclonal antibody specific for ACTH 1-17. When samples are incubated with the tracer and mouse anti-goat coated polystyrene beads, the ACTH-specific antibodies will bind to the ACTH present in the sample. Only ACTH 1-39 present in the sample will bind to both antibodies to form an antibody complex. The unbound radioactivity is washed away and the radioactivity bound to the bead is measured. Concentrations of ACTH present in the sample are directly proportional to the radioactivity measured. Results are calculated by comparing the CPM's of each sample to the CPM's for the ACTH calibrators included with the kit.

4. REAGENTS PROVIDED IN THE KIT

ACTH IRMA Calibrator 0	1 vial/ 5.0 mL	1 vial/ 5.0 mL
ACTH IRMA Calibrators (1-5)	5 vials/ 2.0 mL	5 vials/ 2.0 mL
ACTH IRMA Beads	1 container/ 65 beads	2 container/ 65 beads
ACTH IRMA Tracer	2 vials/ 2 mL	4 vials/ 2 mL
ACTH IRMA Wash Solution	1 vial/ 50 mL	1 vial/ 50 mL
ACTH IRMA Controls (Levels 1-3)	3 vials/ 2.0 mL	3 vials/ 2.0 mL
Number of tests	65	130

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store all reagents at 2-8°, unless mentioned otherwise below, until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

When reconstituting the contents of the vials, mix gently to avoid foaming. Store all reconstituted reagents as recommended below. Reagents from different batches must not be mixed.

4.1 ACTH IRMA Calibrators (0-5): lyophilized reagent

Set of six synthetic ACTH (1-39) calibrators, at concentrations ranging from 0-1500 pg/mL, in ACTH-free human plasma with stabilizers and 0.1% sodium azide added as preservative. Exact concentration values are assigned with each lot. Calibrators are calibrated against synthetic human ACTH (1-39). The HPLC-purified synthetic human ACTH was tested and found to react equivalently to the synthetic human ACTH (1-39) obtained from the National Hormone and Pituitary Program (University of Maryland). The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended. Reconstitute 0 calibrator with 5.0 mL purified water. Reconstitute calibrators 1-5 with 2.0 mL of purified water. At room temperature, rotate calibrators for 25 - 30 minutes or let stand for 60 minutes with frequent mixing until contents are completely dissolved. Mix thoroughly before using. Store the reconstituted calibrators as follows:

For periods up to 5 days: Store at 2-8°C.

For prolonged storage (kit expiration) between use: Store at -15°C or lower.

The calibrators (including the 0 calibrator) can be frozen and thawed up to three times. Calibrators may be held at room temperature for a maximum of four 1 hour intervals during their storage at 2-8°C.

4.2 ACTH IRMA Beads: ready to use reagent

Polystyrene beads are coated with polyclonal mouse anti-goat antibody.

4.3 ¹²⁵I ACTH IRMA Tracer: ready to use reagent

Each vial contains 2 mL of polyclonal goat anti-ACTH and ¹²⁵I labeled-mouse monoclonal anti-ACTH diluted in buffered serum containing red dye and 0.1% sodium azide.

4.4 ACTH IRMA Wash Solution Concentrate: liquid solution

Contains a concentrated buffered surfactant. Prepare a working wash solution by diluting the entire vial contents to a final volume of 500 mL with distilled or deionized water.

4.5 ACTH IRMA Control, Levels 1-3: lyophilized reagent

Human plasma is spiked with the appropriate amount of synthetic human ACTH. 0.1% sodium azide and other stabilizers are added. Reconstitute each vial with 2.0 mL of purified water. At room temperature, rotate controls for 25-30 minutes or let stand with frequent mixing for 60 minutes until contents are completely dissolved. Mix thoroughly before using. Treat each control as an unknown sample. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. Store the reconstituted controls as follows:

For periods up to 5 days: Store at 2-8°C.

For prolonged storage (kit expiration) between use: Store at -15°C or lower.

The controls can be frozen and thawed up to three times. Controls may be held at room temperature for a maximum of four 1 hour intervals during their storage at 2-8°C.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected.

This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 7.2 μCi (266 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material. For practitioners or institutions receiving radioisotopes under a general license: This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license: The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.2 A decrease in maximum binding.
- 6.3 High nonspecific zero binding.
- 6.4 Poor duplicate values.

7. SPECIMEN COLLECTION AND PREPARATION OF PLASMA

Four hundred microliters of EDTA plasma are required to assay the specimen in duplicate in the ACTH IRMA. Heparinized samples should not be utilized in the ACTH IRMA. Grossly icteric, hemolyzed or lipemic samples should not be used.

Blood should be collected aseptically by venipuncture in a 5 or 10 mL EDTA vacutainer tube. Mix sample by gentle inversion and immediately place on crushed ice. Centrifuge sample within 30 minutes of time collected. Centrifugation should be performed for 15 minutes using approximately 1500 x g*. Immediately remove plasma from cells using a plastic pipette and transfer into a polypropylene/polystyrene tube. Store plasma samples at 2-8°C prior to assay, provided the assay can be performed within 4 hours. If the samples can not be assayed within 4 hours, samples must be stored at -20°C. Samples may be stored at 2-8°C for up to 4 hours or at -20°C for up to 6 months.

NOTE: DiaSorin studies on freshly drawn samples show an average decrease of 20-35% in ACTH concentrations after one freeze thaw. Minimal changes were observed after three subsequent freeze thaws.

A study (n = 50) was done using siliconized versus non-siliconized tubes. The geometric mean for the siliconized tubes was 28.8 pg/mL (2SD range 13.9-59.7). The geometric mean for the non-siliconized tubes was 26.7 pg/mL (2SD range 12.4-57.5).

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable polypropylene or polystyrene tubes, 12 x 75 mm.
Borosilicate glass tubes may not be used.
- 8.2 Test tube rack.
- 8.3 Gamma counter capable of counting iodine-125.
- 8.4 Vortex mixer.
- 8.5 Pipetting devices:
 - a. Micropipettor calibrated to deliver 50 µL and 200 µL.
 - b. Repeating dispensers, calibrated to deliver 50 µL, 200 µL and 1 mL.
 - c. Volumetric pipettes for reconstituting controls and calibrators; 2.0 mL and 5.0 mL.
- 8.6 500 mL graduated cylinder.
- 8.7 Parafilm or equivalent for covering test tubes.
- 8.8 Teflon-coated forceps or other device for dispensing beads.
- 8.9 Aspiration unit.
- 8.10 Purified water for component reconstitution.

9. ASSAY PROCEDURE

Allow all reagents to come to room temperature prior to use.

- 9.1 Reconstitution of ACTH IRMA Calibrators / Controls
 - a. Remove lyophilized ACTH IRMA Calibrators and Controls from the refrigerator.
 - b. Using a volumetric pipette, add 2.0 mL of purified water to ACTH IRMA Calibrators and Controls. Similarly, add 5.0 mL of purified water to the ACTH 0 Calibrator.
 - c. Gently tap the bottles to dislodge the lyophilized pellet and mix gently.
 - d. At room temperature, rotate calibrators and controls for 25 - 30 minutes or let stand with frequent mixing for 60 minutes until contents are completely dissolved. Mix thoroughly before using.
 - e. Reconstitution is complete when all solutions appear homogeneous.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.2** Preparation of ACTH IRMA Wash Solution.
The ACTH IRMA Working Wash Solution Concentrate is diluted 10X in purified water to give a working strength solution.
- a. Empty the 50 mL bottle containing the Concentrate into a 500 mL graduated cylinder.
 - b. Following the addition, rinse the bottle contents with purified water and add the rinse of the concentrate to the 500 mL cylinder.
 - c. Bring the final volume of the 500 mL cylinder to 500 mL with the continued addition of purified water. Stir to ensure complete dissolution of any precipitated material and uniform mixing. Label this material ACTH IRMA Wash Working Solution.
- 9.3** Set up properly labeled 12 x 17 mm polypropylene/polystyrene tube, in duplicate. Add reagents to the tubes as follows:
- a. **Calibrator 0**
200 μ L Calibrator 0
50 μ L 125 I ACTH IRMA Tracer (red)
 - b. **Calibrators (1-5)**
200 μ L Calibrator
50 μ L 125 I ACTH IRMA Tracer (red)
 - c. **Controls and unknown samples**
200 μ L sample
50 μ L 125 I ACTH IRMA Tracer (red)
- 9.4** Vortex the rack containing all tubes three times for 1-2 seconds each. Tube contents should not be allowed to foam.
- 9.5** Dispense one ACTH IRMA bead into each tube (except the Total Count tubes) with Teflon-coated forceps or suitable dispenser for handling beads. (Do not use fingers.)
- 9.6** Cover the tubes with parafilm or equivalent.
Incubate the tubes for 20 (\pm 2) hours at 18-25°C. Record temperature at the beginning and the termination of the assay.
Following the 20 (\pm 2) hour incubation, the tubes are washed as outlined below:
- a. Aspirate the reaction mixture into an appropriate container for radioactive waste.
 - b. Wash each tube with the ACTH IRMA Working Wash Solution.
Wash the bead by vigorously dispensing 1 mL of Working Wash Solution into each test tube with sufficient force to raise the bead from the bottom of the test tube. Aspirate following each wash step. Repeat the wash procedure two times for a total of three washes.
- 9.7** Measure the radioactivity present in each tube using a gamma counter. Tubes should be counted for 1 minute or longer (see Limitations of Procedure section).

10. PROCEDURAL COMMENTS

- 10.1** When adding beads to the tubes, tilt the test tube rack slightly and allow beads to roll into the tubes. This will prevent excessive splashing of the test solution. Do not handle beads with fingers.
- 10.2** Pipette 125 I ACTH IRMA Tracer to the lower one third of the test tube.
- 10.3** Samples reading greater than the highest calibrator must be diluted with 0 Calibrator and re-assayed.
- 10.4** All patient samples MUST be kept on ice until assayed.
- 10.5** To obtain optimal kit performance beads must be thoroughly washed.

- 10.6 No assay drift was observed in assays up to 130 tubes.
- 10.7 In order for a laboratory to completely monitor the performance of an IRMA there are additional factors which may be checked. DiaSorin suggests a check of the following parameters to assure consistent kit performance.
 - a. **Maximum Binding**
Counts per minute (CPM) of "5" Calibrator Tube
 - b. **Nonspecific Binding**
CPM of 0 Calibrator Tube

11. QUALITY CONTROL

Each laboratory must include at least one high and one low control in every assay, to monitor assay performance. Commercially available controls or the 3 reference controls provided with the kit may be utilized. The kit controls contain ACTH at concentrations typical of normal and elevated patient ranges. The kit controls have been assayed using the DiaSorin ACTH IRMA kit. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. The controls used should be treated as unknown specimens and assayed in duplicate. Quality control charts should be maintained by the laboratory to follow performance of the controls. Control results should be analyzed by statistical methods that adequately evaluate trends. Acceptable performance limits should be determined by the individual laboratory for each level of control using statistically based methods designed to detect both random and systematic errors. Control results must meet the laboratory's criteria for acceptability prior to reporting patient test results.¹⁷⁻¹⁹

12. CALCULATION OF RESULTS

To determine the concentration of ACTH found in unknown and control samples, a calibrator curve is prepared with every assay using the calibrator concentrations stated on the vial labels. Values for unknowns are generated as follows:

- 12.1 Calculate the average CPM for each calibrator, control, and unknown sample.
- 12.2 Using log-linear graph paper, plot the CPM of each calibrator level on the ordinate (Y-axis) against the calibrator concentration on the abscissa (X-axis).

NOTE: Automated data reduction programs may also be utilized in the analysis of data. DiaSorin utilizes Multi-Calc (Pharmacia) with a log-CPM smooth-SPLINE fit program. Other data reduction methods must be validated before incorporating for regular use.
- 12.3 Interpolate the levels of ACTH in the samples from the plot.
- 12.4 If any unknown sample reads greater than the highest calibrator, it should be diluted with 0 calibrator and reassayed.
- 12.5 If an unknown sample has been diluted, correct for the appropriate dilution factor.
- 12.6 Values less than Calibrator 1 but greater than 1.5 pg/mL (assay sensitivity) may be estimated using the formula:

$$\text{Values of unknown} = \frac{\text{CPM (unknown)}}{\text{CPM (Ca1. 1)}} \times \text{Value of Cal. 1}$$

- 12.7 The reportable range of the assay falls within the concentration of the highest calibrator and assay sensitivity.

Sample Data

Typical sample data and a calibrator curve for ACTH IRMA are shown in TABLE I and Figure 1; this information is for reference only and should not be used for the calculation of any value. This data was generated with fresh tracer.

Data Reduction

The DiaSorin QC lab uses a smoothed spline curve fit.

TABLE I
DiaSorin ACTH IRMA Sample Data

Tube	Duplicate CPM	Average CPM	Conc. (pg/mL)
Calibrator 0	563 591	577	
Calibrators (pg/mL)			
1 (9.9)	1,562 1,692	1,627	
2 (23.3)	3,348 3,385	3,367	
3 (96.1)	11,573 11,464	11,519	
4 (328)	35,676 35,426	35,551	
5 (1330)	104,713 99,583	102,148	
Unknown Samples			
1	3,999 4,264	4,132	30
2	13,668 11,993	12,831	109
3	41,418 41,557	41,488	387

ACTH IRMA SAMPLE CALIBRATOR CURVE


FIGURE 1

13. LIMITATIONS OF THE PROCEDURE

13.1 Counting times should be sufficient to reduce statistical error (for example, accumulation of 2,000 counts will yield 5% error; 10,000 counts will yield 1% error).

13.2 All ACTH results must be interpreted carefully with the overall clinical presentations and other supportive diagnostic tests.

14. EXPECTED VALUES

DiaSorin Inc. recommends that each testing laboratory establish its own reference (normal) range.

Normal Range

Each laboratory should establish its own normal range. A normal range has been established by DiaSorin using 100 frozen EDTA plasma samples from apparently healthy volunteers (50 females, 50 males). Samples were drawn in non-siliconized tubes between 7-10:00 a.m. from fasting (overnight) donors. The geometric mean was found to be 18.5 pg/mL with a 2 S.D. range of 6.0 to 56.7 pg/mL.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Precision

Within-run and total precision were determined using CLSI Document EP5-T2. Each of three sites assayed four controls in duplicate over twelve operating days.

Within - run Precision (values = pg/mL).

Sample Number	Mean Value	S.D.	% C.V.
1	32.7	1.37	4.20
2	58.0	2.05	3.53
3	194	4.90	2.53
4	773	41.97	5.43

Total Assay Precision (values = pg/mL).

The data was combined from three sites.

Sample Number	Mean Value	S.D.	% C.V.
1	32.7	1.88	5.8
2	62.0	6.15	9.9
3	194	9.04	4.7
4	817	87.6	10.7

Low end precision was further evaluated at DiaSorin as follows:

Between - run Precision (values = pg/mL).

Each sample was assayed in 10 separate assays.

Sample Number	Mean Value	S.D.	% C.V.	n
1	8.65	0.496	5.74	10
2	10.67	0.527	4.94	10
3	257.0	8.126	3.16	10

Within - Assay Precision (values = pg/mL).

Each sample was run as 10 replicates in a single assay.

Sample Number	Mean Value	S.D.	% C.V.	n
1	8.6	0.411	4.78	10
2	10.39	0.437	4.21	10
3	260.7	9.22	3.54	10

15.2 Trueness: The assay trueness has been checked by the linearity test and the recovery test.

Linearity

Two EDTA spiked samples and one EDTA spiked pool were diluted with 0 calibrator. Results are summarized below, (values = pg/mL).

Sample Number	Dilution	Measured	Corrected for dilution
1	Undiluted	>E Cal.	
	1:2	755	1510
	1:4	40	1360
	1:8	158	1264
	1:16	84	1344
	1:32	41	1312
	1:67	21	1344
	1:128	11	1408
2	Undiluted	1239	
	1:2	644	1287
	1:4	336	1344
	1:8	158	1264
	1:16	80	1280
	1:32	37	1184
	1:67	20	1280
	1:128	9	1152
Pool	Undiluted	>E Cal.	
	1:2	713	1426
	1:4	342	1368
	1:8	155	1240
	1:16	78	1248
	1:32	39	1248
	1:67	20	1280
	1:128	10	1380

Accuracy

Representative data from three pooled samples were spiked with ACTH 1-39 to determine the recoveries, (values = pg/mL).

Background Value	ACTH Added	Expected Value	Measured Value	Percent Recovery
Pool 1 17.4	5	22.4	20.4	91.1
	25	42.4	37.6	88.7
	50	67.4	59.8	88.7
	75	92.4	81.1	87.8
	125	142.4	118.5	83.2
	500	517.4	571.7	110.5
			Avg. Recovery	91.7
Pool 2 32.1	5	37.1	36.2	97.6
	25	57.1	52.3	91.6
	50	82.1	71.1	86.6
	75	107.1	92.8	86.6
	125	157.1	128.3	81.7
	500	532.1	583.0	109.6
			Avg. Recovery	92.3
Pool 3 31.2	5	36.2	35.5	98.1
	25	56.2	48.9	87.0
	50	81.2	69.0	85.0
	75	106.2	89.6	84.4
	125	156.2	129.7	83.0
	500	531.2	520.6	98.0
			Avg. Recovery	87.5

15.3 Analytical Sensitivity

Sensitivity, defined as the ACTH concentration measured two standard deviations above the 0 calibrator counts, at minimum binding is 1.5 pg/mL or less, with fresh tracer.

15.4 Specificity and Cross Reactivity

The cross-reactivity of the ACTH IRMA was made with the following peptide concentrations added to both a low and medium calibrator.

Peptide	Dose (pg/mL)	Change in medium Calibrator (pg/mL)	% change	Change in low Calibrator (pg/mL)	% change
ACTH 1 - 26	100,000	-346.7	-98.5%	-24.5	-98.8%
	1,000	-181.2	-51.5%	-13.3	-53.6%
	500	-116.8	-33.2%	-7.4	-29.8%
	200	-47.6	-13.5%	-2.5	-9.9%
ACTH 18 - 24	100,000	0.0	0.0%	0.0	0.0%
	50,000	0.0	0.0%	0.0	0.0%
ACTH 26 - 39	100,000	-349.6	-99.3%	-24.3	-98.0%
	5,000	-318.9	-90.6%	-23.1	-93.1%
	2,000	-287.8	-81.8%	-20.3	-81.9%
	1,000	-244.6	-69.5%	-16.7	-67.3%
ACTH 1 - 10	100,000	0.0	0.0%	0.0	0.0%
	50,000	0.0	0.0%	0.0	0.0%
	1,000	0.0	0.0%	0.0	0.0%
	500	0.0	0.0%	0.0	0.0%
α -MSH	100,000	0.0	0.0%	0.0	0.0%
	1,000	0.0	0.0%	0.0	0.0%
	500	0.0	0.0%	0.0	0.0%
	100,000	0.0	0.0%	0.0	0.0%
β -MSH	100,000	0.0	0.0%	0.0	0.0%
	50,000	0.0	0.0%	0.0	0.0%
β -Endorphin	100,000	0.0	0.0%	0.0	0.0%
	50,000	0.0	0.0%	0.0	0.0%

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
1. Sayers, M.A., G. Sayers and L.A. Woodbury, "The Assay of Adrenocorticotrophic Hormone by the Adrenal Ascorbic Acid-Depletion Method," **Endocrinology**, 42:379-393, (1948).
2. Findings, J.W., W.C. Engeland and H. Raff, "The Use of Immunoradiometric Assay for the Measurement of ACTH in Human Plasma," **TEM**, Elsevier Science Publishing Co., Inc., (1990).
3. **Harrison's Principles of Internal Medicine**, 11th edition, ACTH Physiology.
4. Mains, R.E., B.A. Eipper and N. Ling, "Common Precursor to Corticotropins and Endorphins," **Proceedings of the National Academy of Sciences**, U.S.A., 74: 3014 - 3018, (1977).
5. Orth, D.N., "Adrenocorticotrophic Hormone," **Methods of Hormone Radioimmunoassay**, Second Edition, (1979).
6. Reader, S.C.J., J. Alagband-Zadeh, R. Daly and W.R. Robertson, "Negative, Rate-sensitive Feedback Effects on Adrenocorticotropin Secretion by Cortisol in Normal Subjects," **The Journal of Endocrinology**, 92:443, (1982).
7. Tietz, N.W., eds, **Clinical Guide to Laboratory Tests**, W.B. Saunders Company (Philadelphia), (1990).
8. Krieger, D.T., "Factors Influencing the Circadian Periodicity of ACTH and Corticosteroids," **The Medical Clinics of North America**, 62:251, (1978).
9. Kuhn, J.M., M.F. Proeschel, D.J. Seurin, et. al., "Comparative Assessment of ACTH and Lipotropin Plasma Levels in the Diagnosis and Follow-up of Patients with Cushing's Syndrome: A Study of 210 Cases," **The American Journal of Medicine**, 86:678-684, (1989).
10. Hermus, A.R., G.F. Pieters, A.G. Smals, et. al., "Transition from Pituitary-Dependent to Adrenal-Dependent Cushing's Syndrome," **New England Journal of Medicine**, 318:966, (1988).
11. Landolt, A., A. Valavanis, J. Girard and An. Eberle, "Corticotrophin-Releasing Factor Test Used with Bilateral, Simultaneous Inferior Petrosal Sinus Blood Sample for the Diagnosis of Pituitary-Dependent Cushing's Disease," **Clinical Endocrinology**, 25:687, (1986).
12. Findling, J.W. and J.B. Tyrrell, "Occult Ectopic Secretion of Corticotropin," **Archives of Internal Medicine**, 146:929, (1986).
13. Singer, W., N. Kovacs and E. Horvath, "Ectopic ACTH Syndrome: Clinicopathological Correlations," **The Journal of the Association of Clinical Pathologists**, 31:591, (1978).
14. Yalow, R.S., "Ectopic ACTH in Carcinoma of the Lung," **Progress in Cancer Research and Therapy**, eds F.M. Muggia and M. Rozenzweig, Raven Press (New York), 11:209, (1979).
15. Nelson, D.H., J.W. Meakin and G.W. Thorn, "ACTH-Producing Tumors Following Adrenalectomy for Cushing's Syndrome," **Annals of Internal Medicine**, 52:560 - 569, (1960).
16. Krieger, D.T. and M. Luria "Plasma ACTH and Cortisol Responses to TRF, Vasopressin or Hypoglycemia in Cushing's Disease and Nelson's Syndrome," **The Journal of Clinical Endocrinology Metabolism**, 44:361-368, (1977).
17. Tonks, D.B., "A Study of the Accuracy and Precision of Clinical Chemistry Determinations in 170 Canadian Laboratories," **Clinical Chemistry**, 9: 217, (1963).
18. Rodbard, D., et. al., "Statistical Quality Control of Radioimmunoassays," **The Journal of Clinical Endocrinology and Metabolism**, 28:1412, (1968).
19. **National Committee for Clinical Laboratory Standards**, Internal Quality Control Testing: Principles and Definitions, Approved Guideline, CLSI document C24-A (ISBN 1-56238-112-1), NCCLS, 771 East Lancaster Avenue, Villanova, PA, 19085, (1991).

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Solid phase (Coated beads.)	Phase solide (Billes enduites).	Festphase (Beschichtete Kügelchen).	Fase sólida (Perlas recubiertas).	Fase solida (Perle con rivestimento).
	Wash buffer	Tampon de lavage	Waschpuffer	Tampón de lavado	Tampone di lavaggio
	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Trazador: antígeno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo	Nocivo


SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
+44(0) 1344 401 430
FAX: +44(0) 7884 050812

12715

34690
7/10

PRINTED IN U.S.A.