ALDOCTK-2
(P2714)

DiaSorin

English.....	p. 1
Italiano.....	p. 11
Français.....	p. 21
Deutsch.....	p. 32
Español.....	p. 42
Português.....	p. 52
Svenska.....	s. 62
Dansk.....	s. 71
Magyar.....	s. 81
Česky.....	p. 91
Ελληνικά.....	p. 101


ALDOSTERONE RADIOIMMUNOASSAY KIT

Procedure for quantitative determination of aldosterone in human serum, plasma or urine samples

For in vitro use only

1. INTRODUCTION

Aldosterone is a steroid hormone of molecular weight 360.4 daltons, having the following chemical formula:

Aldosterone is secreted by the zona glomerulosa of the adrenal cortex. It is by far the most potent regulator of electrolyte excretion.

Aldosterone binds to specific protein receptors in the responsive renal cells, where it probably induces the DNA-directed synthesis of mRNA with subsequent formation of new protein. The role of this protein remains controversial, but it may either facilitate entry of sodium into the transporting cells or provide these cells with an increased energy supply for active sodium transport. In either case, its net effect is to enhance sodium transport (reabsorption).

In the distal nephron, the action of aldosterone is manifested by increased reabsorption of sodium and chloride ions from the tubular lumen and increased secretion of potassium and hydrogen ions into the tubular lumen. The reabsorption of sodium and chloride tends to increase the osmolality of extracellular fluid. An increase in extracellular osmolality stimulates the secretion of ADH and ADH facilitates the renal conservation of water. Sodium, chloride, and water are commonly retained or excreted together, and aldosterone can be said to promote, indirectly, the renal tubular reabsorption of water.

The mechanism regulating the secretion of aldosterone responds to changes in body hydration. The renin-angiotensin system is the principal mediator of the adrenocortical response to changes in body hydration. Conditions that reduce extracellular fluid volume, and conditions that result in sequestration of venous blood, all tend to stimulate the production of renin, thereby increasing the formation of angiotensin. Angiotensin not only acts directly on the vascular system as a pressor agent, but also stimulates the adrenal cortex to secrete aldosterone.

Aldosterone secretion is also influenced by:

1. Potassium concentration in body fluids: potassium depletion decreases, while potassium retention increases aldosterone secretion.
2. ACTH secretion.
3. Sodium concentration in body fluids: sodium depletion increases, while sodium retention decreases aldosterone secretion.

The assay of aldosterone may be carried out both in serum or plasma and in 24-hour urine. From a clinical point of view it must be remembered that, while the 24-hour urinary measurement of 18-oxo-conjugate of aldosterone is an integrated reflection of the daily aldosterone secretion, the plasma values may reflect only a single point in time and since plasma aldosterone shows typical bursts, which follow a circadian rhythm, it might not be safe to draw valid conclusions from a single determination. Therefore, while plasma aldosterone measurement is more suitable for acute studies (like circadian rhythms, postural changes, acute effect of drugs), the 24-hour urinary measurement is the method of choice for clinical investigation.

2. PRINCIPLE OF THE ASSAY

The principle of the assay is based on the competition between labelled aldosterone and aldosterone contained in calibrators or samples to be assayed for a fixed and limited number of antibody binding sites. After the incubation, the amount of labelled aldosterone bound to the antibody on the tube walls is inversely related to the concentration of unlabelled aldosterone present in calibrators or samples. The method adopted for B/F separation is based on the use of antibody-coated tubes, where the antibody is coated on the tube walls.

3. REAGENTS PROVIDED IN THE KIT

Coated tubes	100
¹²⁵ I-labelled aldosterone	1 bottle
Aldosterone calibrators	6 vials
Control serum	1 vial
Number of tubes	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 2 assay runs if reagents are stored as the manufacturer recommends.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming.

Reagents from different batches must not be mixed.

3.1. Coated tubes

The inner surface of each tube is coated with aldosterone antiserum raised in rabbits.

Before use, bring the coated tubes to room temperature prior to opening the box, to avoid condensation of humidity.

Securely reseal the box containing unused tubes. Do not mix different batches of coated tubes.

3.2. ¹²⁵I-labelled aldosterone (red): ready-to-use reagent

The bottle contains 52 mL aldosterone labelled with ¹²⁵I, steroid-free human serum, phosphate buffer, preservatives and an inert red dye. Radioactivity is 85 kBq (2.3 µCi) or less on the calibration date.

3.3. Aldosterone calibrators

The vials contain increasing amounts of aldosterone, steroid-free human serum and preservatives. The zero calibrator vial contains 3 mL solution.

Reconstitute the contents of the 1-5 calibrator vials with 1 mL distilled water. The resulting solutions contain 50 - 100 - 250 - 500 - 1000 pg/mL aldosterone (0.138 - 0.277 - 0.692 - 1.385 - 2.770 nmol/L) and are stable for one week at 2-8°C or stored in deep-frozen aliquots (-20°C or below) for extended storage.

As no international standard preparation is currently available, the kit calibrators are referenced to an internal reference preparation (99% purity by HPLC). *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends.*

3.4. Control serum: lyophilized reagent

The vial contains aldosterone, human serum and preservatives. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Reconstitute the vial contents with 1 mL distilled water. The resulting solution is stable for one week at 2-8°C or stored in deep-frozen aliquots (-20°C or below) for extended storage.

4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Laboratory reagents (HCl, boric acid) for urine assay.
- Glassware.
- Disposable glass tubes.
- Test tube rack.
- Micropipettes with disposable tips (10, 200, 500 µL) (10 µL: trueness ± 3%, precision 2%; 200, 500 µL: trueness ± 2%, precision 1%).
- Vortex mixer.
- Thermostatically-controlled water bath or heating block capable of maintaining 37° ± 1°C.
- Device for aspiration of incubation mixture.
- Gamma counter suitable for counting ¹²⁵I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

5. SPECIMEN COLLECTION AND PREPARATION

Careful standardization of the patient's preparation and sampling conditions is strongly recommended.

Serum or plasma samples

The anticoagulants potassium EDTA and lithium heparin have been tested and may be used with this assay (see §9.1). Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

If aldosterone levels greater than 1000 pg/mL are expected, the samples should be diluted with the zero calibrator.

Urine samples

Collect 24-hour urine, measure and record the volume. Mix well before withdrawing an aliquot to be assayed. Add 1 g boric acid/100 mL and store at 2-8°C for one day or in deep-frozen aliquots (-20°C or below) for extended storage.

Urine hydrolysis

Aldosterone can be assayed in urine samples after acid hydrolysis of aldosterone 18-glucuronide. In the conditions recommended for the assay, hydrolysis is complete.

- **Mix 100 µL urine and 1 mL 0.1N HCl.** Use glass tubes.
- Cap tubes and **incubate overnight (18-22 hours) at 30 ± 2°C.**

6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each calibrator and sample.

- Dispense reagents *in the bottom of coated tubes*. Operate according to the following scheme:

reagents \ tubes	Zero calibrator	Calibrators 1-5	Serum samples	Urine samples
Zero calibrator	200 µL	–	–	200 µL
Calibrators 1-5	–	200 µL	–	–
Serum samples	–	–	200 µL	–
Urine samples	–	–	–	10 µL*
Tracer	500 µL	500 µL	500 µL	500 µL

* 20 µL hydrolyzed urine should be used for urine samples that are expected to contain low aldosterone levels.

- **Mix** the contents of tubes with a Vortex and **incubate overnight** (18-22 hours) **at room temperature**.
- Carefully **aspirate** the incubation mixture. *Be sure that the aspirator tip touches the bottom of the coated tube so that all the liquid is removed. Failure to remove adhering solution adequately may result in poor reproducibility and spurious results. No trace of dye should still be visible.*
- **Measure the radioactivity** of tubes.

7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Compute the B/Bo ratio for each calibrator and unknown sample as follows:

$$B/Bo\% = \frac{\text{calibrator or sample mean counts}}{\text{zero calibrator mean counts}} \times 100$$

Plot in semilog or linear-linear coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of aldosterone concentration expressed as pg/mL or nmol/L on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1).

Serum or plasma samples

Read the aldosterone concentration of each sample expressed as pg/mL or nmol/L directly from the calibration curve. If the sample was diluted, the aldosterone concentration value derived from the diluted specimen must be multiplied by the dilution factor.

Urine samples

Daily aldosterone excretion expressed as µg or nmol/24 hours is calculated by the following procedure.

- *For a 10-µL sample volume of hydrolyzed urine*
Multiply the value read from the calibration curve by 220 to obtain urinary aldosterone concentration in pg/mL or nmol/L. The coefficient of 220 derives from the dilution factor of hydrolyzed urine (11), multiplied by the calibrator-to-sample ratio (20).
- *For a 20-µL sample volume of hydrolyzed urine*
Multiply the value read from the calibration curve by 110 to obtain urinary aldosterone concentration in pg/mL or nmol/L. The coefficient of 110 derives from the dilution factor of hydrolyzed urine (11), multiplied by the calibrator-to-sample ratio (10).

To obtain the daily aldosterone excretion concentration, multiply the aldosterone concentration by the volume of urine excreted (expressed in mL) using the following formulae:

$$\begin{aligned} \mu\text{g}/24 \text{ hours} &= \text{pg/mL} \times \text{mL urine excreted} \times 10^{-6} \\ \text{nmol}/24 \text{ hours} &= \text{nmol/L} \times \text{mL urine excreted} \times 10^{-3} \end{aligned}$$

Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/Bo x 100
Zero calibrator	10,220	100
50 pg/mL - 0.138 nmol/L	8,534	83.5
100 pg/mL - 0.277 nmol/L	7,409	72.5
250 pg/mL - 0.692 nmol/L	5,243	51.3
500 pg/mL - 1.385 nmol/L	3,352	32.8
1,000 pg/mL - 2.770 nmol/L	1,942	19.0
Serum sample	6,805	66.8
Urine sample	8,262	81.0


Serum sample

By interpolation from the calibration curve, the sample is found to contain 140 pg/mL (0.388 nmol/L) aldosterone.

Urine sample

By interpolation from the calibration curve, the sample is found to contain 70 pg/mL (0.194 nmol/L) aldosterone. Taking into consideration the factor 220 and the volume of urine excreted during 24 hours (e.g. 1100 mL), the daily aldosterone excretion is calculated as follows:

$$\begin{aligned} \mu\text{g}/24 \text{ hours} &= 70 \times 220 \times 1100 \times 10^{-6} = 15400 \times 1100 \times 10^{-6} = 16.94 \\ \text{nmol}/24 \text{ hours} &= 0.194 \times 220 \times 1100 \times 10^{-3} = 42.68 \times 1100 \times 10^{-3} = 46.95 \end{aligned}$$

ALDOSTERONE

Fig. 1

8. EXPECTED VALUES

The ranges given below are merely indicative for a normal population with a sodium intake of 100-150 mEq/24 hours: each laboratory should establish its own reference ranges.

Normal subjects	PLASMA, pg/mL	URINE, µg/24 hours
Supine position	7.5 - 150	2.8 - 30.0
Upright position	35 - 300	

Conversion to nmol aldosterone/L is possible using the following formula:

$$\text{nmol aldosterone/L} = \text{pg aldosterone/mL} \times 2.77 \times 10^{-3}$$

9. SPECIFIC PERFORMANCE CHARACTERISTICS

9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by lipaemia (up to 500 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin). Aldosterone values determined in samples containing 200 mg/dL haemoglobin are about 10% higher than those determined in normal samples.

Matched serum and plasma samples (the latter obtained with lithium heparin and potassium EDTA) were drawn from 41 healthy volunteers to set up an anticoagulant interference study. All samples were tested by the standard method and the following results were observed.

#	Sample	Serum	Heparin	EDTA
41	mean (pg/mL)	181.6 (100%)	183.9 (101.3%)	197.1 (108.5%)

The results of those tests showed that aldosterone values determined in potassium EDTA plasma samples were about 10% higher than those determined in serum or lithium heparin plasma samples.

Cross-reactions. The percentage of cross-reactions, calculated according to Abraham, shows the specificity of the antiserum used.

- Aldosterone	100%
- 3 beta, 5 alpha-tetrahydroaldosterone	4.2%
- 3 alpha, 5 beta-tetrahydroaldosterone	0.3%
- Cortisolone	$6.4 \times 10^{-3}\%$
- 18-OH-deoxycorticosterone	$5.0 \times 10^{-3}\%$
- Corticosterone	$4.0 \times 10^{-3}\%$
- Cortisone	$2.5 \times 10^{-3}\%$
- Deoxycorticosterone	$1.0 \times 10^{-3}\%$
- Testosterone	$8.0 \times 10^{-4}\%$
- Prednisolone	$6.4 \times 10^{-4}\%$
- Cortisol	$1.0 \times 10^{-4}\%$
- Progesterone	$8.0 \times 10^{-5}\%$
- Spironolactone	$6.5 \times 10^{-5}\%$
- Dehydroepiandrosterone	$6.0 \times 10^{-5}\%$
- Pregnenolone	$5.0 \times 10^{-5}\%$
- Pregnanetriol	$3.3 \times 10^{-5}\%$
- Dexamethasone	$3.3 \times 10^{-5}\%$
- Estradiol	$1.0 \times 10^{-6}\%$

9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 10 pg/mL (0.0275 nmol/L) at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator, that is, two standard deviations below zero.

9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	10	10	10
Mean (pg/mL)	86.4	249.9	495.6
Standard deviation	4.60	9.50	8.60
Coefficient of variation (%)	5.3	3.8	1.7

Reproducibility	A	B	C
Number of determinations	15	15	15
Mean (pg/mL)	91.3	230.3	607.5
Standard deviation	5.72	7.90	26.40
Coefficient of variation (%)	7.0	3.4	4.3

9.4. Trueness

The assay trueness has been checked by the dilution and recovery tests.

Dilution test. Two sera with high aldosterone concentration were tested after serially diluting with the zero calibrator.

Dilution	Expected concentration, pg/mL	Measured concentration, pg/mL	% Recovery
neat	–	706.0	–
1:2	353.0	360.5	102.1
1:4	176.5	167.7	95.0
1:8	88.3	92.7	105.0
1:16	44.1	43.0	97.5
neat	–	355.0	–
1:2	177.5	165.0	93.0
1:4	88.8	92.7	104.5
1:8	44.4	46.0	103.7
1:16	22.2	25.7	115.8

Recovery test. Two sera containing aldosterone were tested as such and after mixing with increasing amounts of aldosterone.

Added concentration, pg/mL	Expected concentration, pg/mL	Measured concentration, pg/mL	% Recovery
–	–	147.0	–
520	659.7	713.0	108.1
300	439.7	469.0	106.7
180	319.7	305.0	95.4
40	179.7	177.0	98.5
–	–	291.0	–
520	796.5	815.0	102.3
300	576.5	573.0	99.4
180	456.5	435.0	95.3
40	316.5	321.0	101.4

10. LIMITATIONS OF THE PROCEDURE

The clinical significance of aldosterone determination is invalidated when the subjects are not kept under controlled conditions of sodium and potassium intake and controlled conditions of posture or when drugs are administered (e.g., diuretics, clonidine, betablocking agents, reserpine, alpha-methyldopa, etc.).

Diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration of incubation mixture
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

11. WARNINGS AND PRECAUTIONS

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. However, as no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

13. BASIC RULES OF RADIATION SAFETY

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 2.17 μCi (81 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

SCHEME OF THE ASSAY

- 1 - DIRECTLY ASSAY SERUM OR PLASMA SAMPLES.
HYDROLYZE URINE SAMPLES WITH 0.1N HCl OVERNIGHT AT 30°C.
- 2 - RECONSTITUTE REAGENTS. IDENTIFY COATED TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS ACCORDING TO THE FOLLOWING SCHEME AND MIX THE INCUBATION MIXTURE:

REAGENTS \ TUBES	0 CAL	CAL 1-5	SERUM SAMPLES	URINE SAMPLES
ZERO CALIBRATOR	200 µL	–	–	200 µL
CALIBRATORS 1-5	–	200 µL	–	–
SERUM SAMPLES	–	–	200 µL	–
URINE SAMPLES	–	–	–	10 µL*
TRACER	500 µL	500 µL	500 µL	500 µL

- 4 - INCUBATE OVERNIGHT AT ROOM TEMPERATURE.
- 5 - ASPIRATE THE REACTION MIXTURE.
- 6 - MEASURE THE RADIOACTIVITY OF TUBES.

* 20 µL in case of low aldosterone levels.

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In: «Handbook of radioimmunoassay», G.E. Abraham ed., M. Dekker Inc., New York, p. 591 (1977).
DiaSorin Inc.
Stillwater, MN 55082-0285 U.S.A.
DiaSorin S.p.A.,