**GammaCoat™ Cortisol**  
**<sup>125</sup>I RIA Kit**

For the quantitative determination  
of cortisol levels in serum, plasma or urine

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### 1. INTENDED USE

FOR IN VITRO DIAGNOSTIC USE.

The GammaCoat™ [<sup>125</sup>I] Cortisol Radioimmunoassay Kit is to be used for the quantitative determination of cortisol (hydrocortisone, Compound F) levels in serum, plasma or urine.

### 2. SUMMARY AND EXPLANATION

Cortisol (hydrocortisone, Compound F) is the major glucocorticoid which is secreted by the adrenal cortex. Among its physiological effects are anti-inflammatory activity, blood pressure maintenance, and synthesis of carbohydrate from protein.

Glucocorticoids are synthesized in response to adrenocorticotrophic hormone (ACTH) which is secreted by the anterior lobe of the pituitary gland. In the normal individual, cortisol participates in a negative feedback loop with ACTH through the hypothalamus-pituitary-adrenal cortex axis (HPA axis). Pituitary ACTH is, in turn, regulated by corticotropin-releasing factor (CRF) which is secreted by the hypothalamus. The hypothalamic CRF is responsive to cortisol levels and to stress. Physical and psychological stress, diurnal variation and low blood sugars will affect the rate of cortisol secretion.<sup>1-6</sup>

Malfunction of any organ in the HPA axis will result in alteration of cortisol levels.<sup>7,8</sup> The evaluation of adrenal gland function by measurement of plasma or serum cortisol levels aids in the diagnosis of normal and abnormal states. Combinations of morning and evening measurements and stimulation and suppression tests can offer strong evidence for diagnosis of specific adrenal-related diseases. The diagnosis of Addison's disease (chronic adrenal insufficiency) and Cushing's syndrome (adrenal overproduction) are two examples for which the specific assay of circulating cortisol levels is helpful in patient evaluation.<sup>9,10</sup>

### 3. PRINCIPLES OF THE ASSAY

The GammaCoat [<sup>125</sup>I] Cortisol Radioimmunoassay Kit procedure is based on the competitive binding principles of radioimmunoassay.<sup>11</sup> Calibrators and unknown samples (serum, plasma or urine) are incubated with cortisol tracer in antibody-coated tubes where the antibody is immobilized onto the lower inner wall of the GammaCoat tube. After incubation, the contents of the tube are aspirated or decanted and the tube is counted. A calibrator curve is prepared with 5 serum calibrators ranging from 1-60 µg/dL. Unknown values are interpolated from the calibrator curve. The entire assay is performed in the coated tube. No separate sample dilution, protein denaturation or extraction steps are required for serum, urine or plasma samples.

In addition to unextracted urinary cortisol assays, the GammaCoat [<sup>125</sup>I] Cortisol RIA Kit procedure may also be used to measure urinary free (unconjugated) cortisol. The cortisol present in the urine must first be extracted into methylene chloride. An aliquot of the extract is evaporated directly in the GammaCoat antibody-coated tube. (Experiments have shown that no deleterious effects occur to the antibody-coated tube if good quality methylene chloride is used.) To insure that the samples contain proper levels of protein, 10 microliters of Cortisol Serum Blank (CA-2367) are added to each GammaCoat tube. The remainder of the assay procedure is the same as that used for serum, plasma or unextracted urine cortisol determinations. The efficiency of the cortisol extraction from urine by methylene chloride should be monitored.

### 4. REAGENTS PROVIDED IN THE KIT

Cortisol Calibrators (A-F)	6 vials /1 mL	6 vials /1 mL
Cortisol Coated Tubes	2 bags /50 ea	10 bags / 50 ea
Cortisol Tracer	1 vial / 10 mL	5 vials/10 mL ea
PBS Buffer	1 vial/ 100 mL	5 vials / 100 mL ea
Number of tests	100	500

STORAGE: All reagents including the tracer-buffer reagent are stable through the expiration date printed on the kit label when stored at 2-8°C. The anti-cortisol coated tubes can be stored at 2-27°C. Do not use reagents beyond the kit expiration date. Reagents from different batches must not be mixed.

**4.1 Cortisol Tracer:** ready to use reagent

Each vial contains approximately 4 µCi of tracer (<1 µg/mL Cortisol) in 10 mL of phosphate buffered saline, ANS (8-Anilino-1-naphthalene sulfonic acid) with 0.1% sodium azide as a preservative.

**4.2 Rabbit Anti-Cortisol Serum Coated Tubes:** ready to use reagent

12 x 75 mm tubes are coated with rabbit anti-cortisol serum (titer<1 µg/tube) and packaged in a plastic bag.

**4.3 PBS Buffer:** ready to use reagent

Each bottle contains 100 mL of phosphate buffered saline with 0.1% sodium azide as a preservative.

**4.4 Cortisol Serum Calibrators:** ready to use reagent

Each vial contains cortisol in 1 mL of processed human serum with 0.1% sodium azide as a preservative. Calibrators are provided at the following concentrations: 0, 1, 3, 10, 25 & 60 µg/dL, respectively. The DiaSorin Cortisol Calibrators have been calibrated against the U.S.P. Cortisol Reference Standard. Any comparison with other products or procedures should be done with this reference standard. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

**4.5 Auxiliary Label**

For tracer-buffer reagent

**4.6 AVAILABLE UPON REQUEST:**

Semi-Logarithmic Graph Paper  
50 sheets per pad

**5. WARNINGS AND PRECAUTIONS**

FOR IN VITRO DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

**REAGENTS CONTAINING HUMAN SOURCE MATERIAL**

**Treat as potentially infectious.**

Each serum/plasma donor unit used in the preparation of this product has been tested by an FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4<sup>TH</sup> ed., May 1999 or current edition.

**Reagents Containing Sodium Azide**

**CAUTION:** Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

**European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

**Reagents Containing Iodine-125**

Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

**6. SPECIMEN COLLECTION AND HANDLING**

**Plasma Samples**

Collect a venous blood sample aseptically in an evacuated glass tube containing heparin or EDTA.

Separate the plasma fraction. The plasma may be stored refrigerated overnight at 2-8°C. For longer storage, samples must be stored frozen at -20°C. Hemolyzed, lipemic or icteric samples should be avoided.

**Serum Samples**

Although plasma samples are usually specified for cortisol assays, comparable results are obtained using serum samples.<sup>12</sup>

Separate the serum fraction. The serum may be stored refrigerated overnight at 2-8°C. For longer storage, samples must be stored frozen at -20°C. The use of hemolyzed, lipemic or icteric samples should be avoided.

**Urine Samples**

Collect a 24-hour urine specimen. Record the total volume of urine. Add 10 grams of boric acid per liter of urine as a preservative, and store a well-mixed aliquot frozen at -20°C. Thaw frozen samples at room temperature and mix well prior to use.

## **7. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED.**

- 7.1 Constant temperature water bath,  $37 \pm 2^\circ\text{C}$ .
- 7.2 Precision pipet (10  $\mu\text{L}$ ).
- 7.3 Semi-automatic pipet (1.0 mL).
- 7.4 Gamma counter.
- 7.5 Test tube rack.
- 7.6 Water aspirator (optional).
- 7.7 Vortex mixer.
- 7.8 Plastic-backed absorbent paper (optional).

Additional supplies for use in urine extraction procedure only:

- 1. Disposable glass test tubes (12 x 75 mm) and caps.
- 2. Methylene chloride (dichloromethane; ACS analytical reagent grade).
- 3. Ice bath.
- 4. Precision pipets (100 & 200  $\mu\text{L}$ ).
- 5. Hood or well-ventilated area.
- 6. Means of evaporating the methylene chloride. This can be done with a stream of dry air or nitrogen with mild warming of the extract to about  $37^\circ\text{C}$  in a hood or well-ventilated area.

## **8. URINE EXTRACTION PROCEDURE**

This procedure includes the methodology for determining the efficiency of methylene chloride extraction. Typical recoveries of cortisol from urine samples averaged  $99 \pm 14\%$  in DiaSorin laboratories.

- 8.1 Place an adequate volume of methylene chloride in an ice bath to chill. Allow more than 1 mL per tube. The ice bath should be large enough to accommodate a beaker of methylene chloride and the rack of glass extraction test tubes.
- 8.2 To determine the extraction efficiency, label two disposable glass test tubes for each sample to be assayed with the patient's identification. Label one tube "S" for Sample and the other tube "E" for Extraction Efficiency.
- 8.3 Add 200 microliters of a well-mixed 24 hour urine sample to the appropriate "S" and "E" tubes.
- 8.4 To each "E" tube only, add 10  $\mu\text{L}$  of the 60  $\mu\text{g}/\text{dL}$  Cortisol Serum Calibrator. Mix reagents by shaking the test tube rack or mixing each tube.
- 8.5 Add 1.0 mL chilled methylene chloride to each "S" and "E" tube. Vortex well to extract the cortisol into the organic layer. Cap or cover and place in an ice bath to avoid excessive evaporation of methylene chloride. To facilitate pipetting the methylene chloride, it is important to rinse the pipet 3 to 4 times with COLD methylene chloride before transferring it to the appropriate tubes.
- 8.6 While waiting for phase separation to occur, label two GammaCoat Cortisol tubes for each "S" tube and two GammaCoat Cortisol tubes for each "E" tube.
- 8.7 After phase separation, carefully remove 100  $\mu\text{L}$  from the bottom (methylene chloride) layer of each tube and add to the appropriately labeled GammaCoat tubes. Rinsing the pipet tip in cold methylene chloride prior to pipetting the sample will minimize seepage of extracted sample from the pipet tip. When transferring the 100  $\mu\text{L}$  aliquot of extraction solution to the coated tube, do not attempt to wipe the pipet tip, because the physical properties of methylene chloride cause leakage and air bubbles in the pipet tip.

To ease the transfer of extraction solution to the coated tube, hold both the extraction tube and the GammaCoat tube in the same hand. Rapidly, but carefully, transfer the 100  $\mu\text{L}$  of extract into the bottom of the GammaCoat tube.

- 8.8 Evaporate the methylene chloride from the GammaCoat tubes with a stream of dry air or nitrogen and mild warming of the extract to about 37°C.  
**CAUTION:** Avoid excessive concentrations of halogenated hydrocarbon fumes in work areas. EVAPORATION SHOULD BE DONE ONLY IN A HOOD OR WELL-VENTILATED AREA.
- 8.9 After the methylene chloride is completely evaporated, add 10 µL of the 0 µg/dL Cortisol Serum Blank, CA-2367, to each GammaCoat patient sample tube ("S" and "E" tubes), to ensure that the samples contain the proper levels of protein.
- 8.10 Continue with the RIA procedure starting with step 1 of Assay Procedure.

## 9. ASSAY PROCEDURE

### Tracer-Buffer Reagent:

Add the entire contents of 1 vial of cortisol tracer to 1 bottle of PBS buffer and mix gently. A final rinse of the cortisol tracer vial with the tracer-buffer reagent will maximize transfer of the tracer. Fill out and attach the auxiliary tracer-buffer reagent label to the container for proper identification. The tracer-buffer reagent must be stored in the dark, or in a foil-wrapped container, at 2-8°C.

The assay procedure includes the preparation of a calibrator curve from which the unknown cortisol content in the sample is interpolated.

- 9.1 Allow all reagents to reach room temperature and mix thoroughly without foaming prior to use.
- 9.2 Label the GammaCoat tubes in duplicate according to the following scheme. Total Count [T1,T2] and B0 tubes [1,2] may be required for certain data reduction and quality control programs, but may be omitted if the calibrator curve is plotted on semi-logarithmic graph paper.<sup>13-15</sup>
- 9.3 Add to the bottom of the appropriate duplicate GammaCoat tubes:
- 10 µL of Cortisol Serum Blank, or each Cortisol Serum Calibrator.
  - FOR SERUM, PLASMA AND UNEXTRACTED URINE SAMPLES, add 10 µL of each sample. Care must be taken to pipet calibrators and samples in the same manner.<sup>16</sup>
  - FOR EXTRACTED URINE SAMPLES, use the GammaCoat tubes containing the extracted samples from step 9 of Urine Extraction Procedure.  
 Serum Blank provide with this kit as diluent. With either unextracted or extracted urine samples, dilute as required with distilled water. For Urinary Free Cortisol, extract a 200 µL aliquot of the diluted urine sample (Urine Extraction Procedure, Step 3). Continue with the routine assay. The dilution factor must be included in the calculations.
- 9.4 Add 1.0 mL of tracer-buffer reagent to each tube. Gently mix each tube on a vortex mixer set at a low speed.
- 9.5 Incubate tubes for 45 minutes in a 37 ± 2°C water bath.
- 9.6 Aspirate or decant all tubes except Total Counts.  
**FAILURE TO REMOVE ADHERING SOLUTION ADEQUATELY MAY RESULT IN POOR REPLICATION AND SPURIOUS VALUES.**  
 If the aspiration technique is used, be sure that the plastic tip of the aspirator touches the bottom of the coated tube and that all the liquid is removed.  
 If the decanting technique is used, allow the tubes to drain in an inverted position for 3-5 minutes. Tap the tubes on absorbent paper to remove any adhering liquid before placing the tubes upright.
- 9.7 Count all tubes in a gamma counter for 1 minute with the window suitably adjusted for iodine-125.
- 9.8 Calculate results. Refer to Results.

## 10. SALIVARY PROCEDURE

This kit also allows assays of saliva samples.

It has been shown that there is a good correlation between the concentration of cortisol in saliva and the concentration of free cortisol in plasma and serum (these concentrations represent 60% to 100% of the concentration of the free fraction in circulation and 4% to 6% of the total concentration in serum or plasma).

Numerous studies have demonstrated the practical and diagnostic advantages of salivary cortisol assays and confirmed their validity and clinical value.

Two important points were taken into account in adapting this technique to salivary cortisol assays:

First, because salivary cortisol represents only 5% of total plasma cortisol, a 20-fold increase in sensitivity is required. To attain this level of sensitivity, the saliva sample was increased to a volume of 200  $\mu$ L.

Second, saliva has a lower concentration of proteins (generally consisting of enzymes and glycoproteins). Only trace amounts of plasma proteins are present, while proteins that bind with corticosteroids are absent. To guarantee that the effects of the protein matrix are similar in the sample and standards, 0 calibrator (without cortisol) is added to each saliva sample (as for assays on extracted urine samples).

Thus, only minor modifications are required in the pipetting of samples and standards to adapt the method to saliva sample assays, using only reagents already provided in the kit.

### 10.1 Collecting and preparing saliva samples

Saliva must be collected 10 minutes after patients rinse their mouths with tap water. Saliva production should not be stimulated and the patient must be told the difference between saliva and mucus. A sufficient volume of 1 mL of saliva can be collected in a few minutes, even from children.

The samples are stored frozen at  $-20^{\circ}\text{C}$ . Before the assay, they must be defrosted and centrifuged; the clear supernate is then pumped out for the assay.

### 10.2 Salivary Procedure


For each **saliva sample**, pipette 10  $\mu$ L of 0 calibrator in pairs into the labeled tubes; add 200  $\mu$ L of the sample and 1 mL of tracer.

For the **standard range**, pipette 10  $\mu$ L of 0 calibrator or the standard in pairs into the labeled tubes; add 200  $\mu$ L of PBS buffer, plus 1 mL of tracer

Gently swirl each tube. Incubate for 45 min in a water bath at  $37^{\circ}\text{C}$ . Pump out or decant the content of the tubes and measure the radioactivity bound to the tube sides for 1 minute.

Standard curve

The actual standard range for salivary cortisol is 0.10 – 0.25 – 0.50 – 1.25 and 3.0 μg/dL. The standard curves for the two assay techniques (serum and salivary) covered by this kit are illustrated in the chart to the left.

Although the incubation volume is somewhat higher in the salivary assay than the serum assay (1,210 μL instead of 1,010 μL), this has no effect on the taggant's binding ability within the standard curve's range.

Figure 3  
Standard curves for assay of cortisol in serum or plasma (●) and in saliva (○) with the GammaCoat™ [<sup>125</sup>I] RIA kit.

Note: This curve was obtained during a study conducted in 1984 with a standard range of 0, 2.0, 5.0, 10, 25 and 60 μg/dL. The standard range has since been modified. Please refer to the new range for salivary cortisol calculations, taking into account a ratio of 1:20.

Normal values

Table 1

Normal morning cortisol levels in saliva and serum (μg/dl)						
Group of subjects*	Saliva			Serum		
	n	Median	Range	n	Median	Range
Normal adults	150	0.43	0.18 – 1.01	103	6.5	4.5 – 21.1
Children (4 – 11 years old)	105	0.40	0.14 – 0.91	---	---	---
Pregnant women (3 <sup>rd</sup> trimester)	84	0.47	0.33 – 0.91	84	15.7	12.9 – 39.9

\* Living in the UK

11. QUALITY CONTROL

Each laboratory should utilize controls at several levels to monitor assay performance. The controls should be treated as unknowns. Quality control charts should be maintained to follow the performance of the controls. Appropriate statistical methods should be used to evaluate trends. Acceptable performance limits should be ascertained for the individual laboratory.

## 12. CALCULATION OF RESULTS

Unknown serum, plasma and urine samples may be assayed simultaneously and are all interpolated from a calibrator curve generated with Cortisol Serum Calibrators. Both unextracted urinary cortisol samples (UC) and extracted urinary free cortisol samples (UFC) are then converted to "µg/24 hours" with the appropriate formula.

- 12.1 Record the counts per minute (CPM) bound for each tube.
- 12.2 Plot CPM bound for cortisol calibrators (vertical axis) versus cortisol concentration (horizontal axis) on semi-logarithmic graph paper.
- 12.3 Draw the best fitting curve. For typical data and graph, refer to TABLE II and FIGURE 1.
- 12.4 Locate CPM bound for each tube corresponding to each sample (controls, serum, plasma, unextracted and extracted urine) on the vertical axis and follow a horizontal line intersecting the calibrator curve. At the point of intersection, read cortisol concentration (µg/dL) on the horizontal axis.

### 12.5 Unextracted Urinary Cortisol

To convert the urinary cortisol (UC) of unextracted urine samples from µg/dL to µg cortisol per 24 hours, use the following equation:

$$UC = \frac{S \times V}{100}$$

where:

- S = unextracted urine sample concentration in µg/dL  
V = total urine volume in mL/24 hours  
100 = conversion factor

### 12.6 Extracted Urinary Free Cortisol

To convert the urinary free cortisol (UFC) of extracted urine samples from µg/dL to µg cortisol per 24 hours, use the following equation:

$$UFC = \frac{S \times V}{33.3 (E - S)}$$

where:

- S = extracted urine sample concentration in µg/dL  
E = efficiency tube concentration for corresponding sample in µg/dL  
V = total urine volume in mL/24 hours

33.3 = combined dilution and extraction efficiency factor

Example:

- E = 13.50 µg/dL  
S = 7.85 µg/dL  
V = 1150 mL/24 hours

$$UFC = \frac{(7.85) (1150)}{33.3 (13.50 - 7.85)} = 48 \text{ } \mu\text{g cortisol/24 hours}$$

For labs wishing to monitor the efficiency of the urine extraction procedure (EE), the following formula can be used to calculate the efficiency as a percentage.

$$EE = \frac{E-S}{6.0} \times 100$$

where:

- S = Extracted urine sample concentration in µg/dL
- E = Efficiency tube concentration for corresponding sample in µg/dL
- 6.0 = Expected difference between the "E" and "S" tube concentrations in µg/dL
- 100 = Conversion factor

### 12.7 S.I. Units

Using a molecular weight of 362.5 daltons, serum calibrators, controls, and patient values (serum and plasma) can be converted from micrograms per deciliter (µg/dL) to nanomoles per liter (nmol/L) of S.I. units by multiplying by 27.59 (Refer to TABLE II for values of cortisol calibrators listed in both units). Similarly, both unextracted and extracted urine controls and patient values can be converted from micrograms per 24 hours (µg/24 hr) to nanomoles per 24 hours (nmol/24 hr) of S.I. units by multiplying by 2.76.


**TABLE II**  
Concentration of Calibrators

Cat. No.	Calibrator Level	S.I. Units
CA-2367	0 µg/dL	0 nmol/L
CA-2358	1 µg/dL	28 nmol/L
CA-2359	3 µg/dL	83 nmol/L
CA-2370	10 µg/dL	276 nmol/L
CA-2371	25 µg/dL	690 nmol/L
CA-2372	60 µg/dL	1655 nmol/L

**TABLE III**  
Recording the Data  
Do not use to calculate unknowns

Tube No.	Contents of Tubes	CPM Bound	Cortisol Concentration (µg/dL)
T <sub>1</sub>	Total Counts (Tracer-Buffer)	44516	
T <sub>2</sub>	Total Counts (Tracer-Buffer)	43651	
1	Cortisol Serum Blank, B <sub>0</sub>	22311	0
2	Cortisol Serum Blank, B <sub>0</sub>	22771	0
3	Cortisol Serum Calibrator	19216	1
4	Cortisol Serum Calibrator	19243	1
5	Cortisol Serum Calibrator	15331	3
6	Cortisol Serum Calibrator	15731	3
7	Cortisol Serum Calibrator	10660	10
8	Cortisol Serum Calibrator	10294	10
9	Cortisol Serum Calibrator	6862	25
10	Cortisol Serum Calibrator	6657	25
11	Cortisol Serum Calibrator	4186	60
12	Cortisol Serum Calibrator	4224	60
13	Patient "A" Sample	18470	1.3
14	Patient "A" Sample	18992	1.1
			Av: 1.2
15	Patient "B" Sample	10377	10.1
16	Patient "B" Sample	10608	9.6
			Av: 9.9
17	Patient "C" Sample	6316	28.8
18	Patient "C" Sample	6318	28.8
			Av: 28.8

**Typical Results**  
Do not use to calculate unknowns

**FIGURE 1**

### 13. LIMITATIONS OF THE PROCEDURE

#### 13.1 Procedural

The tracer-buffer reagent must be protected from light. The coated tubes are ready to use as received and should be stored at 2-27°C. The plastic wrapper should be closed securely when storing unused tubes to prevent condensation.

Incubation conditions should be checked routinely. The water bath temperature must be constant at 37 ± 2°C and the water level kept above that of the solution in the tubes without allowing tubes to float.

The user should note that erroneous results may be caused by improper storage and hand-ling of samples. Hemolyzed, lipemic and icteric samples should be avoided.

Failure to obtain the appropriate cortisol values for controls may indicate imprecise manipulations, improper handling, or deterioration of reagents. Controls should be included in every run. A calibrator curve must be established for every run.

The GammaCoat Cortisol Kit should not be used to qualify cortisol levels in samples from patients treated with prednisolone. See Interpretive Limitations.

#### 13.2 Interpretive

Normal cortisol levels are generally highest from 5 AM to 10 AM and lowest from 8 PM to 4 AM.<sup>5</sup> PM levels are usually 1/2 to 1/3 of AM levels. A mid-day surge of cortisol levels in women, coincident with food intake and markedly attenuated by food deprivation, has been described.<sup>17</sup> Diurnal variation, stimulation and suppression tests, and various types of stressful situations normally alter cortisol levels. For example, ACTH stimulation usually yields values two to three times the reference AM level.<sup>18</sup> Metyrapone and dexamethasone suppression normally yield AM samples with values less than 5 µg/dL.<sup>19</sup> Estrogens, as in oral contraceptives, may cause elevated cortisol levels as a result of altered concentrations of plasma and serum proteins. Thus, each patient's cortisol values are usually evaluated in reference to his/her AM cortisol level.

Cortisol assays of samples from patients under prednisolone treatment are contraindicated because of high cross-reactivity of this kit's cortisol antibody with prednisolone and 6-methylprednisolone. The high cross-reactivity may yield high apparent cortisol levels.

Cortisol determinations generated by the direct urinary protocol will result in clinically valid values when compare with the reported unextracted urinary cortisol normal range. Since it is well known that unidentified substances present in urine can interfere with free cortisolimmunoassays, the values obtained using direct urinary protocols will be higher than those obtained from urine samples that have been purified by techniques such as high pressure liquid chromatography.<sup>19,20</sup>

### 14. EXPECTED VALUES

#### 14.1 Plasma or Serum Levels

Mornings:

7-25 µg/dL or 70-250 ng/mL

[193-690 nmol/L]

Evenings:

2-9 µg/dL or 20-90 ng/mL

[55-248 nmol/L]

These ranges may be used as guides until each laboratory has established its own normal ranges.

#### 14.2 Unextracted Urinary Cortisol

The normal range for urinary cortisol (UC) of unextracted urine samples is 75-270 µg per 24 hours [207-745 nmol/24 hr]. Improper urine specimen collections may result in incorrect patient evaluation. It has been recommended that creatinine be determined to identify grossly inadequate urine collections.<sup>21</sup>

### 14.3 Extracted Urinary Free Cortisol

The normal range of urinary free cortisol is 20-90 µg cortisol per 24 hours [55-248 nmol/24 hr]. Improper urine specimen collection may result in incorrect patient evaluation. It has been recommended that creatinine be determined to identify grossly inadequate urine collections.<sup>21</sup>

## 15. SPECIFIC PERFORMANCE CHARACTERISTICS

### 15.1 Precision

The intra-run precision was determined from the mean of 20 simultaneous assays per serum sample. The inter-run precision was determined from the mean of the average of duplicates for 20 separate runs.

Intra-Run Precision Sample	Number of Assays	Mean (µg/dL)	Standard Deviation (µg/dL)	Coefficient of Variation (%)
Serum Pool A	20	2.9	0.19	6.6
Serum Pool B	20	12.1	0.93	7.7
Serum Pool C	20	47.1	3.18	6.8

Inter-Run Precision Sample	Number of Separate Runs	Mean (µg/dL)	Standard Deviation (µg/dL)	Coefficient of Variation (%)
Serum Pool D	20	3.7	0.33	9.0
Serum Pool E	20	12.1	1.19	9.8
Serum Pool F	20	36.9	3.24	8.8

**15.2 Trueness:** The Trueness of this assay has been checked by the recovery test.

### Recovery from Pooled Sera

Spiked samples were prepared by diluting a known high cortisol sample with a previously analyzed serum pool. Each spiked sample was assayed in quadruplicate. Percent recovery is calculated as (recovered/expected) x 100.

Contribution from Spike (µg/dL)	Contribution from Pool (µg/dL)	Expected Value (µg/dL)	Recovered Value (µg/dL)	Percent (%) Recovery
4.2	1.6	5.8	5.5	95
11.2	1.6	12.7	12.8	101
17.2	1.6	18.8	18.4	98
26.3	1.6	27.9	25.4	91
41.2	1.6	42.8	39.2	92

### 15.3 Analytical Sensitivity

The sensitivity of the calibrator curve is defined as the smallest single value which can be distinguished from zero. A statistical estimation of the minimal detectable concentration (sensitivity) was calculated according to the method of D. Rodbard<sup>22</sup> for 30 replicates at the zero point of the calibrator curve. The calculated sensitivity is 0.21 µg/dL.

### 15.4 Avidity

The calculated affinity constant of the cortisol antiserum is approximately  $2 \times 10^8$  liters/mole.

### 15.5 Analytical Specificity

Data on the cross-reactivity of the antiserum used in this kit are expressed as the ratio of cortisol concentration to the cross-reacting substance concentration at 50% inhibition of maximum binding.

Compound	% Cross-reactivity
Cortisol	100
Prednisolone	25.5
6-Methylprednisolone	14.5
11-Deoxycortisol	9.8
17-Hydroxyprogesterone	0.4
Corticosterone	0.9
Dexamethasone	0.049
Prednisone	3.8
Deoxycorticosterone	0.14
Tetrahydrocortisone	0.2
Aldosterone	<0.1
$\beta$ -Cortol	<0.1
$\beta$ -Cortolone	<0.1
Cortisone	10.3
Dihydrocortisone	0.43
Progesterone	0.015
Spirolactone	<0.1
Tetrahydrocortisol	0.3
6- $\beta$ -Hydrocortisone	4.2

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.