GammaCoat[®] Plasma Renin Activity **¹²⁵I RIA Kit**

For the quantitative determination of plasma renin activity (PRA) by the radioimmunoassay of generated angiotension I

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1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

The GammaCoat [¹²⁵I] Plasma Renin Activity Radioimmunoassay Kit is to be used for the quantitative determination of plasma renin activity (PRA) by the radioimmunoassay of generated angiotensin I.

2. SUMMARY AND EXPLANATION

Renin, a proteolytic enzyme with a molecular weight of approximately 40,000, is released from the juxtaglomerular cells of the kidney. The enzyme acts in general circulation to cleave its substrate, an alpha-2 globulin synthesized by the liver, to produce the deca-peptide angiotensin I.

Angiotensin I is cleaved rapidly by the activity of converting enzyme, primarily in the lungs, to the biologically active octapeptide angiotensin II. In turn, angiotensin II is degraded rapidly to inactive peptide fragments by enzymes present in plasma and tissues, known as angiotensinases. The metabolic pathway of the renin angiotensin system is described below:

renin substrate

↓ renin

angiotensin I

↓ converting enzyme

angiotensin II

↓ angiotensinases

inactive peptide fragments

Angiotensin II has an extremely short *in vivo* half-life, but it is the most potent vaso-pressor known. Angiotensin II plays a key role in several forms of hypertension, as well as in blood pressure regulation. In addition, angiotensin II has been established as the major influence on aldosterone secretion by the adrenal gland.

Technical difficulties associated with the measurement of angiotensin II levels in blood have retarded general acceptance of its assay in the clinical laboratory. Since angiotensin I levels are a direct representation of plasma renin activity, the determination of plasma renin activity has been adopted widely to evaluate the renin-angiotensin system in disease states. Measurement of plasma renin activity in hypertensives is an important aid in the differential diagnosis of primary and secondary aldosteronism. Estimation of renin activity is also valuable in determining the prognosis and most appropriate therapy for persons with essential hypertension.

A radioimmunoassay for angiotensin I and its application as an index of plasma renin activity has been described.¹ The GammaCoat Kit procedure for determining plasma renin activity is adapted from this method.

3. PRINCIPLES OF THE ASSAY

The procedure is based on competitive binding principles of radioimmunoassay.² In the GammaCoat [¹²⁵I] Plasma Renin Activity RIA Kit, the antibody is immobilized onto the lower inner wall of the GammaCoat tube. The PRA determination involves an initial incubation of plasma to generate angiotensin I, followed by quantitation of angiotensin I by radioimmunoassay. The generation of angiotensin I in the human substrate system is governed by:

3.1 pH of the incubation.

3.2 extent of plasma dilution.

3.3 choice of enzyme inhibitors.

3.4 unknown factors in each individual plasma sample such as concentration of the renin substrate level.

3.5 length of incubation.

3.6 temperature of incubation.

Generation pH

The generation rate of angiotensin I *in vitro* is pH dependent. This kit procedure uses a generation pH of 6.0, the usual pH of samples after addition of pH 5.7 buffer. The advantage of incubation at pH 6.0, rather than the physiologic pH 7, is the twofold greater rate of generation of angiotensin I with correspondingly better sensitivity for low renin samples.³⁻¹⁴ As a result of optimization, shorter generation times may be utilized.

Selection of Enzyme Inhibitor

A variety of agents can be employed to block enzymatic conversion of angiotensin I to angiotensin II and to prevent proteolytic degradation during generation of angiotensin I.¹⁵⁻¹⁷ Phenylmethylsulfonyl fluoride (PMSF) is the inhibitor used in this procedure.

Effects of Plasma Dilution

Plasma dilution should be avoided since the renin reaction is substrate-dependent at the concentration of renin substrate which usually exists in plasma. The effects of substrate dilution can be seen at all intervals during generation, and more marked effects are seen at higher dilution.⁴

Length of Generation

The shortest practical time should be employed for determination of enzyme activity to provide the best estimate of initial velocity of the reaction and to minimize substrate consumption. The radioimmunoassay of angiotensin I is sensitive enough to permit PRA determination using optimized conditions of angiotensin I generation at pH 6.0. The GammaCoat Plasma Renin Activity Kit method offers a reduced generation time because the plasma sample is minimally diluted, buffered at an optimized pH, and run with an inhibitor which is effective at the selected pH. It has been suggested that low level samples can be quantified by generating the samples for eighteen hours and then assaying.⁴

Length of Radioimmunoassay

The GammaCoat Plasma Renin Activity Kit method uses a three hour radioimmunoassay incubation. Results are available within the same working day.

4. REAGENTS PROVIDED IN THE KIT

Rabbit anti-Angiotensin I Coated Tubes	1 bag/100 tubes	5 bags/100 tubes
Angiotensin Assay Buffer (11x)	1 vial/10 mL	5 vials/10 mL
Angiotensin Maleate Generation Buffer	1 vial/5 mL	3 vials/5 mL
Angiotensin PMSF	1 vial/1 mL	2 vials/1 mL
Renin Activity Control	1 vial/3 mL	2 vials/3 mL
Angiotensin Calibrators (A-F)	6 vials/2 mL	12 vials/2 mL
Angiotensin Tracer	2 vials/5 mL	10 vials/5 mL
Number of tests	100	500

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent as recommended below. Reagents should not be used past the expiration date on the label. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

Reagents from different batches must not be mixed.

4.1 [¹²⁵I] Angiotensin I Tracer: lyophilized reagent

Reconstitute the [¹²⁵I] Angiotensin I Tracer vial by adding 5 mL purified water. Mix well by inverting or gently vortexing the vial. Add the contents of **two** tracer vials to 100 mL of assay buffer. The tracer-buffer reagent is stable at 2-8°C until the kit expiration date. Upon reconstitution each vial contains approximately 1 µCi tracer (<1 µg/mL Angiotensin I) in 5 mL of phosphate buffered saline with bovine serum albumin with thimerosal added as a preservative.

4.2 Rabbit Anti-Angiotensin I Serum Coated Tubes: ready to use reagent

12 x 75 mm tubes are coated with rabbit anti-angiotensin I serum (titer <1 µg/tube) and packed in a plastic bag. Store at 2-27°C.

4.3 Angiotensin Assay Buffer Concentrate (11X): concentrated reagent

Add the entire contents of the assay buffer concentrate vial to 100 mL purified water and mix well. Each vial contains 10 mL of phosphate buffer with 0.1% sodium azide as a preservative.

4.4 Angiotensin Maleate Generation Buffer: ready to use reagent

Each vial contains 5 mL of maleate buffer, sodium EDTA, neomycin sulfate and inert blue coloring with 0.1% sodium azide as a preservative. If crystal formation occurs during shipping the vial may be warmed at 37°C until crystals dissolve.

4.5 Angiotensin Phenylmethylsulfonyl Fluoride: ready to use reagent

Each vial contains 1 mL of phenylmethylsulfonyl fluoride (PMSF) in ethanol.

4.6 Renin Activity Control: lyophilized reagent

Reconstitute before use by pipetting 3.0 mL purified water at 2-8°C. Keep vial in an ice bath for 10 minutes, and gently vortex or swirl to mix thoroughly. Aliquot 1.0 mL portions into pre-chilled glass or plastic tubes. Seal tightly and store at -20°C. One aliquot may be kept in an ice bath or refrigerated at 2-8°C if the assay is performed at this time.

The reconstituted control must be stored frozen at -20°C and is stable for at least two months. Thawing must be done in an ice bath or in a refrigerator at 2-8°C and the reconstituted control must **not** be refrozen after thawing.

Variations in PRA have been observed upon repeated assay of frozen plasma after various periods of storage.^{19,20} Thus, the use of stored frozen plasma as a control in PRA determinations may lead to unreliable results. The Renin Activity Control is utilized routinely during Angiotensin I **generation** as well as radioimmunoassay, providing a reliable quality control index for the entire assay. Refer to Quality Control section for more information.

Upon reconstitution each vial contains 3 mL of processed human plasma with 0.1% sodium azide.

4.7 Angiotensin I Calibrators (A-F): lyophilized reagent

NOTE: Blank and calibrators are used **only** in the radioimmunoassay portion and **must not** undergo angiotensin I generation.

Reconstitute each vial with 2.0 mL of purified water and mix well before use. Blank and calibrators should then be chilled in an ice bath before assaying. They may be stored at 2-8°C up to two weeks, and at -20°C for longer storage.

Upon reconstitution, each vial contains 2 mL of angiotensin I, BSA in phosphate buffered saline. Calibrators are calibrated at 0, 0.2, 0.8, 3.0, 10.0 and 50.0 ng/mL, respectively. The DiaSorin Plasma Renin Activity calibrators are calibrated using the current production master lot. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus, hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories" 4th ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING THIMEROSAL

Some reagents in this kit contain thimerosal which contains a mercury compound. Disposal of elemental mercury, inorganic mercury, mercury oxides and mercury compounds should be done in strict compliance with all local, state, and federal regulations.

WARNING: This product contains a chemical known to the State of California to cause birth defects or other reproductive harm.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

Reagents Containing Ethanol

R11 - Highly flammable.

S7 - Keep container tightly closed.

S16 - Keep away from sources of ignition — No smoking.

Reagents Containing PMSF

R22 - Harmful if swallowed.

R36/38 - Irritating to eyes and skin.

S45 - In case of accident or if you feel unwell, seek medical advice immediately. (show the label where possible.)

Reagents Containing Potassium Hydroxide

R35 - Causes severe burns.

S26 - In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S37/39 - Wear suitable gloves and eye/face protection.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive materials which do not exceed 2 μCi (74 kBq – kit CA-1533) or 10 μCi (370 kBq – kit CA-1553) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:
The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. COLLECTION AND PREPARATION OF SPECIMEN

Plasma renin activity is highest in the morning because of diurnal variation in renin release, and samples should be collected regularly at that time.^{18,19} A venous blood sample is collected aseptically using an evacuated glass tube containing EDTA. Samples may be collected, transported and centrifuged at room temperature. Blood can sit at room temperature for as long as six hours before centrifugation without significant accumulation of angiotensin I. After centrifugation of the blood, plasma can be stored frozen until assayed. Samples may be frozen at -20°C up to one month. Freeze-thaw cycles, which can occur in self-defrosting freezers, should be avoided. Variations in PRA have been observed both due to storage at 2-8°C, and to prolonged storage at -20°C.^{20, 21}

Prior to assay, frozen samples should be thawed rapidly to room temperature. Hemolyzed, lipemic, icteric or citrated samples should not be used. Do not use heparin as an anti-coagulant. See Procedural Limitations section for additional precautions and information.

7. MATERIALS REQUIRED, BUT NOT SUPPLIED

- 7.1 Purified water.
- 7.2 Constant temperature water bath, 37 ± 2°C.
- 7.3 Ice water bath.
- 7.4 Volumetric pipets (1, 2, 3 and 5 mL).
- 7.5 Graduated cylinder (100 mL).
- 7.6 Precision pipets (10, 100, 500 and 1000 µL).
- 7.7 Gamma counter.
- 7.8 Plastic or glass tubes, 12 x 75 mm (for angiotensin I generation).
- 7.9 Test tube racks.
- 7.10 Vortex mixer (optional).
- 7.11 Water aspirator (optional).
- 7.12 Plastic-backed absorbent paper (optional).

8. ANGIOTENSIN I GENERATION PROCEDURE

- 8.1 Transfer 1.0 mL of plasma to an **uncoated** glass or plastic tube marked with the identifying number and the suffix 37.
- 8.2 Add 10 µL of the PMSF solution and 100 µL of the maleate generation buffer (pH 6.0) to each 1.0 mL aliquot from step 1. Mix well and place in ice bath.
- 8.3 Transfer 500 µL of the plasma from step two to a corresponding chilled tube marked with the suffix 4.
- 8.4 Place the 37 series of tubes into a 37°C water bath. Maintain the 4 series of tubes in an ice bath. Run the maleate buffer generation for 90 minutes.
It has been suggested that patient samples with PRA values less than 1.0 ng/mL/hr may be generated in an eighteen hour incubation and then reassayed.⁴ Calculations must be adjusted for any change in generation ti
- 8.5 At the end of the generation period, transfer the 37 series of samples either to the ice bath if the assay is to be performed immediately, or to the freezer along with the 4 series for future assay.

9. ASSAY PROCEDURE

The assay procedure includes the preparation of a calibrator curve from which the unknown angiotensin I content in both the 37 and 4 series samples is interpolated. The background (4 series) is then subtracted from the corresponding 37 generation samples.

9.1 Allow tracer-buffer reagent to reach ambient temperature and mix well before using. Maintain blank and calibrators, generated control, and generated patient samples in an ice bath and mix well before using.

9.2 Label a set of GammaCoat tubes in duplicate according to the following scheme. Total Counts [T₁,T₂] and B₀ tubes [1, 2] may be required for certain data reduction programs, but may be omitted if the calibrator curve is plotted on semi-logarithmic graph paper.²³⁻²⁵

Tube No	Contents of Tubes		Code Letter	Angiotensin I (ng/0.1 mL)
T ₁ ,T ₂	Total Counts (Tracer-Buffer)			
1, 2	Angiotensin I Blank	0 ng/mL	A	0
3, 4	Angiotensin I Calibrator	0.2 ng/mL	B	0.02
5, 6	Angiotensin I Calibrator	0.8 ng/mL	C	0.08
7, 8	Angiotensin I Calibrator	3.0 ng/mL	D	0.30
9, 10	Angiotensin I Calibrator	10.0 ng/mL	E	1.0
11, 12	Angiotensin I Calibrator	50.0 ng/mL	F	5.0
13, 14	Renin Activity Control, 37°C Generation			
15, 16	Renin Activity Control, 4°C Generation			
17, 18	Patient "X", 37°C Generation			
19, 20	Patient "X", 4°C Generation			

9.3 Pipet into the appropriate duplicate tubes:

- a. 100 µL of each angiotensin I blank or calibrator.
- b. 100 µL of the "37" and "4" generation set renin activity control to the appropriate group of 4 tubes.
- c. 100 µL of each "37" and "4" generation set patient sample to the appropriate group of 4 tubes.

9.4 Immediately add 1.0 mL of tracer-buffer reagent to each tube, including Total Counts. Mix reagents by gently vortexing each tube.

9.5 Incubate all tubes at room temperature (20-27°C) for three hours.

9.6 Aspirate or decant all tubes, except Total Counts.
FAILURE TO REMOVE ADHERING SOLUTION ADEQUATELY MAY RESULT IN POOR REPLICATION AND SPURIOUS VALUES.

If the aspirating technique is used, be sure that the plastic tip of the aspirating tube touches the bottom of the coated tube and that all liquid is removed.

If the decanting technique is used, allow the tubes to drain in an inverted position for 3-5 minutes. Tap the tubes on absorbent paper to remove any adhering liquid before placing them upright.

9.7 Count all tubes in a gamma counter for one minute with the window suitably adjusted for iodine-125.

9.8 Calculate results. Refer to Results section.

10. QUALITY CONTROL

One Renin Activity Control is provided in the kit. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. The control must undergo generation before assaying.

Each laboratory should utilize controls at several levels to monitor assay performance. The controls should be treated as unknowns. Quality control charts should be maintained to follow the performance of the controls. Appropriate statistical methods should be used to evaluate trends. Acceptable performance limits should be ascertained for the individual laboratory.^{22, 23}

11. RESULTS

- 11.1 Record the counts per minute (CPM) bound for each tube.
- 11.2 Plot the CPM bound for angiotensin I calibrators (vertical axis) versus nanograms (ng) of angiotensin I added per tube (horizontal axis) on semi-logarithmic graph paper.
- 11.3 Draw the best fitting curve. For typical data, refer to TABLE I. FIGURE 1 illustrates a typical calibrator curve.
- 11.4 Locate the CPM bound for each tube which corresponds to each sample on the vertical axis and follow a horizontal line intersecting the calibrator curve. At the point of intersection, read ng angiotensin I from the horizontal axis. For each sample and control, determine the average ng "37" and average ng "4".

12. PLASMA RENIN ACTIVITY (PRA)

Plasma renin activity (PRA) is expressed as ng/mL/hr of generated angiotensin I. The following mathematical corrections are made to obtain ng/mL/hr:

- 12.1 Subtract the average "4" value (Background) from the average "37" (Generated):

$$\text{Net ng} = \text{ng "37"} - \text{ng "4"}$$

- 12.2 Size of sample assayed in coated tube (C): 0.1 mL_C

- 12.3 Dilution of sample (S) by generation buffer (B) and inhibitor (I):

$$\frac{(1.0 \text{ mL}_S + 0.1 \text{ mL}_B) + 0.01 \text{ mL}_I}{1.0 \text{ mL}_S} = \frac{1.11 \text{ mL}_{\text{Final}}}{1.00 \text{ mL}_{\text{Initial}}}$$

- 12.4 Time of angiotensin I generation: 1.5 hours

NOTE: If length of generation time is altered, then the corresponding change must be made in this number.

- 12.5 Maleate Generation Calculation:

$$\text{PRA} = \frac{\left(\frac{\text{ng "37"} - \text{ng "4"}}{0.1 \text{ mL}_C} \right) \times \left(\frac{1.11 \text{ mL}}{1.00 \text{ mL}} \right)}{1.5 \text{ hr}}$$

$$\text{PRA} = (7.40/\text{mL/hr}) \times (\text{net ng})$$

Example: Patient "X" Plasma Sample:

$$\text{PRA} = (0.41 - 0.08) \text{ ng} \times (7.40/\text{mL/hr})$$

$$\text{PRA} = 2.44 \text{ ng/mL/hr}$$

TABLE I
Recording the Data
Do not use to calculate unknowns

Tube No.	Contents of Tubes	CPM Bound	Angiotensin I Level (ng/tube)	PRA Final (ng/mL/hr)
T1	Total Counts (Tracer-Buffer)	29242	–	–
T2	Total Counts (Tracer-Buffer)	28992	–	–
1	Angiotensin I Blank 0 ng/tube	13617	–	–
2	Angiotensin I Blank 0 ng/tube	13048	–	–
3	Angiotensin I Calibrator 0.02 ng/tube	13163	–	–
4	Angiotensin I Calibrator 0.02 ng/tube	12625	–	–
5	Angiotensin I Calibrator 0.08 ng/tube	11305	–	–
6	Angiotensin I Calibrator 0.08 ng/tube	11291	–	–
7	Angiotensin I Calibrator 0.3 ng/tube	6991	–	–
8	Angiotensin I Calibrator 0.3 ng/tube	7110	–	–
9	Angiotensin I Calibrator 1.0 ng/tube	3053	–	–
10	Angiotensin I Calibrator 1.0 ng/tube	3056	–	–
11	Angiotensin I Calibrator 5.0 ng/tube	1307	–	–
12	Angiotensin I Calibrator 5.0 ng/tube	1232	–	–
13	Renin Activity Control, 37°C	3149	0.97	
14	Renin Activity Control, 37°C	3053	<u>1.01</u>	
			Av: 0.99	
15	Renin Activity Control, 4°C	8674	0.20	
16	Renin Activity Control, 4°C	8607	<u>0.20</u>	
			Av: 0.20	5.85
17	Patient "X" Plasma Sample, 37°C	4362	0.41	
18	Patient "X" Plasma Sample, 37°C	4459	<u>0.40</u>	
			Av: 0.41	
19	Patient "X" Plasma Sample, 4°C	8546	0.08	
20	Patient "X" Plasma Sample, 4°C	8557	<u>0.08</u>	
			Av: 0.08	2.44

Typical Calibrator Curve
Do not use to calculate unknowns


FIGURE 1

13. LIMITATIONS OF THE PROCEDURE

Procedural

- 13.1 The coated tube plastic bag should be closed securely when storing unused tubes and maintained at 2-27°C.
- 13.2 The water level must be kept above that of the solution in the tubes without allowing the tubes to float.
- 13.3 The user should note that erroneous results can be caused by improper handling of patient samples.
- 13.4 Strict attention must be given to the maintenance of ice water temperature during reagent preparation, angiotensin I generation procedure, and assay procedure.
- 13.5 Frozen samples should be thawed rapidly to room temperature and mixed thoroughly before pipetting.
- 13.6 Heparin has been noted to inhibit the renin reaction at concentrations of 5 IU/mL; and therefore, heparin should not be used as an anticoagulant in blood samples to be assayed for plasma renin activity.⁵
- 13.7 Hemolyzed specimens contain angiotensinases and should not be used.
- 13.8 Published data presented evidence that linearity of angiotensin I generation during a predetermined generation period, even as short as thirty minutes or as long as eighteen hours, cannot be assumed.³ Departure from linearity may be the rule rather than the exception even under controlled buffering at any pH. It may be that factors other than pH can influence the linearity of angiotensin I generation. Until such variables are truly controllable, it seems necessary to anticipate non-linearity characteristics.
- 13.9 Failure to obtain appropriate control values may indicate imprecise manipulations, improper handling or deterioration of reagents. A calibrator curve based on the appropriate values for calibrators must be established for every run.

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It is critical that physiological variables, such as the patient's posture and either sodium intake or 24-hour urinary sodium excretion, be considered for accurate interpretation of assay results. Age, sex, and race are also important factors for consideration.

Changes in PRA can be affected by a spectrum of drugs including diuretics, adrenergic blocking agents, vasodilators, anti-hypertensives, large doses of progestins, angiotensin II antagonists, oral contraceptives or estrogen therapy, and mineralocorticoid anta-gonists. Use of diatrizoate as in renography, can also affect PRA values. Failure to take into account the pharmacological contribution of the patient's drug intake may lead to the misinterpretation of PRA values. Various clinical conditions, including normal pregnancy, affect PRA. These factors have been outlined concisely in review articles published by several investigators.²⁷⁻³⁰

14. EXPECTED VALUES

Clinical evaluation of this PRA kit was performed by an independent investigator, and results are summarized briefly in this section. The data presented are intended as a guide for the user of this kit. Because of variances in patient population, users are advised to establish their own data (for example, PRA vs. Sodium Excretion Nomogram).

14.1 Effect of Salt Intake

Plasma samples and 24-hour urine collections were obtained from normal adult subjects without apparent cardiovascular or renal diseases. Each subject was placed on a diet of known sodium and potassium content. After an equilibration period, plasma samples for PRA and 24-hour urine samples were obtained. The results of these studies are contained in TABLE II. The data were divided into four subgroups based on the range of sodium excretion observed: *ad libitum* salt intake (75-150 mEq/24 hrs.); salt supplemented diet (150 mEq/24 hrs.); and two salt restricted diets. It is feasible in routine clinical circumstances to obtain modest reductions in sodium excretion, to 30-75 mEq/24 hrs. The lower range of sodium excretion, 0-30 mEq/24 hrs., is unusual and requires severe restriction of sodium intake. The predicted inverse relationship between PRA and urine sodium excretion was observed. Since dietary salt content varies regionally, it is suggested that PRA vs sodium excretion data be developed by the individual laboratory.

During the course of clinical evaluation of this kit, samples were obtained from ambulatory or hospitalized adults without apparent cardiovascular or renal disease. These subjects were on *ad libitum* salt intakes and no effort was made to control their activities. All samples were obtained between 8 AM and 10 AM. The data from these observations are contained in TABLE III. They fall within the range of PRA observed for subjects on normal salt intake for this population. The higher values of hospitalized subjects may be ascribed to dietary salt alterations imposed by institutional food, stress induced by hospital environment, or other undetermined factors.

TABLE II
PRA with Restricted Salt Intake

Na ⁺ Excretion mEq/24 hr	Mean ±1 S.D. (ng/mL/hr)
0 - 30	16.34 ± 7.52
30 - 75	5.91 ± 1.82
75 -150	2.12 ± 0.68
>150	0.85 ± 0.46

TABLE III
PRA with *ad libitum* Salt Intake

	Mean ±1 S.D. (ng/mL/hr)
Random Ambulatory Normals	1.67 ± 0.83
Non-Ambulatory Normals	3.30 ± 1.85

14.2 Positional Effect

Accurate assessment of PRA values depends on the clinician's knowledge of the subject's posture, since PRA values are affected by position. The effect of assuming an erect position after a night of supine sleep is indicated in TABLE IV.

TABLE IV
Positional Effect

	Mean \pm 1 S.D. (ng/mL/hr)
Supine	1.24 \pm 1.09
Erect	2.63 \pm 1.32

14.3 Furosemide Stimulation

Clinical assessment of the renin-angiotensin system may require evaluation of PRA change after stimulation. The effect of furosemide on the renin-angiotensin system, as detected by change in PRA, was determined as shown by a study summarized in TABLE V. There was a significant rise in PRA values with furosemide stimulation.

TABLE V
Furosemide Stimulation*

	Mean \pm 1 S.D. (ng/mL/hr)
Pre-stimulation	2.36 \pm 1.23
Post-stimulation	6.92 \pm 2.76

NOTE: A monograph detailing the entire evaluation is available upon request by contacting Customer Service. The monograph contains the results of the following studies:

1. urine aldosterone vs. PRA.
2. renal vein renin studies in patients with renovascular hypertension.
3. PRA determination in anephric patients.
4. PRA results in patients on hemodialysis.
5. renin profiling of patients with essential hypertension.

* five hrs. after furosemide, 40 mg p.o.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Analytical Sensitivity

The sensitivity of the calibrator curve is defined as the smallest single value which can be distinguished from zero. A statistical estimation of the minimal detectable concentration (sensitivity) was calculated according to the method of D. Rodbard, (1978),³¹ for thirty replicates at the zero point of the calibrator curve. The calculated sensitivity is 0.018 ng/tube.

15.2 Precision

The intra-run precision was determined from the mean of twenty simultaneous assays per sample. The inter-run precision was determined from the mean of the average of duplicates for twenty separate runs.

Intra-Run Precision Sample	Number of Assays	Mean (ng/mL/hr)	Standard Deviation (ng/mL/hr)	Coefficient of Variation
Plasma Pool A	20	1.6	0.16	10.0
Plasma Pool B	20	6.2	0.28	4.6
Plasma Pool C	20	17.9	1.68	9.4

Inter-Run Precision Sample	Number of Assays	Mean (ng/mL/hr)	Standard Deviation (ng/mL/hr)	Coefficient of Variation
Plasma Pool A	20	1.6	0.09	5.6
Plasma Pool B	20	10.7	0.82	7.6
Plasma Pool C	20	15.2	1.04	6.8

15.3 Avidity

The calculated affinity constant of this kit's antiserum is approximately 3×10^{10} liters/ mole.

15.4 Analytical Specificity

Data on the cross-reactivity of the antiserum used in this kit are expressed as the ratio of angiotensin I concentration to the cross-reacting substance concentration at 50% inhibition of maximum binding.

Compound	% Cross Reactivity
Angiotensin I	100
Tetradecapeptide*	0.02
Angiotensin II	<0.03
Angiotensin III	<0.03

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.