C-terminal PTH

¹²⁵I RIA Kit

For the quantitative determination of
C-terminal Human Parathyroid Hormone

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1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

This kit contains instructions and materials for the quantitative determination of C-terminal human parathyroid hormone (hPTH) in serum or plasma by radioimmunoassay (RIA).

2. SUMMARY AND EXPLANATION

Native parathyroid hormone is an eighty-four amino acid peptide which acts to increase serum calcium by specific effects in soft tissue and bone metabolism. Parathyroid peptides are secreted from small glands in response to a decrease of serum calcium. High levels of serum calcium generally inhibit PTH secretion, except in the case of PTH secreting tumors. PTH circulates at a relatively low concentration; therefore, sensitive radiometric techniques are needed to measure PTH in blood.

Clinically, the measurement of PTH has proven to be an extremely helpful and efficient tool in the differential diagnosis and management of hypercalcemia; more specifically, PTH measurement can help the physician diagnose tumors and hyperplasia of the parathyroid glands.¹ Radioimmunoassay detection of high blood levels of PTH can often establish that excess PTH secretion is causing hypercalcemia.² PTH can also be useful in localizing hyperfunctioning parathyroid tissue by RIA of samples obtained via venous catheterization.³ The PTH radioimmunoassay is also useful for the diagnosis and management of hypocalcemia. With rare exceptions, hypocalcemic persons with deficient PTH (hypoparathyroidism) have subnormal concentrations of PTH.⁴ Measurement of serum PTH is frequently used to assess the status of renal osteodystrophy in renal failure patients on chronic dialysis.⁵ Serum PTH levels are often markedly elevated in patients with renal failure due to chronic negative calcium balance, hyperphosphatemia, and hypocalcemia with subsequent secondary hyperparathyroidism.

3. PRINCIPLE OF THE ASSAY

The DiaSorin C-terminal RIA is a disequilibrium procedure using delayed tracer addition to increase sensitivity. The antibody was produced against intact beef PTH and is reactive to the 65-84 sequence of human PTH. The antibody has an affinity constant of 1×10^{10} L/mole. In this RIA, sample and first antibody are combined and incubated for 16-24 hours (Option A) or 4 hours (Option B) at 2-8°C. Tracer is then added, followed by a second incubation for 16-24 hours at 2-8°C. Separation is accomplished by use of pre-precipitated double antibody which is accelerated with PEG. The pre-precipitated antibody is added as a single reagent with a subsequent incubation of two hours. The calibrators are expressed in terms of the mass of human PTH 1-84 in ng/mL. Quality control will be performed at DiaSorin using only Option B.

4. REAGENTS PROVIDED IN THE KIT

C-Term PTH Calibrator 0	1 vial/20 mL	1 vial/20 mL
C-Term PTH Calibrators 1-6	6 vials/1.0 mL	6 vials/1.0 mL
C- Term PTH Antiserum	1 vial/14 mL	2 vials/14 mL
¹²⁵ I C-Term PTH Tracer	1 vial/14 mL	2 vials/14 mL
C-Term PTH Controls	2 vials/ 1.0 mL	2 vials/ 1.0 mL
C-Term PTH Precipitating Complex	1 vial/ 35 mL	2 vials/ 35 mL
Number of tests	65	130

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After reconstitution store each reagent at -15°C or lower until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 C-terminal PTH Calibrator 0: lyophilized reagent

BSA-borate buffer with sodium azide (0.2%) and other stabilizers added. Reconstitute the vial with 20 mL of purified water and allow it to stand for 15-20 minutes until completely dissolved; mix thoroughly before using.

4.2 C-terminal hPTH Calibrators (1-6): lyophilized reagent

Synthetic hPTH 65-84 calibrators, at nominal concentrations ranging from 0.312 to 10 ng/mL (31.2-1000 pmol/L), are pre-diluted in Calibrator 0 containing sodium azide (0.1%) and other stabilizers. Exact concentration values are assigned according to each lot. Reconstitute each vial with 1.0 mL of purified water and allow it to stand for 15-20 minutes until completely dissolved; mix thoroughly before using. The DiaSorin synthetic hPTH 65-84 calibrator has been calibrated against the World Health Organization (W.H.O.) preparation 79/500. Any comparison with other products or procedures should be done with this reference standard. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

4.3 C-terminal PTH Antiserum: lyophilized reagent

Chicken anti-human parathyroid hormone antibody diluted in BSA-borate buffer with sodium azide (0.1%) added. Reconstitute the vial with 14 mL of purified water and allow it to stand for 15-20 minutes until completely dissolved; mix thoroughly before using.

4.4 ¹²⁵I C-terminal PTH: lyophilized reagent

Synthetic human C-terminal PTH, labeled with iodine-125 and diluted in BSA-borate-EDTA buffer with sodium azide (0.4%) added. Reconstitute the vial with 14 mL of purified water and allow it to stand for 15-20 minutes until completely dissolved; mix thoroughly before using.

4.5 C-terminal PTH Precipitating Complex: lyophilized reagent

Normal chicken serum pre-precipitated with rabbit or goat anti-chicken serum and polyethylene glycol (PEG), diluted in BSA-borate buffer with sodium azide (0.3%) added. Reconstitute the vial with 35 mL of purified water; mix thoroughly until the suspension appears homogeneous and then allow the vial to stand for a minimum of 30 minutes at room temperature with occasional mixing.

4.6 C-terminal PTH Controls, Level 1 and 2: lyophilized reagent

Human serum spiked, if necessary, with the appropriate amount of synthetic hPTH 65-84 to obtain a concentration within a specified range. Sodium azide (0.1%) and other stabilizers are added. Reconstitute the vial with 1.0 mL of purified water and allow it to stand for 15 - 20 minutes until completely dissolved; mix thoroughly and treat as an unknown sample.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg and antibody to HIV. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4TH ed., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control, Atlanta, GA, U.S. 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

The kits contain radioactive material not exceeding 1.5 μCi (55.5 kBq) of iodine-125 for kit No. 13065 or 3 μCi (111 kBq) of iodine-125 for kit No. 13130. Appropriate precautions and good laboratory practices should be used in the storage, handling and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of particulate matter in any of the reagents.
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

7. COLLECTION AND PREPARATION OF SERUM AND PLASMA

One hundred microliters of serum or EDTA plasma, in duplicate, are required for use in the C-terminal PTH assay.

For serum, collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. Allow the blood to clot at room temperature (15°C-25°C). For plasma, collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. EDTA (7.2 mg/5 mL blood) is used as an anticoagulant. Centrifuge for 15 minutes using approximately 760 x g* to obtain hemolysis free sera. No additives or preservatives are required to maintain integrity of the sample. All plastics, glassware or other material coming into contact with the specimen should be entirely free of any contamination. A fasting specimen is recommended but not required. Serum or plasma samples should be stored frozen for best results. Studies indicate that C-terminal PTH specimens are stable for at least 3-6 months at -20°C. Specimens should not be repeatedly frozen and thawed.

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED.

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Temperature controlled centrifuge to accommodate 12 x 75 mm tubes.
- 8.3 Gamma scintillation counter capable of counting ¹²⁵I.
- 8.4 Vortex.
- 8.5 Pipetting devices:
 - a. Micropipettors calibrated to deliver 100 µL, 200 µL, and 500 µL.
 - b. Repeating dispensers, calibrated to deliver 200 µL and 500 µL.
- 8.6 Purified water.

9. ASSAY PROCEDURE – OPTION A

- 9.1 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 20-25°C. Mix gently before using.
- 9.2 Set up labeled 12 x 75 mm glass tubes in duplicate according to the Scheme of the Assay, on the back page.
- 9.3 Add reagents to the tubes as follows:
 - a. **Total Count Tubes**
Set aside until step number 5
 - b. **Nonspecific Binding (NSB)**
100 µL of calibrator 0
 - c. **Calibrator 0**
100 µL of calibrator 0
200 µL of C-terminal PTH antiserum
 - d. **Calibrators (1-6)**
100 µL of PTH calibrator
200 µL of C-terminal PTH antiserum
 - e. **Controls and Unknown Samples**
100 µL of serum
200 µL of C-terminal PTH antiserum
- 9.4 Vortex gently without foaming and incubate at 2-8°C for 16-24 hours.
- 9.5 Add 200 µL of ¹²⁵I PTH to all tubes.
- 9.6 Vortex gently without foaming and incubate at 2-8°C for 16-24 hours.
- 9.7 Add 500 µL of precipitating complex to all tubes except the total count tubes.

$$* g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 9.8 Vortex gently without foaming and incubate for two hours (± 15 minutes) at 2-8°C.
- 9.9 Centrifuge the tubes using a minimum of 760 x g* for 20 minutes at 2-8°C or 20-25°C.
- 9.10 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a maximum of two minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 9.11 In a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure section.)

10. ASSAY PROCEDURE – OPTION B

- 10.1 Reconstitute the lyophilized reagents and allow any frozen reagents to thaw completely. Do not allow reagents to reach temperatures above 20-25°C. Mix gently before using.
- 10.2 Set up labeled 12 x 75 mm glass tubes in duplicate according to the protocol in TABLE II.
- 10.3 Add reagents to the tubes as follows:
 - a. **Total Count Tubes**
Set aside until step number 5
 - b. **Nonspecific Binding (NSB)**
100 μ L of calibrator 0
 - c. **calibrator 0**
100 μ L of calibrator 0
200 μ L of C-terminal PTH antiserum
 - d. **Calibrators (1-6)**
100 μ L of PTH calibrator
200 μ L of C-terminal PTH antiserum
 - e. **Controls and Unknown Samples**
100 μ L of serum
200 μ L of C-terminal PTH antiserum
- 10.4 Vortex gently without foaming and incubate at 2-8°C for 4 hours (± 15 minutes).
- 10.5 Add 200 μ L of 125 I PTH to all tubes.
- 10.6 Vortex gently without foaming and incubate at 2-8°C for 16-24 hours.
- 10.7 Add 500 μ L of precipitating complex to all tubes except the total count tubes.
- 10.8 Vortex gently without foaming and incubate for two hours (± 15 minutes) at 2-8°C.
- 10.9 Centrifuge the tubes using a minimum of 760 x g* for 20 minutes at 2-8°C or 20-25°C.
- 10.10 Immediately decant the supernatant from all the tubes except the total count tubes by inverting them for a maximum of two minutes. Blot the tubes on absorbent paper to remove any drops of supernatant that may be remaining on the rims before turning the tubes upright.
- 10.11 In a gamma scintillation counter, count the precipitate of each tube and the total count tubes for a sufficient time to achieve statistical accuracy. (See Limitations of the Procedure section.)

* $g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$

11. PROCEDURAL COMMENTS

- 11.1 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 11.2 Some manufacturers' disposable borosilicate glass tubes yield elevated nonspecific bindings.
- 11.3 If you choose to aspirate supernatant from the precipitate, be careful not to disturb the precipitate pellet.
- 11.4 To completely monitor the consistent performance of an RIA there are additional factors which may be checked. DiaSorin suggests a check of the following parameters to assure consistent kit performance:
 - a. **Total Counts**
 - b. **Maximum Binding**
Average counts per minute (CPM) of Calibrator 0 Tubes/Average CPM of Total Count Tubes.
 - c. **Nonspecific Binding**
Average CPM of NSB Tubes/Average CPM of Total Count Tubes.

12. QUALITY CONTROL

Each laboratory should include at least two controls in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included in this kit using Option B. These ranges are printed on the control vials.

13. CALCULATION OF RESULTS

There are many methods in existence for calculating the results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either linear or logarithmic scale. Each of these methods give essentially the same values for controls and samples; although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is % B/B₀ versus log concentration.

- 13.1 Calculate the average CPM for each calibrator, control, and unknown sample.
- 13.2 Subtract the average CPM of the NSB tubes from all counts.
- 13.3 Divide the corrected CPM of each calibrator, control, or unknown sample by the corrected CPM of the Calibrator 0.

$$B/B_0 (\%) = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of Calibrator 0} - \text{CPM of NSB}} \times 100$$

- 13.4 Using 2-cycle semi-log or log-logit graph paper, plot percent B/B₀ for the PTH calibrators (vertical axis) versus the concentration (horizontal axis).
- 13.5 Draw a best-fit line through the points.
- 13.6 Interpolate the levels of PTH in the samples from the plot.
- 13.7 If any unknown sample reads greater than the highest calibrator, it should be diluted appropriately with Calibrator 0 and assayed again.
- 13.8 If an unknown sample has been diluted, correct for the appropriate dilution factor.
- 13.9 Calculate maximum binding by dividing CPM of Calibrator 0 by the average total counts obtained in the total count tubes.

REDUCTION DATA

The DiaSorin QC lab uses a smoothed spline curve fit.

TABLE III
DiaSorin C-terminal PTH RIA Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Percent (B/T)	Percent Bound (B/B ₀)	Conc. (ng/mL)
Total Count	15,724 15,738	15,731				
NSB	302 360	331		2.1		
Calibrator 0	7,262 7,124	7,193	6,862	45.7	100	
Calibrators (ng/mL)						
1 (0.32)	6,734	6,702 6,669	6,371		92.8	
2 (0.62)	6,358	6,392 6,425	6,061		88.3	
3 (1.44)	5,566 5,656	5,611	5,280		76.9	
4 (2.93)	4,139 4,074	4,105	3,774		55.0	
5 (6.03)	2,173	2,195 2,217	1,864		27.2	
6 (11.6)	753 745	749	418		6.1	
Unknown Samples						
1	4,463 4,522	4,493	4,162		60.7	2.48
2	6,382 6,466	6,424	6,093		88.8	0.64

Typical sample data and a calibrator curve for Option B are shown in TABLE III and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

C-TERMINAL PTH SAMPLE CALIBRATOR CURVE


FIGURE 1

14. LIMITATIONS OF THE PROCEDURE

- 14.1 If the initial concentration of an unknown sample is greater than the highest calibrator, dilute with Calibrator 0 only.
- 14.2 Counting times should be sufficient to prevent statistical error, for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error.
- 14.3. Recent studies have shown that if the second antibody incubation is performed at 2-8°C instead of 20-25°C, a substantial increase in zero binding will result. If centrifugation is performed at 20-25°C, the temperature must be controlled so that it does not exceed 25°C.

15. EXPECTED VALUES

Normal Range

Each laboratory should establish its own normal range. Normal values observed in the DiaSorin Laboratory were 0.44 ± 0.22 ng/mL (44.4 ± 22.0 pmol/L). This would lead to a normal range (two std. dev.) of nondetectable (ND) - 0.88 ng/mL (ND - 88 pmol/L).

Interpretation of Clinical Values

Values for C-terminal PTH fragments in primary hyperparathyroidism are generally above normal, although, as with all other C-terminal PTH assays, there can be overlap of pathological values with normal values. Uremic sera will give vastly elevated values for C-terminal PTH.

TABLE IV shows the data for 16 proven primary hyperparathyroid patients and two parathyroidectomized individuals on vitamin D and calcium. In 87% of these cases, there was clear separation from normal, and the hypoparathyroid patients were in the undetectable range.

In hypercalcemia of malignancy C-terminal PTH is not generally elevated out of the normal range in contrast to primary hyperparathyroidism, although with calcium levels less than 11.5 mg/dL there is overlap of these two disease states.

TABLE IV
 C-terminal PTH Values in Sera from
 Surgically Proven Parathyroid Adenomas
Units
 Calcium: mg/dL
 PTH: ng/mL (1-84 equivalents)
Parathyroid Adenoma

Sample No.	Calcium	PTH	Sample No.	Calcium	PTH
1	12.4	2.6	9	12.4	5.9
2	10.7	2.1	10	11.0	2.0
3	11.2	1.3	11	11.3	2.2
4	10.9	3.0	12	10.8	0.8
5	11.3	3.8	13	11.2	2.0
6	13.0	3.4	14	10.8	2.4
7	11.1	1.9	15	11.0	2.9
8	12.2	1.0	16	12.8	2.8
			Mean	11.5	2.5
			Std. Dev.	0.78	1.2

Hypoparathyroid (treated)

Sample No.	Calcium	PTH
1	8.4	<0.3
2	9.2	<0.3

15. PERFORMANCE DATA

15.1 Precision

Within Assay Variation (values = ng/mL)

	Mean Value	S.D.	%C.V.
LOW	0.4	0.02	5.0
MEDIUM	2.8	0.13	4.6
HIGH	5.6	0.28	5.0

Between Assay Variation (values = ng/mL)

	Mean Value	S.D.	%C.V.
LOW	0.7	0.10	14.3
HIGH	5.1	0.41	8.0

15.2 Trueness: The assay trueness has been checked by the linearity test and the recovery test.

Linearity (Parallelism)

Serial Dilution Study of 4 Patient Samples (values = ng/mL)

Sample No.	Undiluted	1/2	1/4
1	5.6	5.4	5.2
2	7.3	6.8	6.8
3	10.0	9.0	8.8
4	>10	10.8	10.4

Recovery

Recovery Study (values = ng/mL)

Background	Calibrator Added	Expected Value	Measured Value	Percent Recovery
Set No. 1				
0.4	1.3	1.7	1.7	100
0.4	2.5	2.9	3.1	107
0.4	5.0	5.4	4.9	91
Set No. 2				
0.3	1.3	1.6	1.7	106
0.3	2.5	2.8	2.7	96
0.3	5.0	5.3	5.3	100
Set No. 3				
0.4	1.3	1.7	1.9	112
0.4	2.5	2.9	2.9	100
0.4	5.0	5.4	5.4	100

15.3 Analytical Sensitivity

When defined as the apparent concentration at three standard deviations from the counts at maximum binding, the minimum detectable amount is 0.01 ng/tube (0.1 ng/mL) if the first incubation is 24 hours, and 0.02 ng/tube (0.2 ng/mL) if the first incubation is 4 hours.

15.4 Analytical Specificity

Comparison of the cross-reactivity of C-terminal PTH antibody was made with the following peptides at 1,000 times the top calibrator for C-terminal PTH:

Peptide	% Cross-reactivity
bPTH	100
hPTH 65-84	100
hCalcitonin	<0.01
hProlactin	<0.01
hFSH	<0.01
hLH	<0.01
hTSH	<0.01
Thyroxine (T ₄)	<0.01

15.5 Interference

Interference studies were performed using NCCLS EP7-A as a guideline to assess the effects of common endogenous interferents. Values obtained demonstrated interference due to the addition of Tryglycerides, cholesterol and hemoglobin.

SEE LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

Option A

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow any reagents to reach above 20-25°C.
2. Identify tubes in duplicate.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrator 0	-	100 µL	100 µL	-
Calibrators (1-5)	-		100 µL	-
Controls	-	-	-	100 µL
Unknown	-	-	-	100 µL
Samples				
Antiserum			200 µL	200 µL

4. Mix well; incubate for 16-24 hours at 2-8°C.
5. Dispense 200 µL of Tracer into all wells.
6. Mix well; incubate for 16-24 hours at 2-8°C.
7. Dispense 500 µL Precipitating Complex into all wells except the Total Count tubes.
8. Mix well, incubate for 2 hours (+/- 15 minutes) at 2-8°C.
9. Centrifuge using 760 x g* for 20 minutes.
10. Decant the supernatants.
11. Count each tube in a gamma counter for 60 seconds or longer.

SCHEME OF THE ASSAY

Option B

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Do not allow any reagents to reach above 20-25°C.
2. Identify tubes in duplicate.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	NSB	Cal 0-5	Controls and unknown samples
Calibrator 0	-	100 µL	100 µL	-
Calibrators (1-5)	-		100 µL	-
Controls	-	-	-	100 µL
Unknown	-	-	-	100 µL
Samples				
Antiserum			200 µL	200 µL

4. Mix well; incubate for 4 hours (+/- 15 minutes) at 2-8°C.
5. Dispense 200 µL of Tracer into all wells.
6. Mix well; incubate for 16-24 hours at 2-8°C.
7. Dispense 500 µL Precipitating Complex into all wells except the Total Count tubes.
8. Mix well, incubate for 2 hours (+/- 15 minutes) at 2-8°C.
9. Centrifuge using 760 x g* for 20 minutes.
10. Decant the supernatants.
11. Count each tube in a gamma counter for 60 seconds or longer.

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




*g = $(1118 \times 10^{-8}) (\text{raggio in cm}) (\text{rpm})^2$

REFERENCES/BIBLIOGRAFIA/BIBLIOGRAFIA

1. Roos, B.A., A.W. Lindall, D.C. Aron, J.W. Orf, M. Yoon, M.B. Huber, J. Pensky, J. Ells and P.W. Lambert, "Detection and Characterization of Small Mid-region Parathyroid Hormone Fragment(s) in Normal and Hyperparathyroid Glands and Sera by Immuno-extraction and Region-specific Radioimmunoassay," **Journal of Clinical Endocrinology and Metabolism**, 53:709, (1981).
2. Lindall, A.W., B.A. Roos and M. Cecchetti, "Potential Clinical Usefulness of New Glandular and Circulating Parathyroid Peptides Illuminated by Sequence Specific Radioimmunoassay," Presented at the 3rd International Symposium on Calcitropic Hormones Methods and Clinical Applications at Gardone, Italy, April 1981, published in *Monoclonal Antibodies and Developments in Immunoassay*, 217-229, A. Albertini and R. Ekins, eds., **Elsevier/North-Holland Biomedical Press**, (1981).
3. Lindall, A.W. and E.T. Wong, "Localization of Substernal Parathyroid Adenomas," **Minnesota Medicine**, 57:87, (1974).
4. Mallette, L.E., S.N. Tuma, R.E. Berger and J.L. Kirkland, "Radioimmunoassay for the Middle Region of Human Parathyroid Hormone Using an Homologous Antiserum with a Carboxy-terminal Fragment of Bovine Parathyroid Hormone as Radioligand," **Journal of Clinical Endocrinology and Metabolism**, 54:1017, (1982).
5. Lindall, A.W., J. Ells, J. Elting and B. Roos, "Estimation of Biologically Active Intact Parathyroid Hormone in Normal and Hyperparathyroid Sera by Sequential N-terminal Immunoextraction and Mid-region Radioimmunoassay," **Journal of Clinical Endocrinology and Metabolism**, 57:1007, (1983).

SYMBOLS USED WITH DEVICES

	English	Español	Italiano
	European Conformity	Conformidad europea	Conformità europea
	Expiration Date	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricante	Fabbricante
	Consult Instructions for Use	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	<i>In vitro</i> diagnostic.	Diagnóstico in vitro.	Diagnostica in vitro.
	Maximum Temperature	Temperatura máxima	Temperatura massima
	Reconstitute with X mL	Reconstitución con X mL	Ricostituire con X mL
	Lot No.	Número de lote	Lotto n°
	Antiserum	Antisuero	Antisiero
	Precipitating reagent	Reactivo precipitante	Reagente precipitante
	Tracer: antigen labelled with ¹²⁵ I	Trazador: antígeno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
	Calibrator	Calibrador	Calibratore
	Control serum	Suero de control	Siero di controllo
	Radioactive	Radioactivo	Radioattivo
	Harmful	Nocivo	Nocivo

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the U.S. and Canada Call Toll Free: 800-328-1482	
In the United Kingdom Call: 44 118 9364200 FAX: 44 118 9792061	
10391	27830 1/04

PRINTED IN U.S.A.