**N-tact[®] PTH SP
IRMA Kit**

For the quantitative determination of biologically active intact
hPTH 1-84 in serum or plasma by immunoradiometric assay

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1. INTENDED USE

FOR IN VITRO DIAGNOSTIC USE.

This kit contains instructions and materials for the quantitative determination of biologically active intact human parathyroid hormone (hPTH 1-84) in serum or plasma by immunoradiometric assay (IRMA).

2. SUMMARY AND EXPLANATION

The maintenance of proper calcium levels in cells and extracellular fluids is of fundamental importance for many biological processes, and parathyroid hormone (PTH) is one of the most important regulators of calcium homeostasis.

Clinically, the measurement of PTH is used in conjunction with calcium determinations to assess disorders in calcium metabolism.¹⁻³ The measurement of PTH has proven to be difficult due to the heterogeneity of peptides that occur in both the gland and circulation.⁴

PTH is synthesized in the parathyroid glands as a 115 amino acid precursor referred to as preproparathyroid hormone. Preproparathyroid hormone is converted to a 90 amino acid intermediate form, proparathyroid hormone.⁵ Additional proteolytic action converts the intermediate hormone to biologically active PTH 1-84.⁶ In addition to the formation of fragments within the gland, there is further proteolytic breakdown of PTH once it is released into the circulation.^{7, 8}

Scientific studies indicate only intact PTH and PTH fragments containing the 1-34 amino acid sequence contain biological activity. Insignificant concentrations of the N-terminal fragment(s) contribute to the circulating pool of immunoreactive parathyroid hormone. Therefore, when using PTH levels to evaluate patients, the measurement of intact PTH correlates best with their state of calcium metabolism.⁹⁻¹³

Parathyroid hormone and its fragments are cleared from the circulation by both the kidneys and liver.¹⁴ The clearance of C-terminal fragments is slower than the clearance of the intact hormone and is more dependent upon renal mechanisms.¹⁵ In patients with severe or end-stage renal failure, C-terminal fragments accumulate to very high levels. When using a C-terminal specific PTH assay, elevated PTH concentrations may be found simply because of diminished clearance of these biologically inactive fragments, complicating the differentiation of hypercalcemia.^{8, 16-18}

3. PRINCIPLE OF THE ASSAY

The DiaSorin intact PTH SP immunoradiometric assay (IRMA) utilizes 2 different polyclonal antibodies that have been purified using affinity chromatography. These purified antibodies are specific for 2 different regions of the PTH molecule.

The first antibody, specific for PTH 39-84, is bound to a solid phase (polystyrene beads). The second antibody is specific for PTH 1-34 and is labeled with iodine-125. Samples are incubated simultaneously with both antibodies. Intact PTH 1-84 contains both the 1-34 and the 39-84 amino acid sequences and is the only form of PTH that will be bound by both the antibody on the bead and the antibody labeled with iodine-125.

Since the antibody coupled to the solid phase is specific for C-terminal and mid-region fragments as well as intact PTH, the capacity of the solid phase has been designed to accommodate very high levels of PTH. This prevents interference by extremely elevated C-terminal and mid-region PTH fragments in unknown samples.

Following the incubation period, each bead is washed to remove any unbound labeled antibody. The radioactivity present in the remaining bound labeled antibody is then measured using a gamma counter.

Concentrations of intact PTH present in the samples are directly proportional to the radioactivity measured.

4. REAGENTS PROVIDED IN THE KIT

N-Tact PTH SP Calibrator 0	1 vial/15 mL
N-Tact PTH SP Calibrators 1-5	5 vials/2.0 mL
N-Tact PTH SP Beads	1 container/100 beads
¹²⁵ I N-Tact PTH SP Antibody	2 vials/5 mL
N-Tact PTH SP Wash Solution Conc.	1 vial/50 mL
N-Tact PTH SP Controls	2 vials/2.0 mL
Number of tests	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent at 2-8°C until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

When reconstituting the contents of the vials, mix gently to avoid foaming. Store all reconstituted reagents at -15°C or lower immediately following use. Reagents from different batches must not be mixed.

4.1 N-tact PTH SP Calibrator 0: ready to use reagent

Contains human serum with 0.1% sodium azide added as a preservative.

4.2 N-tact PTH SP Calibrators 1-5: lyophilized reagent

Set of 5 human PTH 1-84 calibrators, at nominal concentrations ranging from 15-2000 pg/mL, contains 0.1% sodium azide and other stabilizers. Exact concentration values are assigned with each lot. Reconstitute each vial of lyophilized calibrator with 2.0 mL of the 0 calibrator provided, mix by swirling and allow them to stand for 5 minutes at 2-8°C or until the contents are completely dissolved. The kit calibrators are calibrated using Synthetic Human PTH. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

4.3 N-tact PTH SP Beads: ready to use reagent

Polystyrene beads are coated with affinity-purified goat antibody specific for the 39-84 sequence of PTH.

4.4 ¹²⁵I N-tact PTH SP Antibody: ready to use reagent

Each vial contains 5.0 mL of affinity purified goat antibody specific for PTH 1-34, labeled with iodine-125. The labeled antibody is diluted in buffered serum containing red dye and 0.1% sodium azide.

4.5 N-tact PTH SP Wash Solution Concentrate: reagent in solution

Contains a concentrated buffered surfactant. Prepare a working wash solution by diluting the entire vial contents with 450 mL of distilled or deionized water. Store the working wash solution at room temperature.

4.6 N-tact PTH SP Control, Level 1 and Level 2: lyophilized reagent

Human serum is spiked with hPTH 1-84 to obtain a concentration within a specified range. 0.1% sodium azide is added. Reconstitute each lyophilized control with 2.0 mL of distilled or deionized water, mix and allow to stand at 2-8°C until the contents are completely dissolved. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

5. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 20 μ Ci (740kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:
The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.2 A decrease in maximum binding.
- 6.3 High nonspecific zero binding.
- 6.4 Poor duplicate values.

7. SPECIMEN COLLECTION AND PREPARATION

Two hundred microliters of serum or EDTA plasma, in duplicate, are required for use in the N-tact PTH SP assay.

Either human serum or plasma may be used in this kit. The anticoagulants EDTA may be used with this assay. A fasting specimen is recommended, but not required. Blood should be collected aseptically by venipuncture in a 5 or 10 mL evacuated glass tube. Allow the blood to clot at room temperature (15-25°C). Centrifuge for 15 minutes using approximately 760 x g* to obtain hemolysis free sera. No additives or preservatives are required to maintain integrity of the sample. All plastics, glassware or other material coming into contact with the specimen should be entirely free of any contamination. Store serum or plasma samples at -20°C or lower. Specimens may be stored in glass or plastic vials, as long as the vials are tightly sealed to prevent desiccation of the sample.

Normal EDTA plasma and serum samples were compared by DiaSorin laboratories. No significant differences in values were seen.

8. ASSAY PROCEDURE

- 8.1 Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Keep reagents and thawed specimens on ice while setting up assay.
- 8.2 Set up labeled 12 x 75 mm borosilicate glass tubes in duplicate according to the Scheme of the Assay.
- 8.3 Add reagents to the tubes as follows:
 - a. **Total count tubes**
100 µL ¹²⁵I N-tact PTH SP Antibody (red)
 - b. **Calibrator 0**
200 µL calibrator 0
100 µL ¹²⁵I N-tact PTH SP Antibody (red)
 - c. **Calibrators (1-5)**
200 µL calibrator
100 µL ¹²⁵I N-tact PTH SP Antibody (red)
 - d. **Controls and unknown samples**
200 µL sample
100 µL ¹²⁵I N-tact PTH SP Antibody (red)
- 8.4 Vortex all tubes.

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

- 8.5 Dispense one bead into each tube, except the Total Counts, with teflon-coated forceps or suitable dispenser for handling beads. (Do not use fingers.)
- 8.6 Cover the tubes with parafilm or equivalent.
- 8.7 Incubate the tubes for 22 (± 2) hours at 20-25°C.
- 8.8 Aspirate the reaction mixture from each tube.
- 8.9 Wash the beads by vigorously dispensing 1 mL of wash solution into each tube with sufficient force to raise the bead from the bottom of the test tube. Aspirate wash solution. Repeat wash procedure 3 times.
- 8.10 Measure the radioactivity present in each tube using a gamma counter.

9. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 9.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 9.2 Test tube rack.
- 9.3 Gamma counter capable of counting iodine-125.
- 9.4 Vortex mixer.
- 9.5 Pipetting devices:
 - a. Micropipettor calibrated to deliver 200 μ L.
 - b. Repeating dispensers, calibrated to deliver 100 μ L and 1 mL.
 - c. Volumetric pipettes for reconstituting controls and calibrators.
- 9.6 Parafilm or equivalent for covering test tubes.
- 9.7 Teflon-coated forceps or other kit-approved device for dispensing beads.
- 9.8 Aspiration unit to aspirate 125 I antibody and wash solutions.

10. PROCEDURAL COMMENTS

- 10.1 When adding beads to the tubes, tilt the test tube rack slightly and allow beads to roll into the tubes. This will prevent excessive splashing of the test solution. Do not handle beads with fingers.
- 10.2 Pipette 125 I antibody to lower one third of the test tube.
- 10.3 In order for a laboratory to completely monitor the consistent performance of an RIA assay there are additional factors which must be checked. DiaSorin suggests a check for every assay of the following parameters to assure consistent kit performance.
 - a. **Total Counts**
 - b. **Maximum Binding**
Counts per minute (CPM) of 2000 Calibrator Tube
 - c. **Nonspecific Binding**
CPM of 0 Calibrator Tube

11. QUALITY CONTROL

Each laboratory should include both kit control sera in every assay to ensure the validity of each assay's results. A mean and standard deviation should then be determined for each control using a minimum of ten (10) assays. An acceptable range of values may then be obtained for these controls using ± 2 standard deviations of the values previously determined. The DiaSorin Quality Control Laboratory has determined a range for the controls included in this kit.

12. CALCULATION OF RESULTS

To determine the concentration of intact PTH found in unknown and control samples, a calibrator curve is prepared using the calibrator concentrations stated on the vial labels. Values are generated as follows:

- 12.1 Calculate the average CPM for each calibrator, control, and unknown sample.
- 12.2 Subtract the average CPM of the 0 calibrator tubes from all other average counts to obtain corrected CPMs.
- 12.3 Using log-log graph paper, plot the corrected CPM of each calibrator level on the ordinate against the calibrator concentration on the abscissa.
- 12.4 Interpolate the levels of PTH in the samples from the plot.
- 12.5 If any unknown sample reads greater than the highest calibrator, it should be diluted appropriately with 0 calibrator and assayed again.
- 12.6 If an unknown sample has been diluted, correct for the appropriate dilution factor.
- 12.7 Values less than Calibrator 1 but greater than 0.7 pg/mL (assay sensitivity) may be calculated using the formula:

$$\text{Values of unknown} = \frac{\text{Corrected CPM (unknown)}}{\text{Corrected CPM (Cal. 1)}} \times \text{Value of Cal. 1}$$

TABLE I
N-tact PTH SP Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Conc. (pg/mL)	B/T
Total Count	260,242 253,536	256,889			
0 Calibrator	384 374	379			
Calibrators (pg/mL)					
1 (15)	2,324 2,351	2,338	1,959		
2 (50)	6,469 6,630	6,550	6,171		
3 (150)	16,390 16,289	16,340	15,961		
4 (450)	34,605 35,246	34,926	34,547		
5 (2,000)	85,678 84,050	84,864	84,485		33%
Unknown Samples					
1	6,494 6,450	6,472	6,093	49	
2	34,144 34,932	34,538	34,159	442	

Typical sample data and a calibrator curve for intact PTH SP are shown in TABLE I and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.


FIGURE 1

13. LIMITATIONS OF THE PROCEDURE

- 13.1 If the initial concentration of an unknown sample is greater than the highest calibrator, dilute with Calibrator 0 only.
- 13.2 Counting times should be sufficient to reduce statistical error (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).

REDUCTION DATA

The DiaSorin QC lab uses a smoothed spline curve fit.

14. EXPECTED VALUES

NORMAL RANGE

Each laboratory should establish its own normal range. A normal range has been established by DiaSorin using serum from 129 apparently healthy, fasting adults. The geometric mean was calculated to be 26 pg/mL with a 2 S.D. range of 13-54 pg/mL.

CLINICAL SPECIMENS

1. Serum samples from forty-six (46) patients with primary hyperparathyroidism were analyzed in the N-tact PTH SP assay system. Intact PTH values ranged from 43.6-686.1 pg/mL.
2. Fourteen (14) samples from patients diagnosed as having hypoparathyroidism were analyzed in the N-tact PTH SP assay system. Values for intact PTH ranged from 0-20.6 pg/mL.
3. Forty-eight patients (48) with chronic renal disease were analyzed in the N-tact PTH SP assay system. The intact PTH values ranged from 10.3-1042 pg/mL.
4. Four (4) samples for patients with familial hyperparathyroidism were analyzed in the N-tact PTH SP assay system. The intact PTH values ranged from 59.5-100.5 pg/mL.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Precision

Intra assay Variation (values = pg/mL). Three serum controls were run in replicates of 20 in a single assay to determine within assay variation.

Sample Number	Mean Value	S.D	%C.V.
1	26	0.9	3.6
2	276	5.4	2.0
3	1,125	27.5	2.4

Inter assay Variation (values = pg/mL). Three serum controls were run by 4 technicians in 20 separate assays to determine between assay variation.

Sample Number	Mean Value	S.D	%C.V.
1	49	1.6	3.4
2	285	13.1	4.6
3	1,111	54.8	4.9

15.2 Trueness: The assay trueness has been checked by the dilution test and the recovery test.

Parallelism

Serial Dilution Study of 3 Unknown Samples (values = pg/mL). (Samples diluted with 0 calibrator.)

Sample Number	Dilution	Measured	Corrected for dilution	% Obtained
1	Undiluted	164	-----	100%
	1:2	77	154	94%
	1:4	36	144	87%
	1:8	19	152	93%
2	Undiluted	286	-----	100%
	1:2	139	278	97%
	1:4	72	288	100%
	1:8	36	288	101%
3	Undiluted	750	-----	100%
	1:2	357	714	95%
	1:4	184	736	98%
	1:8	100	800	106%

ACCURACY

Recovery Study (values = pg/mL)

Background Value	Calibrator Added	Expected Value	Measured Value	Percent Recovery
Set No. 1				
24.3	50	74.3	77.9	107
24.3	150	174.3	186.3	108
24.3	450	474.3	464	98
Set No. 2				
17.7	50	67.7	69.5	104
17.7	150	167.7	177.6	107
17.7	450	467.7	471.7	101
Set No. 3				
19.5	50	69.5	72.4	106
19.5	150	169.5	179.4	107
19.5	450	469.5	461.2	98
Set No. 4				
10.6	50	60.6	60.4	100
10.6	150	160.6	163.9	102
10.6	450	460.6	435.8	94

15.3 Analytical Sensitivity

When defined as the apparent concentration at 2 standard deviations from the counts at minimum binding, the minimum detectable amount is 0.7 pg/mL.

15.4 Analytical Specificity

The cross-reactivity of the N-tact PTH SP system was made with the following human PTH fragments at 200,000 pg/mL.

Fragment	% Cross-reactivity
hPTH 39-84	<0.1
hPTH 53-84	<0.1
hPTH 39-68	<0.1
hPTH 44-68	<0.1
hPTH 1-34	<0.1
hPTH 13-34	<0.1
hPTH 1-84	100.0

REFER TO LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Keep reagents and thawed specimens on ice while setting up the assay.
2. Identify tubes in duplicate.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	Cal 0-5	Control and unknown samples
Calibrators	-	200 μ L	-
Controls	-	-	200 μ L
Unknown samples	-	-	200 μ L
Tracer	100 μ L	100 μ L	100 μ L

4. Mix well.
5. Dispense 1 bead to all tubes, except the Total Counts tubes.
6. Cover the tubes with parafilm, incubate for 22 hours (+/- 2 hours) at 20 – 25°C.
7. Aspirate the reaction mixture from each tube.
8. Wash the beads by vigorously dispensing 1 mL of wash solution into each tube with sufficient force to raise the bead from the bottom of the test tube. Aspirate wash solution. Repeat wash procedure 3 times.
9. Count each tube in a gamma counter.

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1. Pribor, H.C., and L.F. Fabisinski, "Artificial Intelligence in Laboratory Medicine III: The Parathyroid Hormone Profile," **MLO**, (1983).
2. Habener, J.F. and G.V. Segre, "Parathyroid Hormone Radioimmunoassay," **Annals of Internal Medicine**, 91(5):782, (1979).
3. Kao, P.C., N. Jiang, et. al., "Development and Validation of a New Radioimmunoassay for Parathyrin (PTH)," **Clinical Chemistry**, 28(1):69, (1982).
4. Lindall, A.W. and E.T. Wong, "Column Chromatography of Human Serum Para-thyroid Immunoreactive Peptides," **Proceedings of the Society for Experimental Biology and Medicine**, 148:799, (1975).
5. Wong, E.T. and A.W. Lindall, "Preliminary Evidence for a Microsomal Precursor to Human Parathyroid Hormone," **Proceedings of the Society for Experimental Biology and Medicine**, 148:387, (1975).
6. Habener, J.F. and J.T. Potts, Jr., "Biosynthesis of Parathyroid Hormone," **The New England Journal of Medicine**, 299 (11 & 12):580 and 635, (1978).
7. Segre, G.V., J.F. Habener, et. al., "Parathyroid Hormone in Human Plasma," **The Journal of Clinical Investigation**, 51:3163, (1972).
8. Martin, K.J., K.A. Hruska, et. al., "The Peripheral Metabolism of Parathyroid Hormone," **New England Journal of Medicine**, 301:1092, (1979).
9. Dambacher, M.A., J.A. Fischer, et. al., "Distribution of Circulating Immunoreactive Components of Parathyroid Hormone in Normal Subjects and in Patients with Primary and Secondary Hyperparathyroidism: The Role of the Kidney and of the Serum Calcium Concentration," **Clinical Science**, 57:435, (1979).
10. Chu, S.Y. and A.K. Chu, "Intact vs. C-terminal/Mid-region Parathyrin (PTH) Assay in Diagnosis of Hyperparathyroidism - a Clinical Evaluation," **Clinical Chemistry**, 32(12): 2206, (1986).
11. Goltzman, D., B. Henderson, et. al., "Cytochemical Bioassay of Parathyroid Hormone, Characteristics of the Assay and Analysis of Circulating Hormone Forms," **Journal of Clinical Investigation**, 63:1309, (1980).
12. Lindall, A.W., J. Elting, et. al., "Estimation of Biologically Active Parathyroid Hormone in Normal and Hyperparathyroid Sera by Sequential N-terminal Immuno-extraction and Mid-region Radioimmunoassay," **Journal of Clinical Endocrinology and Metabolism**, 57(5):1007, (1983).
13. Endres, D.B. and R. Villanueva, "Immunochemiluminometric and Immunoradio-metric Determinations of Intact and Total Immunoreactive Parathyrin: Performance in the Differential Diagnosis of Hypercalcemia and Hypoparathyroidism," **Clinical Chemistry**, 37(2):162, (1991).
14. Potts, J.T., and J.F. Habener, "Fundamental Considerations in the Physiology, Biology, and Biochemistry of Parathyroid Hormone," **Metabolic Bone Disease and Clinically Related Disorders**, Avioli and Krane Eds., W.B. Saunders, Second Edition: 69, (1990).
15. Freigtag, J., K.J. Martin, et. al., "Impaired Parathyroid Hormone Metabolism in Patients with Chronic Renal Failure," **New England Journal of Medicine**, 298(1):29, (1978).
16. Schmidt-Gayk, H., M. Schmitt-Fiebig, et. al., "Two Homologous Radioimmuno-assays for Parathyrin Compared and Applied to Disorders of Calcium Metabolism," **Clinical Chemistry**, 32:57, (1986).
17. Parnham, A.J., F.C.R. Wijeyesingh, et. al., "Intact Parathyroid Hormone assay Offer Increased Sensitivity Over C-terminal Assays in The Study of Parathyroid Graft Function," **Nephron**, 44:8, (1986).
18. Hackeng, W.H.L., P. Lips, et. al., "Clinical Implications of Estimation of Intact Parathyroid Hormone (PTH) Versus Total Immunoreactive PTH in Normal Subjects and Hyperparathyroid Patients," **Journal of Clinical Endocrinology and Metabolism**, 63(2):447, (1986).

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
IVD	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
LOT	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
Ab	Antiserum	Antisérum	Antiserum	Antisuero	Antisiero
Ab PEG	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reactivo precipitante	Reagente precipitante
Ag ¹²⁵ I	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Trazador: antigéno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
BUF	Buffer	Tampon	Puffer	Tampón	Tampone
CAL	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
CONTROL	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
REAG	Reagent	Réactif	Reagenz	Reactivo	Reagente
DNR	Dry Natural Rubber	Caoutchouc naturel sec	Trockener Naturkautschuk	Goma natural seca	Gomma naturale secca
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo
	Highly flammable	Facilement inflammable	Leichtentzündlich	Fácilment inflamable	Facilment infiammabile
	Harmful	Nocif	Gesundheitsschädlich	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
+44(0) 1344 401 430
FAX: +44(0) 7884 050812

12621

34654
7/10

PRINTED IN U.S.A.

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DiaSorin Inc.