## 25-Hydroxyvitamin D <sup>125</sup>I RIA Kit

For the quantitative determination of 25-OH-D and other hydroxylated metabolites in serum or plasma

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### 1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE.

This kit is intended for quantitative determination of 25-hydroxyvitamin D (25-OH-D) and other hydroxylated vitamin D metabolites in human serum or plasma to be used in the assessment of vitamin D sufficiency. Assay results should be used in conjunction with other clinical and laboratory data to assist the clinician in making individual patient management decisions in an adult population.

### 2. SUMMARY AND EXPLANATION

As calciferol (Vitamin D) enters the circulation, it is metabolized to several forms, the major of these being 25-hydroxycalciferol (25-OH-D). The first step in the metabolism of vitamin D, 25 hydroxylation, occurs mainly in the liver.<sup>1</sup> Only a small amount of 25-OH-D is metabolized in the kidney to other dihydroxyvitamin D metabolites in man.<sup>2,3</sup> Since 25-OH-D is the predominant circulating form of vitamin D in the normal population, it is considered to be the most reliable index of vitamin D status.<sup>4</sup>

The two principal forms of 25-OH-D are cholecalciferol (vitamin D<sub>3</sub>) and ergocalciferol (vitamin D<sub>2</sub>).<sup>3,5</sup> Vitamin D<sub>3</sub> is derived mainly from actions of ultraviolet light on the skin while D<sub>2</sub> is derived solely from dietary sources. Since these two parent compounds provide various contributions to the overall vitamin D status of the individual, it is important that both forms are measured equally.<sup>6,7</sup> A great deal of research has provided information about the circulating levels of 25-OH-D metabolism and their physiological significance.

The measurement of 25-OH-D is becoming increasingly important in the management of patients with various disorders of calcium metabolism associated with rickets, neonatal hypocalcemia, pregnancy, nutritional and renal osteodystrophy, hypoparathyroidism, and postmenopausal osteoporosis.<sup>6,8-10</sup> NOTE: The performance characteristics of the DiaSorin 25-OH-D assay have not been established in a pediatric population.

### 3. PRINCIPLE OF THE ASSAY

The DiaSorin 25-OH-D assay consists of a two-step procedure. The first procedure involves a rapid extraction of 25-OH-D and other hydroxylated metabolites from serum or plasma with acetonitrile. Following extraction, the treated sample is then assayed using an equilibrium RIA procedure. The RIA method is based on an antibody with specificity to 25-OH-D. The sample, antibody and tracer are incubated for 90 minutes at 20-25°C. Phase separation is accomplished after a 20 minute incubation at 20-25°C with a second antibody precipitating complex. A NSB/Addition buffer is added after this incubation prior to centrifugation to aid in reducing non-specific binding.

### 4. REAGENTS PROVIDED IN THE KIT

|                                   |                 |
|-----------------------------------|-----------------|
| 25-OH-D Calibrators               | 6 vials/1 mL    |
| 25-OH-D NSB/Addition Buffer       | 1 bottle/70 mL  |
| 25-OH-D Antiserum                 | 1 bottle/105 mL |
| 25-OH-D Tracer                    | 1 vial/6 mL     |
| 25-OH-D DAG Precipitating Complex | 2 vials/30 mL   |
| 25-OH-D Controls                  | 2 vials/1 mL    |
| 25-OH-D Acetonitrile              | 2 vials/15 mL   |
| Number of tests                   | 100             |

STORAGE: Upon receipt, the kit should be stored at 2-8°C. After opening, store each reagent at 2-8°C until the expiration date on the label. Reagents should not be used past the expiration date. The expiration date of the kit is reported on the external label and corresponds to the expiration date of the tracer.

Reagents from different batches must not be mixed.

**4.1 25-OH-D NSB/Addition Buffer:** ready to use reagent

Phosphate-gelatin buffer containing 0.1% sodium azide. This component has a dual function. It will serve both as the NSB buffer and as the Addition buffer.

**4.2 25-OH-D<sub>3</sub> Calibrators:** ready to use reagent

Six (25-OH-D<sub>3</sub>) calibrators at concentrations ranging from 0-100 ng/mL are prediluted in processed human serum containing 0.1% sodium azide. Handle these reagents with care (see Warnings and Precautions). The kit calibrators are calibrated using UV quantitation and have been verified by HPLC analysis. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as recommended.

**NOTE:** If desired, an "Optional" calibrator (2.5ng/mL) can be created by diluting the kit calibrator "1" 1:2 in the kit zero calibrator (i.e. 100 uL "1" + 100 uL "0"). This dilution should be performed prior to extraction and then extracted as per product insert.

**4.3 25-OH-D Antiserum:** ready to use reagent.

Goat anti-25-OH-D serum is diluted in phosphate-gelatin buffer containing 0.1% sodium azide.

**4.4 25-OH-D<sub>3</sub> Tracer:** ready to use reagent.

An iodinated [<sup>125</sup>I] analog of 25-OH-D<sub>3</sub> is diluted in an ethanol-phosphate buffer. Handle this reagent with care (see Warnings and Precautions).

**4.5 Donkey Anti-Goat (DAG) Precipitating Complex:** ready to use reagent.

Donkey anti-goat serum, normal goat serum, and polyethylene glycol are diluted in a BSA-borate buffer containing gentamycin sulfate and 0.1% sodium azide. Mix for 5-10 minutes before and during use to ensure that a homogeneous suspension is achieved.

**4.6 25-OH-D Controls:** ready to use reagent.

Human serum containing 0.1% sodium azide is spiked with the appropriate amounts of 25-OH-D<sub>3</sub> to obtain control concentrations within specified ranges. Mix thoroughly and treat as an unknown samples. Control 1 represents a low-normal range and Control 2 represents a high-normal range. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. Handle these reagents with care (see Warnings and Precautions).

**4.7 Acetonitrile:** ready to use reagent.

Acetonitrile. Qualified for use with this kit. Handle this reagent with care (see Warnings and Precautions). This product contains dry natural rubber.

## 5. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

**CAUTION:** This device is to be received, acquired, possessed, and used only by physicians, clinical laboratories, hospitals, veterinaries or research facilities. Testing is to be performed only by properly qualified and trained laboratory personnel.

### REAGENTS CONTAINING HUMAN SOURCE MATERIAL

**Treat as potentially infectious.**

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV),

hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

#### **REAGENTS CONTAINING SODIUM AZIDE**

**CAUTION:** Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up (Council Directive 1999/45/EC). For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA., U.S. 1976.

#### **European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)**

R 20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R 32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

#### **REAGENTS CONTAINING ACETONITRILE**

#### **European Communities Hazardous Substance Risk Phrases (Council Directive 88/379/EEC)**

R11 - Highly flammable.

R23/24/25 - Toxic by inhalation, in contact with skin and if swallowed.

S16 - Keep away from sources of ignition — No smoking.

S36 - Wear suitable protective clothing.

S43 - In case of fire, use dry chemical or carbon dioxide.

S45 - In case of accident or if you feel unwell, seek medical advice immediately.

S51 - Use only in well-ventilated areas.

#### **REAGENTS CONTAINING ETHANOL**

#### **WARNING - FLAMMABLE**

#### **European Communities Hazardous Substance Risk Phrases (Council Directive 88/379/EEC)**

R11 - Highly flammable.

S16 - Keep away from sources of ignition — No smoking.

S43 - In case of fire, use dry chemical or carbon dioxide.

#### **REAGENTS CONTAINING IODINE-125**

This kit contains radioactive material which does not exceed 4  $\mu\text{Ci}$  (148 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

## 6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 The presence of abnormal particulate matter in any of the reagents. (The NSB buffer, antiserum and DAG precipitating reagents may be cloudy.)
- 6.2 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.3 A decrease in maximum binding.
- 6.4 A high nonspecific binding.

## 7. SPECIMEN COLLECTION AND PREPARATION

Fifty (50) microliters of serum or plasma (EDTA or Heparin) is required for the 25-OH-D extraction; a volume of 150 microliters will permit repeat analysis and provide adequate pipetting volume as well.

Either human serum or plasma may be used in this kit. The anticoagulants EDTA or heparin may be used with this assay. A fasting specimen is recommended, but not required. Blood should be collected aseptically by venipuncture in a 5 or 10 mL evacuated glass tube. Allow the blood to clot at room temperature (15-25°C). Centrifuge for 15 minutes using approximately 760 x g\* to obtain hemolysis free sera. No additives or preservatives are required to maintain integrity of the sample. All plastics, glassware or other material coming into contact with the specimen should be entirely free of any contamination. Store serum or plasma samples at -20°C or lower. Specimens may be stored in glass or plastic vials, as long as the vials are tightly sealed to prevent desiccation of the sample.

A study performed by DiaSorin on a limited number of patient samples showed samples to be stable for up to 9 weeks when stored at -20°C. No significant change in values was observed following 3 freeze-thaw cycles of the samples. Repeated freeze-thaw cycles should be avoided. Transported specimens should generally arrive frozen.

## 8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Disposable borosilicate glass tubes, 12 x 75 mm.  
**NOTE:** Plastic tubes are not suitable for use with this kit.
- Temperature controlled centrifuge to accommodate 12 x 75 mm tubes.
- Gamma counter capable of counting 125-iodine.
- Vortex mixer.
- Pipetting devices
  - a. Micropipettors calibrated to deliver 25 µL and 50 µL (imprecision less than or equal to ± 1%).
  - b. Repeating dispensers using tips that are calibrated to deliver 50 µL, 500 µL and 1.0 mL (imprecision less than or equal to ± 1%).

**NOTE:** For extraction and assay procedures below, refer to "Equipment and Materials Required but not Supplied" for pipette specifications.

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

## 9. EXTRACTION PROCEDURE

- 9.1 Set up labeled 12 x 75 mm disposable glass tubes for each calibrator, control and patient sample.
- 9.2 Add 500  $\mu$ L of acetonitrile to each tube.
- 9.3 Place pipet tip containing 50  $\mu$ L of calibrator, control or patient sample below the surface of the acetonitrile and SLOWLY add into the acetonitrile.
- 9.4 Vortex for 10 seconds.
- 9.5 Centrifuge using 1200 x g\* for 10 minutes at 20-25°C.
- 9.6 Pipette duplicate 25  $\mu$ L aliquots from the supernatant into separate appropriately labeled 12 x 75 mm tubes. **CAUTION:** Take care not to disturb pellet.
- 9.7 Assay supernatants according to the assay procedure.

## 10. ASSAY PROCEDURE

- 10.1 Allow all reagents and samples to equilibrate to room temperature. Do not allow reagents to reach temperatures above 25°C.
- 10.2 Set up labeled 12 x 75 mm disposable glass tubes in duplicate according to the Scheme of the Assay on the last page.
- 10.3 Add reagents as follows:
  - a. **Total count tubes**  
50  $\mu$ L of  $^{125}\text{I}$  25-OH-D  
1.0 mL of NSB/Addition buffer
  - b. **Nonspecific binding tubes (NSB)**  
25  $\mu$ L of 0 calibrator (extracted)  
50  $\mu$ L of  $^{125}\text{I}$  25-OH-D  
1.0 mL of NSB/Addition buffer
  - c. **Calibrators, controls, and unknown samples**  
25  $\mu$ L of calibrator, control, or unknown sample (extracted)  
50  $\mu$ L of  $^{125}\text{I}$  25-OH-D  
1.0 mL of 25-OH-D antiserum
- 10.4 Vortex gently without foaming and incubate for 90 (+/- 10) minutes at 20-25°C.
- 10.5 Add 500  $\mu$ L of DAG precipitating complex (DAG precipitating complex should be mixed thoroughly before and during use) to all tubes except the total count tubes.
- 10.6 Mix tubes well and incubate for 20-25 minutes at 20-25°C.
- 10.7 Add 500  $\mu$ L of NSB/Addition buffer to all tubes except the total count tubes. Vortex gently to mix tubes well. Use caution when performing this step to avoid splashing due to high liquid volume in tube.
- 10.8 Centrifuge all tubes for 20 minutes at 20-25°C at 1800 x g\*, except the total counts.
- 10.9 Decant the supernatants, except the total count tubes, using a foam rack tube holder or equivalent by inverting the rack into an appropriate waste container. Place the inverted rack onto absorbent paper for 2-3 minutes. Blot the tubes gently to ensure all liquid is removed.
- 10.10 In a gamma scintillation counter, count each tube for a minimum of 1 minute. Each tube should be counted for a sufficient time to achieve statistical accuracy. (See results section: Limitation of the Procedure.)

$$*g = (1118 \times 10^{-8}) (\text{radius in cm}) (\text{rpm})^2$$

## 11. PROCEDURAL COMMENTS

- 11.1 Add each aliquot of reagent to the lower third of the assay tube to ensure complete mixture of reagents.
- 11.2 To completely monitor the consistent performance of an RIA assay there are additional factors which may be checked. DiaSorin suggests a regular check of the following parameters to assure consistent kit performance.
- Total Counts**
  - Maximum Binding**  
Average counts per minute (CPM) of 0 Calibrator Tube/Average CPM of Total Count Tubes.
  - Nonspecific Binding**  
Average CPM of NSB Tube/Average CPM of Total Count Tubes.
  - Slope of Calibrator Curve**  
For example, monitor the 80 and 50% points of the calibrator line.

## 12. QUALITY CONTROL

Each laboratory should include at least two controls (one at low-normal level and one at high-normal level) in every assay to monitor assay performance. Commercially available controls or the two reference controls provided with the kit may be utilized. The kit controls have been evaluated by DiaSorin using DiaSorin 25-Hydroxyvitamin D <sup>125</sup>I RIA kit. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

The controls should be treated as unknown specimens and assayed in duplicate. Quality control charts should be maintained to follow control performance. Acceptable performance limits should be determined by each individual laboratory of each level of control using statistically based methods designed to detect both systematic and random errors. Control results must meet the laboratory's criteria for acceptability prior to reporting patient test results.<sup>18, 19, 20</sup>

## 13. CALCULATION OF RESULTS

There are many methods in existence for calculating results of RIAs. Each is based on obtaining a calibration curve by plotting the extent of binding against stated concentrations of the calibration calibrators. This graph may be either a linear or logarithmic scale. Each of these methods gives essentially the same values for controls and samples, although certain assays may "fit" better into one particular method versus another. The calculation method for the DiaSorin Quality Control Laboratory is % B/B<sub>0</sub> versus log concentration.

- 13.1 Calculate the average CPM for each calibrator, control, and unknown sample.
- 13.2 Subtract the average CPM of the NSB tubes from all counts.
- 13.3 Divide the corrected CPM of each calibrator, control, or unknown sample by the corrected CPM of the 0 calibrator.

$$B/B_0 (\%) = \frac{\text{CPM of Calibrator or Unknown Sample} - \text{CPM of NSB}}{\text{CPM of 0 Calibrator} - \text{CPM of NSB}} \times 100$$

- 13.4 Using 2 cycle semi-log or log-logit graph paper, plot percent B/B<sub>0</sub> for 25-OH-D calibrators (vertical axis) versus the concentration (horizontal axis).
- 13.5 Draw a best-fit line through the points.
- 13.6 Interpolate the levels of 25-OH-D in the unknown samples from the plot.
- 13.7 If an unknown sample has been diluted, correct for the appropriate dilution factor.

- 13.8** Calculate maximum binding by dividing CPM of 0 calibrator by the average total counts obtained in the total count tubes.
- 13.9** Automated data reduction programs may also be utilized to analyze data. DiaSorin utilizes RIACalc (Pharmacia) with a %B/B<sub>0</sub> versus log concentration, smooth-SPLINE fit program. Other data reduction methods must be validated before incorporating for regular use.

**TABLE I**  
DiaSorin 25-OH-D Vitamin D RIA Sample Data

| Tube            | Duplicate CPM                    | Average CPM | Corrected CPM | Percent Bound (B/T) | Percent (B/B <sub>0</sub> ) | Graph Conc. (ng/mL) |
|-----------------|----------------------------------|-------------|---------------|---------------------|-----------------------------|---------------------|
| Total Count     | 41878<br>39924                   | 40901       |               |                     |                             |                     |
| NSB             | 555<br>508                       | 531         |               | 1.3                 |                             |                     |
| 0 Calibrator    | 19585<br>19673<br>20314<br>19718 | 19822       | 19291         | 48.5                |                             |                     |
| Calibrators     |                                  |             |               |                     |                             |                     |
| 1               | 14503<br>14442                   | 14472       | 13941         |                     | 72.3                        | 5                   |
| 2               | 11005<br>10670                   | 10838       | 10307         |                     | 53.4                        | 12                  |
| 3               | 7182<br>7642                     | 7411        | 6880          |                     | 35.7                        | 20                  |
| 4               | 4687<br>4622                     | 4654        | 4123          |                     | 21.4                        | 40                  |
| 5               | 2768<br>2985                     | 2877        | 2346          |                     | 12.2                        | 100                 |
| Unknown Samples |                                  |             |               |                     |                             |                     |
| 1               | 10365<br>10447                   | 10406       |               |                     | 53.9                        | 12.8                |
| 2               | 4219<br>4720                     | 4634        |               |                     | 24.0                        | 40.4                |

Typical sample data and a calibrator curve are shown in TABLE I and FIGURE 1; this information is for reference only and should not be used for the calculation of any value.

#### 25-HYDROXYVITAMIN D SAMPLE STANDARD CURVE


Figure 1

#### 14. LIMITATIONS OF THE PROCEDURE

- 14.1 Counting times should be sufficient to prevent statistical error (for example, accumulation of 2,000 CPM will yield 5% error; 10,000 CPM will yield 1% error).
- 14.2 Assay results should be used in conjunction with other clinical and laboratory data to assist the clinician in making individualized patient management decisions.
- 14.3 Effects due to hemolysis or lipemia present in patient samples have not been evaluated for this assay; however, the acetonitrile extraction is expected to minimize any interference due to these substances.
- 14.4 The greatest source of imprecision in this method is likely to be the extraction step. Other sources of procedural imprecision may be imprecise micropipettors or out-dated reagents.
- 14.5 The kit antibody will demonstrate some cross-reactivity with all forms of dihydroxyvitamin D2 and D3 steroids; however, in humans, these compounds are naturally present only in picomolar concentrations.

**NOTE:** The performance characteristics of this assay have not been established in a pediatric population.

## 15. EXPECTED VALUES

### Reference Range

It is important for each laboratory to establish its own reference range representative of its typical population. Factors such as UV exposure<sup>9,10</sup> season,<sup>11,12</sup> race,<sup>7</sup> and dietary intake<sup>8</sup> are all known to affect levels of 25-OH-D in humans. A decrease in vitamin D levels with age has been demonstrated in the literature; however, the research conducted for most of these studies has been conducted in the Northern European countries, where there is less sunlight and less pervasive food fortification.<sup>8,15,16</sup> This, coupled with highly variable differences in health criteria and medication status, makes it difficult to generalize the effect of age alone on vitamin D status.<sup>12</sup> Results from studies conducted in the U.S. have yielded equivocal results, with some showing occurrence and others the absence of a decline in vitamin D status with age. Anticonvulsants have been shown to cause a decrease in 25-OH-D levels as have long term use of sunscreen contaminating p-aminobenzoic acid.<sup>17</sup> A pediatric reference range has not been established by DiaSorin using this kit.

DiaSorin evaluated serum collected from 20 male and 24 female, apparently healthy, predominately Caucasian volunteers from the midwestern U.S. using the DiaSorin 25-Hydroxyvitamin D <sup>125</sup>I RIA kit. The age of the volunteers fell within the range of 23 - 67 years. Samples were collected during the month of October. Only a slight difference in 25-OH-D levels was seen between males (mean = 21.7 ng/mL) and females (mean = 24.1 ng/mL) in this reference group. The mean for the entire population (n = 44) was 23.0 ng/mL with a 2 S.D. range of 9.0 – 37.6 ng/mL.

A high prevalence of subclinical 25 OH Vitamin D deficiency among normal apparently healthy populations has been observed in many countries, particularly in winter months.<sup>13,14</sup> Recent literature has suggested the following ranges for the classification of 25 OH Vitamin D status: Deficiency: 0 to 5 ng/mL (0 - 12.5 nmol/L); Insufficiency: 5 to 20 ng/mL (12.5 - 50 nmol/L); Hypovitaminosis D: 20 to 40 ng/mL (50 - 100 nmol/L); Sufficiency: 40 to 100 ng/mL (100 - 250 nmol/L); and Toxicity: greater than 100 ng/mL (> 250 nmol/L).<sup>24</sup>

## 16. SPECIFIC PERFORMANCE CHARACTERISTICS

### 16.1 Precision

Assay precision was evaluated at DiaSorin by testing four control levels spanning the curve over 23 operating days in 40 assays. Each assay included the four controls x 2 extractions (4 replications). The study design was consistent with CLSI Document EP5-T2 guidelines.<sup>22</sup> Data calculations were made using SAS PROC VARCOMP. The statistical model fitted is a random effects linear model.<sup>23</sup>

| Mean (ng/mL) | Within-run (% C.V.) | Total Imprecision (% C.V.) |
|--------------|---------------------|----------------------------|
| 8.6          | 11.7                | 9.4                        |
| 22.7         | 10.5                | 8.2                        |
| 33.0         | 8.6                 | 9.1                        |
| 49.0         | 12.5                | 11.0                       |

### 16.2 Trueness

The assay trueness has been checked by the dilution and recovery test .

#### LINEARITY (PARALLELISM)

Serial Dilution Study of 4 Patient Samples (values = ng/mL)

| Sample | Undiluted | 1/2  | 1/4  |
|--------|-----------|------|------|
| 1      | 74.6      | 79.4 | 74.8 |
| 2      | 87.8      | 96.4 | 95.2 |
| 3      | 85.6      | 88.2 | 82.4 |
| 4      | 55.2      | 57.8 | 60.4 |

### Recovery

Known amounts of 25-OH-D<sub>3</sub> were added to patient serum samples and percent recovery was determined.

| Sample Number | Initial Conc. (ng/mL) | Amount spiked (ng/mL) | *Expected Conc. (ng/mL) | Measured Conc. (ng/mL) | Recovery |
|---------------|-----------------------|-----------------------|-------------------------|------------------------|----------|
| 1             | 15.8                  | 5.0                   | 25.5                    | 24.8                   | 97%      |
|               |                       | 10.0                  | 35.1                    | 41.8                   | 119%     |
|               |                       | 25.0                  | 62.7                    | 61.9                   | 99%      |
| 2             | 51.2                  | 5.0                   | 60.6                    | 62.0                   | 102%     |
|               |                       | 10.0                  | 69.8                    | 73.8                   | 106%     |
|               |                       | 25.0                  | 96.4                    | 99.4                   | 103%     |
| 3             | 26.6                  | 5.0                   | 36.2                    | 33.5                   | 92%      |
|               |                       | 10.0                  | 45.7                    | 54.0                   | 118%     |
|               |                       | 25.0                  | 73.0                    | 82.8                   | 113%     |
| 4             | 57.3                  | 5.0                   | 66.6                    | 63.8                   | 96%      |
|               |                       | 10.0                  | 75.8                    | 78.1                   | 103%     |
|               |                       | 25.0                  | 102**                   | 104**                  | 102%     |
| 5             | 41.3                  | 5.0                   | 50.8                    | 48.5                   | 95%      |
|               |                       | 10.0                  | 60.1                    | 56.1                   | 93%      |
|               |                       | 25.0                  | 87.0                    | 93.2                   | 107%     |

\* Expected 25-OH-D concentration calculation includes dilution factor introduced by spiking solution.

\*\* These are extrapolated values, the measured concentration is > 100 ng/mL (top calibrator)

### 16.3 Analytical Sensitivity (Limits of Detection)

The sensitivity of this assay, when defined as the lowest quantity differentiated from zero at 2 standard deviations below the mean cpm's of the zero calibrator (n = 20), has been shown to be at or below 1.5 ng/mL. This sensitivity level applies both with the use of the calibrator "1" and the "Optional" calibrator.

### 16.4 Analytical Specificity

Data on the cross-reactivity of the antiserum used in this kit is expressed as the ratio of 25-OH-D concentration to the cross reacting substance concentration at 50% inhibition of maximum binding.

| Steroid        | % Cross-reactivity |
|----------------|--------------------|
| Vitamin D2     | 0.8                |
| Vitamin D3     | 0.8                |
| 25-OH-D2       | 100                |
| 25-OH-D3       | 100                |
| 24,25-(OH)2-D2 | 100                |
| 24,25-(OH)2-D3 | 100                |
| 25,26-(OH)2-D2 | 100                |
| 25,26-(OH)2-D3 | 100                |
| 1,25-(OH)2-D2  | 11.0               |
| 1,25-(OH)2-D3  | 11.0               |

**REFER TO LAST PAGE FOR REFERENCES**

## SCHEME OF THE ASSAY

1. Extract samples: Dispense 500 µL of acetonitrile into a glass tube and slowly add 50 µL of calibrator, control or patient sample below the surface of the acetonitrile. Vortex for 10 seconds and centrifuge using 1200 x g\* for 10 minutes at 20-25°C.
2. Dispense reagents according to the following scheme:

| Tubes/Reagents            | Total Counts | NSB    | Cal 0-6 | Controls and unknown samples |
|---------------------------|--------------|--------|---------|------------------------------|
| Extracted 0 Calibrators   | -            | 25 µL  | -       | -                            |
| Extracted Calibrators     | -            | -      | 25 µL   | -                            |
| Extracted Controls        | -            | -      | -       | 25 µL                        |
| Extracted Unknown Samples | -            | -      | -       | 25 µL                        |
| Tracer                    | 50 µL        | 50 µL  | 50 µL   | 50 µL                        |
| NSB/Addition Buffer       | 1.0 mL       | 1.0 mL | -       | -                            |
| 25-OH-D Antiserum         | -            | -      | 1.0 mL  | 1.0 mL                       |

3. Mix well; incubate for 90 minutes (+/- 10 minutes) at 20-25°C.
4. Dispense 500 uL of DAG precipitating complex into all wells, except the total count tubes.
5. Mix well; incubate for 20-25 minutes at 20-25°C.
6. Dispense 500 uL of NSB/Addition Buffer into all well, except the Total Count tubes.
7. Centrifuge using 1800 x g\* for 20 minutes.
8. Decant the supernatants.
9. Count each tube in a gamma counter for 60 seconds or longer.

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


1. Chesney. R.W.. "Current Clinical Applications of Vitamin D Metabolite Research." **Clinical Orthopaedics and Related Research**. 161:285. (1981).
2. DeLuca. H.F. and H.K. Schnoes. "Metabolism and Mechanism of Action of Vitamin D." **Annual Reviews Biochemistry**. 45:631. (1976).
3. Clemens. T.L.. "Vitamin D: Recent Advances in Basic Research and Clinical Assay Methodology." **Journal of Clinical Immunoassay**. 9(4):183. (1986).
4. Preece. M.A.. S. Tomlinson. C.A. Ribot. J. Pietrek. H.T. Korn. D.M. Davies. J.A. Ford. M.G. Dunnigan and J.L.H. O'Riordan. "Studies of Vitamin D Deficiency in Man." **Q.J. Med.**. 44:575. (1975).
5. Kruger. D.W.. J. Lee. A. Raines and E.D. Lyne. "Vitamin D Metabolism and Skeletal Disease in Pediatric Populations." **Journal of Clinical Immunoassay**. 9(4):200. (1986).
6. Hollis. B.W. and W.B. Pittard III. "Relative Concentrations of 25-Hydroxyvitamin D<sub>2</sub>/D<sub>3</sub> and 1.25-Dihydroxyvitamin D<sub>2</sub>/D<sub>3</sub> in Maternal Plasma at Delivery." **Nutrition Research**. 4:27. (1984).
7. Hollis. B.W. and W.B. Pittard. "Evaluation of the Total Fetomaternal Vitamin D Relationships at Term: Evidence for Racial Differences." **Journal of Clinical Endocrinology and Metabolism**. 59(4):652. (1984).
8. Parfitt. A.M.. M.B. BChir. J.C. Gallagher. R.P. Heaney. C.C. Johnston. R. Neer and G.D. Whedon. "Vitamin D and Bone Health in the Elderly." **American Journal of Clinical Nutrition**. 36:1014. (1982).
9. Haddad. J.G. and K.J. Chyu. "Competitive Protein-Binding Radioassay for 25-Hydroxycholecalciferol." **Journal of Clinical Endocrinology and Metabolism**. 33: 992. (1971).
10. Adams. J.S.. T.L. Clemens. J.A. Parrish and M.F. Holick. "Vitamin D Synthesis and Metabolism After Ultraviolet Irradiation of Normal and Vitamin-D-Deficient Subjects." **New England Journal of Medicine**. 306(12):722. (1982).
11. Beadle. P.C.. J.L. Burton and J.F. Leach. "Correlation of Seasonal Variation of 25-Hydroxycalciferol with UV Radiation Dose." **British Journal of Dermatology**. 102: 289. (1980).
12. Sherman. S.S.. B.W. Hollis and J.D. Tobin. "Vitamin D Status and Related Parameters in a Healthy Population: The Effects of Age. Sex. and Season." **Journal of Clinical Endocrinology and Metabolism**. 71(2):405. (1990).
13. McLaughlin. M. A. Fairney. E. Lester. P.R. Raggatt. O.J. Brown and M.R. Willis. "Seasonal Variations in Serum 25-Hydroxycholecalciferol in Man." **Lancet**. 1:536. (1974).
14. Gloth. F.M.. J.D. Tobin. S.S. Sherman and B.W. Hollis. "Is the Recommended Daily Allowance for Vitamin D Too Low for the Homebound Elderly?" **Journal of the American Geriatric Society**. 39(2):137. (1991).
15. Sowers. M.R.. R.B. Wallace. B.W. Hollis and J.H. Lemke. "Parameters Related for 25-OH-D Levels in a Population-based Study of Women." **The American Journal of Nutrition**. 43:621. (1986).
16. Belsey. R.E.. H.F. DeLuca and J.T. Potts. "A Rapid Assay for 25-OH-Vitamin D<sub>3</sub> Without Preparation Without Preparative Chromatography." **Journal of Clinical Endocrinology and Metabolism**. 38:1046. (1974).
17. Effects of Drugs on Clinical Laboratory Tests. 1991 Supplement to the third edition. **AACC Press**. Washington. D.C.

18. Tonks, D.B.. "Study of the Accuracy and Precision of Clinical Chemistry Determinations in 170 Canadian Laboratories." **Clinical Chemistry**, 9:217. (1963).
19. Rodbard, D.. et. al.. "Statistical Quality Control of Radioimmunoassay." **The Journal of Clinical Endocrinology and Metabolism**, 28:1412. (1968).
20. National Committee for Clinical Laboratory Standards. Internal Quality Control Testing: Principles and Definitions. Approved Guideline. NCCLS document C24-A (ISBN 1-56238-112-1). **NCCLS**, 771 East Lancaster Avenue, Villanova, PA 19085. (1991).
21. U.S. Department of Health, Education, and Welfare. National Center for Health Statistics. "HANES I Hematology and Clinical Chemistry Procedure Developed or Utilized by the Center for Disease Control Bureau of Laboratories, 1971-1975." Part 16. (1976).
22. National Committee for Clinical Laboratory Standards. "Evaluation of Precision Performance of Clinical Chemistry Devices." NCCLS document EP5-T2, 12(4). Tentative Guideline (1992).
23. Montgomery, D.C.. "Design and Analysis of Experiments." Wiley, NewYork (1991).
24. Zitterman, A. 2003. Vitamin D in preventative medicine: Are we ignoring the evidence? **British Journal of Nutrition** 89:552-572.

**SYMBOLS USED WITH DEVICES**

|   | English  | Français  | Deutsch                                     | Español  | Italiano   |
|---|--|---|---|--|--|
|    | European Conformity                            | Conformité aux normes européennes                 | CE – Konformitätskennzeichnung              | Conformidad europea                                | Conformità europea                                     |
|    | Expiration Date                                | Date limite d'utilisation                         | Mindesthaltbarkeitsdatum                    | Fecha de caducidad                                 | Data di scadenza                                       |
|    | Manufacturer                                   | Fabricant   | Hersteller                                  | Fabricante   | Fabbricante  |
|    | Consult Instructions for Use                   | Consulter les instructions d'utilisation          | Gebrauchsanweisung beachten                 | Consulte las instrucciones de uso                  | Consultare le istruzioni per l'uso                     |
|    | In vitro diagnostic.                           | Diagnostic in vitro.                              | In-vitro-Diagnostikum.                      | Diagnóstico in vitro.                              | Diagnostica in vitro.                                  |
|    | Lot No.  | No. de lot  | Chargen-Nr.                                 | Número de lote                                     | Lotto n°   |
|    | Temperature limitation.                        | Limitation de température.                        | Temperaturbereich                           | Limitación de temperatura                          | Limite della temperatura                               |
|    | Antiserum                                      | Antisérum   | Antiserum                                   | Antisero   | Antisiero  |
|    | Precipitating reagent                          | Réactif précipitant                               | Fällungsreagenz                             | Reactivo precipitante                              | Reagente precipitante                                  |
|    | Tracer: antigen labelled with <sup>125</sup> I | Traceur : antigène marqué à l'iode <sup>125</sup> | Tracer: <sup>125</sup> I-markiertes Antigen | Trazador: antígeno etiquetado con <sup>125</sup> I | Tracciatore: antigene etichettato con <sup>125</sup> I |
|  | Buffer   | Tampon  | Puffer                                      | Tampón   | Tampone  |
|  | Calibrator                                     | Étalon  | Kalibrator                                  | Calibrador   | Calibratore  |
|  | Control serum                                  | Sérum de contrôle                                 | Kontrollserum                               | Suero de control                                   | Siero di controllo                                     |
|  | Reagent  | Réactif   | Reagenz                                     | Reactivo   | Reagente   |
|  | Dry Natural Rubber                             | Caoutchouc naturel sec                            | Trockener Naturkautschuk                    | Goma natural seca                                  | Gomma naturale secca                                   |
|  | Radioactive                                    | Radioactif  | Radioaktiv                                  | Radioactivo  | Radioattivo  |
|  | Highly flammable                               | Facilement inflammable                            | Leichtentzündlich                           | Fácilment inflamable                               | Facilment infiammabile                                 |
|  | Harmful  | Nocif   | Gesundheitsschädlich                        | Nocivo   | Nocivo   |

**SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.**





CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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
SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.  
13040 Saluggia (VC) Italy

Phone: 651-439-9710  
FAX: 651-351-5669

|  |               |
|--|---------------|
| For Customer Service in the<br>U.S. and Canada Call Toll Free:<br>800-328-1482 |               |
| In the United Kingdom Call:<br>+44(0) 1344 401 430<br>FAX: +44(0) 7884 050812  |               |
| 12654E   | 34659<br>7/10 |

PRINTED IN U.S.A.

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DiaSorin Inc.