N-tact[®] Osteo SP
Osteocalcin IRMA Kit

For the quantitative determination of
intact human osteocalcin 1-49 in serum

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1. INTENDED USE

FOR *IN VITRO* DIAGNOSTIC USE

This kit contains instructions and materials for the quantitative determination of intact human osteocalcin 1-49 in serum by immunoradiometric assay (IRMA).

2. SUMMARY AND EXPLANATION

Osteocalcin is a 5900 dalton vitamin K dependent protein which has been described as one of the most abundant non-collagenous proteins in bone. Osteocalcin may account for up to 3% of total bone protein and contains two or three amino acid residues of gamma carboxyglutamic acid (GLA), hence the designation bone GLA protein (BGP).¹ BGP or osteocalcin is found only in bone tissue and is known to be a product of the osteoblast.²

It is believed to play a role in the mineralization process and is under the influence of the other calcium regulating hormones: calcitonin,³ parathyroid hormone and vitamin D.^{4,5,6} Because the levels of osteocalcin are a direct reflection of bone formation, its measurement closely correlates with the status of the bone metabolism of the patient.

While such measurements as alkaline phosphatase and hydroxyproline have been widely used as biochemical markers of bone metabolism, there have been significant limitations to their usefulness because they tend to be nonspecific. Generally, osteocalcin correlates with alkaline phosphatase, but occasionally, much of the alkaline phosphatase in the blood may be of non-bone origin. This limits its usefulness as a specific marker of bone turnover. Osteocalcin, however, is completely specific for bone. It has been found to be elevated in the following disorders: Osteomalacia,⁷ Paget's disease of the bone,⁸ hyperthyroidism,⁹ primary hyperparathyroidism and renal osteodystrophy.¹⁰ Osteocalcin is variably elevated in postmenopausal osteoporosis depending on whether there is high or low bone turnover.³ Reduced levels of osteocalcin are found in hypoparathyroidism,^{11,12} or long term corticosteroid therapy.^{9,13}

3. PRINCIPLE OF THE ASSAY

The DiaSorin Osteocalcin immunoradiometric assay (IRMA) measures intact human osteocalcin 1-49 with no cross-reactivity to the 1-43 fragment. This assay utilizes human osteocalcin 1-49 for calibrators and controls and two polyclonal antibodies that have been purified using affinity chromatography. These purified antibodies are specific for two different regions of the osteocalcin molecule. The first antibody, specific for the amino terminus of human osteocalcin, is bound to a solid phase (polystyrene beads). The second antibody is specific for the carboxy terminus of human osteocalcin and is labeled with ¹²⁵I. Samples are incubated on an orbital agitator at room temperature for 2 hours. Intact human osteocalcin 1-49 is the only form of osteocalcin that will be bound by both the antibody on the bead and the ¹²⁵I labeled antibody. Following the incubation period, each bead is washed to remove the unbound labeled antibody. The radioactivity present in the remaining bound labeled antibody is then measured using a gamma counter. The concentration of osteocalcin 1-49 present in the sample is directly proportional to the radioactivity measured.

4. REAGENTS PROVIDED IN THE KIT

N-Tact Osteo SP Calibrators 0-5	6 vials/1.0 mL
N-Tact Osteo SP Beads	1 container/100 beads
¹²⁵ I N-Tact Osteo SP Antibody	3 vials/11 mL
N-Tact Osteo SP Wash Solution Conc.	1 vial/ 50 mL
N-Tact Osteo SP Controls	2 vials/1.0 mL
Number of tests	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Store all reconstituted reagents at -15°C, or lower, until the expiration date on the label. Reagents should not be used past the expiration date.

When reconstituting the contents of the vials, mix gently to avoid foaming. Reagents from different batches must not be mixed.

4.1 N-tact® Osteo SP Calibrators (0-5): lyophilized reagent

Six human osteocalcin 1-49 calibrators, at nominal concentrations ranging from 0 to 200 ng/mL (0.0 nmol/L - 33.9 nmol/L). The calibrators must be reconstituted with 1.0 mL of purified water before use. Reconstituted calibrators should not be subjected to more than two freeze thaws. The DiaSorin N-tact® Osteo SP Kit is standardized by quantitative amino acid analysis. The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this *in vitro* diagnostic test as recommended.

4.2 N-tact® Osteo SP Controls, Level 1 and Level 2: lyophilized reagent

Buffered clarified human serum containing Trasylol at 400 KIU/mL and 0.1% sodium azide spiked with human osteocalcin 1-49. The controls must be reconstituted with 1.0 mL of purified water before use. Reconstituted controls should not be subjected to more than two freeze thaws.

4.3 N-tact® Osteo SP Beads: ready to use reagent

Polystyrene beads coated with affinity purified goat antibody specific for the amino terminus of human osteocalcin.

4.4 ¹²⁵I N-tact® Osteo SP Tracer: ready to use reagent

Affinity purified goat antibody, specific for the carboxy terminus of human osteocalcin, labeled with iodine-125 and diluted in buffered goat serum containing 0.2% sodium azide and 0.2% gentamicin sulfate.

4.5 10X Wash Solution: concentrated reagent

10X concentrate of phosphate buffer and Tween 20. Dilute ten fold with purified water prior to use.

5. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE.

Not for internal or external use in humans or animals.

REAGENTS CONTAINING HUMAN SOURCE MATERIAL

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV), or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

REAGENTS CONTAINING SODIUM AZIDE

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Council Directive 1999/45/EC)

R20/21/22 - Harmful by inhalation, in contact with skin and if swallowed.

R32 - Contact with acids liberates very toxic gas.

S28 - After contact with skin, wash immediately with plenty of water.

REAGENTS CONTAINING IODINE-125

This kit contains radioactive material which does not exceed 30 μCi (1110 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

6. INDICATIONS OF POSSIBLE DETERIORATION OF KIT REAGENTS

- 6.1 A shift in the slope or position of the calibrator curve from what is normally obtained.
- 6.2 A decrease in maximum binding.
- 6.3 High nonspecific zero binding.
- 6.4 Poor duplicate values.

7. SPECIMEN COLLECTION AND PREPARATION

Twenty microliters of serum, in duplicate, are required for use in the N-tact® Osteo SP Kit. Osteocalcin has been shown to exhibit a marked diurnal variation with a maximum level at night and a minimum in late morning or early afternoon.¹⁴ Because of this, time of sampling must be taken into account when interpreting results. Collect blood by venipuncture in a 5 or 10 mL evacuated glass tube. Allow the blood to clot at room temperature. Spin blood samples for fifteen minutes in a refrigerated centrifuge to obtain hemolysis free serum. Promptly separate the serum from the cells and freeze immediately. All plastics, glassware or other material coming into contact with the specimen should be free of any contamination. Samples may be stored for up to six months at -20°C. If longer storage is needed, samples should be kept at -70°C or colder.

8. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- 8.1 Disposable borosilicate glass tubes, 12 x 75 mm.
- 8.2 Test tube rack.
- 8.3 Calibrated gamma counter capable of counting 125 iodine.
- 8.4 Vortex mixer.
- 8.5 Pipetting devices:
 - a. Micropipettor calibrated to deliver 20 µL.
 - b. Repeating dispensers calibrated to deliver 300 µL and 1 mL.
 - c. Volumetric pipettes for reconstituting controls and calibrators.
- 8.6 Parafilm or equivalent for covering test tubes.
- 8.7 Teflon-coated forceps or other kit-approved device for dispensing beads.
- 8.8 Aspiration unit to aspirate ¹²⁵I antibody and wash solution.
- 8.9 Orbital agitator capable of 190 ± 10 rpm with an approximate orbital diameter of 2 cm.
- 8.10 Purified water.

9. ASSAY PROCEDURE

- 9.1 Reconstitute the lyophilized reagents and dilute the 10X Wash Solution (see Reagents Section). Mix the reconstituted calibrators and controls thoroughly. Allow the tracer and beads to equilibrate at room temperature. Allow any frozen samples to thaw completely. Keep calibrators, controls, and thawed samples on ice while setting up assay.
- 9.2 Set up labeled 12 x 75 mm borosilicate glass in duplicate according to the scheme of the assay, on the back page.
- 9.3 Add reagents to the tubes as follows:
 - a. **0 Calibrator**
20 µL of 0 Calibrator
300 µL of Tracer (red)
 - b. **Calibrators (A-E)**
20 µL of Calibrator
300 µL of Tracer (red)
 - c. **Controls and unknown samples**
20 µL of control or unknown sample
300 µL of Tracer (red)
 - d. **Total Count tubes**
300 µL of Tracer (red)
- 9.4 Vortex all tubes.
- 9.5 Dispense one bead into each tube (EXCEPT the Total Count tubes) with teflon-coated forceps or suitable dispenser for handling beads. (Do not use fingers.)

- 9.6 Cover the tubes with parafilm or equivalent.
- 9.7 Agitate the tubes on an orbital agitator at 190 ± 10 rpm for 2 hours \pm 5 minutes at 20 - 25°C.
- 9.8 Aspirate the reaction mixture from each tube (EXCEPT the Total Count tubes).
- 9.9 Wash the beads by vigorously dispensing 1 mL of diluted Wash Solution into each tube (EXCEPT the Total Count tubes) with sufficient force to raise the bead from the bottom of the test tube. Aspirate Wash Solution. Repeat wash procedure two more times.
- 9.10 Measure the radioactivity present in each tube using a gamma counter.

10. PROCEDURAL COMMENTS

- 10.1 Assay all samples in duplicate to ensure confidence in values obtained.
- 10.2 Pipette samples and tracer into the lower one third of the test tube.
- 10.3 When adding beads to the tubes, tilt the test tube rack slightly and allow beads to roll into the tubes. This will prevent excessive splashing of the test solution. Do not handle beads with fingers.
- 10.4 For a shaker with an orbital diameter less than 2 cm, adjust the speed accordingly to allow vigorous shaking.
- 10.5 To completely monitor the consistent performance of an IRMA assay, there are additional factors which may be checked. DiaSorin suggests a check of the following parameters.
 - a. **Total Counts**
 - b. **% Maximum Binding**

$$\frac{(\text{CPM of E calibrator} - \text{CPM of 0 calibrator})}{(\text{CPM of Total Counts} - \text{CPM of 0 calibrator})} \times 100$$
 - c. **Nonspecific Binding**
CPM of 0 Calibrator Tubes
- 10.6 To convert ng/mL to nmol/L; $\text{ng/mL} \times 0.1695 = \text{nmol/L}$.

11. QUALITY CONTROL

Each laboratory must include at least one high and one low control in every assay to monitor assay performance. The controls provided with the kit may be utilized for this purpose. The kit controls contain osteocalcin at both high and low concentrations and have been assayed using the DiaSorin N-tact® Osteo SP Kit to establish expected ranges. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs. The controls should be treated as unknown samples and assayed in duplicate. Assays should be considered valid only if results meet the laboratory's criteria for acceptability.

12. CALCULATION OF RESULTS

To determine the concentration of Osteocalcin found in control and unknown samples, a calibrator curve is prepared using the calibrator concentrations stated on the vial labels. Values are generated as follows:

- 12.1 Calculate the average CPM for each calibrator, control and unknown sample.
- 12.2 Subtract the average CPM of the 0 calibrator tubes from all other average counts to obtain corrected CPM.
- 12.3 Counting times should be sufficient to reduce statistical error (for example, accumulation of 2,000 counts will yield 5% error; 10,000 counts will yield 1% error).

- 12.4 Using log-log graph paper, plot the corrected CPM of each calibrator level on the ordinate against the calibrator concentration on the abscissa.
- 12.5 DiaSorin uses the spline smoothed method of data reduction.
- 12.6 Interpolate the levels of Osteocalcin in the samples from the plot.
- 12.7 If any unknown sample reads greater than the highest calibrator, it should be diluted appropriately with 0 calibrator and assayed again.
- 12.8 If an unknown sample has been diluted, correct for the appropriate dilution factor.

13. SAMPLE DATA FOR THE N-TACT® OSTEO SP ASSAY

Typical sample data and a calibrator curve for N-tact® Osteo SP are shown in TABLE I and FIGURE I respectively. This information is for reference only and should not be used for the calculation of any value.

TABLE I
DiaSorin N-tact Osteo SP Sample Data

Tube	Duplicate CPM	Average CPM	Corrected CPM	Conc. (ng/mL)
Total Count	307358 306568	306963		
0 Calibrator	282 240	261		
Calibrators (ng/mL) [nmol/L]				
A (1.5) [0.25]	915 1020	968	707	
B (5.0) [0.85]	2146 1956	2051	1790	
C (20) [3.39]	10623 10226	10425	10164	
D (60) [10.2]	33102 33775	33438	33177	
E (200) [33.9]	78471 77768	78120	77859	
Unknown Samples				
1	1057 1081	1069	808	1.88
2	2401 2205	2303	2042	5.66
3	25772 24402	25087	24826	44.00


FIGURE 1

14. EXPECTED VALUES

Reference Range

Each laboratory should establish its own reference range. In studies performed by DiaSorin, serum samples from forty-seven apparently healthy adults were collected between the hours of 7 and 10 A.M. The males ranged in age from 23 - 51; the premenopausal females ranged in age from 26 - 50. Analysis of these samples gave the following results:

	GEOMETRIC MEAN (ng/mL)	REFERENCE RANGE (ng/mL)
Males (n = 21)	6.25 [1.05 nmol/L]	3.2-12.2 [0.54-2.07 nmol/L]
Females (n = 26)	5.58 [0.94 nmol/L]	2.7-11.5 [0.46-1.95 nmol/L]

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1 Analytical Specificity

The DiaSorin N-tact® Osteo SP assay only measures intact human osteocalcin 1-49. The cross-reactivity to human osteocalcin 1-43 is < 0.01%.

15.2 Precision

TABLE III
Within Assay Variation (values = ng/mL)

Intra-Assay Precision

Sample	Mean Value (n = 20)	S.D.	%C.V.
1	2.21	0.14	6.3
2	4.98	0.27	5.4
3	40.76	1.83	4.5

TABLE IV
Between Assay Variation (values = ng/mL)

Inter-Assay Precision

Sample	Mean Value (n = 10)	S.D.	%C.V.
1	2.21	0.20	9.0
2	5.05	0.36	7.1
3	42.10	4.01	9.5

15.3 Trueness: The assay trueness has been checked by the linearity test and the recovery test.

Linearity (Parallelism)

TABLE V
Serial Dilution Study of three Patient Samples (values = ng/mL)

Sample	Undiluted	1/2	1/4	1/8
1	29.2	15.6	7.6	3.8
2	43.1	23.5	12.2	6.0
3	68.3	33.4	16.6	8.6

Accuracy

TABLE VI
Recovery Study of three Patient Samples Spiked with Osteocalcin 1-49 (values = ng/mL)

Sample Number	Added	Value Measured	Percent Recovery
1	0	2.5	100
	9.6	11.9	98
	29.7	30.9	96
	61.4	50.4	78
2	0	4.8	100
	9.6	13.7	93
	29.7	33.2	96
	61.4	62.8	94
3	0	3.1	100
	9.6	12.6	99
	29.7	33.0	101
	61.4	57.2	88

15.4 Analytical Sensitivity

Three tracers (all at eight days post iodination) were used to determine the sensitivity of the assay. The backfit method of analysis was used in which the CPM at 2 standard deviations from the zero calibrator is calculated and this value is read as an unknown off the calibrator line.

Based on the above studies, the sensitivity of the N-tact® Osteo SP kit has been determined to be 0.2 ng/mL (0.034 nmol/L).

SEE LAST PAGE FOR REFERENCES

SCHEME OF THE ASSAY

1. Reconstitute the lyophilized reagents and allow any frozen specimens to thaw completely. Keep reagents and thawed specimens on ice while setting up the assay.
2. Identify tubes in duplicate.
3. Dispense reagents according to the following scheme.

Tubes/Reagents	Total Counts	Cal 0-5	Control and unknown samples
Calibrators	-	20 µL	-
Controls	-	-	20 µL
Unknown samples	-	-	20 µL
Tracer	300 µL	300 µL	300 µL

4. Mix well.
5. Dispense 1 bead to all tubes (EXCEPT the Total Count tubes).
6. Cover the tubes with parafilm, agitate the tubes at 190 +/- 10 rpm for 2 hours +/- 5 minutes at 20-25°C.
7. Aspirate the reaction mixture from each tube except the total count tubes.
8. Wash the beads by vigorously dispensing 1 mL of diluted wash solution into each tube with sufficient force to raise the bead from the bottom of the test tube. Aspirate wash solution. Repeat wash procedure 2 times.
9. Count each tube in a gamma counter.

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
1. Poser, J.W., F.S. Esch, N.C. Ling and P.A. Price, "Isolation and Sequence of the Vitamin K Dependent Protein from Human Bone: Undercarboxylation of the First Glutamic Acid Residue," **Journal of Biological Chemistry**, 255: 8685-8691 (1980).
2. Price, P.A., M.K. Williamson and J.W. Lothringer, "Origin of the Vitamin K Dependent Bone Protein Found in Plasma and Its Clearance by Kidney and Bone," **Journal of Biological Chemistry**, 256: 12760-12766 (1981).
3. Civitelli, R., S. Gonnelli, F. Zacchei, S. Bigazzi, A. Vattimo, L.V. Avioli and C. Gennari, "Bone Turnover in Postmenopausal Osteoporosis: Effects of Calcitonin Treatment," **Journal of Clinical Investigation**, 82: 1268-1274 (1988).
4. Caniggia, A., R. Nuti, M. Galli, F. Lore, V. Turchetti and G.A. Righi, "Effect of a Long-Term Treatment with 1,25-Dihydroxyvitamin D₃ on Osteocalcin in Postmenopausal Osteoporosis," **Calcified Tissue International**, 38: 328-332 (1986).
5. Zerwekh, J.E., K. Sakhaee and C.Y. Pak, "Short-Term 1,25-Dihydroxyvitamin D₃ Administration Raises Serum Osteocalcin in Patients with Postmenopausal Osteoporosis," **Journal of Clinical Endocrinology and Metabolism**, 60: 615-617 (1985).
6. Kerner, S.A., R.A. Scott and J.W. Pike, "Sequence Elements in the Human Osteocalcin Gene Confer Basal Activation and Inducible Response to Hormonal Vitamin D₃," **Proceedings of the National Academy of Sciences, U.S.A.**, 86: 4455 (1989).
7. Demiaux, B., M.E. Arlot, M. Chapuy, P.J. Meunier and P.D. Delmas, "Serum Osteocalcin is Increased in Patients with Osteomalacia: Correlations with Biochemical and Histomorphometric Findings," **Journal of Clinical Endocrinology and Metabolism**, 74 (5): 1146-1151 (1992).
8. Merle, B. and P.D. Delmas, "Normal Carboxylation of Circulating Osteocalcin (bone gla protein) in Paget's Disease of the Bone," **Bone and Mineral**, 11: 237-245 (1990).
9. Lukert, B.P., J.C. Higgins and M.M. Stoskopf, "Serum Osteocalcin is Increased in Patients with Hyperthyroidism and Decreased in Patients Receiving Glucocorticoids," **Journal of Clinical Endocrinology and Metabolism**, 62: 1056-1058 (1986).
10. Charhon, S.A., P.D. Delmas and L. Malaval, "Serum Bone GLA Protein in Renal Osteodystrophy: Comparison with Bone Histomorphometry," **Journal of Clinical Endocrinology and Metabolism**, 63: 892-897 (1986).
11. Duda, R.J.J., J.F. O'Brien, J.A. Katzmann, J.M. Peterson, K.G. Mann and B.L. Riggs, "Concurrent Assays of Bone Gla Protein and Bone Alkaline Phosphatase: Effects of Sex, Age, and Metabolic Bone Disease," **Journal of Clinical Endocrinology and Metabolism**, 66: 951-957 (1988).
12. Kruse, K., U. Kracht, K. Wohlfart and U. Kruse, "Biochemical Markers of Bone Turnover, Intact Serum Parathyroid Hormone and Renal Excretion in Patients with Pseudohypoparathyroidism and Hypoparathyroidism Before and During Vitamin D Treatment," **European Journal of Pediatrics**, 148: 535-539 (1989).
13. Ekenstam, E.A., S. Ljunghall and R. Hallgren, "Serum Osteocalcin in Rheumatoid Arthritis and other Inflammatory Arthritides: Relation Between Inflammatory Activity and the Effect of Glucocorticoids and Remission Inducing Drugs," **Ann. Rheum. Dis.**, 45: 484-490 (1986).
14. Gundberg, C.M., M.E. Markowitz, M. Mizuchi and J.F. Rosen, "Osteocalcin in Human Serum: A Circadian Rhythm," **Journal of Clinical Endocrinology and Metabolism**, 60: 736-739 (1985).b

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Maximum temperature	Température maximale	Höchsttemperatur	Temperatura máxima	Temperatura massima
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Antiserum	Antisérum	Antiserum	Antisero	Antisiero
	Precipitating reagent	Réactif précipitant	Fällungsreagenz	Reactivo precipitante	Reagente precipitante
	Tracer: antigen labelled with ¹²⁵ I	Traceur : antigène marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markiertes Antigen	Trazador: antigeno etiquetado con ¹²⁵ I	Tracciatore: antigene etichettato con ¹²⁵ I
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Wash buffer	Tampon de lavage	Waschpuffer	Tampón de lavado	Tampone di lavaggio
	Reconstitute with X mL	Reconstituer avec X mL	Mit X mL auflösen	Reconstitución con X mL	Ricostituire con X mL
	Coated beads	Billes enduites	Beschichtete Kügelchen	Perlas recubiertas	Perle con rivestimento
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo
	Harmful	Nocif	Gesundheitsschädlich	Nocivo	Nocivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

For Customer Service in the
U.S. and Canada Call Toll Free:
800-328-1482

In the United Kingdom Call:
+44(0) 1344 401 430
FAX: +44(0) 7884 050812

15744

34688
7/10

PRINTED IN U.S.A.

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DiaSorin Inc.