Prolifigen[®] TPA[®] -M IRMA

Immunoradiometric assay for the quantification of
Tissue Polypeptide Antigen (TPA)

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Intended Use

Prolifigen® TPA®-M IRMA is an immunoradiometric assay for the quantitative measurement of soluble fragments of cytokeratins 8, 18 and 19 in serum or other body fluids.

Principle of the procedure

Prolifigen® TPA®-M IRMA is a two-site monoclonal/polyclonal immunoradiometric (sandwich) assay. The sample is incubated with a monoclonal antibody coated plastic bead and ¹²⁵I-labelled polyclonal antibody to TPA® is added.

Warnings and Precautions for Users

1. General

- This product is for *in vitro* diagnostic use only.
- This kit must not be used after the expiry date printed on the package.
- Do not mix reagents from different master lots.
- The accuracy of the results of this assay is directly related to the degree of care exercised during pipetting, vortexing, aspiration, and adherence to methodology and temperature requirements.
- Sodium azide used as a preservative in this product can cause irritation: avoid contact with the skin and mucosa.
- Avoid microbial contamination of reagents.
- Observe the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in accordance with local guidelines.
- The components of the kit as well as the patient sera should be considered as if they were potentially infectious and should be handled with adequate precautions.
- Material Safety data sheets are available on request.

2. Radioactive material

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 27 µCi (1003 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Components of the kit

The following reagents are supplied with each test kit for 100 assay tubes:

- 1. TPA[®]-M IRMA Antibody-coated beads**
1 bottle containing >100 dry beads coated with monoclonal (mouse) anti-TPA antibodies. Contains bovine serum albumin (BSA).
- 2. TPA[®]-M IRMA Tracer**
1 vial containing 21 mL ¹²⁵Iodine labelled polyclonal anti-TPA antibody. Total radioactivity less than 1003 KBq (27 µCi). Red coloured to aid identification. Contains bovine serum albumin (BSA).
- 3. TPA[®]-M IRMA Calibrators**
6 x 1 mL vials TPA calibrator material; concentrations of 0, 50, 100, 200, 600 and 4000 U/L respectively in phosphate buffer with BSA.
- 4. TPA[®]-M IRMA Diluent**
1 x 2.5 mL vial containing phosphate buffer with BSA.
- 5. TPA[®]-M IRMA Controls 1 and 2**
2 x 1 mL vials containing TPA in phosphate buffer with BSA.
- 6. Quality control report**

Note: Reagents 2, 3, 4 and 5 above contain <0.1% sodium azide as a preservative.

The Prolifigen[®] TPA[®]-M IRMA kit should be stored at 2 – 8 °C prior to use.

Other materials required but not provided

- Polystyrene round bottomed test tubes (12 x 70–75 mm), disposable.
- Test tube rack.
- Micropipette (100 µL) and microdispensers (200 µL).
- Purified water as Wash Fluid.
- Vortex mixer.
- Bead washing system, or 2 mL repeating dispenser for addition of wash fluid and a vacuum aspiration device.
- Horizontal shaker suitable for holding test tube racks.
Optimal speed about 200 rpm.
- Gamma counter set for ¹²⁵Iodine. Efficiency of counter should be >40%.
- Plastic forceps or bead dispenser for adding beads to tubes.
- Facilities for handling and disposal of radioactive material.
- Computer with appropriate program (often integrated in gamma counter), or calculator and graph paper, log-log type, or the RIA-mat[®] 280 for the fully automated performance of the test.

Traceability of Calibrators and Controls

The calibrators and controls of the kit consist of Tissue Polypeptide Antigen (TPA). This TPA originates from purified and well-characterized material. Calibrators and Controls are calibrated against in-house reference TPA material.

Specimen collection, preparation & storage

- Blood samples should be taken before any treatment is given, either initial or repeat.
- Serum is recommended for the quantitative determination of TPA®.
- Enough blood should be collected to be sufficient for 2 x 100 µL serum (duplicate samples).
- Collect samples using standard procedures.
- If the specimen is to be assayed within 24 hours, store at 2 – 8 °C. If the specimen will not be assayed within 24 hours, freeze to below –18 °C. Once frozen, the specimen should only be thawed once.
- Do not use samples which are grossly lipaemic, hemolysed or contaminated.
- Samples which have been found to contain over 4000 U/L can be diluted with TPA®-M IRMA Diluent to get an exact value. The diluted sample should be mixed gently prior to assay and the dilution factor taken into account in the final calculation.
- Frozen sera should be mixed gently after thawing.

Preparation of reagents before use

Before each use, all reagents must be allowed to warm to room temperature (20 – 28 °C) and should be gently mixed. Bacterial contamination of stored reagents must be avoided.

Storage of reagents after use


If all the kit is not used at one time, the reagents should be stored in their original containers at 2 – 8 °C. The reagents, thus stored, are stable until the expiry date printed on the box label.

Assay procedure

Perform each determination in duplicate.

1. Add 100 µL TPA®-M IRMA Calibrator/patient serum/TPA®-M IRMA Control to each appropriately marked tube.
2. Add 200 µL of TPA®-M IRMA Tracer to each tube (the solution has a red color). To estimate total counts (T), add 200 µL Tracer to a tube and count radioactivity without further processing.
3. Add one plastic bead to each tube using plastic forceps or a bead dispenser.
4. Incubate the tubes for 2 hours at room temperature (20 – 28 °C) on a shaker. (Optimal speed about 200 rpm and amplitude 20 mm).
5. Wash each bead 3 times with 2 mL purified water.
6. Count radioactivity in a gamma counter.

Flow Diagram

Processing of results

Use a computer with a program for handling IRMA type data to evaluate the TPA[®] levels in reference sera and patient sera. Alternatively, counts or bound activity (% B/T) may be plotted manually against analyte concentration. For this purpose, log-log diagram paper is recommended.

Quality Control

It is recommended that the laboratory has its own set of QC controls, in addition to the TPA[®]-M IRMA Controls 1 and 2 provided in the kit.

Example of results


This is an **example** only and should not be used in calculations.

Sample	CPM	% B/T	TPA U/L
Total	408 992	–	–
B ₀	845	0.2	0
Calibrator 50	2 014	0.5	50
Calibrator 100	3 364	0.8	100
Calibrator 200	6 220	1.5	200
Calibrator 600	14 811	3.6	600
Calibrator 4000	59 945	14.7	4 000
Patient sample 1	1 927	0.5	46
Patient sample 2	8 110	2.0	274
Patient sample 3	27 056	6.6	1 355

To obtain optimal results (B/T ~12% for calibrator 4000 U/L) the shaking speed should be 200 rpm and amplitude 20 mm.

Example of Calibrator Curve


This is an **example** only and should not be used in calculations.

Expected values and limitations of the procedure

In a study of 400 apparently healthy individuals, 95% in the lognormal distribution were found to have serum TPA® values of less than 75 U/L. The distribution is shown in the figure below.

Distribution of TPA® values in healthy individuals

Test characteristics

Precision

Mean (U/L)	Coefficient of variation	
	CV (%) within run	CV (%) between run
71	4.6	7.0
201	2.6	1.7
349	2.8	4.0
771	2.1	4.1
1 553	2.8	4.1

Sensitivity

Minimum measurable TPA[®] value (B_0 value + 3 SD) <15 U/L.

Measuring range

The measuring range of the TPA[®]-M IRMA kit is up to 4000 U/L.

Values above 4000 U/L can be measured after diluting the sample with TPA[®]-M IRMA Diluent and repeating the assay.

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1. Nathrath WBJ, Heidenkummer P, Björklund V, Björklund B: Distribution of tissue polypeptide antigen (TPA) in normal human tissues: Immunohistochemical study on unfixed, methanol-, ethanol-, and formalin-fixed tissues. *J Histochem Cytochem* 1985; 33:99–109.
2. Kirsch J, Oehr P, Winkler C: Localization of tissue polypeptide antigen in interphase HeLa cells by immunofluorescence microscopy. *Tumordiagn Ther* 1983; 4:189–228.
3. Kirsch J, Oehr P, Winkler C: Localization of tissue polypeptide antigen in HeLa cells during mitosis as revealed by immunofluorescence microscopy. *Tumordiagn Ther* 1984; 5:122–6.
4. Weber K, Osborn M, Moll R, Wiklund B, Lüning B: Tissue polypeptide antigen (TPA) is related to the non-epidermal keratins 8, 18 and 19 typical of simple and non-squamous epithelia: re-evaluation of a human tumour marker. *EMBO J* 1984; 3:2707–14.
5. Lüning B, Nilsson U: Sequence homology between tissue polypeptide antigen (TPA) and intermediate filament (IF) proteins. *Acta Chem Scand B* 1983; 731–3.
6. Mellerick DM, Osborn M, Weber K: On the nature of serological tissue polypeptide antigen (TPA); monoclonal keratin 8, 18 and 19 antibodies react differently with TPA prepared from human cultured carcinoma cells and TPA in human serum. *Oncogene* 1990; 5:1007–1017.
7. Lüning B: New aspects on the specificity of tumour marker TPA. *Lab Med* 1984; 8:445.
8. Oehr P, Vogel J: TPA release from cells: Mechanism for induction of elevated plasma levels in patients with TPA positive tumours. 12 Ann Meet ISOBM, Houston Texas October 1984 Abstr no 106.
9. Barak M, Steiner M, Finkel B, Abrahamson J, Antal S, Gruener N: CA-15.3, TPA and MCA as markers for breast cancer. *Eur J Cancer* 1990; 26:577–580.
10. Maffei S, Rudoni M, Antonini G, Sacchetti G, Krengli M: CA 15-3 associated with CEA and TPA in the follow-up of breast carcinoma (Ita). *Minerva Med* 1990; 81:759–763.
11. Müller MM, Griesmacher A, Hölzel W: Diagnostic relevance of tumour-marker profiles. *Int Congr Ser* 1989; 807:33–50.
12. Söletormos G, Nielsen D, Schiöler V, Skovsgaard T, Dombrowsky P: Does addition of CA549 to CA15.3, CEA and TPA give further information in monitoring metastatic breast cancer? (Meeting abstract). *Proc Annu Meet Am Assoc Cancer Res* 1990; 31:A1008.
13. Nicolini A, Colombini C, Luciani L, Carpi A, Giuliani L: Evaluation of serum CA15-3 determination with CEA and TPA in the post-operative follow-up of breast cancer patients (eng). *Br J Cancer* 1991; 64:154–8.
14. Abbasciano V, Levato F, Zavagli G: Specificity of carcino embryonic antigen, gastrointestinal cancer-associated antigen, tissue polypeptide antigen, fibrinopeptide A and gamma-glutamyltransferase in the diagnosis and follow-up of gastric cancer. *Oncology* 1988; 45:159–61.
15. Barillari P, D Angelo F, Ramacciato G, Ricci M, Santeusano G, De Angelis R, Aurello P, Fegiz G: The role of CEA, TPA and CA 19-9 in the early detection of recurrent gastric cancer. *Cancer J* 1990; 3:139–142.
16. Ståhle E, Glimelius B, Bergström R, Pählman L: Preoperative prediction of late cancer-specific deaths in patients with rectal and rectosigmoid carcinoma. *Int J Color Dis* 1989; 4:182–187.
17. Gronowitz JS, Bergström R, Nöu E, Pählman S, Brodin O, Nilsson S, Källander CF: Clinical and serologic markers of stage and prognosis in small cell lung cancer. A multivariate analysis. *Cancer* 1990; 66:722–732.
18. Haga E: Tissue polypeptide antigen in serum and tissue in patients with lung cancer (Jpn). *Nippon Kyobu Shikkan Gakkai Zasshi* 1990; 28:595–604.
19. Ratto GB, Capponi G, de Grandi R, Augeri C, Secco GB, Fardelli R: Tumour markers. I. Their significance in the preoperative assessment of lung neoplasms. *Minerva Chir* 1990; 45:1265–72.

20. Tizzani A, Casetta G, Piana P, Cavallini A: Tumour-associated antigens in diagnosis and follow-up of superficial bladder carcinomas (Meeting abstract). *European Urol* 1990; 18 (Suppl 1):69.
21. Tizzani A, Casetta G, Cavallini A, Piana P, Piantino P: Blood and urine determinations of tissue polypeptide antigen in patients with bladder carcinoma (Ita). *Minerva Urol Nefrol* 1990; 42:69–71.
22. Oehr P, Fischer L, Kersjes W, Adolphs HD, Biersack HJ, Winkler C: Distribution, sensitivity and specificity of CEA, TPA and (CEAxTPA) marker product values in patients with carcinoma of the bladder. *Nuc Compact* 1984; 15:152–4.
23. Adolphs HD, Oehr P: Significance of plasma tissue polypeptide antigen determination for diagnosis and follow-up of urothelial bladder cancer. *Urol Res* 1984; 12:125–8.
24. Rütger U, Lütgens M, Bäuerle K, Rassweiler J, Bub P, Eisenberger F, Jipp P, Kruse-Jarres JD: Surveillance and M-VEC therapy monitoring of patients with advanced carcinoma of the urinary tract by TPA. *J Tumour Marker Oncol* 1989; 4:303–321.
25. Carbin BE, Ekman P, Eneroth P, Nilsson B: Urine-TPA (tissue polypeptide antigen), flow cytometry and cytology as markers for tumour invasiveness in urinary bladder carcinoma. *Urol Res* 1989; 17:269–272.
26. Yoshida M, Morimoto S, Uekado Y, Yasukawa S, Aoshi H, Yoshida T, Ebisuno S, Ohkawa T: Clinical studies on urinary tissue polypeptide antigen (TPA): problems and improvement in determination (Jpn). *Hinyokika Kyo* 1989; 35:217–223.
27. Lewenhaupt A, Ekman P, Eneroth P, Nilsson B, Nordström L: Tissue polypeptide antigen (TPA) as a prognostic aid in human prostatic carcinoma. *Prostate* 1985; 6:285–91.
28. Kreienberg R, Melchert F: New developments in the area of markers in diagnosis and monitoring of ovarian cancer (Ger). *Onkologie* 1985; 8:253–9.
29. Lahousen M, Stettner H, Pickel H, Urdl W, Pürstner P: The predictive value of a combination of tumour markers in monitoring patients with ovarian cancer. *Cancer* 1987; 60:2228–32.
30. Panza N, Pacilio G, Campanella L, Peluso G, Battista C, Amoriello A, Utech W, Vacca C, Lombardi G: Cancer antigen 125, tissue polypeptide antigen, carcinoembryonic antigen, and b-chain human chorionic gonadotropin as serum markers of epithelial ovarian carcinoma. *Cancer* 1988; 61:76–83.
31. Paulick R, Caffier H, Horner G, Lang K, Filbry E: Clinical significance of different serum tumour markers in gynaecological malignancies. *Cancer Detect Prev* 1985; 8:115–20.
32. Lahousen M, Stettner H, Pürstner P, Pickel H: Tumour marker combination versus second-look operation in ovarian cancer (Ger). *Zentralbl Gynäkol* 1990; 112:561–566.

The data presented in this instruction booklet have been carefully compiled from our records and from the scientific literature and we believe them to be accurate and reliable. We do not, however, give any warranties or make any representations with respect hereto, nor is freedom from any patent to be inferred.

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



SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

	Dansk	Svenska	Magyar	Česky	Ελληνικά
	Europæisk øverensstømmeise	Europeisk øverensstømmeise	European Conformity	Evropská shoda	Ευρωπαϊκή Συμμόρφωση
	Udløbsdato	Utgångsdatum	Lejárati ido	Datum expirace	Ημερομηνία λήξης
	Producent	Tilverkare	Gyártó	Výrobce	Κατασκευαστής
	Se brugsanvisning	Se bruksanvisningen	Felhasználói utmutató	Srovnejte s návodem pro použití	Συμβουλευτείτε τις Οδηγίες Χρήσης
	In vitro diagnostik.	Diagnostik in vitro.	In vitro diagnosztika	In vitro diagnostika	In vitro διαγνωστικό ιατροτεχνολογικό προϊόν.
	Temperaturgrænse	Temperatur- begrænsning.	Homérséklet tartomány	Číslo šarže	Περιορισμοί θερμοκρασίας
	Fast fase (Coatede slanger)	Solid fas (Belagda strångar)	Szilárd fázis, bevonatos cso	Teplotní omezení	Στερεή φάση (/Επικαλυμμένες χάντρες)
	Fortyndingsvæske, prøve	Provutspådningsmedel	Minta hígító	Kalibrátor	Αραιωτικό δείγματος
	Tracer: antistof mærket med ¹²⁵ I	Spårelement: antikropp betecknad med ¹²⁵ I	Tracer: ²⁵ I-tel jelölt antitest	Kontrolní sérum	Ιχνηθέτης: Αντίσωμα σημασμένο με ¹²⁵ I
	Kalibrator	Kalibrator	Kalibrátor	Protílátka značená ¹²⁵ I	Βαθμονομητής
	Kontrolserum	Kontrollserum	Kontroll	Diluent vzorku	Ορός μάρτυρα


DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

Phone: 651-439-9710
FAX: 651-351-5669

13531

30846
1/07

PRINTED IN U.S.A.