Prolifigen® NSE IRMA

Immunoradiometric assay for the quantification of
Neuron Specific Enolase (NSE)

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Intended Use

Prolifigen[®] NSE IRMA is an *in vitro* assay for the quantitative determination of neuron-specific enolase (NSE) in human serum used for the management of patients with suspected or diagnosed neuroblastoma and small cell lung carcinoma.

Principle of the Procedure

The Prolifigen[®] NSE IRMA is a monoclonal two-site single incubation immunoradiometric (sandwich) assay. The sample is incubated with a plastic bead coated with antibody to NSE and an antibody against NSE labelled with ¹²⁵I. After washing unreacted radioactive antibody off the bead, the radioactivity bound to the bead is measured using a gamma counter.

Warnings and Precautions for Users

1. General

- This product is for *in vitro* diagnostic use only.
- This kit must not be used after the expiry date printed on the package.
- Do not mix reagents from different master lots.
- The accuracy of the results of this assay is directly related to the degree of care exercised during pipetting, vortexing, aspiration, and adherence to methodology and temperature requirements.
- Sodium azide used as a preservative in this product can cause irritation, so avoid contact with the skin and mucosa.
- Avoid microbial contamination of reagents.
- Observe the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in accordance with local guidelines.
- The components of the kit as well as the patients' sera should be considered as if they were potentially infectious and should be handled with adequate precautions.
- Material Safety data sheets are available on request.

2. Radioactive material

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 5.1 μ Ci (188 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Components of kit

The following reagents are supplied with each test kit for 50 assay tubes.

1. **Prolifigen® NSE IRMA Antibody-coated beads**
1 bottle containing >50 dry beads coated with monoclonal anti-NSE antibody.
2. **Prolifigen® NSE IRMA Tracer**
1 x 11 mL vial containing ¹²⁵I labelled monoclonal anti-NSE antibody. Total radioactivity less than 188 kBq (5.1 µCi). Red coloured to aid identification.
3. **Prolifigen® NSE IRMA Calibrators**
6 x 0.5 mL vials containing lyophilized NSE material; concentrations of 0, 2, 6, 20, 60 and 200 µg/L respectively.
4. **Prolifigen® NSE IRMA Diluent**
1 x 2.5 mL vial containing TRIS buffer with bovine serum albumin (also used as zero standard).
5. **Prolifigen® NSE IRMA Control, 1 & 2**
2 x 0.5 mL vial containing NSE in buffer, lyophilized.

Note: Reagents 2, 3, 4 & 5 above contain <0.1 % sodium azide as a preservative.

Each kit of Prolifigen® NSE IRMA is sufficient for 17 samples, controls and one calibrator curve in duplicate. If larger numbers of tubes are required, use pooled and mixed reagents from kits bearing the same kit lot number only. Mix gently to avoid foaming.

Prolifigen® NSE IRMA kit should be stored at 2–8°C prior to use.

Other materials required but not provided

- Polystyrene round bottomed test tubes (12x70–75 mm), disposable.
- Test tube rack.
- Micropipette (25 µL) and microdispensers (200 µL, 500 µL, 2.5 mL).
- Purified water of good quality.
- Vortex mixer.
- Bead washing system, or 2 mL repeating dispenser for addition of wash fluid and a vacuum aspiration device.
- Gamma counter set for ¹²⁵Iodine. Efficiency of counter should be >40%.
- Plastic forceps or bead dispenser for adding beads to tubes.
- Facilities for handling and disposal of radioactive material including protective gloves and other personal protective devices.
- Computer with program for evaluation of IRMA type curves or calculator and graph paper, lin-log type.

Traceability of calibrators and controls

The calibrators and controls of the kit consist of neuron-specific enolase (NSE). This NSE originates from purified and well-characterized material. Calibrators and Controls are calibrated against in-house reference NSE material.

Specimen collection, preparation & storage

- Serum only can be used for the quantitative determination of NSE. Plasma is not recommended as NSE may be released from the thrombocytes.
- Blood samples should be taken before any treatment is given, either initial or repeat.
- Collect samples using standard procedures. *Prolonged storage of blood can release NSE from the blood cells.*
- If the specimen is to be assayed within 72 hours, store at 2-8°C. If the specimen will not be assayed within 72 hours, freeze to below -18°C.
- Do not use samples which are visibly haemolytic (the absorbance at 500 nm of the sample should not exceed 0.3), contaminated or grossly lipaemic.
- Samples which have been found to contain >200 µg/L can be diluted with NSE IRMA diluent to get exact values. Diluted samples must be gently but thoroughly mixed prior to assay and the dilution factor taken into account in the final calculation.
- Frozen sera must be gently but thoroughly mixed after thawing.

Preparation of reagents before use

To each of the vials of NSE IRMA Calibrators and Controls, add 500 µL purified water. Mix gently, allow to stand for 15 min at room temperature (18–30°C) and mix gently just prior to use.

Allow kit contents to warm to room temperature before use. Mix each vial well before use.

Storage of reagents after use


If all the kit is not used at one time, the NSE IRMA Tracer and NSE IRMA Antibody-coated beads can be stored at 2–8°C until the expiry date of the kit. The NSE IRMA Calibrators and the NSE IRMA Controls must be used within 72 hours if stored at 2–8°C or else frozen to below -18°C. Repeated freezing and thawing must be avoided. Bacterial contamination of stored reagents must be avoided.

Assay procedure

Perform each determination in duplicate.

1. Add 25 µL NSE IRMA Calibrator, patient serum or NSE IRMA Control to each appropriately marked tube.
2. Add 200 µL NSE IRMA Tracer to each tube. (Solution has a red colour.) To estimate total counts (T), add 200 µL Tracer to a tube. Count the radioactivity without further processing.
3. Add one plastic bead to each tube using plastic forceps or a bead dispenser.
4. Incubate the tubes for 3±0.5 hours at room temperature (18–30°C).
5. Wash each bead 3 times with 2 mL purified water each time.
6. Count radioactivity in a gammacounter.

Flow Diagram of Method

Processing of results

Use a computer with a program for handling IRMA type data to evaluate the NSE levels in reference sera and patient sera. Alternatively, counts or bound activity (% B/T) may be plotted manually against analyte concentration. For this purpose, lin-log diagram paper is recommended.

Quality Control

It is recommended that the laboratory has its own set of quality control calibrators, in addition to the NSE IRMA Controls 1 and 2 provided in the kit.


Example of results

This is an **example** only and should not be used in the calculations.

Sample	CPM	% B/T	NSE µg/L
Total (T)	132 601	–	–
B ₀	720	0.5	–
Calibrator 2 µg/L	1 369	1.0	2.0
Calibrator 6 µg/L	2 655	2.0	6.0
Calibrator 20 µg/L	6 043	4.6	20.0
Calibrator 60 µg/L	16 589	12.5	60.0
Calibrator 200 µg/L	35 488	26.8	200.0
Patient sample 1	3 377	2.6	9.5
Patient sample 2	6 195	4.7	20.5
Patient sample 3	13 930	10.5	48.6

Example of calibrator curve


This is an **example** only and should not be used in the calculations.

Expected Values and Limitations of the procedure

In a study of 100 apparently healthy individuals, >95 % were found to have serum NSE values of less than 12.5 µg/L. The distribution is shown in the figure below.

Distribution of NSE in normal subjects

Test characteristics

Precision

Precision in 10 determinations using 5 different lots of Prolifigen® NSE IRMA.

Sample	Mean µg/L	Coefficient of variation Within assay	% Total
1	7.8	8.2	10.6
2	14.5	9.8	11.6
3	58.7	7.5	11.2
4	112.3	8.4	14.2


Recovery

To a serum sample containing 7.7 µg/L of NSE, varying quantities of NSE were added and the concentration determined:

NSE added µg/L	Expected µg/L	Observed µg/L	Recovery %
170.0	177.7	178.3	100.3
34.0	41.7	41.9	100.5
17.0	24.7	24.7	100.0
3.4	11.1	10.9	98.2
1.7	9.4	9.0	95.7

Dilution curve

Dilution of a serum sample with NSE IRMA diluent.

Sensitivity

Minimum measurable NSE value (B_0 value + 3 SD) 0.5 µg/L.

Measuring range

The measuring range is up to 200 µg/L. Values above 200 µg/L can be measured after diluting the sample with NSE IRMA diluent and repeating the assay.

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1. Ariyoshi Y, Kanefusa KI, Ishiguero Y et al: Evaluation of serum neuron-specific enolase as a marker for carcinoma of the lung. *Gann* 74 (1983) 219–225.
2. Carney DNB, Ihde DC, Cohen MH et al: Serum neuron specific enolase: a marker for disease extent and response to therapy of small-cell lung cancer. *Lancet* 1982;i, 583–585.
3. Cooper EH, Pritchard J, Baily CC, Ninane J: Serum neuron specific enolase in childrens cancer. *Br J Cancer* 56 (1984) 65–67.
4. Esscher T, Steinholtz L, Bergh J et al: Neuron specific enolase: A useful diagnostic serum marker for small cell carcinoma of the lung. *Thorax* 40 (1985) 85–90.
5. Fujita K, Haimoto H, Imaizumi M et al: Evaluation of g-enolase as a tumor marker for lung cancer. *Cancer* 60 (1987) 362–369.
6. Gerbitz K-D, Summer J, Schumacher I et al: Enolase isoenzymes as tumour markers. *J Clin Chem Clin Biochem* 24 (1986) 1009–1016.
7. Johnson DH, Marangos PJ, Forbes JT et al: Potential utility of serum neuronspecific enolase levels in small cell carcinoma of the lung. *Cancer Res* 44 (1984) 5409–5414.
8. Pålman S, Esscher T, Bergh J et al: Neuron specific enolase as marker for neuroblastoma and small cell carcinoma of the lung. *Tumor Biology* 5 (1984) 119–126.
9. Phålman S, Esscher T, Bergvall T et al: Purification and characterization of human neuron-specific enolase: radioimmunoassay development. *Tumor Biol* 5 (1984) 127–139.

The data presented in this instruction booklet has been carefully compiled from our records and from the scientific literature and we believe them to be accurate and reliable. We do not, however, give any warranties or make any representations with respect hereto, nor is freedom from any patent to be inferred.

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





SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

Continued from previous page - Symbols used with devices

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.