**Prolifigen<sup>®</sup> TK-REA**  
Radioenzyme Assay for Thymidine Kinase

**Instruction Manual**

Manuel d'Instructions

Testanleitung

Manual de Instrucciones

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Brugsanvisning

Bruksanvisning

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**REF: 324.250**

**DiaSorin**

**Stillwater, Minnesota 55082-0285, U.S.A.**

## TABLE OF CONTENTS

English .....	Page 1
Français .....	7
Deutsch .....	13
Español .....	19
Italiano.....	25
Dansk .....	31
Svenska .....	37
Česky .....	43

## Intended Use

Prolifigen® TK-REA is an *in vitro* assay for the quantitative determination of Thymidine Kinase in human serum.

## Precautions for Users

### A. General

1. This product is for *in vitro* diagnostic use only.
2. This kit must not be used after the expiry date printed on the package label.
3. Do not mix reagents from different master lots.
4. The accuracy of the results of this assay is directly related to the degree of care exercised during reconstitution, pipetting, vortexing, aspiration, and adherence to methodology and temperature requirements.
5. Sodium azide used as a preservative in this product can cause irritation, so avoid contact with the skin and mucosa.
6. Avoid microbial contamination of reagents.
7. Material Safety data sheets are available on request.

### B. Radioactive material

#### Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 5  $\mu\text{Ci}$  (185 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

**WARNING:** This product contains a chemical known to the State of California to cause cancer.

**ATTENTION:** Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

### Components of the kit

The following components are supplied with each test kit for 50 assay tubes.

1. **TK-REA Assay Buffer**  
1 vial containing buffer, lyophilized
2. **TK-REA Radioactive Substrate**  
2 vials containing <sup>125</sup>I-iododeoxyuridine, total radioactivity less than 88.8 kBq (2.4 µCi) per vial. Total kit contains less than 185 kBq (5 µCi) radioactivity.
3. **TK-REA Sample Diluent, also used as 0-Calibrator**  
1 vial containing diluent, lyophilized
4. **TK-REA Calibrator**  
1 vial containing TK *Calibrator* (40 U/L alt. 80 U/L) in bovine serum, lyophilized
5. **TK-REA Wash Fluid, 30 x concentrated**  
1 vial containing 30 x concentrated wash fluid
6. **TK-REA Separator Tablets**  
1 bottle containing >50 separator tablets
7. **TK-REA Controls, 1 and 2**  
2 vials containing TK in bovine serum, lyophilized.

**Note:** Reagents 4, 5 and 7 above contain <0.1% sodium azide as a preservative.

### Other materials required but not provided

1. Polystyrene round bottomed test tubes (12x70–75 mm), disposable.
2. Micropipette (20 µL), microdispenser (500 µL) and dispenser (1.5 mL).
3. Purified water.
4. Multi-vortex mixer.
5. Waterbath, set at 37° C ± 2° C.
6. Aspiration system for washing, e.g. wash comb.
7. Gamma counter suitable for <sup>125</sup>I.

### Specimen collection, preparation and storage

1. Blood samples should be taken before any treatment is given, either initial or repeat.
2. Serum is recommended for the quantitative determination of TK.
3. Collect samples using standard procedures.
4. If the specimen is to be assayed within 24 hours, store at 2 – 8 ° C.  
If the specimen will not be assayed within 24 hours, freeze to below –18° C.
5. Do not use samples which are grossly lipaemic, haemolysed or contaminated.
6. Samples which have been found to contain >40 U/L can be diluted with the Sample Diluent to get an exact value. The diluted sample must be mixed very well prior to assay and the dilution factor taken into account in the final calculation.
7. Frozen sera must be thoroughly mixed after thawing.
8. It is recommended that the laboratory has its own set of QC controls in addition to the TK-REA Controls, 1 and 2 provided in the kit.

### Preparation of reagents before use

#### 1. TK-REA Radioactive Substrate

Add 27 mL of purified water to the TK-REA Assay Buffer. Mix to ensure total reconstitution.

Before opening the TK-REA Radioactive Substrate, hold it upright and tap it several times on the laboratory work surface to bring the droplets to the bottom of the vial, then open and add 13 mL of TK-REA Assay Buffer to each of the vials of TK-REA Radioactive Substrate. Replace cap securely and mix vial contents thoroughly by repeated inversion.

If all the Substrate is used at the same time, mix the contents of both vials before use.

**2. TK-REA Sample Diluent**

Add 4 mL purified water to the TK-REA Sample Diluent. Mix to ensure total reconstitution.

**3. TK-REA Calibrator, 40 U/L or 80 U/L**

Add 500 µL of purified water, to obtain a stock solution of 40 U/L or 250 µL of purified water to obtain a stock solution of 80 U/L, to the TK-REA calibrator.

Mix well and allow to stand for 15 minutes. Mix well before use. To construct a calibrator curve, a dilution with Sample Diluent of calibrator point 40 U/L (alt. 80 U/L) is needed according to the following:

Calibrator	Dilution
40 U/L	–
10 U/L	100 µL Cal 40+300 µL Sample Diluent
5 U/L	100 µL Cal 10+100 µL Sample Diluent
2.5 U/L	100 µL Cal 5+100 µL Sample Dilluent

A calibrator curve consisting of points 40, 10, 5 and 2.5 U/L is recommended.

**4. TK-REA Wash Fluid**

Add the contents of the TK-REA Wash Fluid to purified water to make a total volume of 600 mL.

**5. TK-REA Controls**

To each of the vials of TK-REA Controls, add 250 µL of purified water. Mix well and allow to stand for 15 minutes. Mix well before use.

**Storage of reagents after use**

If all the kit is not used at one time, the diluted Wash Fluid and the Separator Tablets can be stored in their original containers at 2–8° C until the expiry date printed on the carton label.

The reconstituted Calibrator Stock Solution, Sample Diluent and Assay Buffer are stable 4 weeks when stored at –20° C.

The reconstituted Controls are stable up to three months when stored at –20° C.

The diluted Radioactive Substrate and the diluted calibrator points can **not** be stored and have to be diluted for each assay run.


**Assay procedure**

Perform each determination in duplicate

1. Add 20 µL of TK-REA calibrator, patient sample and TK-REA control to each appropriately marked tube.
2. Add 500 µL TK-REA Radioactive Substrate to each tube.
3. Vortex tubes and incubate for 4 hours at 37° C in a waterbath.
4. Add one TK-REA Separator tablet to each tube using plastic forceps or a small spoon.
5. Vortex tubes for 15 min. on a multi-vortex mixer.  
– or vortex for 1 min. and then let stand for 60 min.

6. Add 1.5 mL diluted wash fluid to each tube and mix well by vortexing.
7. Allow suspension to settle for 1 min. and remove supernatant by aspiration using a wash comb. Take care not to remove any of the sediment.  
**Note:** The supernatant will appear cloudy due to the filler used in the manufacture of the tablets. Removal of this filler will not affect the results. The active ingredient of the tablet will sediment within 30 seconds.
8. Repeat washing steps 6 and 7 a further three times (a total of 4 washes).
9. Count radioactivity in a gamma counter.

**Flow Diagram of Method**

**Processing of results**

Use a computer with a program for handling immunoassay data to evaluate the TK levels in the control sera and the patient sera. Alternatively, counts or bound activity %B/T may be plotted manually against analyte concentration. For this purpose log-log diagram paper is recommended.

The background value (0-cal) should be subtracted from all values to give a corrected cpm-value.


**Example of results**

This is an **example** only and should not be used in calculations.

Sample	CPM (corrected)	B/T %	TK U/L
Total	154 209	–	–
B <sub>0</sub>	656	0.4	0
Calibrator 2.5	2 543	1.6	2.5
Calibrator 5	4 085	2.6	5.0
Calibrator 10	8 116	5.3	10.0
Calibrator 40	20 569	16.2	40.0
Patient sample 1	2 483	1.6	2.4
Patient sample 2	9 234	6.0	11.4
Patient sample 3	13 887	9.0	18.3

### Example of Calibrator curve


This is an **example** only and should not be used in the calculations.

### Reference range

In a study of 160 blood donors, >95 % had a serum TK value of less than 6.1 U/L. The distribution is shown in the figure below.

### Distribution of TK values in blood donors


## Test characteristics

### Precision

Mean U/L	% cv within	% cv between
4.9	5.5	7.7
12.9	3.0	7.0
19.5	4.0	4.1
34.8	3.5	4.6

### Dilution

Dilution of a serum sample with TK-REA Sample Diluent.

### Sensitivity

Minimum measurable TK value ( $B_0$  value + 3 SD) 0.6 U/L.

### Measuring range

The measuring range is up to 40 U/L. Values above 40 U/L can be measured after diluting the sample with TK-REA Sample Diluent and repeating the assay.

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1. Gronowitz, J.S., Källander, C. F. R., Diderholm, H., Hagberg, H. and Pettersson, U. Application of an *in vitro* assay for serum thymidine kinase: results on viral disease and malignancies in humans. *Int. J. Cancer*, **33**, 5-12 (1984)
2. Gronowitz, J.S., Hagberg, H., Källander, C. F. R and Simonsson, B. The use of deoxythymidine kinase as a prognostic marker, and in the monitoring of patients with non-Hodgkin's Lymphoma. *Br. J. Cancer*, **47**, 487-495 (1983).
3. Hagberg, H., Glimelius, B., Gronowitz, S., Killander, A., Källander, C. and Schröder, T. Biochemical markers in non-Hodgkin's Lymphoma stages III and IV and prognosis: A multivariate analysis. *Scand. J. Haematol.*, **33**, 59-67 (1984)
4. Källander, C. F. R, Simonsson, B., Hagberg, H. and Gronowitz, J.S. Serum deoxythymidine kinase gives prognostic information in Chronic Lymphocytic Leukaemia. *Cancer*, **54**, 2450-2455 (1984)
5. Hagberg, H., Gronowitz, S., Killander, A., Källander, C., Simonsson, B., Sundström, C. and Öberg, G. Serum Thymidine kinase in acute leukaemia. *Br.J. Cancer*, **49**, 537-540 (1984)
6. Eriksson, B., Hagberg, H., Glimelius, B., Sundström, C., Gronowitz, S. and Källander, C. Serum thymidine kinase as a prognostic marker in Hodgkin's disease. *Acta Radiol. Oncol.*, **24**, 167-171 (1985)
7. Simonsson, B., Källander, C.F.R., Brenning, B., Killander, A., Åhre, A. and Gronowitz, J.S. Evaluation of serum deoxythymidine kinase as a marker in multiple myeloma. *Br.J. Haematol.*, **61**, 215-244 (1985)
8. Ellims, P.H., Gan, T.E., Medley, G. and Van der Weyden, M.B. Prognostic relevance of Thymidine Kinase Isozymes in Adult Non-Hodgkin's lymphoma. *Blood*, **58**, 926-930 (1981)
9. Gronowitz, J.S., Steinholtz, L., Källander, C.F.R., Hagberg, H. and Bergh, J. Serum deoxythymidine kinase in small-cell carcinoma of the lung: relation to clinical features, prognosis and other biochemical markers. *Cancer*, **58**, 111-118 (1986)
10. Larsson, A., Frödin, L., Tufveson, G., Larsson, E., Källander, C.F.R., and Gronowitz, J.S. Deoxythymidine kinase as a marker monitoring the treatment of cytomegalovirus infection with Foscarnet in a renal transplant patient. *Scand. J. Urol. Nephrol.*, **20**, 75-76 (1986)
11. Gronowitz, J.S., Larson, A., Källander, C.F.R., Claesson, K., Sjöberg, O., Lernestedt, J.-O., Frödin, L. and Tufveson, G. Serum thymidine kinase in transplantation patients: its relation to cytomegalovirus activity, renal transplant rejection and its use for monitoring anti-viral therapy. *Ann. Clin. Research*, **18**, 71-75 (1986)
12. Gronowitz, J.S., Tötterman, T.H., Källander, C.F.R., Hagström, Å. and Tufveson, G. Serum thymidine kinase, CMV-specific antibodies and total immunoglobulins in renal transplant patients during immunosuppression with cyclosporine. *A. Clin. Transplantation*, **2**, 26-35 (1988)
13. Hagberg, H., Gronowitz, J.S., Killander, A., and Källander, C. F. R. Serum thymidine kinase in Vitamin B<sub>12</sub> deficiency. *Scand.J. Haematol.*, **32**, 41-45 (1984)
14. Hallek, M., Schick, H.D., Dusch, R., Strohmeyer, S. at al. Serum Thymidine kinase in low and intermediate grade non-Hodgkin's lymphoma. *Folia Oncologica*, **13**, 87, (1990)
15. Hallek, M., Wanders, L., Strohmeyer, S., Emmerich, B., Thymidine kinase: a tumour marker with prognostic value for non-Hodgkin's lymphoma and a broad range of potential clinical applications. *Ann. Hematol.* **65**, 1-5 (1992)

The data presented in this instruction booklet have been carefully compiled from our records and from the scientific literature and we believe them to be accurate and reliable. We do not, however, give any warranties or make any representations with respect hereto, nor is freedom from any patent to be inferred.

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SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

**CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES**

	Dansk	Svenska	Česky
	Europæisk øverensstämmeise	Europeisk øverensstämmeise	Evropská shoda
	Udløbsdato	Utgångsdatum	Datum expirace
	Producent	Tilverkare	Výrobce
	Se brugsanvisning	Se bruksanvisningen	Srovnejte s návodem pro použití
<b>IVD</b>	In vitro diagnostik.	Diagnostik in vitro	In vitro diagnostika
	Temperaturgrænse	Temperatur- begrænsning	Číslo šarže
<b>SORB</b>	Fast fase (Coatede slanger)	Solid fas (Belagda strångar)	Teplotní omezení
<b>DIL</b>   <b>SPE</b>	Fortyndingsvæske, prøve	Provutspådningsmedel	Kalibrátor
<b>Ab</b>   <sup>125</sup> I	Tracer: antistof mærket med <sup>125</sup> I	Spårelement: antikropp betecknad med <sup>125</sup> I	Kontrolní sérum
<b>CAL</b>	Kalibrator	Kalibrator	Protilátka značená 125I
<b>CONTROL</b>	Kontrolserum	Kontrollserum	Diluent vzorku


DiaSorin Inc.  
1951 Northwestern Avenue  
P.O. Box 285  
Stillwater, MN 55082-0285

DiaSorin S.p.A.,  
13040 Saluggia (VC) Italy

Phone: 651-439-9710  
FAX: 651-351-5669

13533

30845  
1/07

PRINTED IN U.S.A.