CA 125 II™
¹²⁵I IRMA Kit

For the quantitative determination of
CA 125 antigen in serum

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For professional use only!

Intended Use

In vitro test for the quantitative determination of CA 125™ antigen (OC 125 defined antigen) in serum as well as in EDTA plasma and heparin plasma during follow-up of patients with primary invasive ovarian carcinoma.

Summary and Explanation of the Test

The two antigen determinants defined by the monoclonal antibodies OC125(1) and M11(5) are found on a heterogeneous group of macromolecular glycoproteins (molecular weight 200,000–1,000,000).

These glycoproteins can be detected in most nonmucinous epithelial ovarian carcinomas, besides some fetal tissues (amnion, periderm, derivatives of the coelomic epithelium) and in the epithelium of the fallopian tubes, endometrium, apocrine sweat glands, mammary glands and endocervix.

Elevated CA 125II assay values are found in the serum of most patients with an active ovarian carcinoma in early stages of the disease (4,6).

The diagnostic value of the CA 125II determination lies in the monitoring of the therapy and progression of patients with ovarian carcinoma (9).

Principles of the Procedure

Immunoradiometric assay based upon the "sandwich principle". The high affinity monoclonal antibody M11 and the monoclonal antibody OC 125 are used for the coating of the solid phase (coated tube) and for the tracer.

One-Step Method

Coating antibodies and tracer antibodies react in one reaction step with the CA 125 present in the calibrators or patient samples. Unbound material is removed by a washing step. After washing out the excess tracer the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

Two-Step Method

During the first incubation, the CA 125 present in patient samples and calibrators is bound by the antibody immobilized on the test tube wall. Unbound material is removed by a washing step. During the second incubation, the tracer antibody reacts with the CA 125 already bound. After removal of excess tracer by a second washing step, the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

CA125™ is a trademark of Fujirebio Diagnostics Inc.

CONTENTS Material Provided

Determinations	50
¹²⁵ I-anti-CA 125, monoclonal (mouse), red radioactive content (kBq / μ Ci)	11 mL < 185 /5
6 calibrators A-F (1.0 mL) in human serum albumin	1 set
Assay buffer, blue	22 mL
Diluent (0 U/mL) in human serum albumin	6 mL
Test tubes, coated with anti-CA 125, monoclonal (mouse)	50
Control, human, lyophilic	1 mL

Material required but not provided

- Micropipettes (100 μ L, 200 μ L, 300 μ L, 1000 μ L) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively, an appropriate automated analyser system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- Purified water

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 4.83 μ Ci (179 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:
The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly. **(Do not allow foam to form.)**

Open the control carefully and reconstitute with 1.0 mL of purified water. **(Do not allow foam to form.)** Ensure that lyophilised material adhering to the cap is also dissolved.

Storage of Reagents

Reconstituted control: 1 week at 2–8°C or 4 weeks at –20°C.

Store all other reagents at 2–8°C. Use before the expiry date printed on the packaging.

- **Store upright.**

- **Keep away from direct light.**

Sample Collection, Material and Storage

- Collect samples using standard procedures.
- Sample material: serum, EDTA and heparin plasma
- Storage at 2–8°C: 24 hours
- To store for longer periods, freeze samples at –20°C
- Freezing and thawing samples once does not affect the test results.
- Stored samples should be thoroughly mixed prior to use (vortex mixer).
- Do not use sera which are agglutinated, lipemic, hemolytic, icteric or contaminated.

Interfering Substances

No interference with test results has been observed with concentrations of bilirubin < 0.125 mg/mL, hemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The individual components of the kit are carefully matched to one another. If any components from different lots are exchanged or mixed up, the manufacturer does not guarantee reliable results.
- Strictly adhere to the sequence of the pipetting steps.
- Set the measuring time to at least 1 minute.
- Observe the quality control guidelines for medical laboratories.
- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that calibrators and samples are run in duplicate.

If values above the highest calibrator are expected, samples should be further diluted with the diluent (e.g. by dilution factors 10, 100, 1000).

As an alternative to manual performance and evaluation, an appropriate automated analyzer system can also be used at the responsibility of the laboratory.

One-Step Method

1. Pipette 100 μL of calibrator, control or patient sample into the bottom of a coated test tube.
2. Add 200 μL ^{125}I -anti-CA 125, mix (vortex mixer). Do not confuse with the assay buffer (blue).
3. Incubate the tubes for 3 hours (\pm 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 3 times with 2 mL of purified water.
6. Measure radioactivity (CPM) in all tubes (at least 1 min).

100 μL	Pipette calibrator, control or patient sample into the bottom of a test tube
200 μL	Add ^{125}I -anti-CA 125 Do not confuse with assay buffer (blue)! Mix
3 hrs (\pm 5 min)	Incubate at room temperature (18–25°C) on a shaker* Aspirate
3 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

*** Keep agitation conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

Two-Step Method

Alternatively to the one-step protocol, the test may also be performed according to a two-step protocol.

1. Pipette 100 μL each of calibrator, control or patient sample into the bottom of a coated test tube.
2. Add 200 μL assay buffer (blue), mix (vortex mixer).
3. Incubate the tubes for 3h (\pm 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 3 times with 2 mL purified water.
6. Add 100 μL of assay buffer.**
7. Add 200 μL ^{125}I -anti-CA 125, mix (vortex mixer)**

** Assay buffer and ^{125}I -anti-CA 125 solution may also be pooled (1+2, only in glass containers). 300 μL of this solution can then be pipetted.

8. Leave all tubes for 20 – 24 h at room temperature (18–25°C).
9. Aspirate the liquid.
10. Wash all tubes 3 times with 2 mL of purified water.
11. Measure radioactivity (CPM) in all tubes (at least 1 min.).

100 µL	Pipette calibrator, control or patient sample into the bottom of a test tube
200 µL	Add assay buffer (blue)
3 hrs (± 5 min)	Incubate at room temperature (18–25°C) on a shaker* Aspirate
3 x 2 mL	Wash with purified water
100 µL	Add assay buffer (blue)**
200 µL	Add ¹²⁵ I-anti-CA 125**
20-24 hours	Leave to incubate at room temperature (18–25°C) Aspirate
3 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

* Keep agitation conditions constant!

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

** Assay buffer and ¹²⁵I-anti-CA 125 solution may also be pooled (1+2, only in glass containers). 300 µL of this solution can then be pipetted.

Calculation of Results

The calibrator curve can be established manually as follows:

1. Determine the mean CPM for each pair of tubes (double determination).
2. Divide the mean CPM of each calibrator (B) by the mean CPM of the highest calibrator (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding ($\%B/B_{max}$) for each calibrator.
3. On semi-log paper, plot the relative bindings of each calibrator ($\%B/B_{max}$) on the Y-axis versus the corresponding concentrations (U/mL) on the X-axis.
4. Sample concentrations (U/mL) can be read directly off the calibrator curve from their corresponding relative binding ($\%B/B_{max}$).

If the radioactivity measured is above that of the highest calibrator, the samples must be diluted with the diluent and tested again. With diluted samples the actual serum concentration and the appropriate dilution factor have to be established.

The instrumental calculation of radioimmunological measured values is performed using a spline approximation.

Quality Control

Observe the quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory. This control should be simultaneously tested with each test run and treated like patient samples. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Interpretation of the Results

Evaluations of 372 subjects have confirmed the reference range of 35 U/mL for obviously healthy subjects (n=25) given in the literature (12).

Variations from these concentrations have been found, however, depending on age (2) and menstrual cycle (11).

Since CA 125 values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Patients with carcinomas may exhibit CA 125II assay values below the cut-off.

Elevated values may be observed in patients with benign diseases such as pericarditis, severe liver diseases, severe endometriosis, ovarian cysts or with carcinoma of the uterus, pancreas, liver or lung (3,8,10). In pregnancy, CA 125II assay values may also rise above the cut-off (7).

CA 125II levels alone may therefore not be used as proof of the presence or exclusion of a malignant disease and should only be interpreted in conjunction with the clinical picture and other diagnostic procedures.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may in principle lead to falsely elevated or decreased values. HAMA-neutralizing agents are added to the test. However the influencing of results cannot be completely ruled out. These samples should not be used for the CA 125 II IRMA.

Analytical Data

Calibration

The CA 125II IRMA has been calibrated using the CA 125II assay from Fujirebio.

Measuring range

Measuring range is 1 - 500 U/mL.

High-dose hook

With the one step method no high-dose hook effect was observed for concentrations up to 50,000 U/mL, whereas with the two step method no high dose hook was found for a concentration up to 100,000 U/mL.

Precision

Intra-assay variation			Inter-assay variation		
Mean value (U/mL)	CV (%)	n=	Mean value (U/mL)	CV (%)	n=
42	4.0	9	44	4.8	13
97	4.1	9	98	7.7	10
339	5.5	9	309	6.1	5

Analytical sensitivity

The lower detection limit is < 1.0 U CA 125/mL. This detection limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest CA 125II concentration that can be differentiated from zero with statistical significance.

Specificity

Interferences with drugs

No cross-reactivities with cyclophosphamide, doxorubicin, fluorouracil, methotrexate and tamoxifen in therapeutic ranges were found.

Linearity upon Dilution

A patient's serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 298.3 U/mL

Dilution	Measured value (U/mL)	Expected value (U/mL)	Recovery (%)
1 : 1.25	223	219	102
1 : 2.5	114	112	101
1 : 5	64	59	108
1 : 10	31	33	95

Recovery

A patient's serum with low CA 125 content was spiked with different amounts of CA 125.

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1. Bast RC, Feeney M, Lazarus H, Nadler LM, Colvin RB, Knapp RC. Reactivity of a Monoclonal Antibody with Human Ovarian Carcinoma. *J Clin Invest* 1981; **68**: 1331-1337
2. Bonfrère JMG, Korse CM, Verstraeten RA, Van Kamp GJ, Hart GAM, Kenemans P. Clinical evaluation of the Byk LIA-mat CA 125II assay: discussion of a reference value. *Clin Chem* 1997; **43** (3): 491-497
3. Daoud E, Bodor G. CA-125 Concentrations in Malignant and Nonmalignant Disease (Washington University Case Conference). *Clin Chem* 1991; **37** (11): 1968-1974
4. Frasci G, Conforti S, Zullo F, Mastrantonio P, Comella G, Comella P, Persico G, Iaffaioli RV. A Risk Model for Ovarian Carcinoma Patients Using CA 125. *Cancer* 1996; **77**: 1122-1130
5. Hardardóttir H, Parmley TH, Quirk JG, Sanders MM, Miller FC, O'Brien TJ. Distribution of CA 125 in embryonic tissues and adult derivatives of the fetal periderm. *Am J Obstet Gynec* 1990; **163** (6): 1925-1931
6. Hasholzner U, Baumgartner L, Stieber P, Meier W, Reiter W, Pahl H, Fateh-Moghadam A. Clinical Significance of the Tumour Markers CA 125 and CA 72-4 in Ovarian Carcinoma. *Int J Cancer (Pred Oncol)* 1996; **69**: 329-334
7. Jacobs IJ, Fay TM, Yovich J, Stabile I, Frost C, Turner J, Oran DH, Grudzinskas JG. Serum levels of CA 125 during the first trimester of normal outcome, ectopic and anembryonic pregnancies. *Human Reproduction* 1990; **5** (1): 116-122
8. Koninckx PR, Meuleman C, Oosterlynck D, Cornillie FJ. Diagnosis of deep endometriosis by clinical examination during menstruation and plasma CA-125 concentration. *Fertil Steril* 1996; **65**(2): 280-287
9. Rustin GJS, Nelstrop AE, Tuxen MK, Lambert HE. Defining progression of ovarian carcinoma during follow-up according to CA 125: A North Thames Ovary Group study. *Ann Oncol* 1996; **7**: 361-364
10. Sproston ARM, Roberts SA, Davidson SE, Hunter RD, West CML. Serum tumour markers in carcinoma of the uterine cervix and outcome following radiotherapy. *Br J Cancer* 1995; **72**: 1536-1540
11. Zeimet AG, Daxenbichler G, Müller-Holzner E, Dapunt O, Marth C. Tumor marker CA-125 in tissues of the female reproductive tract and in serum during the normal menstrual cycle. *Fertil Steril* 1993; **59** (5): 1028-1035
12. Klug TL, Bast RC, Nilloff JM, Knapp RC, Zurawski VR. Monoclonal antibody immunoradiometric assay for an antigenic determinant (CA 125) associated with human epithelial ovarian carcinomas. *Cancer Res* 1984; **44**: 1084-53

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
IVD	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
LOT	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
CAL	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
CONTROL	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
Ab ¹²⁵ I	Tracer: antibody labelled with ¹²⁵ I	Traceur : anticorps marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markierter Antikörper	Trazador: anticuerpo etiquetado con ¹²⁵ I	Tracciatore: anticorpo etichettato con ¹²⁵ I
DIL SPE	Sample diluent	Diluant d'échantillon	Proben-Diluent	Diluyente de muestras	Diluyente campione
RCNS X mL	Reconstitute with X mL	Reconstituer avec X mL	Mit X mL auflösen	Reconstitución con X mL	Ricostituire con X mL
BUF	Buffer	Tampon	Puffer	Tampón	Tampone
SORB	Solid phase Coated tubes	Phase solide Tubes enduits	Festphase Beschichtete Teströhrchen	Fase sólida Tubos recubiertos	Fase solida Cannule con rivestimento.
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.