CA 15-3[®]
¹²⁵I IRMA Kit

For the quantitative determination of
CA 15-3 antigen in serum

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For professional use only!

Intended Use

In vitro test for the quantitative determination of the CA 15-3[®] antigen (DF3 defined antigen) in human serum during the follow-up of patients with breast carcinoma.

Summary and Explanation of the Test

The tumour-associated antigen CA 15-3 bears the epitopes of two different monoclonal antibodies, Mab 115D8 and Mab DF3. Mab 115D8 was produced against the surface antigen MAM-6 (MW 400,000) which is found on mammary gland epithelial cells and most breast carcinoma cells (7).

Mab DF3 was formed against membrane antigens of mammary carcinoma cells (8).

The monoclonal antibody reacts with an epitope termed DF3 on a glycoprotein with the molecular weight 290,000. Elevated CA 15-3 assay values are found in patients with breast carcinoma (4,6).

Due to its high sensitivity for metastatic breast carcinoma, CA 15-3 determination is particularly suitable for the diagnosis of recurrence and in monitoring recurrence therapy (9, 10, 11).

Principles of the Procedure

Immunoradiometric assay based upon the "sandwich principle". The monoclonal antibody 115D8 is used for the coating of the solid phase (coated tube) and the monoclonal antibody DF3 is used for the tracer.

The coating antibody and the tracer antibody react in one reaction step with the CA 15-3 present in calibrators and patient samples. After excess tracer is removed by a washing step, the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

CONTENTS Material Provided

Determinations	50
¹²⁵ I-anti-CA 15-3, monoclonal (mouse), red radioactive content (kBq / µCi)	6 mL < 385 /10
6 calibrators A-F (1.0 mL) in phosphate buffer (for concentration, see Quality control report)	1 set
Diluent (0 U/mL) in phosphate buffer	110 mL
Test tubes, coated with anti-CA 15-3, monoclonal (mouse)	50
Control, human, lyophilic	0.5 mL

CA 15-3[®] is a registered trademark of Fujirebio Diagnostics Inc.

Material required but not provided

- Micropipettes (20 μ L, 100 μ L, 500 μ L) with disposable plastic tips, Multipipette (2000 μ L)
- Vortex Mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively, an appropriate automated analyzer system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- Purified Water

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 10.1 μ Ci (372 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly. **(Do not allow foam to form.)**

Open the control carefully (and reconstitute with 0.5 mL of purified water. **(Do not allow foam to form.)** Ensure that lyophilized material adhering to the cap is also dissolved.

Storage of Reagents

Reconstituted control: 1 week at 2–8°C or 4 weeks at –20 °C

Store all other reagents at 2–8°C. Use before the expiry date printed on the packaging.

Opened tubes should be stored in the tube rack. They can then be used until the expiry date given on the kit label.

- **Store upright.**

- **Keep away from direct light.**

Sample Collection, Material and Storage

- Collect samples using standard procedures.

- Sample material: serum. Disturbances due to plasma are not known.

- Storage at 2–8°C: 24 hours

- To store for longer periods: freeze at –20°C

- Freezing and thawing samples once does not affect the test results.

- Stored samples should be thoroughly mixed prior to use (vortex mixer).

- Do not use sera which are agglutinated, lipemic, hemolytic, icteric or contaminated.

Interfering Substances

No interference with test results has been observed by concentrations of bilirubin < 0.125 mg/mL, haemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The individual components of each kit are carefully matched to one another. If any components from different lots are exchanged or mixed up, the manufacturer does not guarantee reliable results

- Strictly adhere to the sequence of the pipetting steps

- Set the measuring time for at least 1 minute.

- This testing kit must not be used after the expiry date printed on the packaging.

- Observe the quality control guidelines for medical laboratories.

- Avoid microbial contamination of the reagents.

- **Rack-by-rack washing is not appropriate. The washing water must be aspirated from the tubes immediately after being added.**

Test Procedure

It is recommended that calibrators and samples are double determined.

As an alternative to manual performance and evaluation, an appropriate automated analyser system can be used at the responsibility of the laboratory.

Pre-dilution of samples and controls

Patient sera and controls must be diluted prior to assaying (1:101).

1. Pipette 2000 µL of diluent into uncoated tubes.

2. Add 20 µL of control or sample and mix thoroughly.

If values above the highest calibrator are expected, samples should be further diluted with the diluent (e.g. by factors 10, 100, 1000).

Important:

All samples and controls have to be pre-diluted prior to testing, even if elevated CA 15-3 serum levels (beyond the calibrator range) necessitate additional dilution steps. The calibrators are to be used undiluted.

3. Pipette 100 µL of undiluted calibrator, pre-diluted control or patient sample into the bottom of a coated tube.
4. Add 100 µL ¹²⁵I-anti-CA 15-3, mix (vortex mixer).
5. Incubate the tubes for 2 hours (± 5 min) at room temperature (18–25°C) on a horizontal shaker.*
6. Aspirate the liquid.
7. Wash all tubes 3 times with 2 mL of purified water.
8. Measure radioactivity (CPM) in all tubes (at least 1 min.).

Pre-dilution of patient samples and controls

2000 µL	Pipette diluent into an uncoated tube
20 µL	Pipette control or patient sample
	Mix
100 µL	Pipette undiluted calibrator, pre-diluted control or pre-diluted patient sample into the bottom of a test tube
100 µL	Add ¹²⁵ I-anti-CA 15-3
	Mix
2 hrs (± 5 min)	Incubate at room temperature (18–25°C) on a shaker
	Aspirate
3 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

*** Keep agitation conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm
Amplitude 10 mm = 220 rpm
Amplitude < 8 mm = 300 rpm

Calculation of Results

The calibrator curve can be established manually as follows:

1. Determine the mean CPM for each pair of duplicate tubes.
2. Divide the mean CPM of each calibrator (B) by the mean CPM of the highest calibrator (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding (%B/B_{max}) for each calibrator.
3. On semi-log paper, plot the relative binding of each calibrator (%B/B_{max}) on the Y-axis versus the corresponding concentrations on the X-axis.
4. The sample concentrations (U/mL) can be read directly off the calibrator curve from their corresponding relative binding (%B/B_{max}).

If the radioactivity measured is above that of the highest calibrator, the samples must be diluted with the diluent and tested again. With diluted samples the actual serum concentration and the appropriate dilution factor have to be established.

The instrumental calculation of radio-immunological measured values is performed using a spline approximation.

Quality Control

Observe the quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory.

This control should be simultaneously tested with each test run and treated like patient samples. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Interpretation of the Results

In healthy adults (n=200), CA 15-3 serum concentrations of 15.5 U/mL \pm 12.5 U/mL were found (mean \pm 2 x standard deviation).

Since CA 15-3 values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Patients with carcinomas may exhibit CA 15-3 assay values within the normal range. Elevated levels may be found in benign disorders of the breast, ovaries or the liver (2,3), or in ovarian or endometrial carcinoma (1,5).

Therefore, CA 15-3 assay levels may only be interpreted in conjunction with the clinical picture and other diagnostic procedures.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may in principle lead to falsely elevated or decreased values. HAMA-neutralising agents are added to the test. However the influencing of results cannot be completely ruled out. These samples should not be used for the CA 15-3 IRMA.

Analytical Data

Calibration

The CA 15-3 IRMA has been calibrated using the Fujirebio CA 15-3 assay.

Measuring range

The measuring range is 2 - 300 U/mL.

High-dose hook

No high-dose hook effect was observed for concentrations up to 12,300 U/mL.

Precision

Intra-assay variation			Inter-assay variation		
Range (U/mL)	CV (%)	n=	Range (U/mL)	CV (%)	n=
20 - 50	2.9	17	20 - 50	4.9	3
70 - 130	2.8	53	70 - 130	4.2	3

Analytical sensitivity

The lower detection limit is < 2.0 U CA 15-3/mL. This detection limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest CA 15-3 concentration that can be differentiated from zero with statistical significance.

Specificity

Interferences with drugs

No cross-reactivities with doxorubicin, fluorouracil, methotrexate, cyclophosphamide and tamoxifen in therapeutic ranges were found.

Linearity upon Dilution

A patient serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 241 U/mL

Dilution	Measured value (U/mL)	Expected value (U/mL)	Recovery (%)
1 : 1.25	194	193	101
1 : 2.5	95	97	98
1 : 5	45	48	94
1 : 10	23	24	96

Recovery

A serum with low CA 15-3 content was spiked with different amounts of CA 15-3.

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


1. Babilonti L, Riccardi A, Tateo S, Pavesi F, Lotzniker M, Polatti F. Tumor Antigens CA 125 and CA 15-3 as Markers of Endometrial Adenocarcinoma. *J Nucl Med Allied Sci* 1990; **34 (Suppl. 4)**: 79-83
2. Collazos J, Genollà J, Ruibal A. Breast Cancer-Associated Antigen CA 15.3 in Liver Cirrhosis. *Acta Oncologica* 1992; **31 (7)**: 741-744
3. Colomer R, Ruibal A, Genollà J, Salvador L. Circulating CA 15-3 antigen levels in nonmammary malignancies. *Br J Cancer* 1989; **59**: 283-286
4. Dnistrian AM, Schwartz MK, Greenberg EJ, Smith CA, Schwartz DC. CA 15-3 and carcinoembryonic antigen in the clinical evaluation of breast cancer. *Clinica Chimica Acta* 1991; **200**: 81-94
5. Einhorn N, Knapp RC, Bast RC, Zurawski VR. CA 125 Assay used in Conjunction with CA 15-3 and TAG-72 Assays for Discrimination between Malignant and Non-Malignant Diseases of the Ovary. *Acta Oncologica* 1989; **28 (5)**: 655-657
6. Hayes DF, Zurawski VR, Kufe DW. Comparison of Circulating CA 15-3 and Carcinoembryonic Antigen Levels in Patients With Breast Cancer. *J Clin Oncol* 1986; **4 (10)**: 1542-1550
7. Hilkens J, Buijs F, Hilgers J, Hageman P, Calafat J, Sonnenberg A, Van Der Valk M. Monoclonal Antibodies against Human Milk-Fat Globule Membranes detecting Differentiation Antigens of the Mammary Gland and its Tumors. *Int J Cancer* 1984; **34**:197-206
8. Kufe D, Inghirami G, Abe M, Hayes D, Justi-Wheeler H, Schlom J. Differential Reactivity of a Novel Monoclonal Antibody (DF3) with Human Malignant versus Benign Breast Tumors. *Hybridoma* 1984; **3 (3)**: 223-231
9. Molina R, Zanón G, Filella X, Moreno F, Jo J, Daniels M, Latre ML, Giménez N, Pahisa J, Velasco M, Ballesta AM. Use of serial carcinoembryonic antigen and CA 15.3 assays in detecting relapses in breast cancer patients. *Breast Cancer Res & Treatment* 1995; **36**: 41-48
10. Vizcarra E, Lluch A, Cibrián R, Jarque F, Alberola V, Belloch V, García-Conde J. Value of CA 15.3 in breast cancer and comparison with CEA and TPA: A study of specificity in disease-free follow-up patients and sensitivity in patients at diagnosis of the first metastasis. *Breast Cancer Res & Treatment* 1996; **37**: 209-216
11. Vizcarra E, Lluch A, Cibrián R, Jarque F, García-Conde J. CA 15.3, CEA and TPA Tumor Markers in the Early Diagnosis of Breast Cancer Relapse. *Oncology* 1994; **51**: 491-496

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Tracer: antibody labelled with ¹²⁵ I	Traceur : anticorps marqué à l'iode ¹²⁵ I	Tracer: ¹²⁵ I-markierter Antikörper	Trazador: anticuerpo etiquetado con ¹²⁵ I	Tracciatore: anticorpo etichettato con ¹²⁵ I
	Sample diluent	Diluant d'échantillon	Proben-Diluent	Diluyente de muestras	Diluyente campione
	Reconstitute with X mL	Reconstituer avec X mL	Mit X mL auflösen	Reconstitución con X mL	Ricostituire con X mL
	Buffer	Tampon	Puffer	Tampón	Tampone
	Solid phase Coated tubes	Phase solide Tubes enduits	Festphase Beschichtete Teströhrchen	Fase sólida Tubos recubiertos	Fase solida Cannule con rivestimento.
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.