CA 19-9™
¹²⁵I IRMA Kit

For the quantitative determination of
CA 19-9 in serum

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IVD **For professional use only!**

Intended Use

In vitro test for the quantitative determination of the CA 19-9™ antigen (1116NS — 19-9 defined antigen) in human serum during the follow-up of patients with gastrointestinal tumours.

Summary and Explanation of the Test

The monoclonal antibody 1116NS-19-9 (6) recognizes the epitope sialyl-lacto-N-fucopentose II, which is related in biochemical terms to the Lewis^a blood group antigen. This epitope was found in tumours both on glycolipids and mucin-like glycoproteins.

CA 19-9 assay levels are elevated in patients with pancreatic tumours or with hepatic, gastric and colorectal tumours and are suitable for follow-up and therapeutic monitoring (1, 2, 5, 7, 8, 9).

Principles of the Procedure

Immunoradiometric assay based upon the “sandwich principle”. The same highly specific monoclonal antibody 1116NS-19-9 is used for the coating of the solid phase (coated tube) and the tracer.

During the first incubation, the CA 19-9 present in patient samples and calibrators is bound to the antibody immobilized on the test tube wall. Unbound material is removed by a washing step. During the second incubation, the tracer antibody reacts with the CA 19-9 already bound.

After removal of excess tracer by a second washing step, the radioactivity bound at the tube wall is measured in a gamma scintillation counter.

CONTENTS Material Provided

Determinations	50
¹²⁵ I-anti-CA 19-9, monoclonal (mouse), red radioactive content (kBq / µCi)	17 mL < 129/3.5
6 calibrators A-F (1.0 mL) in human serum albumin	1 set
Assay buffer, blue	12 mL
Diluent (0 U/mL) in human serum albumin	6 mL
Test tubes, coated with anti-CA 19-9, monoclonal (mouse)	50
Control, human, lyophilic	1 mL

CA 19-9™ is a trademark of Fujirebio Diagnostics Inc.

Material required but not provided

- Micropipettes (100 μ L, 200 μ L, 300 μ L, 1000 μ L) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively, an appropriate automated analyzer system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- Purified water

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 10.1 μ Ci (372 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly.

(Do not allow foam to form.)

Open the control carefully and reconstitute with 1.0 mL of purified water. **(Do not allow foam to form.)** Ensure that lyophilized material adhering to the cap is also dissolved.

Storage of Reagents

Reconstituted control: 1 week at 2–8°C or 4 weeks at –20°C

Store all other reagents at 2–8°C. Use before the expiry date printed on the packaging.

- **Store upright.**

- **Keep away from direct light.**

Sample Collection, Material and Storage

- Collect samples using standard procedures.
- Sample material: serum. Disturbances due to plasma are not known.
- Storage at 2–8°C: 24 hours
- To store for longer periods freeze at –20°C.
- Freezing and thawing samples once does not affect the test results.
- Stored samples should be thoroughly mixed prior to use (vortex mixer).
- Do not use sera which are agglutinated, lipemic, hemolytic, icteric or contaminated.

Interfering Substances

No interference with test results has been observed with concentrations of bilirubin < 0.125 mg/mL, haemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The individual components of the kit are carefully matched to one another. If any components from different lots are exchanged or mixed up, the manufacturer does not guarantee reliable results.
- Strictly adhere to the sequence of the pipetting steps.
- Set the measuring time to at least 1 minute.
- This testing kit must not be used after the expiry date printed on the package label.
- Observe the quality control guidelines for medical laboratories.
- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that calibrators and samples are double determined.

If values above the highest calibrators concentration are expected, samples should be further diluted with the diluent (e.g. by factors 10, 100, 1000).

As an alternative to manual performance and evaluation, an appropriate automated analyzer system can be used on the responsibility of the laboratory.

1. Pipette 200 µL assay buffer in all tubes.
2. Pipette 100 µL of calibrator, control or patient sample into the bottom of a coated tube, mix (vortex) and cover with plastic foil.
3. Incubate the tubes for 3 hours (± 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 2 times with 2 mL purified water.

6. Add 300 μL ^{125}I -anti-CA 19-9 and cover with plastic foil.
7. Incubate the tubes for 2 hours (\pm 5 min) at room temperature (18–25°C) on a horizontal shaker.*
8. Aspirate the liquid.
9. Wash all tubes 2 times with 2 mL purified water.
10. Measure radioactivity (CPM) in all tubes (at least 1 min.).

200 μL	Pipette assay buffer
100 μL	Pipette calibrator, control or patient sample into the bottom of a test tube
	Mix
3 hrs (\pm 5 min)	Incubate at room temperature (18–25°C) on a shaker*
	Aspirate
2 x 2 mL	Wash with purified water
300 μL	Add ^{125}I -anti-CA 19-9
	Mix
2 hrs (\pm 5 min)	Incubate at room temperature (18–25 °C) on a shaker*
	Aspirate
2 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

*** Keep agitation conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm
 Amplitude 10 mm = 220 rpm
 Amplitude < 8 mm = 300 rpm

Calculation of Results

The calibrator curve can be established manually as follows:

1. Determine the mean CPM for each pair of tubes (double determination).
2. Divide the mean CPM of each calibrator (B) by the mean CPM of the highest calibrator (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding ($\%B/B_{\text{max}}$) for each calibrator.
3. On semi-log paper, plot the relative binding of each calibrator ($\%B/B_{\text{max}}$) on the Y-axis versus the corresponding concentrations on the X-axis.
4. Sample concentrations (U/mL) can be read directly off the calibrator curve from their corresponding relative binding ($\%B/B_{\text{max}}$).

If the radioactivity measured is above that of the highest calibrator, the samples must be diluted with the diluent and tested again. With diluted samples the actual serum concentration and the appropriate dilution factor have to be established.

The instrumental calculation of radio-immunological measured values is performed using a spline approximation.

Quality Control

Observe the quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory.

This control should be simultaneously tested with each test run and treated like patient samples. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Interpretation of the Results

In sera from 29 evidently healthy blood donors CA 19-9 values of up to approx. 24 U/mL were found. In patients (n=95) with various benign diseases (pancreatitis, infections of the gastrointestinal tract) CA 19-9 levels up to 37 U/mL were normally observed. Since CA 19-9 values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Tumour patients may exhibit CA 19-9 assay values below the cut-off.

CA 19-9 levels above the cut-off are observed in patients with cystic fibrosis and severe benign liver disorders (3, 4).

Therefore, the CA 19-9 assay value must only be interpreted in conjunction with the clinical picture and other diagnostic procedures. In patients who are Lewis^a negative the CA 19-9 antigen cannot be detected either.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may in principle lead to falsely elevated or decreased values. HAMA-neutralising agents are added to the test. However the influencing of results cannot be completely ruled out. These samples should not be used for the CA 19-9 IRMA.

Analytical Data

Calibration

The CA 19-9 IRMA has been calibrated using the Fujirebio CA 19-9 assay.

Measuring range

The measuring range is 3 - 240 U/mL.

High-dose hook

No high-dose hook effect was observed for concentrations up to 30,000 U/mL.

Precision

Intra-assay variation			Inter-assay variation		
Mean value (U/mL)	CV (%)	n=	Mean value (U/mL)	CV (%)	n=
15.7	6.3	9	12.3	9.4	14
31.4	4.9	9	41.3	2.7	17
61.2	3.1	9	57.1	5.8	17
122.1	2.1	9	116.4	3.5	17

Analytical sensitivity

The lower detection limit is < 3.0 U CA 19-9/mL. This detection limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest CA 19-9 concentration that can be differentiated from zero with statistical significance.

Specificity

No cross-reactivities with mitomycin-C, doxorubicin, fluorouracil or methotrexate in therapeutic ranges were found.

Linearity upon Dilution

A patient serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 103.8 U/mL

Dilution	Measured value (U/mL)	Expected value (U/mL)	Recovery (%)
1 : 2	51.4	51.8	99.2
1 : 4	27.8	26.7	104.1
1 : 8	13.6	14.1	96.2
1 : 16	7.7	7.9	97.6

Recovery

A patient serum with low CA 19-9 content was spiked with different concentrations of CA 19-9.

Original concentration: 7.4 U/mL

Measured value (U/mL)	Expected value (U/mL)	Recovery (%)
78.6	78.8	99.7
43.0	42.8	100.6
25.5	24.8	102.9
15.0	15.7	95.1

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


1. Bartoloni C, Guidi L, Baroni R, Tricceri A, Barone C, Garufi C, Gambassi G. Longitudinal Study of CA 19-9 and CEA in Gastrointestinal Cancer Patients at Different Stages. *J Tumor Marker Oncol* 1990; **5** (4): 351-362
2. Diez M, Cerdán FJ, Pollán M, Maestro ML, Ortega MD, Martínez S, Moreno G, Balibrea JL. Prognostic Significance of Preoperative Serum CA 19.9 Assay in Patients with Colorectal Carcinoma. *Anticancer Res* 1994; **14**: 2819-2826
3. Encabo G, Ruibal A. Seric CA 19.9 levels in patients with non tumoral pathologies. Our experience in 892 cases. *Bull Cancer (Paris)* 1986; **73** (3): 256-259
4. Fabris C, Falleti E, Pirisi M, Soardo G, Toniutto P, Vitulli D, Bortolotti N, Gonano F, Bartoli E. Non-specific increase of serum carbohydrate antigen 19-9 in patients with liver disease associated with increased circulating levels of adhesion molecules. *Clinica Chimica Acta* 1995; **243**: 25-33
5. Gattani AM, Mandeli J, Bruckner HW. Tumor Markers in Patients with Pancreatic Carcinoma. *Cancer* 1996; **78**: 57-62
6. Koprowski H, Steplewski Z, Mitchell K, Herlyn M, Herlyn D, Fuhrer P. Colorectal Carcinoma Antigens Detected by Hybridoma Antibodies. *Somatic Cell Genetics*, 1979; **5** (6): 957-972
7. Malesci A, Montorsi M, Mariani A, Santambrogio R, Bonato C, Bissi O, Tacconi M, Wizemann G, Spina G. Clinical Utility of the Serum CA 19-9 Test for Diagnosing Pancreatic Carcinoma in Symptomatic Patients: A Prospective Study. *Pancreas* 1992; **7** (4): 497-502
8. Satake K, Takeuchi T. Comparison of CA 19-9 with Other Tumor Markers in the Diagnosis of Cancer of the Pancreas. *Pancreas* 1994; **9** (6): 720-724
9. Willett CG, Daly WJ, Warshaw AL. CA 19-9 Is an Index of Response to Neoadjuvant Chemoradiation Therapy in Pancreatic Cancer. *Am J Surg* 1996; **172**: 350-352

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Tracer: antibody labelled with ¹²⁵ I	Traceur : anticorps marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markierter Antikörper	Trazador: anticuerpo etiquetado con ¹²⁵ I	Tracciatore: anticorpo etichettato con ¹²⁵ I
	Sample diluent	Diluant d'échantillon	Proben-Diluent	Diluyente de muestras	Diluyente campione
	Reconstitute with X mL	Reconstituer avec X mL	Mit X mL auflösen	Reconstitución con X mL	Ricostituire con X mL
	Buffer	Tampon	Puffer	Tampón	Tampone
	Solid phase Coated tubes	Phase solide Tubes enduits	Festphase Beschichtete Teströhrchen	Fase sólida Tubos recubiertos	Fase solida Cannule con rivestimento.
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.