CA 72-4[®]
¹²⁵I IRMA Kit

For the quantitative determination of
TAG 72 in serum

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For professional use only!

Intended Use

In vitro assay for the quantitative determination of TAG 72[®] (tumor-associated glycoprotein 72) in serum during the follow-up of patients with gastric carcinomas or with mucinous-type ovarian carcinoma.

Summary and Explanation of the Test

In 1981, Colcher and coworkers prepared the monoclonal antibody B 72.3 after immunization of mice with a liver metastasis fraction of a mammary carcinoma. The high-molecular glycoprotein (MW 220,000 - 1,000,000) detectable with B 72.3 and belonging to the group of mucins, was termed TAG 72. The second monoclonal antibody cc49 with a 5- to 10-fold higher affinity was raised against highly purified TAG 72. Both antibodies recognize different oligosaccharide epitopes of TAG 72.

Histological detection of the tumor-associated antigen TAG 72 is possible on neoplastic cells of the breast, lung, colon, stomach and the ovaries. Other malignant and benign tumors as well as healthy, mature tissues do not exhibit any reaction to TAG 72 antibodies. The only exception is the secretory endometrium.

Principles of the Procedure

Two-site immunoradiometric assay (sandwich principle) using two highly specific monoclonal antibodies. The monoclonal antibody cc49 is used for the coating of the solid phase (coated tube) and the monoclonal antibody B 72.3 is used for the tracer.

During the first incubation, the CA 72-4 present in patient samples and calibrators is bound to the antibody immobilized at the test tube wall. Unbound material is removed by a washing step.

During the second incubation, the tracer antibody reacts with the CA 72-4 already bound. After removal of excess tracer by a second washing step, the radioactivity remaining at the tube wall is measured in a gamma scintillation counter.

CONTENTS Material Provided

Determinations	50
¹²⁵ I-anti-CA 72-4, monoclonal (mouse), red	11 mL
radioactive content (kBq / μ Ci)	< 222 /6
5 calibrators A-E (0.8 mL) in human serum albumin	1 set
Assay buffer, blue	6 mL
Diluent (0 U/mL) in human serum albumin	6 mL
Test tubes, coated with anti-CA 72-4, monoclonal (mouse)	50
Control, human, lyophil.	1 mL

CA 72-4[®] is a registered trademark of Fujirebio Diagnostics Inc.

Material Required but not Provided

- Micropipets (100 μ L, 200 μ L, 1000 μ L) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively an appropriate automated analyzer system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- Purified water

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 6.04 μ Ci (224 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly (**Avoid foam formation**).

Control must be opened carefully (vacuum) and reconstituted with 1 mL of purified water. Avoid heavy shaking when dissolving (foaming). Make sure that lyophilized material adhering to the screw cap is also dissolved.

Storage of Reagents

Reconstituted control: 1 week at 2–8°C or 4 weeks at –20°C.

Store all reagents at 2–8°C until the expiry date printed on the package label.

- **Keep upright for storage.**
- **Keep away from direct light.**

Sample Collection, Material and Storage

- Collect samples using standard procedures.
- Sample material: serum; disturbances with plasma are not known.
- Storage at 2–8°C: 24 hr.
- For longer storage periods, freeze to below –20°C.
- Frozen samples should only be thawed once.
- Stored samples should be thoroughly mixed prior to use (Vortex mixer).
- Do not use samples which are agglutinated, lipemic, hemolysed, icteric or contaminated.

Interfering Substances

No interference with test results is seen by concentrations of bilirubin < 0.125 mg/mL, hemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The single components of each kit are carefully matched. In case of exchange or mixture of any components from different lots the manufacturer does not guarantee reliable results.
- Strictly adhere to the pipetting step sequence.
- The measuring time at the gamma counter must be adjusted for counting at least 1 minute.
- This kit must not be used after the expiry date printed on the package label.
- Observe quality control guidelines for medical laboratories.
- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that calibrators and samples are assayed in duplicate.

If values are expected above the highest calibrator concentration, samples should be further diluted with the diluent (e.g. by factors 10, 100, 1000).

Alternatively to the manual performance and evaluation an appropriate automated analyzer system can be used on responsibility of the laboratory.

1. Pipet 100 µL assay buffer in all tubes.
2. Pipet 100 µL each of calibrator, control or patient sample, mix (vortex) and cover with plastic foil.
3. Incubate for 4 hours (± 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 3 x with 2 mL purified water.
6. Add 200 µL ¹²⁵I-anti-CA 72-4, vortex.
7. Incubate for 18-20 hours at 6°C (± 2°).
8. Aspirate the liquid.
9. Wash all tubes 3 x with 2 mL of purified water.
10. Measure radioactivity (CPM) in all tubes (at least 1 min).

* **Keep shaking conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

100 µL	Pipet assay buffer
100 µL	Add calibrator, control or patient sample
	Mix
4 h (± 5 min)	Incubate at room temperature (18–25°C) on a shaker* Aspirate
3 x 2 mL	Wash with purified water
200 µL	Add ¹²⁵ I-anti-CA 72-4
	Mix
18-20 h	Incubate standing at 6°C (± 2°C) Aspirate
3 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

* **Keep shaking conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

Calculation of Results

The calibrator curve is established manually as follows:

1. Determine the mean CPM for each pair of duplicate tubes.
2. Divide the mean CPM of each calibrators (B) by the mean CPM of the highest calibrator (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding ($\%B/B_{max}$) for each calibrator.
3. On semi-log paper, plot the relative binding of each calibrator ($\%B/B_{max}$) on the Y-axis versus the corresponding concentrations on the X-axis.
4. Read sample concentrations (U/mL) directly off the calibrator curve by their corresponding relative binding ($\%B/B_{max}$).

Samples with counts above the highest calibrator concentration must be diluted with the kit diluent and re-assayed. For the final concentration of such samples, the appropriate dilution factor has to be considered.

The instrumental calculation of radio-immunological measuring values is based on a spline approximation.

Quality Control

Observe quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory.

This control should be simultaneously tested with each test run and treated like patient samples. The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Expected Values

Sera of apparently healthy blood donors (n=100) show CA 72-4 values up to approx. 4 U/mL (95th percentile).

Since CA 72-4 values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Patients with malignancies may exhibit CA 72-4 values within the normal range. Elevated CA 72-4 values may be observed in patients with benign diseases such as chronic liver impairment or benign ovarian tumors. Therefore, CA 72-4 serum levels provide no evidence of the presence or absence of malignancies. They may only be interpreted in context with the clinical picture and other diagnostic procedures. The CA 72-4 assay as the only criterion should not be used for cancer screening. Elevated concentrations of CA 72-4 may also occur in patients with non small cell lung cancer.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may give falsely elevated or decreased values. Although HAMA-neutralising agents are added, extremely high HAMA serum concentrations may occasionally influence results. These samples should not be used for the CA 72-4 IRMA assay.

Analytical Data

Calibration

The test has been calibrated using the Fujirebio IRMA CA 72-4® assay.

Measuring range

CA 72-4 IRMA measures concentrations between 3 and 100 U/mL.

High-dose hook

No high-dose hook effect was observed for CA 72-4 concentrations up to 1250 U/mL.

Precision

Intra-assay variation			Inter-assay variation		
Mean value (U/mL)	CV (%)	n=	Mean value (U/mL)	CV (%)	n=
9.1	3.0	3	9.4	14.0	7
24.9	2.5	3	25.3	8.2	7
71.2	1.4	3	68.1	5.4	7

Specificity

The CA 72-4 IRMA assay uses the monoclonal antibodies B 72.3 and cc49. These are only available at Fujirebio, their licensees and representatives. Information about the specificity of the antibodies is described in the respective literature (1-5).

Linearity upon Dilution

A patient's serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 42.7 U/mL

Dilution	Measured value (U/mL)	Expected value (U/mL)	Recovery (%)
1 : 1.25	34.6	34.7	99.6
1 : 1.7	27.2	25.6	102.6
1 : 2.5	18.3	18.4	99.4

Recovery

A patient's serum with low CA 72-4 content was spiked with different amounts of CA 72-4 and then measured.

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


1. Johnson BG, Schlom J, Paterson AJ, Bennett J, Magnani JL, Colcher D. Analysis of a human tumor-associated glycoprotein (TAG-72) identified by monoclonal antibody B72.3. *Cancer Res.* 1986; **46(2)**: 850-857.
2. Lan MS, Bast RC, Jr, Colnaghi MI, Knapp RC, Colcher D, Schlom J, Megzgar RS. Co-expression of human cancer associated epitopes on mucin molecules. *J Cancer* 1987; **39(1)**: 68-72.
3. Colcher D, Paterson A, Sears H, Schlom J. Use of monoclonal antibody B72.3 to detect a circulating tumor-associated glycoprotein (TAG-72) in the serum of colon carcinoma patients (meeting abstract). *Hybridoma* 1985; **4(1)**: 71.
4. Paterson AJ, Schlom J, Sears HF, Bennet J, Colcher D. A radioimmunoassay for the detection of a human tumor-associated glycoprotein (TAG-72) using monoclonal antibody B72.3. *Int. J.Cancer* 1986; **37(5)**: 659-666.
5. Klug TL, Sattler MA, Colcher D, Schlom J. Monoclonal antibody immunoradiometric assays for an antigenic determinant (CA 72-4) on a novel pancarcinoma antigen (TAG-72). *Int. J. Cancer* 1986; **38(5)**: 661-669.

SYMBOLS USED WITH DEVICES

	English	Français	Deutsch	Español	Italiano
	European Conformity	Conformité aux normes européennes	CE – Konformitätskennzeichnung	Conformidad europea	Conformità europea
	Expiration Date	Date limite d'utilisation	Mindesthaltbarkeitsdatum	Fecha de caducidad	Data di scadenza
	Manufacturer	Fabricant	Hersteller	Fabricante	Fabbricante
	Consult Instructions for Use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	In vitro diagnostic.	Diagnostic in vitro.	In-vitro-Diagnostikum.	Diagnóstico in vitro.	Diagnostica in vitro.
	Lot No.	No. de lot	Chargen-Nr.	Número de lote	Lotto n°
	Temperature limitation.	Limitation de température.	Temperaturbereich	Limitación de temperatura	Limite della temperatura
	Calibrator	Étalon	Kalibrator	Calibrador	Calibratore
	Control serum	Sérum de contrôle	Kontrollserum	Suero de control	Siero di controllo
	Tracer: antibody labelled with ¹²⁵ I	Traceur : anticorps marqué à l'iode ¹²⁵	Tracer: ¹²⁵ I-markierter Antikörper	Trazador: anticuerpo etiquetado con ¹²⁵ I	Tracciatore: anticorpo etichettato con ¹²⁵ I
	Sample diluent	Diluant d'échantillon	Proben-Diluent	Diluyente de muestras	Diluyente campione
	Reconstitute with X mL	Reconstituer avec X mL	Mit X mL auflösen	Reconstitución con X mL	Ricostituire con X mL
	Solid phase Coated tubes	Phase solide Tubes enduits	Festphase Beschichtete Teströhrchen	Fase sólida Tubos recubiertos	Fase solida Cannule con rivestimento.
	Radioactive	Radioactif	Radioaktiv	Radioactivo	Radioattivo

SYMBOLS USED WITH DEVICES CONTINUED ON NEXT PAGE.

CONTINUED FROM PREVIOUS PAGE - SYMBOLS USED WITH DEVICES

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DiaSorin Inc.
DiaSorin S.p.A.,
PRINTED IN U.S.A.