DiaSorin Inc.
 1951 Northwestern Avenue, P.O. Box 285,
 Stillwater, MN 55082-0285 U.S.A.
 Phone: 651-439-9710 - Fax: 651-351-5669

IRMA-mat[®] AFP

Instructions for Use

REF 324.471

100 tests

English

IVD

For professional use only!

Intended Use

In vitro test for the quantitative determination of alpha fetoprotein (AFP) in human serum and amniotic fluid for pregnancy monitoring and during the follow-up of tumour patients. The assay performance characteristics have not been established for diagnosis of trisomy 21.

Summary and Explanation of the Test

Alpha fetoprotein is a macromolecular glycoprotein (with a molecular weight of approx. 68,000) consisting of a single polypeptide chain. AFP, which belongs to the group of oncofetal proteins, is produced by the yolk sac and in the foetal liver.

In oncology, AFP is determined in patients with liver-cell carcinoma (1, 5) or germ-cell tumours (non-seminomatous tumours of the testes; endodermal sinus tumour of the ovaries) (3, 6, 7) and in pregnancy monitoring (2, 4). During pregnancy, AFP levels in maternal blood continuously increase. After reaching a maximum between weeks 28 to 32 of gestation, the level again falls until delivery. In the amniotic fluid, the maximum is already achieved between the 13th and 15th week of gestation. Elevated AFP levels in early pregnancy indicate neural tube defects (spina bifida, anencephaly).

The determination of serum AFP during monitoring of therapy and the progression of patients with carcinoma provides valuable information about the success of treatment and whether recidivation occurs.

Principles of the Procedure

Two-site immunoradiometric assay (sandwich principle) using two highly specific monoclonal antibodies for coating of the solid phase (coated tubes) and the tracer. The tracer antibody and the coated antibody react simultaneously with the AFP present in patient samples or standards. Excess tracer is removed by a washing step and the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

CONTENTS

Determinations	100
¹²⁵ I-anti-AFP, monoclonal (mouse), red	33 mL
radioactive content (kBq / μ Ci)	< 424 / 11.5
7 standards A-G (0.6 mL) in human serum albumin (The exact concentration is indicated on each vial label)	1 set
Diluent (0 IU/mL) in phosphate buffer (concentrate)	2 x 11 mL
Test tubes, coated with anti-AFP, monoclonal (mouse)	2 x 50
Control serum, human, lyophilic	0.6 mL
Quality control report	1

Material required but not provided

- Micropipettes (25 µL, 300 µL) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively, an appropriate automated analyser system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- Purified water

Warnings and Precautions for Users

Reagents Containing Human Source Material

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th., May 1999 or current edition.

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 11.5 µCi (424 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts," in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly. (Avoid foam formation).

Open the control carefully and reconstitute with 0.6 mL of purified water. (Avoid foam formation). Make sure that lyophilised material adhering to the cap is also dissolved.

Dilute 20 mL of diluent concentrate with 100 mL of purified water.

Storage of Reagents

Reconstituted controls: 1 week at 2–8°C or 4 weeks at – 20°C.

Store all reagents at 2–8°C until the expiry date printed on the package label.

- **Keep upright for storage.**

- **Keep away from direct light.**

Sample Collection and Storage

- Collect samples using standard procedures.

- Sample material: serum, amniotic fluid; disturbances with plasma are not known.

- Storage at 2 - 8°C: 24 h.

- For longer storage periods: freeze to below -20°C.

- Freezing and thawing samples once does not affect the test results.

- Stored samples should be thoroughly mixed prior to use (vortex mixer).

- Do not use samples which are agglutinated, lipemic, hemolysed, icteric or contaminated.

Interfering Substances

No interference with test results has been observed with concentrations of bilirubin < 0.125 mg/mL, haemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The single components of each kit are carefully matched. In case of exchange or mixture of any components from different lots the manufacturer does not guarantee reliable results. See bottom of the kit for the lot numbers of all components.

- Strictly adhere to the procedures.

- The measuring time at the gamma counter must be adjusted for counting at least 1 minute.

- This kit must not be used after the expiry date printed on the package label.

- Observe quality control guidelines for medical laboratories.

- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that standards and samples are assayed in duplicate.

If values above the highest standard are expected, samples should be diluted (e.g. by dilution factor 10, 100, 1000).

Amniotic fluid must always be diluted. During the diagnostically relevant weeks of pregnancy (between 16th and 19th week) 1:100 dilution is sufficient. In pathological cases a higher dilution may be necessary (parallel testing of 1:100 and 1:1000 dilutions).

As an alternative to manual performance and evaluation, an appropriate automated analyser system can be used at the responsibility of the laboratory.

1. Pipette 25 µL of standard, control or patient sample onto the bottom of a coated tube.
2. Add 300 µL ¹²⁵I-anti-AFP, mix (vortex mixer).
3. Incubate the tubes for 3 hours (± 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 3 times with 2 mL of purified water.
6. Measure radioactivity (CPM) in all tubes (at least 1 min.)

25 µL	Pipette standard, control or patient sample onto the bottom of a coated tube.
300 µL	Add ¹²⁵ I-anti-AFP Mix
3 hrs (± 5 min)	Incubate at room temperature (18–25°C) on a shaker* Aspirate
3 x 2 mL	Wash with purified water
1 min	Measure (gamma scintillation counter)

*** Keep shaking conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

Calculation of Results

The standard curve can be established manually as follows:

1. Determine the mean CPM for each pair of tubes (double determination).
2. Divide the mean CPM of each standard (B) by the mean CPM of the highest standard (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding (%B/B_{max}) for each standard.
3. On semi-log paper, plot the relative binding of each standard (%B/B_{max}) on the Y-axis versus the corresponding concentrations (IU/mL) on the X-axis.
4. Sample concentrations (IU/mL) can be read directly off the standard curve from their corresponding relative binding (%B/B_{max}).

If the radioactivity measured is above that of the highest standard, the samples must be diluted with the diluent and tested again. With diluted samples the actual serum concentration and the appropriate dilution factor have to be established. An example of a standard curve is given in the QC report. This curve must not be used for the calculation of unknown samples.

The instrumental calculation of radio-immunological measured values is performed using a spline approximation.

Quality Control

Observe the quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory. This control should be simultaneously tested with each test run and treated like patient samples.

The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Interpretation of the Results

In healthy men and non-pregnant women (n=63) AFP values do not rise above 6 IU/mL (95th percentile). Elevated levels point to neoplastic activity.

During pregnancy, the AFP concentrations in the maternal serum rise until they reach a maximum between the 28th and 32nd week of pregnancy. After week 32 a fall in values is observed.

AFP levels outside the reference range are indicative of fetal damage. Elevated AFP concentrations indicate a neural tube defect. The serum findings are verified by AFP determination in the amniotic fluid.

Based on healthy, pregnant women, the following values for serum and amniotic fluid** have been evaluated for the IRMA - mat[®] AFP assay:

Normal AFP values in early pregnancy

Week of gestation	AFP in maternal serum		AFP in amniotic fluid**	
	(Median IU/mL)	(n*)	(Median IU/mL)	(n*)
15	26.7	236	17,083	15
16	28.0	236	14,679	107
17	34.1	236	12,532	135
18	38.4	236	10,075	106
19	41.9	236	8,381	37
20	60.0	236	6,877	24

** Values for amniotic fluid have been evaluated with the LIAISON[®] AFP assay. As both assays use the same antibodies, the values determined using the LIAISON[®] AFP can be considered standard values.

Since AFP values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Patients with carcinomas may also exhibit AFP values within the normal range.

AFP concentrations may in comparison be elevated in the case of benign diseases such as cirrhosis of the liver, hepatitis or tyrosinaemia (8). For this reason AFP determination is most suitable for therapeutic monitoring or during follow-up and for confirmation of histological results. The AFP assay should not be used as the only criterion for tumour screening.

Since no data on the diagnosis of Down's syndrome has been ascertained with IRMA-mat[®] AFP, the assay is not recommended for use in triple screening.

All tests, in which antigen is incubated together with labelled antibodies and immobilised antibodies in a liquid phase, bear the risk that undiluted samples containing extremely high concentrations of the antigen, will give measuring values below those of the highest standard. In case of the IRMA-mat[®] AFP, this phenomenon is observed at concentrations exceeding 30,000 IU/mL. If such values are suspected, measurement should be repeated after further dilution (e.g. by factors 10, 100, 1000) of the specimen.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may give falsely elevated or decreased values. Although HAMA-neutralising agents are added, extremely high HAMA serum concentrations may occasionally influence results. These samples should not be used for the IRMA-mat[®] AFP assay.

Analytical Data

Calibration

IRMA-mat® AFP has been calibrated using the Reference Standard MRC 72/225.

1ng AFP = 0.83 IU AFP (MRC 72/225)

1IU AFP (MRC 72/225) = 1.21 ng AFP

Measuring range

Measuring range is 2 - 600 IU/mL.

High-dose hook

No high-dose hook effect was observed for AFP concentrations up to 30,000 IU/mL.

Precision

Intra-assay Mean value (IU/mL)	CV (%)	n=	Inter-assay Mean value (IU/mL)	CV (%)	n=
19	2.6	11	11	6.9	11
99	2.8	11	106	4.5	11
338	1.8	11	170	4.2	11

Analytical sensitivity

The lower detection limit is < 2.0 IU AFP/mL. This detection limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest AFP concentration that can be differentiated from zero with statistical significance.

Specificity

No cross-reactivities with mitomycin-C, doxorubicin or fluorouracil in therapeutic ranges were found.

Dilution

A patient's serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 457.5 IU/mL

Dilution	Measured value (IU/mL)	Expected value (IU/mL)	Recovery (%)
1 : 1.25	368	366	101
1 : 2.5	202	183	110
1 : 5	100	92	109
1 : 10	50	46	109

Recovery

A patient's serum with low AFP content was spiked with different amounts of AFP and then measured.

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- 1 . Aoyagi Y. Carbohydrate-based measurements of alpha-fetoprotein in the early diagnosis of hepatocellular carcinoma. *Glycoconjugate Journal* 1995; 12: 194-199
- 2 . Bock JL. Current Issues in Maternal Serum: Alpha-Fetoprotein Screening (Review Article). *Am J Clin Pathol* 1992; 97 : 541-554
- 3 . Brewer JA, Tank ES. Yolk sac tumors and alpha-fetoprotein in first year of life. *Urology* 1993; 42 (1): 79-80
- 4 . Brizot ML et al. First trimester maternal serum alpha-fetoprotein in fetal trisomies. *Br J Obstet Gynaec* 1995; 102: 31-34
- 5 . Burditt LJ et al. Detection of Hepatocellular Carcinoma-Specific Alpha-Fetoprotein by Isoelectric Focusing. *Cancer* 1994; 74 (1): 25-29
- 6 . Johnson PJ et al. Germ Cell Tumors Express a Specific Alpha-Fetoprotein Variant Detectable by Isoelectric Focusing. *Cancer* 1995; 75 (7): 1663-1668
- 7 . Mora J et al. Alpha-fetoprotein-Concanavalin A Binding as a Marker to Discriminate Between Germ Cell Tumours and Liver Diseases. *Eur J Cancer* 1995; 31A (13/14): 2239-2242
- 8 . Stránsky J et al. Assessment of Alpha-Fetoprotein in chronic HBsAg Carriers and in HBsAg Negative Cirrhosis of the Liver. *Sborník lékařsky* 1993;

SYMBOLS USED WITH IVD DEVICES
SIMBOLI USATI CON I DIAGNOSTICI IN VITRO
SYMBOLES UTILISÉS AVEC LES DIAGNOSTICS IN VITRO
MIT IN-VITRO-DIAGNOSTIKA GEBRAUCHTE SYMBOLE
SÍMBOLOS USADOS CON LOS DIAGNÓSTICOS IN VITRO
SÍMBOLOS UTILIZADOS COM OS DIAGNÓSTICOS IN VITRO
SYMBOLER SOM ANVÄNDS MED IN VITRO-DIAGNOSTIK
SYMBOLER SOM BRUGES VED IN VITRO-DIAGNOSTIK
JELMAGYARÁZAT IVD ESZKÖZÖKHÖZ
ΣΥΜΒΟΛΑ ΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΑ ΜΕ ΤΙΣ ΣΥΣΚΕΥΕΣ IVD

CONT.

Kit contents / Contenido del kit / Contenu de la trousse / Inhalt des Kits
 Inhalt des Kits / Contenido del kit / Conteúdo do dispositivo / Satsinnehåll
 Kittets indhold / Készlet tartalma / Περιεχόμενα συσκευασίας.

Ab ¹²⁵I

Tracer: antibody labelled with ¹²⁵I / Tracciante: anticorpi marcati con ¹²⁵I
 Traceur: anticorps marqués à l'¹²⁵I / Tracer: mit ¹²⁵J markierte Antikörper
 Trazador: anticuerpos marcados con ¹²⁵I
 Marcador: anticorpos marcados com ¹²⁵I
 Spårämne: antikroppar märkta med ¹²⁵I
 Tracer: antistoffer mærket med ¹²⁵I
 Tracer: ¹²⁵I-tel jelölt antitest
 Ιχνηθέτης: αντίσωμα σημασμένο με ¹²⁵I.

SORB

Solid phase (Coated tubes / Coated beads).
 Fase solida (Provette sensibilizzate / Sferette sensibilizzate).
 Phase solide (Tubes revêtues / Billes revêtues).
 Feste Phase (Beschichtete Röhrchen / Beschichtete Kugeln).
 Fase sólida (Tubos recubiertos / Bolas recubiertas).
 Fase sólida (Tubos revestidos / Bolas revestidas).
 Fast stadium (Belagda rör / Belagda kulor).
 Fast stadium (Sensibiliserede rør / sensibiliserede kugler).
 Szilárd fázis (bevonatos cső/ bevonatos gyöngy)
 Στερεά φάση (επικαλυμμένοι δοκιμαστικοί σωλήνες / επικαλυμμένα σφαιρίδια).

DIL


Sample diluent / Diluente campioni / Diluant pour échantillons
 Probenverdünnungslösung / Diluyente de muestras / Diluente das amostras
 Provspädning / Fortyndingsmiddel til prøver / Minta hígító / Διαλύτης δειγμάτων.

CAL

Calibrator / Calibratore / Etalon / Kalibrator / Calibrador / Kalibrator
 Kalibrator / Kalibrátor / Μέσο βαθμονόμησης.

CONTROL


Control serum / Siero di controllo / Sérum de contrôle / Kontrollserum
 Suero de control / Soro de controlo / Kontrollserum / Kontrolserum / Kontroll szérum /
 Ορός ελέγχου.

For XX tests / Per XX dosaggi / Pour XX dosages / Für XX Bestimmungen / Para XX ensayos / Para XX testes / För XX dosering Til XX test / XX teszthez / Περιεχόμενο επαρκές για XX εξετάσεις.

Radioactive / Radioattivo / Radioactif / Radioaktiv / Radioactivo / Radioaktiv Radioaktiv / Radioaktív / Ραδιενεργός.

Reconstitute with X mL / Ricostituire con X mL / Reconstituer avec X mL Mit X mL auflösen / Reconstituya con X mL / Reconstitua com X mL Återställ med X mL / Rekonstruer med X mL / X mL-ben feloldandó / Ανασύσταση με X mL.