DiaSorin Inc.
 1951 Northwestern Avenue, P.O. Box 285,
 Stillwater, MN 55082-0285 U.S.A.
 Phone: 651-439-9710 - Fax: 651-351-5669

IRMA-coat[®] CEA

Instructions for Use

REF 324.321

100 tests

English

IVD

For professional use only!

Intended Use

In vitro assay for the quantitative determination of carcinoembryonic antigen (CEA) in serum during the follow-up of cancer patients.

Summary and Explanation of the Test

CEA was first isolated from colon tumours (4). It belongs to a family of glycoproteins with a molecular weight of approx. 180,000 – 200,000 and a carbohydrate content varying between 50 and 60% (7).

Elevated CEA serum concentrations are found in patients with carcinoma of the lung, colon or breast (1, 2, 3). CEA measurement, alone or in combination with other markers, is intended for the early diagnosis of relapses and for therapy monitoring (5, 8).

Principles of the Procedure

Two-site immunoradiometric assay (sandwich principle) using two highly specific monoclonal antibodies for coating of the solid phase (coated tubes) and the tracer. The tracer antibody and the coated antibody react simultaneously with the CEA present in patient samples or standards. Excess tracer is removed by a washing step and the radioactivity bound to the tube wall is measured in a gamma scintillation counter.

CONTENTS

Determinations	100
¹²⁵ I-anti-CEA, monoclonal (mouse), red	11 mL
Radioactive content (kBq / μ Ci)	< 705 /19
6 standards A-F (1.0 mL) in buffer (The exact concentration is indicated on each vial label)	1 set
Diluent (0 ng/mL) in buffer	11 mL
Test tubes, coated with anti-CEA, monoclonal (mouse)	2 x 50
Control serum, A and B, human, lyophilic.	2 x 1mL
Quality control report	1

Material Required but not Provided

- Micropipets (100 μ L) with disposable plastic tips
- Vortex mixer
- Manual or automatic washer with aspiration device
- Horizontal shaker
- Gamma scintillation counter
- Alternatively an appropriate automated analyser system, if available
- Uncoated polystyrene tubes for the dilution of sera and controls
- 0.9 % NaCl solution for washing steps
- Purified water

Warnings and Precautions for Users

Reagents Containing Human Source Material

Treat as potentially infectious.

Each serum/plasma donor unit used in the preparation of this product has been tested by a U.S. FDA approved method and found non-reactive for the presence of HBsAg, antibody to HCV, and antibody to HIV 1/2. While these methods are highly accurate, they do not guarantee that all infected units will be detected. This product may also contain other human source material for which there is no approved test. Because no known test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C virus (HCV), Human Immunodeficiency Virus (HIV), or other infectious agents are absent, all products containing human source material should be handled in accordance with good laboratory practices using appropriate precautions as described in the U.S. Centers for Disease Control and Prevention/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 4th ed., May 1999 or current edition.

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 19 μ Ci (705 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

WARNING: This product contains a chemical known to the State of California to cause cancer.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

Reagents Containing Sodium Azide

CAUTION: Some reagents in this kit contain sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to “Decontamination of Laboratory Sink Drains to Remove Azide Salts,” in the Manual Guide-Safety Management No. CDC-22 issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

Preparation of Reagents

Allow test components to reach room temperature (18–25°C) prior to testing and mix thoroughly (avoid foam formation).

Controls must be opened carefully and reconstituted with 1 mL of purified water. Avoid heavy shaking when dissolving (foaming). Make sure that lyophilised material adherent to the screw cap is also dissolved.

Storage of Reagents

Reconstituted controls: 1 week at 2-8°C or 4 weeks at – 20°C.

Store all reagents at 2-8°C until the expiry date printed on the package label.

- **Keep upright for storage.**
- **Keep away from direct light.**

Sample Collection and Storage

- Collect samples using standard procedures.
- Sample material: serum; disturbances with plasma are not known.
- Storage at 2-8°C: 24 h.
- For longer storage periods: freeze to below -20°C.
- Freezing and thawing samples once does not affect the test results.
- Stored samples should be thoroughly mixed prior to use (vortex mixer).
- Do not use samples which are agglutinated, lipemic, hemolysed, icteric or contaminated.

Interfering Substances

No interference with test results is seen by concentrations of bilirubin < 0.125 mg/mL, haemoglobin < 500 mg/dL or triglycerides < 12.5 mg/mL.

Procedural Notes

- The single components of each kit are carefully matched. In case of exchange or mixture of any components from different lots the manufacturer does not guarantee reliable results. See bottom of the kit for the lot numbers of all components.
- Strictly adhere to the procedures.
- The measuring time at the gamma counter must be adjusted for counting at least 1 minute.
- This kit must not be used after the expiry date printed on the package label.
- Observe quality control guidelines for medical laboratories.
- Avoid microbial contamination of the reagents.

Test Procedure

It is recommended that standards and samples are assayed in duplicate.

If values are expected above the highest standard concentration, samples should be further diluted with the diluent (e.g. by factors 10, 100, 1000).

Alternatively to the manual performance and evaluation an appropriate automated analyser system can be used on responsibility of the laboratory.

1. Pipet 100 μL each of standard, control or patient sample onto the bottom of the corresponding coated tube.
2. Add 100 μL ^{125}I -anti-CEA, mix gently (do not vortex).
3. Incubate for 4 hours (\pm 5 min) at room temperature (18–25°C) on a horizontal shaker.*
4. Aspirate the liquid.
5. Wash all tubes 3 x with 2 mL of 0.9% NaCl solution.
6. Measure radioactivity (CPM) in all tubes (at least 1 min).

100 μL	Pipet standard, control or patient sample
100 μL	Add ^{125}I -anti-CEA
	Mix
4 h (\pm 5 min)	Incubate at room temperature (18–25°C) on a shaker*
	Aspirate
3 x 2 mL	Wash with 0.9% NaCl solution
1 min	Measure (gamma scintillation counter)

*** Keep shaking conditions constant!**

Optimum:

Amplitude 20 mm = 150 rpm

Amplitude 10 mm = 220 rpm

Amplitude < 8 mm = 300 rpm

According to experience an incubation time of 2 h at room temperature is sufficient to perform the test automatically. This should be checked in individual cases.

Calculation of Results

The standard curve is established manually as follows:

1. Determine the mean CPM for each pair of duplicate tubes.
2. Divide the mean CPM of each standards (B) by the mean CPM of the highest standard (B_{max}) and multiply by 100 in order to obtain the percentage of relative binding ($\%B/B_{\text{max}}$) for each standard.
3. On semi-log paper, plot the relative binding of each standard ($\%B/B_{\text{max}}$) on the Y-axis versus the corresponding concentrations (ng/mL) on the X-axis.
4. Read sample concentrations (ng/mL) directly off the standard curve by their corresponding relative binding ($\%B/B_{\text{max}}$).

Samples with counts above the highest standard concentration must be diluted with the kit diluent and re-assayed. For the final concentration of such samples, the appropriate dilution factor has to be considered. An example of a standard curve is given in the QC-certificate. This curve must not be used for the calculation of unknown samples.

The instrumental calculation of radio-immunological measuring values is based on a spline approximation.

Quality Control

Observe quality control guidelines for medical laboratories.

Validity and precision of the results should be checked with control sera or pool sera prepared by the laboratory.

The control included with the kit is well suited for in-house quality control in the laboratory. This control should be simultaneously tested with each test run and treated like patient samples.

The range of concentrations of each control is reported on the certificate of analysis and indicates the limits established by DiaSorin for control values that can be obtained in reliable assay runs.

Expected Values

The reference range of CEA described in the literature is method dependent and varies between 1.5 – 5 ng/mL (9).

The upper reference range of IRMA – coat® CEA in apparently healthy patients (n=60) was found to be 4.4 ng/mL (mean ± 2 standard deviation). A CEA value of 2.1 ng/mL (mean ± 2 standard deviation) was found for non-smokers (n=30) and 5.8 ng/mL was found for smokers (mean ± 2 standard deviation).

Since CEA values may vary depending on the laboratory method used, each laboratory should establish its own reference range.

Limitations of the Procedure

Patients with malignancies may exhibit CEA values within the reference range. Elevated CEA values may also be observed in patients with benign diseases such as liver cirrhosis, viral hepatitis, pancreatic or gastrointestinal disorders. Smoking and alcohol consumption may also lead to an increase in CEA concentrations (6).

Therefore, CEA serum levels may only be interpreted in context with the clinical picture and other diagnostic procedures.

All tests, in which antigen is incubated together with labelled antibodies and immobilised antibodies in a liquid phase, bear the risk that undiluted samples containing extremely high concentrations of the antigen, will give measuring values below those of the highest standard. In case of the IRMA-coat® CEA, this phenomenon is observed at concentrations exceeding 10 300 ng/mL. If such values are suspected, measurement should be repeated after further dilution (e.g. by factors 10, 100, 1000) of the specimen.

HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may give falsely elevated or decreased values. Although HAMA-neutralising agents are added, extremely high HAMA serum concentrations may occasionally influence results. These samples should not be used for the IRMA-coat® CEA assay.

Analytical Data

Calibration

The test has been calibrated using an internal reference standard.

Measuring range

IRMA-coat® CEA allows to measure concentrations between 0.6 and 100 ng/mL.

High-dose hook

A high-dose hook effect was observed for CEA concentrations > 10,300 ng/mL.

Precision

Intra-assay variation			Inter-assay variation		
Mean value (ng/mL)	CV (%)	n=	Mean value (ng/mL)	CV (%)	n=
4.5	4.3	11	4.4	6.6	10
20.2	2.4	11	20.1	4.0	10

Analytical sensitivity

The detection limit of the IRMA-mat® CEA assay is < 0.6 ng CEA/mL. This limit is defined as a value exceeding the zero standard by three standard deviations; it is the lowest CEA concentration that can be differentiated from zero with statistical significance.

Analytical specificity

No cross-reactivity is observed due to the presence of AFP (10000 ng/mL), Ferritin (8000 ng/mL) and PSA (1000 ng/mL).

Dilution

A patient's serum was diluted with diluent and then measured. The measured values were compared with expected values obtained from linear regression.

Original concentration: 70.5 ng/mL

Dilution	Measured value (ng/mL)	Expected value (ng/mL)	Recovery (%)
1 : 1.25	54.9	56.1	98
1 : 2.5	28.8	28.5	101
1 : 5	16.4	14.7	112
1 : 10	7.8	7.8	100

Recovery

A patient's serum with low CEA content was spiked with different amounts of CEA and then measured.

Original concentration: 2.2 ng/mL

Measured value (ng/mL)	Expected value (ng/mL)	Recovery (%)
47.0	47.8	98
25.5	25.0	102
13.4	13.6	98
8.3	7.9	105

DiaSorin Inc.
1951 Northwestern Avenue, P.O. Box 285,
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710 - Fax: 651-351-5669

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1. Chu DZJ, Erickson CA, Russel, MP, Thompson C, Lang NP, Broadwater RJ, Westbrook KC. Prognostic Significance of Carcinoembryonic Antigen in Colorectal Carcinoma. Arch Surg 1991; **126**: 314-316
2. Díez M, Torres A, Maestro mL, Ortega MD, Gómez A, Pollán M, Lopez JA, Picardo A, Hernando F, Balibrea JL. Prediction of survival and recurrence by serum and cytosolic levels of CEA, CA125 and SCC antigens in resectable non-small-cell lung cancer. Br J Cancer 1996; **73**; 1248-1254
3. Dnistrian AM, Schwartz MK, Greenberg EJ, Smith CA, Schwartz DC. CA 15-3 and carcinoembryonic antigen in the clinical evaluation of breast cancer. Clinical Chimica Acta 1991; **200**: 81-94
4. Gold P, Freedman SO. Demonstration of Tumor-Specific Antigens in Human Colonic Carcinomata by Immunological Tolerance and Absorption Techniques. J Exp Med 1965; **121**: 439-462
5. Lamerz R. CEA in Clinical Practice. J Nucl Med Allied Sci 1990; **34 (Suppl. 3)**: 41-48
6. Sikorska HM, Fuks A, Gold P. Clinical Applications of Carcinoembryogenic Antigen. J Nucl Med Allied Sci 1990; **34 (Suppl. 3)**: 7-29
7. Von Kleist S, Zimmermann W, Thompson J. The Carcinoembryonic Gene Family. J Nucl Med Allied Sci 1990; **34 (Suppl. 3)**: 30-33
8. Ward U, Primrose JN, Finan PJ, Perren TJ, Selby P, Purves DA, Cooper EH. The use of tumour markers CEA, CA-195 and CA-242 in evaluating the response to chemotherapy in patients with advanced colorectal cancer. Br J Cancer 1993; **67**: 1132-1135
9. Thomas L. Labor und Diagnose. 5. Auflage.-Frankfurt/Main : TH-Books-Verl.-Ges., 1998: 984

SYMBOLS USED WITH IVD DEVICES
SIMBOLI USATI CON I DIAGNOSTICI IN VITRO
SYMBOLES UTILISÉS AVEC LES DIAGNOSTICS IN VITRO
MIT IN-VITRO-DIAGNOSTIKA GEBRAUCHTE SYMBOLE
SÍMBOLOS USADOS CON LOS DIAGNÓSTICOS IN VITRO
SÍMBOLOS UTILIZADOS COM OS DIAGNÓSTICOS IN VITRO
SYMBOLER SOM ANVÄNDS MED IN VITRO-DIAGNOSTIK
SYMBOLER SOM BRUGES VED IN VITRO-DIAGNOSTIK
JELMAGYARÁZAT IVD ESZKÖZÖKHÖZ
ΣΥΜΒΟΛΑ ΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΑ ΜΕ ΤΙΣ ΣΥΣΚΕΥΕΣ IVD

CONT.

Kit contents / Contenido del kit / Contenu de la trousse / Inhalt des Kits
 Inhalt des Kits / Contenido del kit / Conteúdo do dispositivo / Satsinnehåll
 Kittets indhold / Készlet tartalma / Περιεχόμενα συσκευασίας.

Ab ¹²⁵I

Tracer: antibody labelled with ¹²⁵I / Tracciante: anticorpi marcati con ¹²⁵I
 Traceur: anticorps marqués à l'¹²⁵I / Tracer: mit ¹²⁵J markierte Antikörper
 Trazador: anticuerpos marcados con ¹²⁵I
 Marcador: anticorpos marcados com ¹²⁵I
 Spårämne: antikroppar märkta med ¹²⁵I
 Tracer: antistoffer mærket med ¹²⁵I
 Tracer: ¹²⁵I-tel jelölt antitest
 Ιχνηθέτης: αντίσωμα σημασμένο με ¹²⁵I.

SORB

Solid phase (Coated tubes / Coated beads).
 Fase solida (Provette sensibilizzate / Sferette sensibilizzate).
 Phase solide (Tubes revêtues / Billes revêtues).
 Feste Phase (Beschichtete Röhrchen / Beschichtete Kugeln).
 Fase sólida (Tubos recubiertos / Bolas recubiertas).
 Fase sólida (Tubos revestidos / Bolas revestidas).
 Fast stadium (Belagda rör / Belagda kulor).
 Fast stadium (Sensibiliserede rør / sensibiliserede kugler).
 Szilárd fázis (bevonatos cső/ bevonatos gyöngy)
 Στερεά φάση (επικαλυμμένοι δοκιμαστικοί σωλήνες / επικαλυμμένα σφαιρίδια).

DIL


Sample diluent / Diluente campioni / Diluant pour échantillons
 Probenverdünnungslösung / Diluyente de muestras / Diluente das amostras
 Provspädning / Fortyndingsmiddel til prøver / Minta hígító / Διαλύτης δειγμάτων.

CAL

Calibrator / Calibratore / Etalon / Kalibrator / Calibrador / Kalibrator
 Kalibrator / Kalibrátor / Μέσο βαθμονόμησης.

CONTROL


Control serum / Siero di controllo / Sérum de contrôle / Kontrollserum
 Suero de control / Soro de controllo / Kontrollserum / Kontrolserum / Kontroll szérum /
 Ορός ελέγχου.

For XX tests / Per XX dosaggi / Pour XX dosages / Für XX Bestimmungen / Para XX ensayos / Para XX testes / För XX dosering
Til XX test / XX teszthez / Περιεχόμενο επαρκές για XX εξετάσεις.

Radioactive / Radioattivo / Radioactif / Radioaktiv / Radioactivo / Radioaktiv
Radioaktiv / Radioaktív / Ραδιενεργός.

Reconstitute with X mL / Ricostituire con X mL / Reconstituer avec X mL
Mit X mL auflösen / Reconstituja con X mL / Reconstitua com X mL
Återställ med X mL / Rekonstruer med X mL / X mL-ben felodandó / Ανασύσταση με X mL.