AB-HTGK-3
(P000172)

DiaSorin

English.....	p. 1
Italiano.....	p. 9
Français.....	p. 18
Deutsch.....	p. 27
Español.....	p. 36
Português.....	p. 45
Czesky.....	p. 54
Ελληνικά.....	p. 62

**THYROGLOBULIN ANTIBODIES
IMMUNORADIOMETRIC ASSAY KIT**

**Procedure for quantitative determination of
autoantibodies to human thyroglobulin (IgG to hTg)
in human serum or plasma samples**

For in vitro use only

1. INTRODUCTION

Autoimmune thyroid diseases encompass a wide spectrum of different clinical disorders varying from hypothyroidism (Hashimoto's disease) or hyperthyroidism (Graves-Basedow disease) to subclinical disease in asymptomatic euthyroid or biochemical hypothyroid subjects. The link between these extremes is the presence of serum autoantibodies directed to thyroid antigens.

Three distinct antigen-antibody systems – commonly found in normal thyroid tissue - have been described in human thyroiditis (autoimmune thyroid disease), involving respectively hTg, a colloid and a microsomal antigen (thyroid peroxidase).

Autoantibodies to hTg are immunoglobulins mainly belonging to the IgG class (IgG to hTg). The assessment of their presence is of great clinical usefulness in the diagnosis of autoimmune thyroid diseases including Graves-Basedow disease, Hashimoto's disease and idiopathic myxoedema.

The level of IgG to hTg is not related to the level of either T₃ or T₄; in fact IgG to hTg may be present in the hypothyroid, euthyroid or hyperthyroid status.

Thyroid autoimmunity is more frequent in women, in whom antibody prevalence increases with age, rising from approximately 10% in the early twenties to 30% in the late sixties. Such autoantibody levels can lead to the development of chronic thyroiditis, resulting in hypothyroidism.

In case of idiopathic myxoedema, significant levels of autoantibodies are observed, indicating the end-stage of autoimmune atrophic chronic thyroiditis. In younger patients, a firm goitre found in combination with high levels of autoantibodies is generally an indication of Hashimoto's disease, characterized by progressive reduction of thyroid function, leading to hypothyroidism.

In the hyperthyroid Graves-Basedow disease, toxic goitre is associated with chronic thyroiditis, as confirmed by high serum levels of autoantibodies.

In all cases, monitoring of autoantibody levels is frequently used to follow up therapy and to identify family members at high risk of developing autoimmune disorders.

2. PRINCIPLE OF THE ASSAY

The assay is an immunoradiometric (IRMA) method. In the first assay stage, IgG to hTg contained in calibrators or samples binds to solid-phase hTg. In the second stage, the ¹²⁵I-tracer (¹²⁵I-labelled protein A) is added, which binds to the solid-phase hTg-IgG complex. After a second incubation, the amount of ¹²⁵I-tracer bound to the solid phase is proportional to the concentration of IgG to hTg present in calibrators or samples to be assayed. At the end of each incubation, the unbound material is removed by aspiration and washing. The method adopted for B/F separation is based on the use of coated tubes, where hTg is coated on the tube walls.

The assay of IgG to hTg employs protein A as tracer, which is a bacterial cell wall protein derived from *Staphylococcus aureus*, capable of binding the Fc fragment of IgG molecules.

3. REAGENTS PROVIDED IN THE KIT

Coated tubes	100
¹²⁵ I-tracer	1 vial
Anti-thyroglobulin calibrators	6 vials
Control serum	1 vial
Sample diluent	2 bottles
Incubation buffer	1 bottle
Wash buffer	2 bottles
Number of tubes	100

STORAGE: Upon receipt, the kit should be stored at 2-8°C. Do not freeze. Once opened, the reagents of this kit are stable until the kit expiry date when properly stored. The kit has been designed to perform 4 assay runs when used throughout the day at room temperature and stored overnight at 2-8°C.

Reagents should not be used past the expiry date. The expiry date of the kit is reported on the external label and corresponds to the expiry date of the tracer. The expiry date of each component is reported on the respective vial label.

When reconstituting the contents of the vials, mix gently to avoid foaming.

Reagents from different batches must not be mixed.

3.1. Coated tubes

The inner surface of each tube is coated with highly purified biotinylated hTg.

Before use, bring the coated tubes to room temperature prior to opening the box, to avoid condensation of humidity.

Securely reseal the box containing unused tubes. Do not mix different batches of coated tubes.

3.2. ¹²⁵I-tracer (red): ready-to-use reagent

The vial contains 52 mL protein A labelled with ¹²⁵I, BSA, phosphate-citrate buffer, preservatives and an inert red dye. Radioactivity is 296 kBq (8 µCi) or less on the calibration date.

3.3. Anti-thyroglobulin calibrators: ready-to-use reagent

Each vial contains 0.5 mL of prediluted human serum containing IgG to hTg, calf serum, BSA, PBS buffer and preservatives. The calibrator concentrations are the following: 0 - 100 - 200 - 500 - 1500 - 6500 IU/mL. The kit calibrators are referenced to MRC 65/93 International Standard. *The kit calibrators demonstrate commutability with patient samples when used with reagents and operating procedure of this in vitro diagnostic test as the manufacturer recommends. Because the calibrators have been prediluted 1:51 and values assigned accordingly, the IU/mL concentrations of the unknown specimens can be determined directly by interpolation from the calibration curve. A multiplying factor is necessary only when the specimen dilution is greater than 1:51.*

3.4. Control serum: lyophilized reagent

The vial contains prediluted human serum, anti-hTg antibodies, calf serum, BSA, PBS buffer and preservatives. The reference range is reported on the vial label.

Reconstitute the vial contents with 1 mL distilled water. The resulting solution is stable at 2-8°C until the kit expiry date. The control serum *should not be further diluted*. Read the value directly from the calibration curve.

3.5. Sample diluent: ready-to-use reagent

Each bottle contains 50 mL PBS buffer, calf serum, BSA and preservatives.

The reagent is common to AB-TPOK-3 and AB-HTGK-3 kits.

3.6. Incubation buffer (blue): ready-to-use reagent

The bottle contains 52 mL TRIS buffer solution, BSA, detergents, preservatives and an inert blue dye.

3.7. Wash buffer: reagent in solution (10x)

Each bottle contains 50 mL 0.5% Triton X-100 and saline solution.

Dilute the contents of each bottle to 500 mL with deionized water. The resulting solution is stable at 2-8°C until the kit expiry date. The reagent is used to rinse coated tubes.

4. EQUIPMENT AND MATERIALS REQUIRED, BUT NOT SUPPLIED

- Distilled or deionized water.
- Glassware.
- Disposable polystyrene tubes.
- Micropipettes with disposable tips (10, 50, 500, 1000 µL) (10, 50 µL: trueness ± 3%, precision 2%; 500, 1000 µL: trueness ± 2%, precision 1%).
- Test tube rack.
- Vortex mixer.
- Horizontal shaker capable of achieving a shaking speed of 300-350 rpm.
- Device for dispensing and aspiration of wash buffer capable of delivering 2-3 mL per wash cycle for two wash cycles.
- Gamma counter suitable for counting ¹²⁵I (counter window setting: 15-80 keV - counter efficiency: 70% - counting time: 1 min). If counter efficiency is below 60%, counting time should be prolonged to 2 min.

5. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulants citrate, EDTA and heparin have been tested and may be used with this assay. Blood should be collected aseptically by venipuncture, allowed to clot, and the serum separated from the clot as soon as possible. Samples having particulate matter, turbidity, lipaemia, or erythrocyte debris may require clarification by filtration or centrifugation before testing. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. If the assay is performed within 24 hours of sample collection, the samples should be kept at 2-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freeze-thaw cycles.

Sample dilution

All samples should be diluted 1:51 with sample diluent (3.5) before assaying. Dispense 10 µL sample and 500 µL sample diluent into polystyrene tubes and mix with a Vortex. If high levels of IgG to hTg are expected, further dilution with the sample diluent supplied in the kit should be performed. *Calibrators and control serum are ready to use and must not be diluted.*

6. ASSAY PROCEDURE

Bring all reagents to room temperature (20-25°C) before assaying. Perform the assay at least in duplicate. Calibrators must be run with each series of patient specimens. Calibrators and samples should be subjected to the same process and incubation time. Perform all assay steps in the order given and without any appreciable delays between the steps. A clean, disposable tip should be used for dispensing each calibrator and sample.

- Dispense reagents *in the bottom of coated tubes*. Operate according to the following scheme:

	tubes	Calibrators 0-5	Diluted samples
reagents			
Calibrators		50 µL	-
Diluted samples		-	50 µL
Incubation buffer		500 µL	500 µL

- **Mix** the contents of tubes and **incubate for 2 hours at room temperature** while continuously shaking (300-350 rpm).

- Carefully **aspirate** the incubation mixture and **wash** twice with 2 mL wash buffer. *Be sure that the aspirator tip touches the bottom of the coated tube so that all the liquid is removed. Failure to remove adhering solution adequately may result in poor reproducibility and spurious results. No trace of dye should still be visible.*
- **Dispense 500 μ L tracer** into all tubes. Prepare two non-coated tubes for total activity computation containing only 500 μ L tracer and set them aside until counting.
- **Incubate for 2 hours at room temperature** while continuously shaking (300-350 rpm).
- Carefully **aspirate** tracer and **wash** twice with 2 mL wash buffer. Operate as above.
- **Measure the radioactivity** of tubes.

7. CALCULATION OF RESULTS

Compute the mean net counts for each group of tubes. Compute the B/T ratio for each calibrator and unknown sample as follows:

$$B/T\% = \frac{\text{calibrator or sample mean counts}}{\text{total activity mean counts}} \times 100$$

Plot in log-log coordinates the mean percent value for each calibrator on the ordinate (y axis) as a function of IgG to hTg concentration expressed as IU/mL on the abscissa (x axis). A calibration curve is thus obtained (Fig. 1). Directly from the calibration curve, read the IgG to hTg concentration of each sample expressed as IU/mL. If the sample was diluted beyond the initial 1:51 factor, the antibody concentration value derived from the diluted specimen must be multiplied by the additional dilution factor. Keep in mind that the calibrators are already diluted 1:51.

Calculation example

The following data must only be considered an example and should not be employed instead of the data obtained by the user.

Description	cpm	B/T x 100
Total activity	78,581	—
Zero calibrator/sample diluent	85	—
100 IU/mL	1,313	1.6
200 IU/mL	2,663	3.3
500 IU/mL	6,445	8.2
1500 IU/mL	19,955	25.3
6500 IU/mL	40,869	52.0
Sample	2,977	3.7

By interpolation from the calibration curve, the sample is found to contain 230 IU/mL IgG to hTg.


Fig. 1

Interpretation of results

The first calibrator (100 IU/mL) represents the system's cut-off.

Samples with IgG to hTg values greater than or equal to the cut-off value should be graded positive. Samples with IgG to hTg values less than the cut-off value should be graded negative.

However, samples whose IgG to hTg values range within $\pm 20\%$ of the cut-off value should be graded disputable.

Each laboratory should establish its own reference ranges.

8. EXPECTED VALUES

During the clinical investigation 211 subjects were taken into consideration. In this population the subjects were graded negative, positive and equivocal with a reference test. 96% of the subjects graded negative with DiaSorin kit showed values of IgG to hTg below 100 IU/mL and 100% of the subjects graded positive showed values of IgG to hTg above 100 IU/mL.

However, the values reported above are merely indicative: each laboratory should establish its own reference ranges.

9. SPECIFIC PERFORMANCE CHARACTERISTICS

9.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive analytes.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by anticoagulants (citrate, EDTA, heparin), slight haemolysis (up to 20 mg/dL haemoglobin), lipaemia (up to 10000 mg/dL triglycerides), bilirubinaemia (up to 50 mg/dL bilirubin).

Cross-reactions. The use of highly purified hTg fixed on the solid phase curbs the possibility of interference from non-specific proteins.

No cross-reactivity is observed due to the presence of TPO autoantibodies up to a concentration of 1000 U/mL.

9.2. Analytical sensitivity

Analytical sensitivity may also be expressed as the limit of detection, which is the minimal amount of specific analyte detectable by the assay. The limit of detection is 6 IU/mL at 95% confidence limit. This was calculated as the apparent concentration of analyte which was distinguishable from the zero calibrator (sample diluent), that is, two standard deviations above zero.

9.3. Precision

Different sample pools, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability).

Repeatability	A	B	C
Number of determinations	10	10	10
Mean (IU/mL)	115	951	2924
Standard deviation	3.3	42.7	105.0
Coefficient of variation (%)	2.9	4.4	3.6

Reproducibility	A	D	E
Number of determinations	10	10	10
Mean (IU/mL)	119	793	2361
Standard deviation	13.9	53.0	177.0
Coefficient of variation (%)	11.6	6.7	7.5

9.4. Trueness

The assay trueness has been checked by the dilution test.

Dilution test. Two sera with high IgG to hTg concentration were tested after serially diluting with sample diluent.

Dilution	Expected concentration, IU/mL	Measured concentration, IU/mL	% Recovery
neat	–	3089.00	–
1:2	1544.50	1640.00	106.2
1:4	772.30	776.70	100.6
1:8	386.10	381.80	98.9
1:16	193.10	204.50	105.9
1:32	96.50	107.60	111.5
neat	–	2457.00	–
1:2	1228.50	1246.50	101.5
1:4	614.25	609.4	99.2
1:8	307.13	321.80	104.8
1:16	153.56	157.50	102.6

9.5. High-dose saturation effect

Whenever samples containing extremely high antibody concentrations are tested in a two-step sandwich method, the saturation effect can mimic concentrations lower than real. However, a well-optimized method excludes grossly underestimated results (as it may occur in one-step sandwich methods).

Concentrations of IgG to hTg up to 100,000 IU/mL were tested in this assay and it was observed that such levels produce analytical signals that are consistently above the top calibrator concentration (saturation curve). For correct quantification, samples containing IgG to hTg levels greater than that of the top calibrator should be diluted with sample diluent and retested. The results must then be multiplied by the additional dilution factor beyond the initial 1:51 to obtain the IgG to hTg levels of the neat specimens.

10. LIMITATIONS OF THE PROCEDURE

The use of radioimmunoassay has strongly reduced false positive results. However, diagnosis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement.

Bacterial contamination or repeated freeze-thaw cycles of the specimens may affect the assay results.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results. In particular, precise pipetting and accurate aspiration and washing are essential.

Non-reproducible results may arise from methodological factors, such as:

- cross-exchange of vial caps
- use of the same tip when withdrawing from different vials or dispensing different samples
- leaving the vials open for long
- exposure of reagents or samples to intense heat or heavy sources of bacterial contamination
- inadequate aspiration of incubation mixture and rinsing of tubes
- contamination of tube rims by tracer or samples
- casual oscillations or inadequate handling of the gamma counter
- use of reagents from different master batches.

11. WARNINGS AND PRECAUTIONS

Test components contain sodium azide as a preservative. Because sodium azide may form explosive lead or copper azide in plumbing, it is recommended that drains be thoroughly flushed with water after disposal of solutions containing sodium azide (Council Directive 99/45/EC).

R 22 – Harmful if swallowed.

R 31 – Contact with acids liberates toxic gas.

S 28 – After contact with skin, wash immediately with plenty of water.

S 45 – In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, and anti-HIV-1/2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

12. SAFETY PRECAUTIONS

- Do not eat, drink, smoke or apply cosmetics in the assay laboratory.
- Do not pipette solutions by mouth.
- Avoid direct contact with all potentially infectious materials by using protective clothing such as lab coats, protective glasses and disposable gloves. Wash hands thoroughly at the end of each assay.
- Avoid splashing or forming an aerosol. Any reagent spills should be washed with a 5% sodium hypochlorite solution and disposed of as though potentially infectious.
- All samples, biological reagents and materials used in the assay must be considered potentially able to transmit -infectious agents. They should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction over the laboratory, and the regulations of each Country. Disposable materials must be incinerated; liquid waste must be decontaminated with sodium hypochlorite at a final concentration of 5% for at least half an hour. Any materials to be reused must be autoclaved using an *overkill* approach (USP 24, 2000, p. 2143). A minimum of one hour at 121°C is usually considered adequate, though the users must check the effectiveness of their decontamination cycle by initially validating it and routinely using biological indicators.

13. BASIC RULES OF RADIATION SAFETY

Reagents Containing Iodine-125

This kit contains radioactive material which does not exceed 8.12 µCi (300 kBq) of iodine-125. Appropriate precautions and good laboratory practices should be used in the storage, handling, and disposal of this material.

For practitioners or institutions receiving radioisotopes under a general license:

This radioactive material may be received, acquired, possessed, and used only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories, or hospitals, and only for in vitro clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and the general license of the U.S. Nuclear Regulatory Commission or of the state with which the Commission has entered into an agreement for the exercise of regulatory authority.

1. Storage of radioactive material should be limited to a specifically designated area.
2. Access to radioactive materials must be limited to authorized personnel only.
3. Do not pipette radioactive material by mouth.
4. Do not eat or drink within designated radioactive work areas.
5. Areas where spills may occur should be wiped up, then washed with an alkali detergent or radiological decontamination solution. Any glassware used must be rinsed completely with water before washing with other laboratory glassware.

For practitioners or institutions receiving radioisotopes under a specific license:

The receipt, use, transfer, and disposal of radioactive materials are subject to the regulations and conditions of your specific license.

ATTENTION: Radioactivity printed in the package insert may be slightly different from the radioactivity printed on the box label and on the tracer vial label. The box label and the tracer vial label indicate the actual amount of radioactivity at the calibration date where the package insert indicates the theoretical radioactivity of the kit.

SCHEME OF THE ASSAY

- 1 - PREPARE REAGENTS. DILUTE SAMPLES 1:51.
- 2 - IDENTIFY COATED TUBES IN DUPLICATE.
- 3 - DISPENSE REAGENTS BY THE FOLLOWING SCHEME AND MIX INCUBATION MIXTURE:

REAGENTS \ TUBES	T	CAL 0-5	DILUTED SAMPLES
CALIBRATORS	-	50 μ L	-
DILUTED SAMPLES	-	-	50 μ L
INCUBATION BUFFER	-	500 μ L	500 μ L

- 4 - INCUBATE FOR 2 HOURS AT ROOM TEMPERATURE WHILE SHAKING.
- 5 - ASPIRATE INCUBATION MIXTURE AND WASH TWICE WITH 2 mL WASH BUFFER.
- 6 - DISPENSE 500 μ L TRACER INTO ALL TUBES.
- 7 - INCUBATE FOR 2 HOURS AT ROOM TEMPERATURE WHILE SHAKING.
- 8 - ASPIRATE TRACER AND WASH TWICE WITH 2 mL WASH BUFFER.
- 9 - MEASURE THE RADIOACTIVITY OF TUBES.

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P.J. BANGA P.S. BARNETT, A.M. McGREGOR

Review – Immunological and molecular characterization of the thyroid peroxidase autoantigen, *Autoimmunity*, **8** : 335-343 (1991).

M.J. BAYER, J.P. KRISS

A solid phase, sandwich-type radioimmunoassay for antithyroglobulin: elimination of false positive results and semiquantitative measurement of antithyroglobulin in the presence of elevated thyroglobulin. *J. Clin. Endocrinol. Metab.*, **49** : 565 (1979)

L.C.P. DE CARVALHO, I.M. ROITT, G. WICK

A three layer immunoradiometric assay for antibodies in different immunoglobulin classes and its application to the detection of chicken thyroglobulin autoantibodies and of antibodies to sheep erythrocytes. *J. Immunol. Methods*, **39** : 15 (1980)

H. ENGLER et al.

Assessment of anti-thyroglobulin and anti-microsomal autoantibodies in patients with autoimmune thyroid diseases: comparison of hemagglutination assay, enzyme-linked immunoassay and radioligand assay. *Chim. Clin Acta*, **179** : 251 (1989).

M. GREINER et al.

A modified ROC analysis for the selection of cut-off values and the definition of intermediate results of serodiagnostic tests. *J. Immunol. Meth.*, **185** : 123 (1995)

S.H. INGBAR, K. A. WOEBER

The thyroid gland.

In: *Textbook of Endocrinology*, R.H. Williams ed., W.B. Saunders, Philadelphia, p. **117** (1981)

R. LINDMARK, C. BIRIELL, J. SJÖQUIST

Quantitation of specific IgG antibodies in rabbits by a solid-phase radioimmunoassay with ¹²⁵I protein A from *Staphylococcus aureus*. *Scand. J. Immunol.*, **14** : 409 (1981).

M. LUDGATE, G. VASSART

Review – The molecular genetics of three thyroid autoantigens: thyroglobulin, thyroid peroxidase and the thyrotropin receptor. *Autoimmunity*, **7** : 201-217 (1990).

T. MORI, J. P. KRISS

Measurement by competitive binding radioassay of serum anti-microsomal and anti-thyroglobulin antibodies in Graves' disease and other thyroid disorders. *J. Clin. Endocrinol. Metab.*, **33** : 688 (1971).

R.L. PEAKE et al.

Radioimmunoassay for antithyroglobulin antibodies.

J. Lab. Clin. Med., **86** : 907 (1974).

A. PINCHERA et al.

Significance of thyroid autoantibodies in autoimmune thyroid diseases.

In: *Autoimmunity and the Thyroid*, P. G. Wolfish, R. J. Wall, R. Volpi eds., Academic Press, Orlando, p. 139 (1985).

S.H. ROMAN et al.

Enzyme-linked immunosorbent microassay and hemagglutination compared for detection of thyroglobulin and thyroid microsomal autoantibodies.

Clin. Chem., **30** (2) : 246 (1984).

G.B. SALABÉ et al.
Radioimmunoassay for human antithyroglobulin antibodies. II: Determination of antigen-binding capacity.
J. Clin. Endocrinol Metab., **39** : 1125 (1974)

S. SHULMAN
Thyroid antigens and autoimmunity.
Adv. Immunol., **14** : 85 (1971)


G. TORRIGIANI, I.M. ROITT, D. DONIACH
Quantitative distribution of human thyroglobulin autoantibodies in different immunoglobulin classes.
Clin. Exp. Immunol., **3** : 621 (1968).

R. VOLPI
Autoimmune thyroid diseases. A perspective.
Mol. Biol. Med., **3** : 25 (1986).

A.P. WEETMAN, A.M. MCGREGOR
Autoimmune thyroid disease: development in our understanding.
Endocrinol. Rev., **5** : 309 (1984).

DiaSorin Inc.
1951 Northwestern Avenue
P.O. Box 285
Stillwater, MN 55082-0285 U.S.A.
Phone: 651-439-9710
FAX: 651-351-5669

DiaSorin S.p.A.,
13040 Saluggia (VC) Italy

APM0129

32790